
BAUSCH + STRÖBEL

AFV 1010 / ~~2010~~

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AFV 1010 / 2010

FILLING AND CLOSING MACHINES FOR
AMPOULES AND VIALS

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These combined filling and closing machines are designed for processing ampoules and vials fully automatically. Available either for single or double position processing these fully automatic machines are designed for an output of 3,000 or 6,000 pcs/h respectively and are therefore suitable both for laboratory use and for production in the medium capacity range. With additional equipment either ampoules or vials can be selected for processing. Brief change-over times and fast, simple size changes enhance the versatility of these machines for processing different containers and products.

These intermittent-motion, rotary machines conform to all the requirements of the pharmaceutical industry.

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The overall design of the machine and the arrangement of the equipment optimize the laminar air flow, enabling trouble-free operation in clean-room conditions. In this way the GMP and FDA guidelines and all worldwide standards applicable to the processing of high-quality pharmaceutical products are met.

Both machines are equally suitable as stand-alone machines or for integration into a production line. Virtually all parts in the main work area are made of non-corrosive materials such as stainless steel and high-grade plastics. The design and surface finish of all the parts allow easy cleaning and care.

All the machine operations are PLC-controlled and adhere to the principle of „quality before quantity“. Production data can be called up on the central operating panel at the press of a button and are displayed in plain text.

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

Major features:

- Versatility (for ampoules and vials)
- Streamlined design for laminar flow
- Non-abrasive container transport cushioned with ball bearings
- Container rotation during filling
- Precise dosing with valveless rotary piston pumps
- Fast and simple size change

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Fully automatic production line for processing ampoules

Washing with AWU 5000
Sterilizing with DHT 2531
Filling and sealing with AFV 2010

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■ Work sequence


- Manual container feeding from trays or automatic feeding by a sterilizing tunnel connected upstream
- Intermittent transfer to the star wheel
- Container registration by photo-sensors for control of the following work steps (no container – no dosing)
- Melting open (if closed ampoules are processed)
- Nitrogen flushing before filling
- Precise dose metering with valveless rotary piston pumps
- Nitrogen flushing during filling
- Nitrogen flushing after filling
- For ampoules:
Nitrogen blanket during sealing
Sealing by melting the stem and pulling off the ampoule tip
- For vials:
Closing by insertion of a rubber stopper
- Removal in trays or magazines, or transfer to a machine downstream

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Dosing

The dosing station is equipped with either one or two valveless rotary piston pumps. These high-quality, precision pumps comprising only two parts ① guarantee utmost dosing accuracy. These pumps are simple to remove, install, clean and sterilize as are all the parts in contact with the product.

Depending on the product to be filled, the pumps can be supplied in stainless steel, glass or ceramic.

The dosing system is also available as a CIP/SIP version (clean in place – sterilize in place). All the parts in contact with the product can then be cleaned and sterilized without having to be dismantled.

Both the central dose adjustment and the fine adjustment for each pump are made mechanically with setting spindles. They are equipped with digital position displays enabling the dosing quantities to be repeated, so that they can be reproduced with the same accuracy at any time.

An optional drip retraction device prevents problematic products from dripping at the end of the filling procedure.

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CIP-SIP

„Clean in place – sterilize in place“ is a time and energy-saving system of cleaning or sterilizing the parts in contact with the product. Without removing them from the machine, these parts are cleaned with suitable agents such as hot water or steam, etc. and rinsed thoroughly.

The pumps are moved to the required positions for the CIP/SIP procedure either manually or program-controlled. The filling needles are inserted manually into the needle receptacle.

When the dosing system has been sterilized, aseptic production can be resumed.

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

■ Opening station

Closed ampoules can be processed with the addition of an opening station. This optional extra opens the domed closed ampoule stems immediately before filling.

Opening is carried out in two steps:

First a hole is melted in the dome to compensate pressure.

In the second stage this hole is widened to a circular filling opening while the ampoule is constantly rotated.

Condensation in the ampoule stem arising from the heat is dried by a burner system.

This two-stage process ensures ideal opening of the ampoules for filling without any detriment to the interior of the ampoule.

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■ Sealing station for ampoules

Ampoules are closed by melting and pulling off the tip of the ampoule stem ①.

The flames of the burners are fed with a mixture of fuel gas and oxygen to ensure that they burn steadily. The flow can be regulated at each individual work station separately ②.

An automatic monitor switches off the fuel gas supply if the flames go out.

The ampoule stems are heated up in stages and then pulled off and collected in a bin. The heat arising during the sealing process is withdrawn from the work area by an extraction unit.

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■ Processing vials

With additional equipment both ampoules and vials can be filled and closed on the same machine.

The closures are fed to the closing station by a sorting bowl and a track vibrator.

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

■ Size change

Changing over the machine for different containers is quick and simple as complete size part sets are exchanged. The machine can then be adjusted to any container height.

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Size part set for processing ampoules

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■ Optional extras and alternatives

- Special layout of the infeed and discharge areas
- Automatic burner ignition
- Station to open closed ampoules
- Nitrogen flushing before, during and after filling and during sealing
- Flame heat extraction
- CIP/SIP version (clean in place – sterilize in place)
- Centering device for the ampoule stem
- Drip retraction
- Mobile version
- Program-controlled dose adjustment
- Various supply vessels
- IPC control outlet
- Equipment for processing vials
- Laminar flow unit
- Clean-room cabin
- Encapsulation of the work area by an isolator
- A range of checking and alarm devices and counters
- Operational data logging system for the documentation of all production data. With a suitable interface, this system can be linked to mainframe computers.

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