



Pre-Install Sign-Off For MasteRad VX Series System


Dear Customer;

The checklist below is required before scheduling the installation of your new digital X-Ray System. Please read and make sure that your X-ray room is ready. Please check off completed items, sign, and fax to 949-679-2882.

Please note that we cannot schedule your installation until we receive your complete and signed form.

Please Check off boxes for completed items

Structural Requirements	
	Minimum Size Requirements- 7.7 'w x 6' x 7.2' h
	Recommended Room Requirements – 9'w x 8'd x 8' h
	Room is complete with dry wall and paint.
	All Door openings into X-Ray Room are 30" wide.
	Acoustic Ceiling tiles installed.
	Floor Covering complete.
	Cabinet, molding and fixtures installed.
	Radiation protection such as lead shielding (if required) installed.



Electrical Requirements

Consult your electrician on these requirements

	Safety disconnect is required and should be placed above the generator.
	Safety disconnect should be 5' from the ground.
	Safety disconnect should have a true earth grounding wire going back to the main panel Additional empty ground lug for 6 AWG wire required
	Appropriate wire size used from main to safety disconnect: ☐ For Runs 0' to 50' wire size 2 AWG recommended ☐ For Runs 51' to 100' wire size 1/0 AWG recommended Run Length are from the building main electrical panel to room mounted safety disconnect. For longer runs please call for recommendations
	All wire sizes stated above are for copper wire
	Single phase requires 100 Amp safety disconnect (208 VAC to 240 VAC)
	Available 110V Power Outlets (2 minimum) Elevating Table Requires 110V / 10A



Pre-Install Sign-Off for MasteRad VX Series Systems

Dear Customer;

The checklist below is required before scheduling the installation of your new digital X-Ray System. Please read and make sure that your X-ray room is ready. Please check off completed items, sign, and fax to 949-679-2882.

Please note that we cannot schedule your installation until we receive your complete and signed form.

Network and Telecommunications Requirements

	Network drop in X-Ray Room
	Internet access in X-Ray Room
	Telephone in X-Ray Room

General Notes

- | | |
|--|--|
| <ul style="list-style-type: none">☐ All systems come with single phase generator unless otherwise ordered.☐ For shielding requirement in your state please visit http://www.crcpd.org/Map/default.aspx☐ For rooms smaller in size than the minimum please call for options.☐ Please inspect all delivered materials for any damage and note it on the delivery receipt. | <ul style="list-style-type: none">☐ If during the delivery, damage occurred to your building, you must document the damage and sign the delivery receipt.☐ To avoid charges, please make sure all requirements are completed before the scheduled day of installation.☐ For question or to contact Medicatech USA please Call: 800-817-5030 |
|--|--|

- **Practice Name:** _____
- **Printed Name:** _____
- **Title:** _____
- **Practice Phone Number:** _____ - _____ - _____
- **Date:** ____/____/____
- **Signature:** _____

After completing this form, Please Sign and Fax to (949) 679-2882