### > Computed Radiography



# CR 75.0™ DIGITIZER

for Computed Radiography

Maximizing productivity for the complete range of clinical applications



- No waiting times for improved patient care
- Input/output buffer for maximized productivity
- For a broad range of applications

- > CR 75.0 IS A MULTI-USER DIGITIZER FEATURING A UNIQUE DROP-AND-GO BUFFER THAT ELIMINATES WAITING TIMES AND MAXIMIZES PRODUCTIVITY
- > CR 75.0 IS A MULTI-APPLICATION DIGITIZER, BENEFITING FROM THREE DIFFERENT IMAGE RESOLUTION MODES

### **Highest productivity**

The cassette buffer eliminates waiting times and allows for a continuous workflow within the department. Zero-button operation with automated cassette handling makes CR 75.0 a highly productive and user-friendly system with a throughput of up to 115 plates an hour, depending on size and application.

### No waiting

The CR 75.0 digitizer requires no manual interaction and all the user has to do is to deposit the cassettes in the input buffer (up to 10 cassettes). The digitizer automatically takes cassettes from the input buffer and reads the demographic data from the memory on the cassette. It then scans the imaging plate, digitizes the



# CR 75.0 DIGITIZER

for Computed Radiography



Integrated CR User Station for time-saving identification and optimized workflow





image and returns the cassette to the output buffer for new exposures.

### Full data

CR 75.0 reads imaging plates at a standard resolution of 6 pixels/mm. 10 pixels/mm high resolution capability is available for all image plate sizes. 20 pixels/mm resolution will be available for dedicated 18 x 24 cm and 24 x 30 cm extremities cassettes and plates.

### Compact footprint & optimal accessibility

CR 75.0 occupies a very small floorspace and at the same time provides unhindered access to several users, both at the input and the output buffer, resulting in a smooth flow of operations. This concept makes CR 75.0 the state-of-the-art solution for centralized CR environments.

### Universal CR User Station

Optionally, a fully integrated CR User Station is available. The CR User Station is suitable for all CR environments:

- Decentralized CR (Surgery, Intensive Care Unit, Emergency Room,...)
- Personal CR
- In-room CR solutions.

Its modular and ergonomic design includes:

- Cassette identification functions
- Space for:
  - Workstation for image handling, processing and dispatching
  - Monitor, network switches and UPS
  - · Cassette storage

### An economical way to go digital

CR is compatible with all existing X-ray systems allowing X-ray departments to go digital without significant additional investments and workflow adaptations.

## > CASSETTE SIZES

ACCEPTED CASSETTE SIZES	SPATIAL RESOLUTION	PIXEL MATRIX SIZE
Standard resolution		
35 x 43 cm (14 x 17 in)	6 pixels / mm	2320 x 2826
5 x 35 cm (14 x 14 in)	6 pixels / mm	2320 x 2320
ligh resolution		
35 x 43 cm (14 x 17 in)	10 pixels / mm (option)	3480 x 4240
35 x 35 cm (14 x 14 in)	10 pixels / mm (option)	3480 x 3480
35 x 43 cm automatic collimation to 21 x 43 cm)	10 pixels / mm	2020 x 4240
24 x 30 cm	10 pixels / mm	2320 x 2920
8 x 24 cm	10 pixels / mm	1720 x 2320
L5 x 30 cm	10 pixels / mm	1420 x 2920
3 x 10 in	10 pixels / mm	1950 x 2460
10 x 12 in	10 pixels / mm	2460 x 2970
Extremities		
24 x 30 cm	20 pixels / mm	4760 x 5840
18 x 24 cm	20 pixels / mm	3560 x 4640



## > SAFETY

REGION	REGULATION	X-RAY	LASER
Europe	EN 60601-1: 1990 + A1:	Regulation: 1987	EN 60825 - 1:2001
	1993 + A2: 1995		
	EN 60601-1-2: 2001		
USA	UL 2601	DHHS/FDA 21 CFR	DHHS/FDA 21
	21CFR part 820: good manu-	part 1002, subchapter B	CFR parts 1040, 10
	facturing practice for medical		and 1040, 11
	devices		
 Canada	CSA22.2 No.601.1 No.601.1.2		

### > TECHNICAL SPECIFICATIONS

### General

LCD display

### Cassette buffer capacity and performance

• 10 cassettes of mixed sizes, both in input and output buffer

• Throughput: up to 115 plates/h (depending on size and application)

Machine status and error conditions

### **Greyscale resolution**

• Data acquisition: 12 bits/pixel · Output to processor: 12 bits/pixel

### **Dimensions and weight**

• W x D x H: 84 x 115 x 142 cm (33 x 45 x 56 in)

• At foot: 84 cm (33 in) • At buffer: 142 cm (56 in)

• Weight: Approx. 320 kg (750.5 lbs)

### Power

50/60 Hz single phase 240V +10%, max. fuse 16A

230V  $\pm 10\%$ , max. fuse 16A

208V  $\pm 10\%$ , max. fuse 15A (e.g. USA)

200V ±10%, max. fuse 15A (e.g. Japan)

### **Environmental conditions**

• Temperature: 20 - 30 °C (68 - 86 °F)

• Humidity: 10 - 80% RH

• Magnetic fields: max. 12.60 μT

• Rate of change of temperature: 0.5 °C/minute

### **Environmental effects**

• Noise level: max. 65 dB (A) · Heat dissipation: standby 350 W, max. 2000 W

### Safety

### **Approvals**

TüV, UL, cUL, CE

### **Transport details**

• Temperature: -25 to +55  $^{\circ}$ C (-4 to 131  $^{\circ}$ F), -25 °C for max. 72 hours, +55 °C for max. 96 hours

• Humidity: 5 - 95% RH



Agfa-Gevaert has been approved by Lloyd's Register Quality Assurance limited to the following Quality Management System Standards: ISO 9001:1994; EN ISO 9001:1994, and ANSI/ASQC Q9001-1994. The Quality Management System is applicable to the development, production and distribution of Agfa Medical Films.

Agfa-Gevaert has been awarded the Approval of Conformity Certificate by Lloyd's Register Quality Assurance. It certifies that the Quality Management System for our X-Ray films conforms to the requirements of Annex V of the EEC Directive 93/42 and Medical Devices Regulation 1994:3017.

Agfa HealthCare has been approved by Lloyd's Register Quality Assurance limited to the following Quality Management System Standards: ISO 9001:1994; EN ISO 9001:1994, and ANSI/ASQC Q9001-1994.

The Quality Management System is applicable to: Selling, Servicing, Distribution, and Design of Marketing of Agfa's Product Assortment of Equipement and Sensitive Materials Used in Medical Diagnositic Imaging, Non-destructive Testing, Microfilm and Motion Picture Applications.



Agfa-Gevaert had been awarded the ISO 9001 certificate by TÜV Zertifizierungsgemeinschaft e.V. This is applicable to Agfa's Quality Management System for design, production and servicing of Agfa Medical Equipment.

Products distributed in North America are manufactured by/for Agfa Corporation, 10 South Academy St., Greenville, South Carolina 29601.

Agfa, the Agfa rhombus, Point of Knowledge and See More. Do More. are trademarks of Agfa-Gevaert N.V., Belgium or its affiliates. All other trademarks are held by their respective owners and are used in an editorial fashion with no intention of infringementt.

The data in this publication are for illustration purposes only and do not necessarily represent standards or specifications which must be met by Agfa. All information contained herein is intended for guidance purposes only, and characteristics of the products described in this publication can be changed at any time without notice. Products may not be available for your local area. Please contact your local sales representative for availability information. Agfa diligently strives to provide as accurate information as possible, but shall not be responsible for any typographical error.

© Copyright 2004 Agfa-Gevaert N.V. All rights reserved.

MI 00761 10/04

### Agfa Corporation

HealthCare Headquarters 10 South Academy Street, Greenville, SC 29601 Tel: 1-864-421-1600 Fax: 864-421-1622

### Canada

Agfa, Inc. 77 Belfield Road, Etobicoke, M9W 1G6 Tel.: 416-241-1110 Fax: 416-240-7359

### Mexico

Agfa de México, S.A. de C.V. Benjamin Franklin #98, Colonia Escandon, C.P. 11800, México, D.F. Tel.: +52 (55) 5276 7611 Fax: +52 (55) 5515 8719

