Supplement

Scio Four modules

en/	Supplement Scio Four modules,	da	Supplement Scio Four-moduler,
enUS	page 2		side 140
de	Ergänzung Scio Four-Module,	no	Tillegg Scio Four-moduler, side 158
	Seite 20	sv	Tillägg Scio Four-moduler, sida 176
fr	Supplément Modules Scio Four,	fi	Lisälehti Scio Four -moduulit,
	page 40		sivu 194
es	Suplemento Módulos Scio Four,	zh	增补 Scio Four 模块 , 第 212 页
	página 60		
it	Supplemento per i moduli Scio Four,		
	pagina 80		
ptBR	Suplemento Módulos Scio Four,		
	página 100		
nl	Aanvulling Scio Four-modules,		
	pagina 120		

WARNING

To properly use this medical device,
read and comply with the instructions
for use of the patient monitor and this
supplement.

Supplement for patient monitors

Information to the supplement	3	Technical data	13
Typographical conventions	3	Operation characteristic values	13
Safety information definitions	3	Ambient conditions	14
User group requirements	3	Device outlets	14
Duties of the operating organization	3	Gas measurement	14
User groups	3	EMC declaration for IEC 60601-1-2, Edition 3.2	15
For your safety and that of your patients	4	EMC declaration for IEC 60601-1-2, Edition 4.1	18
General safety information	4	Relevant standards	18
Product-specific safety information	6	List of accessories	19
Further information	6		
Essential performance	6		
Application	7		
Intended use	7		
Indications	7		
Contraindications	7		
Environment of use	7		
Overview	7		
Symbols	8		
Abbreviations	8		
Functional scope	8		
Functional principles	8		
Assembly and preparation	9		
Mounting description	9		
Connecting the mains power supply	9		
Establishing potential equalization	9		
Connecting the gas analyzer to a patient monitor	9		
Water trap	10		
Connecting the sample line	10		
Establishing gas scavenging	10		
Getting started	10		
Switching on the gas analyzer	10		
Operation	11		
Monitoring	11		
Reprocessing	11		
Cleaning the fan filter mat	11		
Service	11		
Definition of service terminology	11		
Inspection	11		
Maintenance	12		
Calibration	12		
Disposal	12		
Disposal of the medical device	12		
Disposal of accessories	13		

Information to the supplement

This document supplements the instructions for use of the following products:

- Infinity Acute Care System
- Infinity M540 patient monitor
- Infinity Delta Series
- Infinity Vista XL
- Vista family of patient monitoring devices

This supplement provides information about Scio Four modules with part number 6871801, 6871802, 6871803, and 6871804.

Availability of versions of the listed products and software versions compatible with the Scio Four Modules is approval and country specific. Please consult Dräger for compatible product and software versions for a specific country.

Product names in use

In this supplement all Scio Four modules (Scio Four, Scio Four Oxi, Scio Four plus, and Scio Four Oxi plus) are referred to as "gas analyzer" or "medical device".

Typographical conventions

- 1 Consecutive numbers indicate steps of action, with the numbering restarting with "1" for each new sequence of actions.
- Bullet points indicate individual actions or different options for action.
- Dashes indicate the listing of data, options, or objects.
- (A) Letters in parentheses refer to elements in the related illustration.
- A Letters in illustrations denote elements referred to in the text.

Safety information definitions

WARNING

A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or in damage to the medical device or other property.

NOTE

A NOTE provides additional information intended to avoid inconvenience during operation.

User group requirements

The term "user group" describes the personnel responsible who have been assigned by the operating organization to perform a particular task on a product.

Duties of the operating organization

The operating organization must ensure the following:

- Every user group has the required qualifications (e.g., has undergone specialist training or acquired specialist knowledge through experience).
- Every user group has been trained to perform the task.
- Every user group has read and understood the relevant chapters in this document.
- The product is used only by authorized personnel and is protected against unauthorized access.

User groups

Clinical users

This user group operates the product in accordance with the intended use.

Clinical users have medical specialist knowledge in the field of anesthesia and intensive care and knowledge of device monitoring and perioperative care.

Reprocessing personnel

This user group carries out the necessary activities to reprocess the product.

Reprocessing personnel has specialist knowledge in the reprocessing of medical devices.

Service personnel

This user group performs service activities.

Service personnel have specialist knowledge in electrical and mechanical engineering and experience in the servicing of medical devices.

Where product-specific knowledge or tools are required, the service activities must be performed by specialized service personnel. The specialized service personnel have been trained by Dräger for these service activities on this product.

For your safety and that of your patients

General safety information

The following WARNING and CAUTION statements apply to general operation of the medical device.

WARNING and CAUTION statements specific to subsystems or particular features of the medical device appear in the respective sections of these instructions for use or in the instructions for use of another product being used with this medical device.

Strictly follow these instructions for use

WARNING

Risk of incorrect operation and of incorrect use
Any use of the medical device requires full understanding and strict observation of all sections of this supplement and the relevant instructions for use. The medical device must only be used for the purpose specified under "Intended use".

Strictly observe all WARNING and CAUTION statements throughout this supplement and the relevant instructions for use and all statements on medical device labels. Failure to observe these safety information statements constitutes a use of the medical device that is inconsistent with its intended use.

Service

WARNING

Risk of medical device failure and of patient injury
The medical device must be inspected and serviced regularly by service personnel. Repair and complex maintenance carried out on the medical device must be performed by specialized service personnel.
If the above is not complied with, medical device failure and patient injury may occur. Observe chapter "Service".

Monitoring

WARNING

Risk due to gas measurement failure
If the gas measurement fails, the patient can no longer be adequately monitored.
Ensure corresponding substitute monitoring

WARNING

Do not make therapeutic decisions based solely on individually measured values and parameters. Therapeutic decisions must be made solely by the user. For further information see "User groups" on page 3.

Accessories

WARNING

Risk due to incompatible accessories
The use of faulty or incompatible accessories may compromise the functional integrity of the product. Personal injury and property damage may occur as a consequence.

- Use only compatible accessories. The accessories that are compatible with this product are listed in the list of accessories supplied with the product.
- Use only accessories that are in good working order.

WARNING

Risk of operating errors and incorrect use
Strictly observe the instructions for use of all accessories.

Connecting to other devices

WARNING

Risk of electric shock and of device failure
Any connected devices or device combinations not complying with the requirements mentioned in these instructions for use can compromise the correct functioning of the medical device and lead to an electric shock. Before operating the medical device, strictly comply with the instructions for use of all connected devices or device combinations.

WARNING

Risk of electric shock
Connecting devices to the RS232 or MEDIBUS interfaces can lead to an increased leakage current. If the protective ground of one of these devices fails, the leakage current may rise above the permissible values. Only connect with the approval of the respective device manufacturer. Have the leakage current checked by service personnel. If the permissible value is exceeded, disconnect the devices from the RS232 or MEDIBUS interface.

WARNING

Risk of operating errors and incorrect use
Strictly observe assembly instructions and instructions for use for each connected device.

WARNING

Risk of electric shock
Electrical connections to equipment not listed in these instructions for use or these assembly instructions must only be made when approved by each respective manufacturer.
Alternatively the operating organization must ensure that the device combination complies with the applicable editions of the relevant standards for medical devices.

WARNING

Risk of device failure
The use of faulty or incompatible cabling and mounting equipment may compromise the functional integrity of the product. Personal injury and property damage may occur as a consequence.

- Use only compatible cables and mounting equipment. The accessories that are compatible with this product are listed in the list of accessories supplied with the product.
- Use only accessories that are in good working order.

WARNING

Risk of device failure
Usage of this equipment adjacent or stacked with other equipment should be avoided. It could result in improper operation.
If such use is necessary, this device and the other equipment must be observed to verify that they are operating normally.

This device can be operated in combination with other Dräger devices or with devices from other manufacturers. Follow the accompanying documents of the individual devices.

If a device combination is not approved by Dräger, the safety and the functional integrity of the individual devices can be compromised. The operating organization must ensure that the device combination complies with the applicable editions of the relevant standards for medical devices.

Device combinations approved by Dräger meet the requirements of the following standards:

- IEC 60601-1 (Edition 3.2)
Medical electrical equipment
Part 1: General requirements for basic safety and essential performance
- IEC 60601-1-2
Medical electrical equipment
Part 1-2: General requirements for safety and essential performance
Collateral standard: Electromagnetic compatibility – Requirements and tests
- IEC 60601-1 (Edition 2)
Medical electrical equipment
Part 1: General requirements for safety
- IEC 60601-1-1
Medical electrical equipment
Part 1-1: General requirements for safety
Collateral standard: Safety requirements for medical electrical systems
- IEC 60601-1-2
Medical electrical equipment
Part 1-2: General requirements for safety and essential performance
Collateral standard: Electromagnetic compatibility - Requirements and tests
- IEC 60601-1-4
Medical electrical equipment
Part 1-4: General requirements for safety
Collateral standard: Programmable electrical medical systems

If Dräger devices are connected to other Dräger devices or third-party devices and the resulting combination is not approved by Dräger, the correct functioning of the devices may be compromised. The operating organization is responsible for ensuring that the resulting system meets the requirements of the applicable standards.

Electromagnetic compatibility (EMC)

Medical electrical equipment is subject to special precautionary measures concerning electromagnetic compatibility (EMC) and must be installed and put into operation in accordance with the EMC information provided in the instructions for use of the patient monitor.

WARNING

Risk of device failure

Electromagnetic fields e.g., those generated by radio frequency communication equipment such as mobile phones, high-frequency electrical surgery equipment, defibrillators or shortwave therapy devices can disrupt the function of the medical device.

Only operate radio frequency devices at a sufficient safety clearance.

WARNING

Risk due to electromagnetic disturbance

Wireless communication devices (e.g., cellular phones) and medical electrical equipment (e.g., defibrillators, electro-surgical devices) and their peripherals (e.g. antenna cables and external antennas) emit electromagnetic radiation.

When such devices are operated too close to this device or its cables, the functional integrity of this device may be compromised by electromagnetic disturbances. As a result, the patient could be put at risk.

- **Maintain a distance of at least 0.3 m (1.0 ft) between this device and wireless communication devices, to ensure that the essential performance of this device is fulfilled.**
- **Maintain an adequate distance between this device and other medical electrical equipment.**

Connections to IT networks

During operation, this device can exchange information with other devices by means of IT networks. An IT network can be any data interface (e.g., printer interface, USB interface, RS 232 interface) that is described in standards and conventions.

The exchange of data using wired technologies enables the following functions in the network:

- Display of waveforms and parameter data
- Recording, storing, and printing
- Service mode, access to logbooks

Connecting this device to a network that incorporates other devices or making subsequent changes to that network can lead to new risks for patients, users, and third parties. Before the device is connected to the network or the network is changed, these risks must be identified, analyzed, and evaluated, and appropriate measures taken.

Examples of subsequent changes to the network:

- Changing the network configuration
- Removing devices from the network
- Adding new devices to the network
- Performing upgrades or updates on devices that are connected to the network

Information on connecting to the network

Prerequisites

This device must only be connected to the network by service personnel. The IT representative of the hospital must be consulted in advance.

The following documents must be followed:

- Accompanying documents of this device
- Description of the network interface
- Description of the network-based alarm systems

Dräger recommends complying with IEC 80001-1 (risk management for IT networks with medical devices).

Serial interfaces

The following interfaces are supported:

- RS232 interfaces conforming to EIA RS-232 (CCITT V.24/V.28) for the following applications:
 - MEDIBUS
 - Connections to medical devices from other manufacturers
 - Service Interface

Using the MEDIBUS protocol

WARNING

Risk of patient injury

All data transferred through the MEDIBUS interface is for information only and must not be used as the sole basis for diagnostic or therapeutic decisions.

WARNING

Risk of patient injury

Data transmitted over the MEDIBUS interface is not intended for a distributed alarm system (in the sense of remote monitoring) according to the applicable standards ("Relevant standards" on page 18). They are for information only and must not be used as the sole basis for diagnostic or therapeutic decisions.

Carry out further diagnostic measures.

NOTE

The RS-232 port is used to establish a point-to-point connection. Unauthorized access to the device can impair the device function.

- Unauthorized persons must not have access to the data sent via the RS-232 port.
- All device interfaces must be protected against malware and computer viruses.
- If commercially available RS-232-to-Ethernet converters are used, the operating organization is responsible for the use and correctness of the data.

Information on the security level as per IEC 60601-4-5 standard

This device has been tested in accordance with the requirements of the IEC 60601-4-5 standard.

Further information on the security level achieved by this device and on measures that can increase the security level are available on the following web page:

<https://www.draeger.com/productsecurity>

Storing the instructions for use

WARNING

Risk of incorrect use

Instructions for use must be kept accessible to the user.

Product-specific safety information

WARNING

Risk of fire

To avoid a fire hazard, do not introduce drugs or other substances based on flammable solvents, such as alcohol, into the device and do not use explosive anesthetic agents such as ether or cyclopropane or spray them into the device. If highly flammable substances are used for disinfection, adequate ventilation is required.

WARNING

Risk of malfunction

Modifications to the medical device may lead to malfunctions.

WARNING

Risk of insufficient oxygen supply

The sample flow diverted by the gas monitor may reduce the breathing system volume in case of low-flow anesthesia.

Compensate either by increasing the fresh-gas flow of the anesthesia machine accordingly or by returning the sample gas to the breathing system.

In some anesthesia systems the sample flow may influence the measurement of the expiratory minute volume.

WARNING

Risk due to inaccurate gas measurement values

The use of a sample line other than specified or provided by the manufacturer and with an incorrect length and/or diameter may lead to erroneous agent concentration readings and waveforms or watertap/sample line alarms.

Never use standard pressure-sensor tubing or IV lines because these can absorb anesthetic agents, which are released later (degassing) resulting in erroneous agent concentration readings.

CAUTION

Risk of wrong measured oxygen value

Do not use oxygen concentrators in minimal or low flow anesthesia.

CAUTION

Risk of gas measurement failure

To prevent damage to the gas measurement, especially during cleaning and disinfection, always keep the water trap installed.

Further information

Mandatory reporting of incidents

Serious incidents with this product must be reported to Dräger and the responsible authorities.

Training

Training for users is available via the Dräger organization responsible (see www.draeger.com).

Essential performance

The essential performance consists in:

- Accuracy of gas measurement

If the accuracy of gas measurement is not fulfilled, the device transmits status information with the purpose of generating a technical alarm on the patient monitor.

Application

Intended use

The Scio Four module samples gas from the breathing gas of pediatric patients and adults. It continuously measures the concentration of CO₂, N₂O, and anesthetic agents (Halothane, Isoflurane, Enflurane, Sevoflurane, and Desflurane) in the breathing gas as well as the O₂ concentration (optional). All measured values as well as derived values are communicated to a patient monitor.

Indications

The intended medical indication for the Scio Four module is the need to continuously monitor respiratory gas concentrations according to the intended use.

Contraindications

No known contraindications.

Environment of use

The gas analyzer is designed for use inside professional health-care facilities in rooms in which therapeutic or diagnostic interventions are carried out.

WARNING

Risk of explosion or fire

This medical device must not be used in areas where oxygen concentrations exceed 25 %, or where combustible or explosive gas mixtures are likely to occur.

WARNING

Risk of device failure and/or patient injury and user injury
Magnetic fields can negatively influence the correct functioning of the medical device and therefore endanger the patient or user.
Do not use the medical device near magnetic resonance imagers (MRI, NMR).

WARNING

Risk of device failure

Do not obstruct or close off the fan on the medical device. Air must be able to enter freely. Otherwise the medical device may become too hot.

WARNING

Risk of fire

Place the medical device at least 25 cm (10 in) from any ignition sources such as electrosurgery and laser surgery.

WARNING

Risk of infection

Clean and disinfect the device before usage (including first usage), maintenance, and disposal.

WARNING

Risk due to falling device

To avoid injury, place the medical device on a secure and level platform or use optional mounting accessories available from Dräger.

WARNING

Risk of strangulation

Negligent placement of sample line, cables, and similar device components can endanger the patient.

Use particular diligence when establishing connections to the patient.

CAUTION

Risk of tripping


Do not place the power supply unit, cables, and sample line on the ground. Place them carefully and securely.

Furthermore, do not use the gas analyzer outside of solid buildings, during patient transport or in ground vehicles or aircraft.

Observe the environmental specifications like temperature, humidity, and atmospheric pressure.

Overview

Front panel


A Water trap

B Model name (example only)

C Power indicator

Back panel

A Cooling fan

B On/Off switch

C Gas exhaust

D Potential equalization connector

E DC power inlet

F MEDIBUS connector

G RS232-Connector to patient monitor - X2

Symbols

Symbol	Explanation
	Follow instructions for use
	Power
	Exhaust for sample gas for recirculation or scavenging
	On/Off
	Manufacturer
	Date of manufacture
	Part number
	Serial number
	Storage temperature
	Relative humidity
	Atmospheric pressure
	WEEE label, Directive 2002/96/EC
	Keep away from rain
	Nominal weight
	Applied part, protection class BF (Body Floating)
	Potential equalization
	Do not use this device in the vicinity of MRI scanners.
	The product is a medical device (CE conformity assessment procedure)
Rx only.	CAUTION: Federal law restricts this device to sale by or on the order of a physician.

Abbreviations

Abbreviation	Explanation
ATPS	Ambient temperature and pressure, saturated
BTPS	Body temperature and pressure, saturated
COM	Serial interface
EMC	Electromagnetic compatibility

Abbreviation	Explanation
HME	Heat and moisture exchanger
HMEF	HME filter
I:E	Ratio of inspiratory time to expiratory time
ID	Identification
MEDIBUS	Dräger communication protocol for medical devices with uniform data definition for all devices
MRI	Magnetic resonance imaging
NMI	Nuclear magnetic imaging
NMR	Magnetic resonance
NTPD	Normal Temperature Pressure Dry, 20 °C (68 °F), 1013 hPa, dry
SW	Software

Functional scope

The gas analyzer samples inhaled and exhaled gases from adult and pediatric patients in non-, partial-, and total rebreathing systems. It measures inspiratory and expiratory gases and communicates in real time to the patient monitor.

The gas analyzer is available in 4 variants with different functions as listed below:

	O ₂	CO ₂ , N ₂ O	Agent	Agent ID	Mixtures
Scio Four	No	Yes	1 out of 5	No	No
Scio Four Oxi	Yes	Yes	1 out of 5	No	No
Scio Four plus	No	Yes	2 out of 5	Yes	Yes
Scio Four Oxi plus	Yes	Yes	2 out of 5	Yes	Yes

Functional principles

The gas analyzer utilizes IR spectroscopy for carbon dioxide, nitrous oxide, and volatile anesthetics measurement. Each gas, with its own specific absorption characteristic, is transported through a sample cell. Multiple optical IR filters select specific bands of IR light, which then are analyzed with a detector. The amount of IR light received by the detector is a measure of gas concentration. This means the higher the gas concentration, the less IR light will be received by the detector.

For oxygen measurement, the paramagnetic characteristic of oxygen molecules is used. The gas is transported through a sample cell, where the paramagnetic effect of oxygen is measured. The higher the oxygen concentration, the greater is the measured effect.

WARNING

Risk due to inaccurate gas measurement values.
For the variants Scio Four and Scio Four Oxi only.
Use care when setting agent ID manually on the patient monitor. Measurements will be inaccurate if the wrong agent ID is selected.

WARNING

Risk of incorrect gas measurement values
If Scio Four or Scio Four Oxi are used with anesthetic gas mixtures, inaccurate gas measurements will result.
Do not use Scio Four or Scio Four Oxi with anesthetic gas mixtures.

CAUTION

Risk due to inaccurate gas measurement values
 Due to the response time of the sensors and the gas sample flow rate, the stated accuracy of O₂, CO₂, N₂O and anesthetic agents is limited by the respiratory rate and by the inspiratory to expiratory ratio (I:E). Observe the technical data of the medical device.

Assembly and preparation

Mounting description

The power supply unit shall be mounted in one of the following mounting options:

- Inlet for mains power downwards or
- Inlet for mains power horizontally with the power supply label downwards.

The power supply unit must not be mounted with the inlet for mains power pointing upwards.

Connecting the mains power supply

The mains voltage must correspond to the voltage range indicated by the rating plate on the rear of the device.

The power supply supplies the Scio Four module with 12 VDC. The metal collar of the 12 VDC supply connector on the back of the Scio four establishes a functional ground connection of the Scio housing to earth.

WARNING


Risk of electric shock
If the power supply of Scio Four modules is connected to a power socket with incorrect mains voltage or a power socket without a protective ground, an electric shock may occur. Connecting the device to power socket strips can lead to an increased leakage current. This may result in an electric shock and the user and patient may be put at risk. Connect the device only to power sockets with correct mains voltage and a protective ground.

WARNING

Risk due to penetrating liquid
Penetrating liquid may cause the following:

- Damage to the device
- Electric shock
- Device malfunctions

Mount the power supply unit only in the orientation described in the chapter "Mounting description".
Ensure that no liquid e.g. disinfectant penetrates the device.
Do not place any containers containing liquids above or on the device.

- 1 Connect the power supply with the "X16" connector to the rear of the gas analyzer.
- 2 Connect the power cable to the rear of the power supply.
- 3 Plug the power cable into the mains power socket on the wall.

NOTE

Keep the power cable clear, so that it can be disconnected at any time.

NOTE

Unplug the power cable from the mains power socket to disconnect the gas analyzer and its power supply from mains (e. g. for maintenance or emergency reasons).

Establishing potential equalization

Differences in electrical potential between devices can be reduced by potential equalization.

Potential equalization does not replace the protective ground connection. During operation, the potential equalization connectors must be readily accessible and the connection must be able to be disconnected without the use of tools.


Connecting the potential equalization cable

- 1 Connect the potential equalization cable to the potential equalization pin on the device (see page 7).
- 2 Connect the potential equalization cable to a potential equalization connector of the hospital (e.g., wall, ceiling supply unit).

Connecting the gas analyzer to a patient monitor

For information on patient monitors compatible to the Scio Four modules, see "Information to the supplement" on page 3.

- 1 Connect and secure data interface cable with the "Pat. Mon. - X2" connector (A) on the back panel of the gas analyzer.

- 2 Connect the data interface cable with the connector at the Infinity Docking Station or directly with the patient monitor.

Water trap

Place the water trap on the water trap holder.

WARNING

Risk due to inaccurate gas measurement values and gas measurement failure
Silicon can get into the measuring cuvette and influence the gas measurement permanently.
Do not spray the O-rings of the water trap holder with silicon spray.

WARNING

Risk due to gas measurement failure
If the gas measurement fails, the patient can no longer be adequately monitored.
Check the water trap for damage or blockages.

CAUTION

Risk of device failure
Do not operate the medical device without water trap.

Connecting the sample line

WARNING

Risk due to leakage from improperly connected or damaged sample line
Check the sample line for damage and connect it properly, otherwise faulty gas measurements may result.

WARNING

Risk due to gas measurement failure
If gas measurement fails, the patient can no longer be adequately monitored.
Check sample line for blockage and resolve these as needed.

WARNING

Risk of patient cross infection
To avoid patient cross infection, if the sampled gas is returned to the breathing system, use the sample gas return kit including its bacterial filter. Replace the bacterial filter according to the schedule in the chapter "Preventive maintenance".

WARNING

Risk of medical device failure and of patient injury
Solely connect the sample line to a sample point of the breathing circuit of the patient. Do not connect the sample line to the patient, infusion equipment (e.g., infusion pumps etc.) or any other equipment near the patient.

- 1 Connect the sample line to the Y-piece or HME filter, or hose adapter with Luer Lock.
- 2 Connect the sample line to the water trap.

Establishing gas scavenging

NOTE

As a result of automated zeroing room air is added to the breathing system when using the sample gas return kit.

WARNING

Risk of inaccurate gas measurement values and adverse effects on bystanders due to contamination of the ambient air
Always connect the gas exhaust of the medical device and anesthesia machine to the scavenging system or return the medical device's sample gas to the breathing system of the anesthesia machine.
Ensure proper ventilation of the place where the medical device is located.

- 1 Use the sample gas return kit to connect the exhaust port on the gas analyzer with the sample gas return port on the anesthesia workstation.
- 2 If no sample gas recirculation is possible, connect the gas analyzer exhaust port to the anesthetic gas scavenging system using the sample gas scavenging hose.

Getting started

Switching on the gas analyzer

WARNING

Risk of explosion and fire
Do not place the medical device into operation if oxygen leakage is suspected in the medical device or its vicinity.
Stop all oxygen supplies and contact service personnel.

WARNING

Risk of inaccurate gas measurement values
During warm-up, reported CO₂ measurement values may not be accurate. Wait until the medical device has completed warm-up.

CAUTION

Risk of gas measurement failure
To prevent condensation and resulting failure of electrical components, do not switch on the medical device after significant temperature changes for 1 to 2 hours (e.g., after storage in unheated rooms).

- 1 Connect the mains plug if necessary and set the On/Off switch on the rear of the device to "ON".
- 2 Check the Power indicator on the front of the gas analyzer. The Scio Four modules are ready to use, when the indicator light on the front of the gas analyzer is constantly on.

When starting, the gas analyzer passes through an initialization and warm-up period. During this time, concentrations for certain gases may not be available and the anesthetic agent may not be identified.

NOTE

In case of a power supply failure, the gas analyzer shuts down and when power is restored, the gas analyzer powers back up.

Operation

WARNING

Risk of inaccurate gas measurement values and gas measurement failure
To avoid temporary influence on the gas measurement and prevent damage to the water trap and measuring system do not use nebulizers/aerosols in the breathing system, when the medical device is connected.
Disconnect the sample line temporarily when using nebulizers in the breathing system.

CAUTION

Risk of airway pressure drop
The medical device will continue to draw air from the patient's ventilation system when the mechanical ventilation is stopped during cardiopulmonary bypass surgery.
Disconnect the sample line during cardiopulmonary bypass surgery. In case ventilation is to be ceased, switch the medical device to standby mode or disconnect the sample line and seal the sampling site to avoid leakage.

Monitoring

Set up and use the patient monitor in accordance with its instructions for use.

Reprocessing

Clean and disinfect the gas analyzer before patient change according to the document "Scio Four modules Reprocessing Supplement".

CAUTION

Risk of damaging the pneumatic system
Do not apply excessive pressure (e.g., syringe, or compressed air) to the inlet, exhaust port, or the water trap of the medical device, e.g., during cleaning and disinfection.

Cleaning the fan filter mat

The fan filter mat should be cleaned once every month.

- 1 Locate the fan at the rear of the gas analyzer.
- 2 Remove the fan filter mat from its holder.
- 3 Vacuum up any accumulation of dust on the filter mat as well as on and inside the fan port.
- 4 Insert the cleaned filter mat or a new filter mat.

Service

This chapter describes the maintenance measures required to maintain the proper functioning of the medical device. Maintenance measures must be performed by the personnel responsible.

WARNING

Risk of infection

Users, reprocessing personnel, service personnel, and specialized service personnel may become infected by disease-causing germs.

Clean and disinfect device or device components before each maintenance step and before returning from repair.

WARNING

Risk of electric shock

There are conducting components under the housing cover.

- Do not remove the housing cover.
- Maintenance measures must be performed by service personnel.

WARNING

Only use replacement parts that are qualified to Dräger standards.

Definition of service terminology

Concept	Definition
Service	All measures (inspection, maintenance, repair) intended to maintain and restore the functional condition of a medical device
Inspection	Measures intended to determine and assess the current state of a product
Maintenance	Regular specified measures intended to maintain the functional integrity of a product
Repair	Measures intended to restore the functional integrity of a product after a failure

Inspection

Inspections must be carried out regularly according to the following guidelines and within the specified intervals. Technical documentation is available on request.

Measure	Interval	User group
Inspection and safety checks	Every 12 months	Service personnel

Safety checks

Safety checks are not a substitute for maintenance, which includes the preventive replacement of wearing parts as specified by the manufacturer.

If safety checks are not performed regularly, the proper operation of the medical device can be compromised. Perform safety checks at the indicated intervals.

NOTE

The gas analyzer needs to be in operation mode for more than 10 min before checking the accuracy of the gas measurement.

- 1 Check accompanying documents:
 - Instructions for use are available in the correct version
- 2 Verify that the gas analyzer is in good condition:
 - Labels are complete and legible
 - No visible damage on:
 - Housing
 - External power supply and power cable
 - Water trap holder, especially o-rings
 - Gas exhaust
 - Power indicator: Green steady light
- 3 Check the electrical safety according to IEC 62353.

CAUTION

Risk of device damage
If a test current of 10 A is used in the protective conductor resistant test, the medical device may be damaged.
Configure the test equipment prior testing that the test will be performed only with 200 mA (DC).

NOTE

Only the appliance inlet of the power supply unit is connected to the protective ground.

NOTE

To dispose of the anesthetic gas, connect the device to the scavenging system.

- 4 Check the accuracy of the gas measurement based on a certified test gas concentration:
 - Anesthetic gas measurement
 - For Scio Four and Scio Four Oxi only:
2 Vol% Desflurane
Accuracy ± 0.50 Vol%
 - For Scio Four plus and Scio Four Oxi plus only:
1 Vol% Isoflurane
1 Vol% Sevoflurane
Accuracy ± 0.35 Vol%
 - N₂O measurement
 - For Scio Four and Scio Four Oxi only:
60 % N₂O
Accuracy ± 6.8 Vol%
 - For Scio Four plus and Scio Four Oxi plus only:
70 % N₂O
Accuracy ± 7.6 Vol%
 - CO₂ measurement, 5 Vol%
Accuracy ± 0.83 Vol%
- 5 Check the accuracy of the O₂ measurement:
 - Ambient air 21 Vol%
Accuracy ± 3 Vol%
- 6 Check the sampling rate of the gas analyzer:
 - Accuracy 200 ± 20 mL/min
- 7 Check the gas analyzer for leakage:

- Leakage at -200 hPa (cmH₂O)
- <20 hPa/min (cmH₂O/min)

Maintenance

WARNING

Risk of faulty components
Device failure is possible due to wear or material fatigue of the components.
To maintain the proper operation of all components, this device must undergo inspection and preventive maintenance at specified intervals.

WARNING

Risk of electric shock
Disconnect all electrical connectors from power supply, before opening the enclosure.

The following table shows the preventive maintenance intervals:

Component	Interval	Task	Responsible person
Fan filter mat	Every year	Exchange	Service personnel
Bacterial filter of sample gas return kit	Every year	Exchange	Service personnel
O-rings on the water trap holder	Every 2 years	Exchange	Service personnel

Calibration

The gas analyzer is self-zeroing and does not need calibration. A yearly check of the calibration components shall be performed by service personnel.

For further information, see chapter "Inspection".

Disposal

Disposal of the medical device

WARNING

Risk due to inappropriately reprocessed products
The product may be contaminated with infectious agents.
Before disposal, reprocess the product in accordance with the chapter "Reprocessing".

When disposing of the medical device:

- Consult the relevant waste disposal company for appropriate disposal.
- Observe the applicable laws and regulations.

For countries subject to the EU Directive 2002/96/EC

This device is subject to EU Directive 2002/96/EC (WEEE). In order to comply with its registration according to this directive, this device may not be disposed of at municipal collection points for waste electrical and electronic equipment. Dräger has authorized a company to collect and dispose of this device. To initiate collection or for further information, visit Dräger on the internet at

www.draeger.com. Use the search function with the keyword "WEEE" to find the relevant information. If access to Dräger's web-site is not possible, contact the local Dräger Organization.

- Filter, HME, HMEF
- Water trap
- Sample line
- Sample gas return kit
- Scavenging hose

Disposal of accessories

When disposing of the following accessories, observe the hospital hygiene regulations. Observe the instructions for use of the respective accessories (if provided):

Technical data

Operation characteristic values

Dimensions with water trap

Height	115 mm (4.8 in)
Width	190 mm (7.5 in)
Depth	270 mm (10.5 in)

Weight

Scio Four	2.9 kg + 0.1 kg (6.4 lbs + 0.22 lbs)
Scio Four Oxi	3.2 kg + 0.1 kg (7.0 lbs + 0.22 lbs)
Scio Four plus	2.9 kg + 0.1 kg (6.4 lbs + 0.22 lbs)
Scio Four Oxi plus	3.2 kg + 0.1 kg (7.0 lbs + 0.22 lbs)

Mode of operation Continuous

Device noise emissions ≤45 db(A)

Cooling Fan, temperature-regulated

Dimensions of the external power supply

Height	33 mm (1.3 in)
Width	68 mm (2.7 in)
Depth	109 mm (4.3 in)

Mains power supply 100 to 240 V~ 50/60 Hz

Rated current 1.5 A (1.5 A - 0.7 A)

Inrush current 40 A, cold start at 264 V~ input

Power consumption ≤60 W

Housing leakage current ≤300 µA (per UL 60601-1:2003)
≤500 µA (per UL 60601-1:1988 incl. corrections)

The power supply supplies the Scio Four module with 12 VDC.

The metal collar of the 12 VDC supply connector on the back of the Scio four establishes a functional ground connection of the Scio housing to earth.

Protection classes

Device Class I

Electric shock protection Type BF

Scio Four modules are intended to connect to a patient by means of a WaterLock2 water trap, a sample gas line and patient connectors.

These parts and accessories are subjected to the requirements for a 'type BF applied part' (IEC 60601-1 & ISO 80601-2-55).

Penetration of liquids IP21

Penetration of liquids (power supply) IP22

Classification according to European Medical Device Regulation IIb

UMDNS Code Universal Medical Device 17-445

Nomenclature System - nomenclature system for medical devices

Use of latex No latex (natural rubber) is used in parts with indirect or direct contact to patients or operators.

Ambient conditions

During operation	
Temperature	10 to 40 °C (50 to 104 °F)
Atmospheric pressure	620 to 1100 hPa (9.0 to 15.9 psi), with automatic barometric pressure compensation
Relative humidity	20 to 95 %, without condensation
CO2 concentration	300 to 1000 ppm
O2 concentration	20.95 Vol% ±0.2 Vol% abs. (dry air)
During storage and transport	
Temperature	−20 to +70 °C (−4 to 158 °F)
Atmospheric pressure	500 to 1060 hPa (7.3 to 15.4 psi)
Relative humidity	5 to 95 %, without condensation

Device outlets

Serial interfaces	RS232 and MEDIBUS Only connect devices that meet the requirements of IEC 60950-1 for ungrounded SELV circuits or the requirements of IEC 60601-1 (as of the 2nd edition) for touchable secondary circuits with maximum 24 VDC nominal voltage.
Protocol	STORM / MEDIBUS
Connector	9-pole Sub-D, galvanically isolated 1.5 kV against internal electronics, 0.5 kV against housing.
Baud rate	1200, 2400, 4800, 9600, 19200, 38400 baud
Data bits	8
Parity	Even
Stop bits	1
Pin assignment	
Pin 1, 4, 6, 7, 8, 9	Open
Pin 2	TXD
Pin 3	RXD
Pin 5	GND
Housing	SHLD

Gas measurement

Method	Sidestream gas measurement All values are measured under ATPS conditions, sampling rate calibrated under NTPD conditions. The endtidal gas concentration is the value measured at the time of maximum CO2 concentration in the expiratory phase.
Sampling rate	200 mL/min ±20 mL/min
Maximum time to empty water trap	41 hrs (sample gas under BTPS conditions, ambient air 23 °C)
Gas leakage	<10 mL/min
Data sample rate	<50 ms
Delay time	<2.4 s
Total system response time	<3 s Sum of delay time and gas-specific rise time (T10...90)
Cross sensitivity	<0.2 Vol% for volatile agents, CO2, and N2O <1 Vol% for O2 against ethanol (<3000 ppm), isopropanol (<10000 ppm), acetone (<1000 ppm), methane (<1000 ppm), saturated humidity at 37 °C
Drift	Compensated by automated cyclical zeroing
O2 ^{1) 2)}	
Range	0 to 100 Vol%
Accuracy	±(2.5 Vol% + 2.5 % rel.)

Rise time (T10...90)	<500 ms
CO₂ ^{1) 2)}	
Range	0 to 10 Vol%
Accuracy	±(0.43 Vol% + 8 % rel.)
Rise time (T10...90)	<300 ms
N₂O ^{1) 2)}	
Range	0 to 100 Vol%
Accuracy	±(2 Vol% + 8 % rel.)
Rise time (T10...90)	<300 ms
Anesthetic agents ^{1) 2)}	
Range	
Halothane	0 to 8.5 Vol%
Isoflurane	0 to 8.5 Vol%
Enflurane	0 to 10 Vol%
Sevoflurane	0 to 10 Vol%
Desflurane	0 to 20 Vol%
Accuracy	±(0.20 Vol% + 15 % rel.)
Rise time (T10...90)	<450 ms
Automatic detection	
Primary gas	At the latest at 0.3 Vol%
Secondary gas	At the latest at 0.4 Vol%
	With a Desflurane concentration greater than 4 Vol%, mixture detection occurs at the latest when the concentration of the second anesthetic agent rises above 10 % of the Desflurane concentration.

The mentioned identification thresholds (primary and secondary) apply to rising agent concentrations (e.g., beginning of an anesthetic case).

When anesthetic agent concentration decreases (e.g., end of an anesthetic case), the gas analyzer measures agent concentrations down to 0.05 Vol% based on the last identified agent. Below that value the gas analyzer indicates loss of agent identification and a value of 0 Vol%.

Respiratory rate (RR)

Range	0 to 100 /min (Respiratory rate is derived from CO ₂ real-time data)
Accuracy	0 to 60 /min: ±1 /min with an I:E ratio of 1:1 >60 /min: not specified 0 to 80 /min: ±1 /min with an I:E ratio of 1:2 >80 /min: not specified
Resolution	1 /min

- 1) The accuracy of the measurement of O₂, CO₂, N₂O and anesthetic agents depends on the respiratory rate and on the I:E ratio. The accuracies of the expiratory values specified in the technical data apply at a maximum respiratory rate of 60 /min with an I:E ratio of 1:1. The influence of respiratory rate and I:E ratio on the accuracy has been determined for the measured gas concentrations in a simulated breathing system with step changes in concentration.
- 2) The specified accuracy is obtained at the latest 450 seconds after start-up.

EMC declaration for IEC 60601-1-2, Edition 3.2

General information

The device was tested for electromagnetic compatibility using accessories from the list of accessories. Other accessories may only be used if they do not compromise the electromagnetic compatibility. The use of non-compliant accessories may result in increased electromagnetic emissions or decreased electromagnetic immunity of the device.

Electromagnetic environment

The medical device is intended for use in an electromagnetic environment as specified below. The user must ensure its use in such an environment.

The device may be used in the direct vicinity of other devices only if Dräger has authorized this device arrangement. If no authorization has been given by Dräger, it must be ensured that the device functions correctly in the desired arrangement before using it. The instructions for use for the other devices must be followed.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The medical device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
	Class A	The medical device is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the medical device or shielding the location.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Electromagnetic immunity

The medical device is intended for use in an electromagnetic environment as specified below. The user must ensure that the medical device is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact	±6 kV	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
	±8 kV air	±8 kV	
Electrical fast transients/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV	Mains power quality should be that of a typical commercial or hospital environment.
	±1 kV for input/output lines	±1 kV	
Surge IEC 61000-4-5	±1 kV line(s) to line(s)	±1 kV	Mains power quality should be that of a typical commercial or hospital environment.
	±2 kV line(s) to earth	±2 kV	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycles	>95%, 0.5 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the medical device requires continued operation during power mains interruptions, it is recommended that the medical device is powered from an uninterruptible power supply or a battery.
	40% UT (60% dip in UT) for 0.5 cycles	60%, 5 cycles	
	70% UT (30% dip in UT) for 25 cycles	30%, 25 cycles	
	<5% UT (>95% dip in UT) for 5 s	>95%, 5 s	

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the medical device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3.5}{3} \right] \sqrt{P}$ $d = \left[\frac{3.5}{3} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{3} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey¹⁾, should be less than the compliance level in each frequency range.²⁾ Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	3 V/m	

- 1) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the medical device is used exceeds the applicable RF compliance level above, the medical device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the medical device.
- 2) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances from mobile high-frequency communication equipment

The medical device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the medical device can help prevent elec-

tromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the medical device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.8	3.8	7.4
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

EMC declaration for IEC 60601-1-2, Edition 4.1

General information

This device was tested for electromagnetic compatibility using accessories from the list of accessories. Other accessories may only be used if they do not compromise the electromagnetic compatibility. The use of non-compliant accessories may result in increased electromagnetic emissions or decreased electromagnetic immunity of the device.

This device may be used in the direct vicinity of other devices only if Dräger has approved this device arrangement. If no approval has been given by Dräger, it must be ensured that this device functions correctly in the desired arrangement before use. The instructions for use for the other devices must be followed.

Electromagnetic environment

This device may only be used in environments specified in section "Environment of use" on page 7.

Emissions	Compliance
Radiated emissions	Class A, group 1 (30 MHz to 1 GHz)
Conducted emissions	Class A, group 1 (150 kHz to 30 MHz)

NOTE

The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required), this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment. The radiated emissions testing was conducted according to CISPR 16-2-3 as an alternate test method.

Immunity against	Test level and required electromagnetic environment
Electrostatic discharge (ESD) (IEC 61000-4-2)	Contact discharge: ± 8 kV Air discharge: ± 15 kV
Fast transient electrical disturbances (bursts) (IEC 61000-4-4)	Power cable: ± 2 kV Longer signal input lines/output lines longer than 3 m (10 ft): ± 1 kV
Impulse voltages (surges) (IEC 61000-4-5)	Voltage, line – line: ± 1 kV Voltage, line – ground: ± 2 kV
Magnetic fields at mains frequency (IEC 61000-4-8)	50 Hz: 30 A/m
Voltage dips and short interruptions in the supply voltage (IEC 61000-4-11)	Voltage dips at various phase angles according to the IEC 60601-1-2 standard
Radiated high-frequency disturbances (IEC 61000-4-3)	80 MHz to 2.7 GHz: 3 V/m
Conducted high-frequency disturbances (IEC 61000-4-6)	150 kHz to 80 MHz: 3 V, ISM bands: 6 V

Immunity against	Test level and required electromagnetic environment
Electromagnetic fields in the vicinity of wireless communication devices	Various frequencies were tested according to table 9 of IEC 60601-1-2
Proximity magnetic fields (IEC 61000-4-39)	134.2 kHz: 65 A/m 13.56 MHz: 7.5 A/m

Recommended separation distances from wireless communication devices

To ensure that the functional integrity of this device is maintained, there must be a separation distance of at least 1.0 m (3.3 ft) between this device and wireless communication devices.

The gas analyzer SCIO Four modules complies with the requirements of IEC TR 60601-4-2.

Relevant standards

This medical device meets other standards in addition to those listed here. Additional information is available on request.

- IEC 60601-1 (Edition 3.2)
Medical electrical equipment
Part 1: General requirements for safety
- IEC 60601-1-2 (3rd edition)
Medical electrical equipment
Part 1-2: General requirements for safety, collateral standard: Electromagnetic compatibility; Requirements and tests
- IEC 60601-1-2 (Edition 4.1)
Medical electrical equipment
Part 1-2: General requirements for safety, collateral standard: Electromagnetic disturbances; Requirements and tests
- ISO 80601-2-55:2018
Medical electrical equipment
Part 2-55: Particular requirements for basic safety and essential performance of respiratory gas monitors

List of accessories

Accessory	Part number
Power supply, 100 to 240 VAC, 50/60 Hz	6874204
WaterLock2	6872130
Sample line (set of 10)	8290286
Sample gas return kit	M32692
Sample gas scavenging hose	6873252
Connecting cables	
SCIO/M500 cable, 0.8 m	MS33986
SCIO/M500 cable, 2.2 m	MS33987
Infinity Acute Care System (M540), 0.8 m	MS33168
Infinity Acute Care System (M540), 2.2 m	MS33167
Infinity Gamma XL with interface plate or Infinity Docking Station (X5 connector), 0.3 m	MS13503
Infinity Gamma XL with interface plate or Infinity Docking Station (X5 connector), 1.5 m	MS13864
Infinity Gamma XL (X8 connector), 2.1 m	MS15849 and MS23211
Infinity Delta series/Gamma XXL/Vista XL with Infinity Docking Station (X3 connector), 0.3 m	MS13504
Infinity Delta series/Gamma XXL/Vista XL with Infinity Docking Station (X3 connector), 1.5 m	MS13865
Infinity Delta series/Gamma XXL/Vista XL/Kappa XLT (X8 connector), 1.8 m	6871581
Y-cable for X8 connector	MS24539
Vista 120 (RS232 connector), 2 m	8601473

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© Drägerwerk AG & Co. KGaA
**Drägerwerk AG & Co. KGaA**