

Technical Service Manual

Fabius Tiro® Anesthesia System



Revision B
5/25/04
6020.002
4118302-002

Copyright by Dräger Medical AG & Co. KGaA, Lübeck, Germany.

No reproduction allowed for commercial purposes.

Read and understand the Instructions for Use/Operator's Manual.

This Technical Documentation does not replace the Instructions for Use/Operator's Manual.

The warranty and liability conditions of the general terms and conditions for business transactions of Dräger Medical AG & Co. KGaA are not extended by this Technical Documentation.

Observe all applicable technical laws and regulations.

Insofar as reference is made to laws, regulations or standards, these are based on the legal system of the Federal Republic of Germany. Observe the laws and regulations applicable in your country.



A Dräger and Siemens Company

Fabius Tiro Service Manual Table of Contents

[What's New in Rev. B](#)

| DESCRIPTION | PAGE |
|---|----------|
| General | 1 |
| Recommendations | 3 |
| How To Use This Manual | 3 |
| General Troubleshooting Guidelines..... | 3 |
| Copyright..... | 6 |
| Trademark Notices..... | 6 |
| Disclaimer | 6 |
| Function Description | 7 |
| General Information about the Fabius Tiro | 9 |
| Fabius Tiro Function Diagram..... | 14 |
| Battery Backup..... | 16 |
| Fabius Tiro Piping Diagram..... | 17 |
| Function Description of Gas Box..... | 19 |
| SORC (Sensitive Oxygen Ratio Controller) | 20 |
| Compact Breathing System, Cosy II | 22 |
| Manual Ventilation..... | 25 |
| Spontaneous Breathing..... | 28 |
| Volume/Pressure Mode Ventilation | 32 |
| Cosy II Absorber | 37 |
| Ventilator | 37 |
| High Pressure Safety Valve | 39 |
| Negative Pressure Relief Valve..... | 40 |
| Pneumatic System | 40 |
| PEEP/Pmax Valve Control..... | 41 |
| APL Bypass Valve Control | 41 |
| Electronic Block Diagram..... | 43 |
| Control PCB | 44 |
| Flowmeter | 44 |
| Patient Interface..... | 44 |
| Ventilator | 44 |
| Front Panel Functions..... | 44 |
| Pneumatic Assembly Control..... | 44 |
| Serial Port Interface | 44 |
| Battery..... | 44 |
| Control Panel Assembly..... | 45 |
| FiO2 Measurement | 47 |
| Respiratory Flow Measurement | 48 |
| Gas Flow Rate Measurement | 49 |

| | |
|-----------------------------|--------------------|
| CONTENTS (continued) | FABIUS TIRO |
|-----------------------------|--------------------|

| DESCRIPTION | PAGE |
|---|-------------|
| Vaporizer | 50 |
| ELECTROMAGNETIC TESTING AND COMPLIANCE | 51 |
| Electromagnetic Testing and Compliance | 53 |
| Guidance and manufacturer's declaration-electromagnetic emissions | 53 |
| Guidance and manufacturer's declaration-electromagnetic immunity | 54 |
| Guidance and manufacturer's declaration-electromagnetic immunity | 55 |
| Recommended separation distances between portable and mobile RF communications equipment and the Fabius Tiro | 57 |
| TROUBLESHOOTING GUIDE | 61 |
| Troubleshooting | 63 |
| Power Supply and Voltage Distribution | 63 |
| Battery | 68 |
| Troubleshooting Guides | 68 |
| DIAGNOSTICS | 81 |
| Diagnostics | 83 |
| System Service Screen | 87 |
| Main Service Screen | 89 |
| Service Log | 91 |
| Service Log Entries | 91 |
| Real Time Values | 97 |
| Preventive Maintenance | 101 |
| Activate Preventive Maintenance Date | 101 |
| General | 103 |
| Serial Number | 103 |
| Reset Hours Run | 104 |
| Reset Last Service Date | 105 |
| Select Market Kit | 106 |
| Calibration | 109 |
| Fresh Gas Flow Calibration | 109 |
| Pressure Calibration | 110 |
| O2 Zero Calibration | 111 |
| PEEP Valve Calibration | 112 |
| Configure | 115 |
| System Settings (US and Non-US) | 115 |
| Model Type | 116 |
| Standard Options | 117 |
| Flowmeter | 118 |
| O2 Position (Virtual Flowtubes) | 119 |
| Gas Selection | 120 |
| Flowtube Resolution | 121 |
| O2 Whistle | 122 |

| DESCRIPTION | PAGE |
|--|-------------|
| Alarms | 123 |
| No Fresh Gas | 124 |
| Fresh Gas Low Alarm | 125 |
| Threshold Low Alarm | 126 |
| Pressure | 127 |
| Ambient Pressure | 128 |
| Plateau-Mean Display Screen | 128 |
| Secure Options | 129 |
| Serial Port | 131 |
| Serial Port Parameters | 131 |
| Pump On/Off | 133 |
| | |
| REPLACEMENT PROCEDURES | 135 |
| Replacement Procedures | 137 |
| Core Module Inversion | 139 |
| Core Module Inversion Disassembly | 139 |
| Core Module Inversion Reassembly | 139 |
| Wall Mount Service Access Procedure (If Applicable) | 143 |
| Cylinder Yokes (Fixed) and Regulators | 145 |
| Cylinder Pressure Gauges | 149 |
| Auxiliary Oxygen Flow Meter | 151 |
| Vaporizers | 153 |
| O2 Flush Valve | 155 |
| Caster | 157 |
| Ventilator | 159 |
| Gas Inlet Assembly (including O2 supply pressure switch) | 171 |
| Battery | 173 |
| Power Supply | 175 |
| Control PCB Assembly | 177 |
| Pneumatic (PEEP Control) Assembly | 183 |
| Flow Meter Bezel Assembly | 185 |
| Fresh Gas Flow Meter | 187 |
| Flow Control Valves | 189 |
| Fresh Gas Flow Sensors and Filter Assembly | 191 |
| Pipeline Pressure Gauges | 193 |
| SORC (Sensitive Oxygen Ratio Controller) | 195 |
| Fresh Gas Display PCB | 197 |
| Monitor Bezel Assembly | 199 |
| Spirolog Sensor and Cable | 201 |
| | |
| ADJUSTMENT AND CALIBRATION PROCEDURES | 203 |
| Cylinder Pressure Regulator Adjustment | 205 |
| Gas Inlet Regulator Output Adjustment | 207 |
| Oxygen Supply Pressure Alarm Switch Adjustment | 209 |
| Sensitive Oxygen Ratio Controller (SORC) Adjustment | 211 |

| | |
|-----------------------------|--------------------|
| CONTENTS (continued) | FABIUS TIRO |
|-----------------------------|--------------------|

| DESCRIPTION | PAGE |
|---|-------------|
| Oxygen Sensor Calibration | 214 |
| Pressure Calibration..... | 216 |
| Fresh Gas Flow Calibration..... | 217 |
| PEEP Valve Calibration..... | 218 |
| Vacuum Adjustment | 219 |
| Vacuum Pressure Pump Assembly (Software Version 2.N)..... | 221 |
| Pump Calibration - If Applicable..... | 221 |
| SOFTWARE UPDATE PROCEDURE | 225 |
| Software Update Procedure..... | 229 |
| Requirements..... | 229 |
| Boot Strap Download Procedure..... | 229 |
| PMS PROCEDURE | 241 |
| SPARE AND REPLACEMENT PARTS LIST | 307 |

General

1 Recommendations

Because of the sophisticated nature of Draeger Medical, Inc. anesthesia equipment and its critical importance in the operating room setting, it is highly recommended that only appropriately trained and experienced professionals be permitted to service and maintain this equipment. Please contact DrägerService® at (800) 543-5047 for service of this equipment in North America. For service in Europe, call 49 (451) 882-4222. For service in other countries, call (215) 721-5402.

Draeger Medical, Inc. recommends that the Fabius Tiro be serviced at six month intervals. Periodic Manufacturer's Certification agreements are available for equipment manufactured by Draeger Medical, Inc. Please contact us for further information concerning these agreements.

Draeger Medical, Inc. products/material in need of factory repair shall be sent to:

| | |
|--|---|
| For North America and other countries except Europe: | For service in Europe: |
| DrägerService® 3124 Commerce Drive Telford, PA 18969 U.S.A. (Include RMA Number) | Dräger Medical AG & Co. KGaA Moislinger Allee 53-55 Reparaturannahme 23542 Lübeck Germany |

2 How To Use This Manual

The manual is divided into several sections. The DIAGNOSTICS section describes self-test and service diagnostics for checking the system functions. An understanding of the on-board service capabilities is necessary before any attempt is made to troubleshoot the unit. The TROUBLESHOOTING section lists error codes and provides troubleshooting guides to assist the Technical Service Representative in locating the source of a problem. The REPLACEMENT PROCEDURES section contains instructions for removal and replacement of assemblies that are considered field-replaceable. The ADJUSTMENT AND CALIBRATION PROCEDURES section contains the field procedures needed to restore original system specifications. The Periodic Manufacturer's Service (PMS) PROCEDURE section outlines the steps required to verify the electrical, mechanical, and pneumatic safety of the unit and also, identifies components requiring periodic replacement. The SPARE PARTS section is provided for use as a reference only to obtain part numbers and descriptions for parts and assemblies for replacement purposes. For items not shown in the Spare Parts section, contact DrägerService®.

3 General Troubleshooting Guidelines

Troubleshooting the Fabius Tiro should always begin by communicating with those who observed or experienced a problem with the unit. This may eliminate unnecessary troubleshooting steps. Once a general problem is identified, refer to the flow charts in the Troubleshooting Section to determine the proper corrective action to be taken.

After any component is replaced, verify that the unit is operating properly by running the appropriate diagnostic procedure. The PMS PROCEDURE must also be performed after any component is replaced.

There are two possible mounting configurations for the Fabius Tiro, a wall mounted unit and a trolley mounted unit. A general arrangement of both mounting configurations of the Fabius Tiro anesthesia system are shown in Figure 1 and Figure 2.

SV88002




Figure 1 Fabius Tiro - Wall Mount Configuration




Figure 2 Fabius Tiro - Trolley Mount Configuration

- 4 Copyright** Copyright© 2004 by Draeger Medical, Inc. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, or stored in a retrieval system in any form or by any means, electronic or mechanical, including photocopying and recording, without written permission of Draeger Medical, Inc.
- 5 Trademark Notices** Datagrip, DrägerService, ORM, Quality Service For Life, Respitone, Narkomed, Vigilance Audit, Vitalert, Vitalink, Narkomed GS, Fabius, Fabius Tiro, and Fabius GS are registered trademarks of Draeger Medical, Inc. All other products or name brands are trademarks of their respective owners.
- 6 Disclaimer** The content of this manual is furnished for informational use only and is subject to change without notice. Draeger Medical, Inc. assumes no responsibility or liability for any errors or inaccuracies that may appear in this manual.

Function Description

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

Function Description

Fabius Tiro

1 General Information about the Fabius Tiro

The Fabius Tiro comprises the following assemblies:

- Bezel assembly: Display and Control Panel
- Flowmeter assembly
- Gas Box: Gas Inlet Assembly and related items
- Breathing system
- Pneumatic Assembly
- Ventilator
- Anesthetic Vaporizer(s)
- Trolley Mount and Wall Mount

SV88002




Figure 1 Front View of Fabius Tiro Anesthesia System - Wall Mount Configuration

Fabius Tiro

Function Description




Figure 2 Front View of Fabius Tiro Anesthesia System - Trolley Mount Configuration

F6020002_Function_Description_USA.fm 17.06.04
All rights reserved. Copyright reserved.

Function Description

Fabius Tiro

SV88003




Figure 3 Back View - Interface Panel and Power Entry

SV88057





Figure 4 Side View - Gas Pipeline Inlets and Tethered Yoke - Wall Mount Configuration

F6020002_Function_Description_USA.fm 17.06.04
All rights reserved. Copyright reserved.



SIDE VIEW





Figure 5 Back/Side View - Gas Pipeline and Cylinder Hose Connections - Trolley Mount Configuration

2 Fabius Tiro Function Diagram



SV88038

Figure 6 Functional Diagram of Fabius Tiro - Typical for Breathing System P/N 4116398 or 4117529

SV88078




Figure 7 Function Diagram of Fabius Tiro - Typical for Breathing System P/N 4118378 or 4118379

3 Battery Backup

Fabius Tiro backup power is provided by two series-connected 12 V rechargeable batteries. These batteries remain on charge as long as the machine is plugged into an active AC outlet. Should power to the main controller PCB assembly fail while the machine is in operation, the batteries will allow the machine to continue operating for a minimum of 45 minutes, provided the batteries are fully charged.

The batteries are located within the controller housing and are accessible by opening the ventilator compartment. The 3.15A battery fuse is located on the power entry panel at the back of the machine.




Figure 8 Battery Backup Arrangement

AB602002001

Function Description

Fabius Tiro

**4 Fabius Tiro Piping
Diagram**





Figure 9 Fabius Tiro 3-Gas Piping Diagram

F6020002_Function_Description_USA.fm 17.06.04
 All rights reserved. Copyright reserved.

5 Function Description of Gas Box

The supply gases flow through the filters and non-return valves in the gas inlet assembly. Pipeline supply pressures are indicated on gauges located on the flowmeter assembly. Cylinder pressure gauges are located on the base machine assembly. The pressures of O₂ and N₂O delivered to the flowmeter assembly are set by regulators on the gas inlet assembly.

Should the O₂ supply fail or its pressure decrease below a certain limit, the O₂ supply pressure alarm switch signals an alarm.



SV88039

Figure 10 Gas Box Functional Diagram (Part 1 of 2)

If the O₂ flush button is pressed, oxygen is delivered to the fresh gas outlet. The fresh-gas ejector prevents the fresh gas from flowing back into the anesthetic vaporizer. This avoids an increase in the anesthetic gas concentration.




Figure 11 Gas Box Functional Diagram (Part 2 of 2)

6 SORC (Sensitive Oxygen Ratio Controller)

The SORC is a control element that functions like an N₂O shut-off device and ensures a vital O₂ concentration in the fresh gas. In the event of an O₂ shortage, the SORC limits the N₂O flow such that the O₂ concentration in the fresh gas does not decrease below 23 vol.%.

If the O₂ flow control valve is closed or if the O₂ flow is lower than or equal to 200 mL/min, the SORC interrupts the N₂O flow.

N₂O can be added when the O₂ flow is approx. 300 mL/min. In this case, the SORC also prevents O₂ concentrations below 23 vol.%.

The SORC bypass allows O₂ to bypass the restrictor in the SORC when O₂ flows above 10 L/min. are needed.

F6020002_Function_Description_USA.fm 17.06.04
All rights reserved. Copyright reserved.




Figure 12 SORC Functional Diagram (Part 1 of 2)

The flow control valves are used to adjust the O2 and N2O flows.

Restrictors located at the outlets of the SORC generate back-pressures. These back-pressures exert a force on the control diaphragms of the SORC. The O2 back-pressure opens the SORC. The N2O back-pressure closes the SORC. The pressure ratio at the control diaphragm affects the N2O flow.

The restrictors and the spring tension are dimensioned such that a minimum concentration of 23 vol.% O2 is always ensured. The maximum O2 flow is approx. 12 L/min.




Figure 13SORC Functional Diagram (Part 2 of 2)

7 Compact Breathing System, Cosy II

The Cosy II compact breathing system allows various modes of patient ventilation: manual and spontaneous breathing, volume controlled, and pressure controlled.


In the "MAN" position, the compact breathing system is closed to atmosphere. This position is used for manual ventilation of the patient. The APL valve opening pressure can be adjusted from 5 to 70 cmH₂O.

In the "SPONT" position the APL valve is open to atmosphere. This position is used for spontaneous patient breathing.

The pressure limit (P_{max}) can also be adjusted (through front panel interface) during volume control from 15 cmH₂O to 70 cmH₂O using the control box and the PEEP/P_{max} valve.

Function Description

Fabius Tiro



SV88061

Figure 14 Compact Breathing System, Cosy II - Typical for P/N 4116398 or 4117529

Table 1 Key, Figure 14

| | |
|----|---|
| 1 | PEEP/Pmax valve |
| 2 | Flow sensor (Spirolog) |
| 3 | MAN/SPONT-APL Valve |
| 4 | Fresh-gas decoupling valve |
| 5 | Breathing bag hook |
| 6 | Inspiratory port |
| 7 | Inspiratory valve and O ₂ sensor port |
| 8 | Breathing bag terminal and Y-piece plug |
| 9 | Expiratory valve |
| 10 | Expiratory port |
| 11 | APL Bypass valve |
| 12 | Fresh-gas port |
| 13 | Ventilator port |
| 14 | Anesthetic gas scavenging port |
| 15 | Absorber |
| 16 | Pressure sensor connector |
| 17 | Exhaust port, anesthetic monitor return (non-U.S. systems only) |




Figure 15 Compact Breathing System, Cosy II - Typical for P/N 4118378 or 4118379

Table 2 Key, Figure 15

| | |
|----|---|
| 1 | PEEP/Pmax valve |
| 2 | Flow sensor (Spirolog) |
| 3 | MAN/SPONT-APL Valve |
| 4 | Fresh-gas decoupling valve |
| 5 | Breathing bag hook |
| 6 | Inspiratory port |
| 7 | Inspiratory valve and O2 sensor port |
| 8 | Breathing bag terminal and Y-piece plug |
| 9 | Expiratory valve |
| 10 | Expiratory port |
| 11 | APL Bypass valve |
| 12 | Fresh-gas port |
| 13 | Ventilator port |
| 14 | Anesthetic gas scavenging port |
| 15 | Absorber |
| 16 | Pressure sensor connector |
| 17 | Exhaust port, anesthetic monitor return (non-U.S. systems only) |

7.1 Manual Ventilation

Manual Ventilation: General

During manual ventilation, the APL valve is set to the "MAN" position. The patient system safety valve is activated. The piston of the ventilator is in the upper end position in order to reduce the volume of the ventilator.

Manual Ventilation: Inspiration

During inspiration, expiratory valve remains closed. When the clinician compresses the breathing bag the gas mixture (expiratory gas and fresh gas) flows through the fresh-gas decoupling valve, the inspiratory valve, the O₂ sensor, the inspiratory hose, and the Y-piece into the patient's lung. The pressure sensor measures the airway pressure. The APL valve limits the ventilation pressure. Any excess amount of the gas mixture flows through the APL valve and the non-return valve to the anesthetic gas scavenging system.




Figure 16 Manual Ventilation (Inspiration) - Typical for Breathing System P/N 4116398 or 4117529




Figure 17 Manual Ventilation (Inspiration) - Typical for Breathing System P/N 4118378 or 4118379

Manual Ventilation: Expiration - Typical for Breathing System P/N 4116398 or 4117529

During expiration, the inspiratory valve remains closed and thus prevents the expiratory gas from flowing back into the inspiratory branch.

After releasing the breathing bag, the expiratory gas from the lung flows through the expiratory hose, the flow sensor, the PEEP/Pmax valve, the expiratory valve, and through the absorber into the breathing bag. At the same time, new fresh gas flows into the breathing bag.




Figure 18 Manual Ventilation (Expiration) - Typical for Breathing System P/N 4116398 or 4117529

Manual Ventilation: Expiration - Typical for Breathing System P/N 4118378 or 4118379

During expiration, the inspiratory valve remains closed and thus prevents the expiratory gas from flowing back into the inspiratory branch.

After releasing the breathing bag, the expiratory gas from the lung flows through the expiratory hose, the flow sensor, the PEEP Pmax valve, the expiratory valve, onto the breathing bag and through the absorber. At the same time, new fresh gas flows into the breathing bag.




Figure 19 Manual Ventilation (Expiration) - Typical for Breathing System P/N 4118378 or 4118379

7.2 Spontaneous Breathing

Spontaneous Breathing: General

A prerequisite for spontaneous breathing is that the patient is supplied with a sufficient amount of fresh gas. The APL valve selector must be set to the "SPONT" position. No pressure builds up in the compact breathing system.

Spontaneous Breathing: Inspiration

During inspiration, the expiratory valve remains closed thus preventing rebreathing of expiratory gas containing CO₂.

Function Description

Fabius Tiro

The patient inhales the gas mixture (expiratory gas and fresh gas) from the breathing bag. The gas mixture flows through the fresh-gas decoupling valve, the inspiratory valve, the O₂ sensor, the inspiratory hose, and through the Y-piece into the lung. The pressure sensor measures the airway pressure.




Figure 20 Spontaneous Breathing (Inspiration) - Typical for Breathing System P/N 4116398 or 4117529




Figure 21 Spontaneous Breathing (Inspiration) - Typical for Breathing System P/N 4118378 or 4118379

Spontaneous Breathing: Expiration - Typical for Breathing System P/N 4116398 or 4117529

During expiration, the inspiratory valve remains closed thus preventing the expiratory gas from flowing back into the inspiratory branch.

The APL valve is open, regardless of its pressure setting.

The expiratory gas flows from the lung through the expiratory hose, the flow sensor, the PEEP control valve, the expiratory valve, and through the absorber into the breathing bag. At the same time, new fresh gas flows into the breathing bag.

Function Description

Fabius Tiro

When the breathing bag is full, any excess gas mixture flows through the non-return valve into the anesthetic gas scavenging system.




Figure 22 Spontaneous Breathing (Expiration) - Typical for Breathing System P/N 4116398 or 4117529

The CO₂ is scrubbed from the expiratory gas by the soda lime contained in the absorber. The fresh gas replaces the anesthetic and oxygen taken up by the patient.

Spontaneous Breathing: Expiration - Typical for Breathing System P/N 4118378 or 4118379

During expiration, the inspiratory valve remains closed thus preventing the expiratory gas from flowing back into the inspiratory branch.

The APL valve is open, regardless of its pressure setting.

The expiratory gas flows from the lung through the expiratory hose, the flow sensor, the PEEP control valve, the expiratory valve, the breathing bag, and through the absorber. At the same time, new fresh gas flows into the breathing bag.

When the breathing bag is full, any excess gas mixture flows through the non-return valve into the anesthetic gas scavenging system.




Figure 23 Spontaneous Breathing (Expiration) - Typical for Breathing System P/N 4118378 or 4118379

7.3 Volume/Pressure Mode Ventilation

Volume Control Mode: General

A prerequisite for volume control is that the patient is supplied with a sufficient amount of fresh gas.

The APL bypass valve opens in volume mode, allowing excess gas to be vented to the scavenging system regardless of the MAN-SPONT valve setting.

Function Description

Fabius Tiro

The safety valve of the patient system makes sure that no pressures greater than 75 cmH₂O build up in the system.

During ventilation, the pressure limit (P_{max}) can be adjusted on the control box.

Volume/Pressure Control Mode: Inspiration

During inspiration, the PEEP/P_{max} valve remains closed. The control pressure present at the PEEP/P_{max} valve varies with the set pressure limit (P_{max}).

The pressure generated by the ventilator's piston closes the fresh-gas decoupling valve. The gas mixture (expiratory gas and fresh gas) flows through the inspiratory valve, the O₂ sensor, the inspiratory hose, and through the Y-piece into the lung. The pressure sensor measures the airway pressure. The ventilation pressure cannot exceed the pressure limit (P_{max}) set on the control box because the PEEP/P_{max} valve opens. The fresh gas then fills the breathing bag.

Any excess fresh gas flows through the open APL bypass valve, and through the non-return valve into the anesthetic gas scavenging system.




Figure 24 Volume Control Ventilation (Inspiration) - Typical for Breathing System P/N 4116398 or 4117529




Figure 25 Volume Control Ventilation (Inspiration) - Typical for Breathing System P/N 4118378 or 4118379

Volume/Pressure Control Mode: Expiration

During expiration, the inspiratory valve remains closed thus preventing rebreathing into the inspiratory branch.

The expiratory gas from the lung flows through the expiratory hose, the flow sensor, the PEEP/Pmax valve, the expiratory valve, and through the absorber back into the breathing bag mixing with fresh gas also flowing into the breathing bag.

The ventilator's piston moves back drawing the gas mixture needed for the next inspiration into the piston space.

Fabius Tiro

Function Description

Any excess fresh-gas flows through the APL bypass valve, and through the non-return valve into the anesthetic gas scavenging system.




Figure 26 Volume Control Ventilation (Expiration) - Typical for Breathing System P/N 4116398 or 4117529




Figure 27 Volume Control Ventilation (Expiration) - Typical for Breathing System P/N 4118378 or 4118379

7.4 Cosy II Absorber

The absorber canister is filled with fresh soda lime. The soda lime scrubs CO₂ from the respiratory expired gas.

Expired soda lime changes its color. The soda lime must be replaced when two thirds of the soda lime in the absorber canister is discolored.

8 Ventilator

The ventilator is located in a swing-out compartment at the left side of the Fabius Tiro. A hose terminal is provided on the left side of the compartment for connection to the breathing system. Fresh gas is delivered to the patient by a piston that is driven by a motor and ball-screw arrangement. A sight window on the compartment allows the operator to verify movement of the piston.

Fabius Tiro

Function Description

Two diaphragms (upper and lower) comprise a bag-type rolling seal that surrounds the piston. Vacuum from the pneumatic assembly (described in a later paragraph) is provided between the outside of the seal and the cylinder, to ensure proper operation of the seal during piston movement.


During inspiration the ventilator delivers fresh gas at a given volume, pressure and frequency. These parameters are set at the control panel. Refer to the Operator's Manual for details on ventilator settings, displays and controls. During expiration, the bag-type rolling seal fills with expired gas from the patient and with fresh gas stored in the breathing bag.

Power for the ventilator motor is distributed from the control PCB. A position sensor on the ventilator signals the control PCB when the piston reaches its lower limit. An incremental encoder on the motor shaft determines the number of revolutions and provides piston travel information to the control PCB.

Ventilator pressure is monitored by a transducer on the control PCB. Should the negative pressure relief valve on the patient assembly open, a Fresh Gas Low alarm is displayed if enabled via service mode.

The ventilator pressure transducer is the same type as the one used for measuring airway pressure. A hose connects the transducer's positive pressure port to a hose barb located on the top cover of the ventilator. The purpose of this transducer is to allow the software to sense when a condition exists that would cause the ventilator negative pressure relief valve to open. The threshold used by the software for this condition is -8 mbar. In normal use the primary cause for this condition is an insufficient amount of reserve gas in the breathing bag. The operator is alerted when this condition exists, with a medium priority FRESH GAS LOW alarm. This alarm may be disabled via service mode.

The ventilator assembly is illustrated on [Figure 28](#).



SV00513

Figure 28 Ventilator (piston shown in down position)

The top of the ventilator assembly (patient system) contains two valves:

8.1 High Pressure Safety Valve

If the pressure limit control fails, the patient system high pressure safety valve limits the gas pressure. This valve is set to open at approximately 75 cmH₂O.




Figure 29 Sectional View of the Safety Valve

8.2 Negative Pressure Relief Valve

The negative pressure relief valve allows the patient to spontaneously breathe ambient air should the medical gas supply and/or Fabius Tiro fail. The opening pressure of this valve is -8 mbar.




Figure 30 Sectional View of the Negative Pressure Relief Valve

9 Pneumatic System

The pneumatic assembly provides pressure for the PEEP valve control, and also provides vacuum for the ventilator bag-type rolling seals and the APL bypass valve control.

F6020002_Function_Description_USA.fm 17.06.04
 All rights reserved. Copyright reserved.




Figure 31 Pneumatic Control System Schematic

9.1 PEEP/Pmax Valve Control

When the Fabius Tiro is operating in the automatic mode, the pump on the pneumatic assembly is running, and the electronic PEEP valve receives a signal from the main control PCB. The amount of current supplied to the coil of the electronic PEEP valve is proportional to the PEEP value set by the operator, and controls the position of the diaphragm within the electronic PEEP valve. This then determines the control pressure applied to the proportional PEEP valve in the breathing system, which maintains the desired amount of PEEP during patient expiration. The V1 reservoir smooths out pressure variations caused by the pump. See [Figure 29](#).

9.2 APL Bypass Valve Control

When the Fabius Tiro is operating in the automatic mode, the pneumatic assembly provides a vacuum signal to hold open the APL bypass valve in the breathing system. The V2 reservoir and filter provide noise damping, and the variable restrictor is used to set the vacuum level in the range of -150 to -240 cmH2O.

Fabius Tiro

Function Description

When the machine is operating in the Manual mode, the pump on the pneumatic assembly (and the ventilator) is stopped, and the spring-loaded APL bypass valve in the breathing system closes, directing exhaled gas through the APL valve.

10 Electronic Block Diagram




Figure 32 Electronic Block Diagram

F6020002_Function_Description_USA.fm 17.06.04
All rights reserved. Copyright reserved.

| Fabius Tiro | Function Description |
|--|---|
| 11 Control PCB | The control PCB provides the following functions in the Fabius Tiro: |
| 11.1 Flowmeter | Receives flow rate data from the N ₂ O, Air, and O ₂ flow sensors and provides fresh gas display information for these gases. |
| 11.2 Patient Interface | Receives information from the oxygen sensor and processes it for display. Receives flow information from the Spirolog sensor and processes it for display. Converts airway breathing pressure to an electrical signal and processes it for display. |
| 11.3 Ventilator | Provides ventilator motor drive, receives ventilator piston position and movement information. Provides ventilator pressure information. |
| 11.4 Front Panel Functions | Provides power for the display panel and LED table lamp. Provides video signals to the display. Receives information from the keypad and rotary encoder for making display and operating selections, setting alarm limits, and service functions. |
| 11.5 Pneumatic Assembly Control | Provides power to the pump on the pneumatic assembly. Provides control signal to the electronic PEEP valve on the pneumatic assembly in response to operator setting. Monitors pump vacuum to stabilize PEEP control. |
| 11.6 Serial Port Interface | Provides isolated port for connecting external monitors and downloading software. |
| 11.7 Battery | Provides charging current for battery and monitors state of battery charge. |




Figure 33 Controller Functional Block Diagram

12 Control Panel Assembly

The control panel assembly consists of a 320 x 240 pixel graphical display, a table lamp with six LEDs, a membrane keypad, rotary encoder and speaker.

Data and power for the display comes from the main controller PCB via a 20-conductor ribbon cable. The keypad interface is connected to the main controller PCB by a 30-conductor ribbon cable. A block diagram of the control panel assembly is shown in [Figure 34](#).




Figure 34 Control Panel Block Diagram

An illustration of the control panel is shown in [Figure 35](#).




Figure 35 Fabius Tiro Control Panel (Main Service Screen Illustrated)

Descriptions of the numbered items in [Figure 35](#) are given in the [Table 3](#).

Table 3 Fabius Tiro Control Panel Key for [Figure 35](#)

| ITEM | FUNCTION |
|------|--|
| 1 | Selects volume controlled ventilation mode. Refer to Operator's Manual |
| 2 | Select pressure controlled ventilator mode. Refer to Operator's Manual |
| 3 | Reserved for future use |
| 4 | Reserved for future use |
| 5 | Controls table lamp: Off/On |
| 6 | Places ventilator in Man Spont mode. Refer to Operator's Manual |
| 7 | Soft Keys: active the corresponding function that appears on screen above the key |
| 8 | For setting alarm limits. Refer to Operator's Manual |
| 9 | Setup Key: activates sub-screens for monitoring functions. Refer to Operator's Manual |
| 10 | Home key: returns display to main screen shown before standby |
| 11 | Rotary control: moves the cursor on the screen; confirms selection when pressed |
| 12 | Alarm Status Indicators: Flashing Red: Warning; Flashing Yellow: Caution; Solid Yellow: Advisory |
| 13 | Alarm Silence key: silences all active alarms for two minutes |
| 14 | Power ON indicator: lighted when machine is plugged into an active AC outlet |
| 15 | Returns unit to Standby mode |

13 FiO2 Measurement

The O2 sensor measures the fraction of inspired O2 (FiO2) in the respiratory gas.

The O2 sensor contains an alkaline electrolyte, a lead anode, two gold cathodes, and a Teflon membrane. The spatial separation of the two gold cathodes allows a voltage comparison to be made as explained below.

The O2 sensor is an electrochemical cell that generates a voltage from the ion current.




Figure 36 O₂ Sensor

The O₂ to be measured diffuses through the Teflon membrane, reacts at the gold cathodes (negative polarity) and forms lead oxide and water at the lead anode (positive polarity). During this chemical process, a voltage is generated which is proportional to the O₂ partial pressure.

The internal resistance of the cell is determined by the surface area of the gold cathodes, the O₂ diffusion velocity, the distance between the gold cathodes and the lead anode. This resistance is approximately 700 ohms.

The chemical process is temperature-sensitive. Therefore temperature-sensitive resistors are connected in parallel with the O₂ sensor. These resistors and the internal resistance of the O₂ sensor correct the measuring voltage. Because two cathodes used in the O₂ sensor cell, two different voltages are generated. These voltages are compared with each other. If their difference exceeds a certain value, the machine prompts the operator to check the cell.

If the O₂ sensor fails, the control box will indicate an error on the graphics display.

14 Respiratory Flow Measurement

The flow sensor functions according to the constant temperature hot-wire anemometer principle. Respiratory gas flows past a thin platinum wire. This platinum wire (A) is located in a measuring tube and is electrically heated. The platinum wire is held at a constant temperature. Gas flow removes heat from the hot wire. The higher the gas flow rate, the greater the heat removal. The amount of electrical current needed to maintain a constant platinum wire temperature is thus proportional to the gas flow rate.

Function Description

Fabius Tiro

A second platinum wire (B) in the measuring tube is used to compensate for interferences from different gases present in the respiratory gas. The heat removed from the second platinum wire is measured during inspiration when the gas flow is zero.

The different gases present in the respiratory gas have a different thermal conductivity. The amount of heat removed from the second platinum wire is thus an indicator of respiratory gas composition.

Internal calibration tables for O₂/N₂O mixtures, Air and 100% O₂ are used to linearize the measured flow.




Figure 37 Respiratory Flow Sensor

15 Gas Flow Rate Measurement

The gas flow sensors operate on the principle of specific heat for individual gases. In each sensor, as the gas flows through a heated chamber the gas molecules carry away a certain amount of heat relative to the specific heat index for that gas.

A known amount of electrical current is required to maintain the temperature in the heated chamber. The higher the gas flow rate, the more heat is removed from the chamber and more current is required to maintain the temperature in the chamber. This current is then scaled and displayed as liters per minute flow rate for each gas.

Fabius Tiro

Function Description




Figure 38 Flow Sensor Details




Figure 39 Flow of Gases through Sensors

16 Vaporizer

Refer to separate technical documentation of the anesthetic vaporizer.

F6020002_Function_Description_USA.fm 17.06.04
All rights reserved. Copyright reserved.

ELECTROMAGNETIC TESTING AND COMPLIANCE

1 Electromagnetic Testing and Compliance

The following CAUTIONS and WARNINGS are applicable to the Fabius Tiro Anesthesia Machine:



CAUTION

Do not use Fabius Tiro in the environment of NMR tomography equipment. Malfunctions may result, thereby endangering the patient.



CAUTION

The use of portable and mobile radio frequency communications equipment can affect medical electrical equipment. Do not use mobile phones within a distance of 10 meters from the machine. Mobile phones can cause malfunctions in electrical medical equipment, thereby endangering the patient and the operator.



WARNING

No third-party components shall be attached to the anesthesia machine, ventilator, or breathing system (except for certain approved exceptions). For more information, contact your local Authorized Service Organization or DrägerService at:

DrägerService
Draeger Medical, Inc.
3122 Commerce Drive
Telford, PA 18969
Tel:(215) 721-5402
(800) 543-5047
Fax:(215) 721-5784



CAUTION

Although the Fabius Tiro is designed to minimize the effects of ambient radio-frequency interference, machine functions may be adversely affected by the operation of electrosurgical equipment or short wave or microwave diathermy equipment in the vicinity.



CAUTION

Communications with external equipment may be temporarily affected by electromagnetic interference due to the use of electrosurgical equipment.

Tables 1 through 4 contained in this section are provided for informational purposes and general guidance, as required by IEC 60601-1-2:2001, Medical Electrical Equipment, Part 1-2, General Requirements for Safety - Collateral Standard: Electromagnetic Compatibility.

2 Guidance and manufacturer's declaration-electromagnetic emissions

The Fabius Tiro is intended for use in the electromagnetic environment specified below. The customer or the user of the Fabius Tiro should assure that it is used in such an environment.

Table 1 Electromagnetic Emissions

| Emissions test | Compliance level | Recommendations |
|---|------------------|---|
| RF emissions CISPR 11 | Group 1 | The Fabius Tiro uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | The Fabius Tiro is suitable for use in all establishments, except NMR environments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 Class A | Complies | |
| Voltage fluctuations/flicker emissions IEC 61000-3-3 | Complies | |

3 Guidance and manufacturer's declaration-electromagnetic immunity

The Fabius Tiro is intended for use in the electromagnetic environment specified below. The customer or the user of the Fabius Tiro should assure that it is used in such an environment.

Table 2 Electromagnetic Immunity

| Immunity test | IEC 60601 test level | Compliance level | Recommendations |
|---|--|------------------|--|
| Electrostatic discharge (ESD) IEC61000-4-2 | +/-6 kV contact +/-8 kV air | Complies | Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast Transient/burst IEC61000-4-4 | +/-2 kV for power supply lines +/-1 kV for input/output lines | Complies | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC61000-4-5 | +/-1 kV differential mode +/-2 kV common mode | Complies | Mains power quality should be that of a typical commercial or hospital environment. |

All rights reserved. Copyright reserved.
L6020002_EMC.fm 17.06.04

| Immunity test | IEC 60601 test level | Compliance level | Recommendations |
|---|--|------------------|---|
| Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11 | <5% U_t (>95% dip in U_t) for 0.5 cycle 40% U_t (60% dip in U_t) for 5 cycles 70% U_t (30% dip in U_t) for 25 cycles <5% VT (>95% dip in U_t) for 5 seconds | Complies | Mains power quality should be that of a typical commercial or hospital environment. The Fabius Tiro provides battery back-up in the event of a power failure. |
| Power frequency (50/60 Hz) magnetic field IEC61000-4-8 | 3 A/m | Complies | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |



NOTE

U_t is the a.c. mains voltage prior to application of the test level.

4 Guidance and manufacturer's declaration-electromagnetic immunity

The Fabius Tiro is intended for use in the electromagnetic environment specified below. The customer or the user of the Fabius Tiro should assure that it is used in such an environment.

Table 3 Electromagnetic Immunity

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment-guidance |
|---------------|----------------------|------------------|--|
| | | | Portable and mobile RF communications equipment should be used no closer to any part of the Fabius Tiro, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. |

Electromagnetic Testing and Compliance

Fabius Tiro

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment-guidance |
|------------------------------|---|----------------------|--|
| | | | Recommended separation distance |
| Conducted RF IEC61000-4-6 | 3 Vrms 150 kHz to 80 MHz outside ISM bands (ref. NOTE A) | [V1] V Complies | $d=[3.5/V1]vP$ |
| | 10 Vrms 150 kHz to 80 MHz in ISM bands (Ref. NOTE A) | [V2] V Complies | $d=[12/V2]vP$ |
| Radiated RF IEC61000-4-3 | 10 V/m 80 MHz to 2.5 GHz | [E1] V/m Complies | $d=[12/E1]vP$ 80 MHz to 800 MHz $d=[23/E1]vP$ 800 MHz to 2.5 GHz |
| | | | where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m) (Ref. NOTE B). |
| | | | Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (ref. NOTE C), should be less than the compliance level in each frequency range (Ref. NOTE D). |
| | | | Interference may occur in the vicinity of equipment marked with the following symbol: ((e»)) |



NOTE

At 80 MHz and 800 MHz, the higher frequency range applies.



NOTE

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.



NOTE A

The ISM (industrial, scientific, and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.



NOTE B

The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.



NOTE C

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Fabius Tiro is used exceeds the applicable RF compliance level above, the Fabius Tiro should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Fabius Tiro.



NOTE D

Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 1 V/m.

5 Recommended separation distances between portable and mobile RF communications equipment and the Fabius Tiro

The Fabius Tiro is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Fabius Tiro can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Fabius Tiro as recommended below, according to the maximum output power of the communications equipment.

Table 4 Distance Recommendations between RF Communications

Equipment and the Fabius Tiro

| Rated maximum output power of transmitter W(atts) | Separation distance according to frequency of transmitter (meters) | | | |
|--|--|--|---|--|
| | 150 kHz to 80 MHz outside ISM bands $d=[3.5/\sqrt{P}] \sqrt{P}$ | 150 kHz to 80 MHz in ISM bands $d=[12/\sqrt{E1}]\sqrt{P}$ | 80 MHz to 800 MHz $d=[12/\sqrt{E1}]\sqrt{P}$ | 800 MHz to 2.5 GHz $d=[23/\sqrt{E1}]\sqrt{P}$ |
| 0.01 | 0.116 | 0.120 | 0.120 | 0.230 |
| 0.1 | 0.368 | 0.379 | 0.379 | 0.727 |
| 1 | 1.166 | 1.200 | 1.200 | 2.300 |
| 10 | 3.689 | 3.794 | 3.794 | 7.273 |
| 100 | 11.66 | 12.000 | 12.000 | 23.000 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.



NOTE

At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.



NOTE

The ISM (industrial, scientific, and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.



NOTE

An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.



NOTE

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

TROUBLESHOOTING GUIDE

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

Troubleshooting

This section contains information to assist the Draeger qualified Service Representative in locating electrical faults affecting operation of the Fabius Tiro. A simplified electrical block diagram is given in [Figure 1](#). Since most troubleshooting efforts begin with verifying power supply voltages, the following paragraph outlines the voltage distribution scheme within the machine along with test points for each of the voltages.

A piping diagram ([Figure 4](#)) and a pneumatic control diagram ([Figure 5](#)) are also included in this section.

1 Power Supply and Voltage Distribution

In the Fabius Tiro, the power supply delivers 28 VDC to the main controller unit PCB assembly. Some of the other voltages derived on this assembly can be measured at the connectors listed in Table 1. Fuse data is listed in Table 2. Controller PCB connector locations are shown in [Figure 2](#).

Table 1 Test Points and Allowable Ranges

| Controller PCB | Function | Acceptable Range |
|----------------|------------------------|-----------------------|
| J6-1, + 28 VDC | Power In | 26.6 VDC to 29.4 VDC |
| J6-4, Common | | |
| J1-1, + 5 VDC | Fresh Gas Display | 4.75 VDC to 5.25 VDC |
| J1-14, Common | | |
| J2-1, +12 VDC | Front Panel Display | 11.5 VDC to 12.5 VDC |
| J2-20, Common | | |
| J16-1, +10 VDC | Fresh Gas Flow Sensors | 9.99 VDC to 10.01 VDC |
| J16-9, Common | | |

Table 2 Fuse Data

| Fuse | Rating | Location | Function |
|----------|----------|----------------|------------|
| External | 2.5 A, T | Power Entry | Mains |
| F4 | 3.15 A | Power Entry | Battery |
| F1 | 4.0 A, T | Controller PCB | Battery |
| F2 | 2.5 A,T | Controller PCB | +28 V dist |
| F3 | 1.6 A,T | Controller PCB | +28 V dist |




Figure 1 Fabius Tiro Block Diagram (Electrical)

R6020002_Troubleshooting_USA.fm 17.06.04
All rights reserved. Copyright reserved.




Figure 2 Main Controller PCB Connector Locations (P/N 4116632)




Figure 3 Main Controller PCB Connector Locations (P/N 4118079)




Figure 4 Fabius Tiro Piping Diagram

R6020002_Troubleshooting_USA.fm 17.06.04
 All rights reserved. Copyright reserved.




Figure 5 Fabius Tiro Pneumatic Control Diagram

2 Battery

The Fabius Tiro backup battery comprises two series-connected 12 V rechargeable batteries. While the machine is operating from AC, the total battery voltage at full charge should be within the range of 27.0 to 29.6 V. Battery voltage can be measured at J5 on the main controller PCB assembly. During battery operation, the low battery cutoff voltage should be within the range of 21.0 to 20.0 V.

3 Troubleshooting Guides

Table 3 lists common failure modes and symptoms (excluding simultaneous multiple faults) for the monitoring and display devices in the Fabius Tiro. Each failure mode or symptom is keyed to a troubleshooting guide flow chart at the back of this section to assist the service representative in locating the cause of a problem. These flow charts assume that the machine is plugged into an AC outlet with the correct voltage, and the machine is not running on its backup battery.

Table 3 Fabius Tiro Failure Mode and Symptom List

| Failure Mode/Symptom | Corrective Action |
|---|---------------------------|
| Loss of O2 Monitor | Figure 6 |
| Loss of Breathing Pressure Monitor | Figure 7 |
| Loss of Respiratory Volume Monitor | Figure 8 |
| No Audio Alarms | Figure 9 |
| Serial Port Communication failure | Figure 10 |
| No Oxygen Supply Pressure Alarms | Figure 11 |
| Display Blank Upon System Power-up | Figure 12 |
| Keypad Inoperative | Figure 13 |
| Ventilator Inoperative | Figure 14 |
| Fresh Gas Display Blank or Not Working Properly | Figure 15 |




Figure 6 Loss of O2 Monitor

R6020002_Troubleshooting_USA.fm 17.06.04
All rights reserved. Copyright reserved.




Figure 7 Loss of Breathing Pressure Monitor




Figure 8 Loss of Respiratory Volume Monitor




Figure 9 No Audio Alarms

R6020002_Troubleshooting_USA.fm 17.06.04
All rights reserved. Copyright reserved.




Figure 10 Serial Port Communication Failure






Figure 11 No O2 Supply Pressure Alarms



SV88070

Figure 12 Main Display Blank Upon System Power-up

R6020002_Troubleshooting_USA.fm 17.06.04
 All rights reserved. Copyright reserved.



R6020002_Troubleshooting_USA.fm 17.06.04
All rights reserved. Copyright reserved.

Figure 13 Keypad Inoperative




Figure 14 Ventilator Inoperative




Figure 15 Fresh Gas Display Blank or Not Working Properly

DIAGNOSTICS

Diagnostics



NOTE

The screen illustrations contained in this section are for reference only and therefore may or may not reflect the software version currently installed.

The Fabius Tiro diagnostic system monitors and records the status of its internal hardware when the machine is turned on. The status of each test is displayed on the Power-up screen as shown in [Figure 1](#). This screen is displayed for several seconds before proceeding to the Standby screen. The power-up screen also displays one of three messages at completion of the diagnostics:

| | |
|--------------------------|---|
| FUNCTIONAL | This message indicates that the Fabius Tiro has passed all power-up tests and is fully functional. The machine will proceed to the Standby screen (Figure 2) after a short delay. |
| CONDITIONALLY FUNCTIONAL | This message indicates that a minor problem has been detected. The Fabius Tiro may be used, but your local authorized service organization or DrägerService should be notified to correct the problem. Press the rotary control to proceed to the Standby screen. |
| NON-FUNCTIONAL | This message indicates that a serious problem has been detected, and the machine will not proceed to the Monitor screen. Do not use the machine. Immediately notify your local Authorized Service Organization or DrägerService to correct the problem. |

The following quick reference table is for use in quickly locating the Service Screens necessary to perform specific functions within the software. The indentions shown below indicates the relationship of the Service Screen/function relative to the location of each Service Screen within the software.

| | SERVICE ACCESS SCREEN | ENTRIES/ACTIONS |
|----------------------------------|-------------------------|---------------------------------------|
| Service Log: | | |
| | Entries | Clear Service Log |
| | Real Time Entries | Displays current data values and data |
| Preventative Maintenance: | | Month/Day/Year |
| General: | Serial Number | 8 Digit Serial Number |
| | Reset Hours Run | Yes/No |
| General: | Reset Last Service Date | Yes/No |

Fabius Tiro

Diagnostics

SERVICE ACCESS SCREEN

ENTRIES/ACTIONS

Select Market Kit U.S./Non-U.S.

Calibration:

Fresh Gas Flow Store zero/Exit
Pressure Store zero/Exit
O2 Store zero/Exit
Peep Calibrate/Exit

Configure:

Model Type Fabius GS/Tiro
Standard Options:

Flowmeter:

O2 Position Left/Right
Gas Selection 2 Gas/3 Gas
Flowtube Resolution Normal/High Res

O2 Whistle:

Change State Enable/Disable

Alarms:

No Fresh Gas Enable/Disable
Fresh Gas Low Enable/Disable
Threshold Low Enable/Disable

Pressure:

Ambient Pressure Adjust
Plateau/Mean Display Change Display


Model Type Enable/Disable

Serial Port:

Parameters:




Baud Rate 1200, 2400, 4800, 9600, 19200, 38400
Parity NONE/EVEN/ODD
Stop Bits 1,2
Data Bits 7,8
Protocol VITALINK/MEDIBUS

The Preventive Maintenance Due message will appear on the screen if the current date exceeds the Periodic Manufacturer's Service due date stored in the machine.

| SYSTEM DIAGNOSTICS | | System Status | |
|----------------------------|------|---|--|
| Watch Dog | Pass |  | |
| System RAM | Pass | | |
| Program Memory | Pass | | |
| Video Test | Pass | | |
| Interrupts | Pass | | |
| A/D Converter | Pass | | |
| NV RAM | Pass | | |
| Serial Port | Pass | | |
| Clock | Pass | | |
| Speaker | Pass | | |
| Main Power | Pass | | |
| Battery | Pass | | |
| Dräger | | | |
| Fabius Tiro SW XXX CRC XXX | | | |

AB6020002011

Figure 1 Power-up Diagnostics Screen

| Standby | | | 10:51 AM | | | | | |
|---|---|---|---|-----------------------------|---------------------------|-------------------------|------------------------|-----------------------------|
|  |  |  | Last system test run on 06/13/03 10:50 AM Sleep Mode will activate in 1 min 40 sec To start operation press one of the keys located to the left of the display SWV 2.02 3EF3 | | | | | |
| O2 | Air | N2O | Run System Test | Calibrate Flow Sensor | Calibrate O2 Sensor | Leak / Compl Test | Access Alarm Log | Restore Site Defaults |

AB6020002012

Figure 2 Standby Screen

Further diagnostic functions are available through service screens that can be called up at the display panel. If no display is present upon system power-up, refer to [Troubleshooting Guide](#) Section for troubleshooting assistance.



NOTE

During display of the Standby screen, a five minute count-down appears on the screen, after which the display changes to a Sleep Mode. Press any key on the panel to return to the Standby screen.

To access the System Service screen, press and hold the Home and Standby keys, and press the rotary knob.

1 System Service Screen

The System Service Screen shown in [Figure 3](#) displays system service data.

Press the rotary knob to go to the Main Service screen.





Figure 3 System Service Screen

2 Main Service Screen

The Main Service Screen shown in [Figure 4](#) lists the service functions available under each of the displayed categories.

Turn the rotary knob to move the scroll box to the desired category, then press the knob to confirm. The Service Log category is illustrated.

If desired, press the Standby key to return to the standby screen.



AB6020002014

Figure 4 Main Service Screen

3 Service Log

When the Service Log category is selected, the screen shown in [Figure 5](#) will appear.

Turn the rotary knob to scroll down to select the desired function, then press the knob to confirm.




Figure 5 Service Log Category Selected

To return to the main service screen, scroll to the return arrow and press the rotary knob to confirm.

3.1 Service Log Entries


When Entries is selected, the screen shown in [Figure 6](#) appears.

Use the rotary knob to scroll through the list of service log entries. The service log holds up to 256 entries. When new entries are made, the oldest entries are replaced with the new entries. The most recent entries appear upon screen call-up.



CAUTION

Selecting **Clear Service Log** and pressing the knob will bring up another screen that allows you to delete all entries. Selecting **YES** will delete all entries.

| | | | |
|---|-------|----------|--------------------|
| Service | | 11:30 | |
| Service Log | | | |
|  | | | |
| Clear Service Log | | | |
| Date | Time | Code | Description |
| 06-13-03 | 15:00 | V043 | MOTOR LIM STOP |
| 06-13-03 | 12:46 | V043 | MOTOR LIM STOP |
| 06-13-03 | 12:46 | V043 | MOTOR LIM STOP |
| 06-13-03 | 11:33 | V044 98 | MEMBRANE PRES LOW |
| 06-13-03 | 11:32 | V045 271 | MEMBRANE PRES HIGH |
| 06-13-03 | 11:32 | V044 98 | MEMBRANE PRES LOW |
| 06-13-03 | 11:31 | V044 97 | MEMBRANE PRES LOW |

AB6020002016

Figure 6 Service Log Entries Screen

Scroll to the return arrow and press the rotary knob to return to the Service Log category screen.

Tables of service log entries, error codes, and their descriptions can be found in [Table 1](#) through [Table 7](#).

Fabius Tiro Service Log Entries:

Table 1 Communication Errors

| Code | Description | Possible Cause |
|------|---------------------------|--|
| C001 | Handler Busy | Bad processor, external device incorrect operation |
| C002 | Incoming Overflow | Bad processor, bad cable or connector, incorrect settings, external device incorrect operation |
| C003 | Outgoing queue overflow | Bad processor, external device incorrect operation |
| C004 | Message not registered | External communications error |
| C005 | Character not transmitted | Bad processor, external device incorrect operation |
| C006 | Interrupt problem | Bad processor |
| C007 | Continuous interrupts | Bad processor |
| C008 | Excessive Breaks | Bad processor, incorrect settings, external device incorrect operation |

R6020002T03_Diagnostics_Service_Log.fm 17.06.04
All rights reserved. Copyright reserved.

| Code | Description | Possible Cause |
|-------------|--------------------|--|
| C009 | Excessive Framing | Bad processor, incorrect settings, external device incorrect operation |
| C010 | Excessive Parity | Bad processor, incorrect settings, external device incorrect operation |
| C011 | Excessive Overrun | Bad processor, incorrect settings, external device incorrect operation |
| C012 | Non-existent UART | Bad processor |

Table 2 Monitor Errors

| Code | Description | Possible Cause |
|-------------|--------------------|---|
| M001 | O2 <10% | O2 cell not working properly |
| M002 | O2 LOW INOP | O2 cell inoperative or disconnected |
| M003 | O2 DIFF HIGH | Difference between Cell A and Cell B >25% |
| M004 | BATTERY <10% | Battery level below 10% of capacity while operating without mains power |
| M005 | BATTERY <20% | Battery level below 20% of capacity while operating without mains power |
| FG01 | N2O CAL ERR | N2O sensor value out of range |
| FG02 | AIR CAL ERR | AIR sensor value out of range |
| FG03 | O2 CAL ERR | O2 sensor value out of range |

Table 3 Diagnostic Errors

| Code | Description | Possible Cause |
|-------------|---------------------|--|
| F001 | Diag Watchdog | Bad hardware |
| F002 | Diag System Ram | System RAM fails on pwr up, bad RAM chip, bad CPU, bad trace, etc. |
| F003 | Diag Program Memory | Bad flash download, bad flash chip |
| F004 | Diag Video Test | Video memory failed, bad video chip |
| F005 | Diag Interrupts | Bad hardware |
| F006 | Diag A/D Converter | Bad A/D hardware |
| F007 | Diag NM Ram | Corrupt entry on pwr up |
| F008 | Diag Serial Port | Failed pwr up loopback test, bad UART chip |
| F009 | Diag Clock | Bad hardware |
| F010 | Diag Speaker | Speaker not connected on pwr up |
| F011 | Diag Main Power | Mains pwr not detected on pwr up |

| Code | Description | Possible Cause |
|------|--------------|-------------------------------------|
| F012 | Diag Battery | Battery voltage on pwr up below 20% |

Table 4 System Errors

| Code | Description | Possible Cause |
|------|--|--|
| I011 | L2 SW Watchdog | Bad processor board |
| I012 | L2 Recovered, w/ Level-Two Error Count | Bad processor board |
| I021 | L3 SW Watchdog | Bad processor board |
| I022 | L3 Recovered, w/ Level-Three Error Count | Bad processor board |
| I031 | L4 SW Watchdog | Bad processor board |
| I032 | L4 Recovered, w/ Level-Four Error Count | Bad processor board |
| WD00 | Watch dog time exp | Watch dog timer fired |
| NVP1 | Pat Check Pri Fail | NVRAM Primary pattern write failed |
| NVP2 | Pat Check Sec Fail | NVRAM Secondary pattern write failed |
| NVF0 | Pri and Sec NV Files Inv | Primary and secondary NVRAM files invalid |
| NVF1 | Pri NV File Inv | Primary NVRAM file invalid |
| NVF2 | Sec NV File Inv | Secondary NVRAM file invalid |
| NVC0 | Non-Matching CRCS | CRCs for primary and secondary NVRAM did not match |
| PS01 | Charger Error | Check Battery, Connectors |
| PS02 | -2.5 V Supply Failure | Bad Processor Board |
| PS03 | +12 V Supply Failure | Bad Processor Board |
| PS04 | -5 V Supply Failure | Bad Processor Board |
| PS05 | Battery Less Than 18 V | Check Battery, Connectors |
| PS06 | +28 V Supply Failure | Check AC Power, Fuse, or 28 V Supply |
| PS07 | Temp Sensor Failure | Bad Processor Board |
| PS08 | 2.5 V Supply Failure | Bad Processor Board |
| PS09 | 3.3 V Supply Failure | Bad Processor Board |
| PS10 | +10 V Supply Failure | Bad Processor Board |
| PS11 | +16 V Supply Failure | Bad Processor Board |

Table 5 Service Log Events

| Code | Description | Possible Cause |
|------|-------------------|---|
| S006 | Service Log Clear | Log cleared by service tech, SW upgrade |
| SS90 | test message | Factory test - no failure |
| I001 | CPU Diag Recovery | CPU recovered from bad initialization |

Table 6 User Interface Events

| Code | Description | Possible Cause |
|------|-----------------------|---|
| U001 | New Device ID Created | Machine detected a new device ID was needed because of NVRAM corruption or previously did not exist |
| U002 | Option Enabled | A secure option was enabled |
| U003 | Option Disabled | A secure option was disabled |
| U007 | Model Changed to FGS | Model-type change |
| U008 | Model Changed to Tiro | Model-type change |
| U009 | Monitoring Disabled | Monitoring-mode change |
| U010 | Monitoring Enabled | Monitoring-mode change |
| U012 | SOS Recovered | A Secure Option Setting was recovered after a partial NVRAM corruption |
| U013 | FGS SOS Invalid | The FGS Secure Options Setting NVRAM file was found corrupt and no Options were recovered |
| U014 | Tiro SOS Invalid | The Tiro Secure Options Setting NVRAM file was found corrupt and no Options were recovered |

Table 7 Ventilator Errors

| Code | Description | Possible Cause |
|------|-------------------|---|
| V001 | VENT PEEP NOT CAL | Auto ventilation was attempted without valid PEEP calibration |
| V002 | ZERO POS INVALID | Crossing of light barrier at bottom vent not detected |
| V037 | EX VALVE LEAK | Greater than 15 mL of Tidal Volume during inspiration, and this Vt is greater than 10% of Vt flow during expiration |

| Code | Description | Possible Cause |
|------|----------------------------|---|
| V042 | Motor Voltage Error | Motor voltage has dropped below minimum required for auto ventilation while in Volume or PCV Ventilation mode |
| V043 | MOTOR LIM STOP | Incident of excessive motor current during auto ventilation |
| V044 | MEMBRANE PRES LOW | Membrane vacuum below normal range: possible leak |
| V045 | MEMBRANE PRES HIGH | Membrane vacuum above normal operating range |
| V049 | MOTOR OVER CURRENT | Persistent excessive motor current during auto ventilation: caused ventilator motor shut off |
| V050 | LT Pres Failed>350 | Leak in path from vent to hoses to peep valve to bag |
| V051 | LT Pres Rose>250 | Leak in path from vent to hoses to peep valve |
| V052 | LT SYS FAIL | Leak in path from vent to hoses to peep valve to bag |
| V053 | LT PC FAIL | Leak in path from vent to hoses to peep valve |
| V060 | Comp Unexp State | Compliance test entered an unexpected test state |
| V061 | Comp Bag Pres Too High | Compliance test start pressure too high |
| V062 | Comp Cant Fully Pres | Compliance test could not attain complete pressurization |
| V063 | Comp Cant Hold Pres | Compliance test could not maintain plat. pressure |
| V064 | COMP Cant DePres | Pressure not dropping during compliance test |
| V065 | COMP Inv Pres | Pressure changing in the wrong direction during compliance test |
| V066 | COMP Inv Pres | Pressure changing in the wrong direction during compliance test |
| V067 | COMP Sys < VentHose | Irregular compliance test results |
| V069 | Comp Unexpected Error | Compliance test entered an unexpected error state |
| V070 | COMP Sys Fail in Range | Unusable compliance value |
| V071 | COMP Sys Fail out of Range | Unusable compliance value |

| Code | Description | Possible Cause |
|------|---------------------------|-------------------------------------|
| V072 | COMP VH Fail in Range | Unusable vent hose compliance value |
| V073 | COMP VH Fail out of Range | Unusable vent hose compliance value |

3.2 Real Time Values

When Real Time Values is selected, the typical screen shown in [Figure 7](#) appears.

Refer to Table 8 for Real Time Value item descriptions and acceptable range data.




Figure 7 Real Time Values Screen Appearance with PCB 4116632

Refer to Table 9 for Battery/Charger Status error message list and Table 10 for PSC (Power Supply Controller) error message list as shown in [Figure 8](#).




Figure 8 Real Time Values Screen Appearance with PCB 4118079

Press the knob to return to the Service Log category screen.

Table 8 Real Time Values Screen Key and Acceptable Range Data

| Label | Description | Acceptable Range |
|-------|--|--|
| Paw | (AD) Pressure Sensor Calib. value | Max 825 (ambient) Min 650 (ambient) |
| F1 | (AD) Flow Sensor Calib. value | 0 - 4095 |
| O2 | O2 Sensor Channel (in mV) Channel 1 + 2 | OFFSET Min 39 - 9,64 mV (21%) / OFFSET Max 541 132,00 mV(100%) |
| CAL | Calibration State of O2 Sensor and Calibration value | Average stored valve between OFFSET Min and OFFSET Max |
| h | Operating Hours (full hours) | N/A |
| U+24 | Supply voltage in V | Max - 29 Volts Min - 21 Volts |
| Pvent | Ventilator Pressure | Max 825 (ambient) Min 650 (ambient) |
| PP | Pump Vacuum Pressure | 0- 120 mB (Pump OFF) 120 - 270 mB (Pump ON) |

| Label | Description | Acceptable Range |
|------------|---|---|
| Pump Speed | Pump Motor Speed in RPM | No MIN/MAX, typically 2400 to 2800 RPM when running |
| Pump Power | Percentage of full power in use (only on 4118079 PCB) | 0 to 100%, typically 6 to 15% running |

Table 9 Battery/Charger Status Error Message List

| Message | Description | Condition |
|---------------------|--|--|
| Charger OFF | The PSC has detected a +28 V supply failure or a Charge error. | Refer to Table 10 for additional conditions. |
| Slow Charging | The charger is trickle charging the battery until voltage rises to 24V. | Normal after the machine has had no AC power for an extended period of time. |
| Bulk Charging | The charger delivers the maximum charging current to the battery. | Normal after the machine has run on battery. |
| Full Charge | The charger keeps the battery voltage at a consistent charging voltage, which is above the rated voltage, for 3 hours. | Normal after the machine has run on battery and been recharged. |
| Holding Full Charge | The charger is maintaining full charge on the battery. | Normal when battery is in good condition and has been charged correctly. |
| No PSC Data | The PSC is not communicating. | If condition persists, check service log, possible replace PSC. |

Table 10 PSC (Power Supply Controller) Error Message List

| Message | Description | Condition |
|------------------------|--------------------------------|---------------------------|
| Charger Error | Charging process has timed out | Check battery, connectors |
| -2.5 V supply failure | <-2.5 or >-2.25 Volts | PCB problem |
| +12 V supply failure | <10.8 or >13.2 Volts | PCB problem |
| +5 V supply failure | <4.5 or >5.5 Volts | PCB problem |
| Battery less than 18 V | <18 Volts | Check battery, connectors |

| Message | Description | Condition |
|----------------------|----------------------|-------------------------------------|
| +28 V supply failure | <25.2 or >30.8 Volts | Check AC power, fuse or 28 V supply |
| Temp sensor failure | <.5 or >1.5 Volts | PCB problem |
| 2.5 V supply failure | <2.25 or >2.75 Volts | PCB problem |
| 3.3 V supply failure | <3.0 or >3.6 Volts | PCB problem |
| +10 V supply failure | <9.0 or >11.0 Volts | PCB problem |
| +16 V supply failure | <14.4 or >17.6 Volts | PCB problem |

4 Preventive Maintenance

When the Preventive Maintenance category is selected, the screen shown in [Figure 9](#) will appear.

Turn the rotary knob to select Activate, then press the knob to confirm. The next preventive maintenance due date can then be set using the screen as shown in [Figure 10](#).




Figure 9 Preventive Maintenance Category Selected

4.1 Activate Preventive Maintenance Date

With Activate selected, the screen will appear as shown in [Figure 10](#). Turn the rotary knob to position the box under Mth and press the knob to highlight the box. Turn the knob to select the next month preventive maintenance is due, and press the knob to enter the month.

With the box under Day, press the knob to highlight the box. Turn the knob to select the desired day, then press the knob to enter the day.

With the box under Year, press the knob to highlight the box. Turn the knob to select the year, and press the knob to enter the year.




Figure 10 PMS Due Date Screen

To return to the main service screen, scroll to the return arrow and press the knob.

5 General

When the General category is selected, the screen shown in [Figure 11](#) will appear.

Turn the rotary knob to scroll down to select the desired function, then press the knob to confirm.



Figure 11 General Category Selected

5.1 Serial Number

When Serial Number is selected, the screen will appear as shown in [Figure 12](#).

Turn the rotary knob to position the box on the first digit of the serial number, then press the knob to highlight the number. Turn the knob to display the desired number, then press the knob to enter the number. The box will then step to the next digit. Press the knob to highlight the number. Turn the knob to select the next digit of the serial number, and press the knob to enter the digit. Continue entering each digit of the serial number in the same manner. The serial number entered shall be the same as on the rear of the machine.




Figure 12 Serial Number Screen

To return to the General category screen, scroll to the return arrow and press the knob.

5.2 Reset Hours Run

When Reset Hours Run is selected, the screen will appear as shown in [Figure 13](#). Turn the rotary knob to select Reset Hours Run, and press the knob to highlight the selection. Then turn the knob to select Yes or No. Selecting Yes will reset the hours run to zero. Press the knob to enter the selection.




Figure 13 Reset Hours Run Screen

To return to the General category screen, scroll to the return arrow and press the knob.

5.3 Reset Last Service Date

When Reset Last Service Date is selected, the screen will appear as shown in [Figure 14](#). Turn the rotary knob to select Reset Last Service Date, and press the knob to highlight the selection. Then turn the knob to select Yes or No. Selecting Yes will reset the last service date to the current date. Press the knob to enter the selection.




Figure 14 Reset Last Service Date Screen

To return to the General category screen, scroll to the return arrow and press the knob.

5.4 Select Market Kit

When Select Market Kit is selected, the screen will appear as shown in [Figure 15](#) or [Figure 16](#). Turn the rotary knob to scroll to Select Market Kit, and press the knob to highlight the selection. Then turn the knob to select U.S. or Non-U.S.

When U.S. is selected, the only language supported is English and the pressure unit of measure is set to cmH₂O. The options to choose a different language or pressure unit are eliminated while in the "U.S." mode.

Selecting Non-U.S. mode enables the option to select language and pressure units of measure from the standby menu.




Figure 15 Select Market Kit Screen (US)





Figure 16 Select Market Kit Screen (Non-US)

To return to the General category screen, scroll to the return arrow and press the knob.

6 Calibration

When the Calibration category is selected, the screen shown in [Figure 17](#) will appear.

Turn the rotary knob to scroll down to select the desired function, then press the knob to confirm.



AB6020002027

Figure 17 Calibration Category Selected

6.1 Fresh Gas Flow Calibration

When Fresh Gas Flow is selected, the screen shown in [Figure 18](#) appears.

Follow the on-screen instructions to perform the calibration. The Store Zero and Exit keys are located just below their legends at the bottom of the screen.




Figure 18 Fresh Gas Flow Calibration Screen

Press the Exit key to return to the Calibration category screen, or the Standby key to return to the Standby screen.

6.2 Pressure Calibration

When Pressure is selected, the screen shown in [Figure 19](#) appears.

Follow the on-screen instructions to perform the calibration. The Store Zero and Exit keys are located just below their legends at the bottom of the screen.




Figure 19 Pressure Calibration Screen

Press the Exit key to return to the Calibration category screen, or the Standby key to return to the Standby screen.

6.3 O2 Zero Calibration

When O2 is selected, the screen shown in [Figure 20](#) appears.

Follow the on-screen instructions to perform the calibration. The Store Zero and Exit keys are located just below their legends at the bottom of the screen.

If the calibration process fails, a "value not stored" message is displayed in place of the Stored zero cell A or cell B values.




Figure 20 O2 Zero Calibration Screen

Press the Exit key to return to the Calibration category screen, or the Standby key to return to the Standby screen.

6.4 PEEP Valve Calibration



NOTE

Before performing the PEEP valve calibration, a valid pressure calibration must be performed. Otherwise, an 'Inv pres.cal' message will appear.

Follow the on-screen instructions to perform the calibration. During calibration, an "In Progress" message is displayed. When the procedure is complete a PASS or FAIL message is displayed. The Calibrate and Exit keys are located just below their legends at the bottom of the screen.




Figure 21 PEEP Calibration Screen

Press the Exit key to return to the Calibration category screen, or the Standby key to return to the Standby screen.

7 Configure

When the Configure category is selected, the screen shown in [Figure 22](#) will appear.

Turn the rotary knob to scroll down to select the desired function, then press the knob to confirm.




Figure 22 Configure Category Selected

7.1 System Settings (US and Non-US)

When the System Settings category is selected, the screen shown in [Figure 23](#) or [Figure 24](#) will appear.

Turn the rotary knob to scroll down to select the desired market kit, then press the knob to confirm.




Figure 23 System Settings Screen

7.1.1 Model Type

When the Model Type category is selected, the screen shown in [Figure 24](#) will appear.

Turn the rotary knob to scroll down to select the model type 'Tiro', then press the knob to confirm.



NOTE

Machines from the factory are configured with the appropriate model type. Other than performing a software download or PCB replacement, the model type should not be changed. Changing the model type will disable some ventilator options previously configured and will require reconfiguration.

To reconfigure the unit, contact DrägerService - Technical Support: Phone 1-800-4-Dräger, Phone: 215-721-5402, or e-mail to techsupport@draegermed.com and provide the following information to receive the necessary release codes:

- Machine Type (Fabius GS or Fabius Tiro)
- Feature Description and Part Number
- Machine Serial Number
- Device ID Number




Figure 24 Typical Model Type Screen

7.2 Standard Options

When the Standard Options category is selected, the screen shown in [Figure 25](#) will appear.

Turn the rotary knob to scroll down to select the desired function, then press the knob to confirm.




Figure 25 Standard Options Screen

To return to the Configure screen, scroll to the return arrow and press the knob.

7.2.1 Flowmeter

When Flowmeter is selected, the screen shown in one of the following illustrations will appear ([Figure 26](#), [Figure 27](#), [Figure 28](#) or [Figure 29](#)). Turn the knob to scroll down to the desired selection, and press the knob to enter that function.




Figure 26 Flowmeter Selection Screen

To return to the Standard Options category screen, scroll to the return arrow and press the knob.

7.2.1.1 O2 Position (Virtual Flowtubes)

When O2 Position for virtual flowtubes is selected, the screen will appear as shown in [Figure 27](#). Turn the rotary knob to select O2 Position, and press the knob to highlight the choices. Turn the knob to select either LEFT or RIGHT, then press the knob to enter that selection.

By entering the desired position (left or right), the virtual O2 flow position is relative to the air virtual flow tube. The air virtual flow tube is always located in the center position.




Figure 27 O2 Flowmeter Position Screen (Virtual Flowtubes)

To return to the Flowmeter selection screen, scroll to the return arrow and press the knob.

7.2.1.2 Gas Selection

When Gas Selection is chosen, the screen will appear as shown in [Figure 28](#). Turn the rotary knob to choose Gas Selection, and press the knob to highlight the choices. Turn the rotary knob to select either 2 or 3, then press the knob to enter that selection.

The gas selection number 2 or 3 shall match the total gas options pneumatically configured on the unit. Otherwise a 'ERR' will occur on the fresh gas digital display, where applicable.




Figure 28 Gas Selection Screen

To return to the Flowmeter selection screen, scroll to the return arrow and press the knob.

7.2.1.3 Flowtube Resolution

When Flowtube Resolution is chosen, the screen will appear as shown in [Figure 29](#). Turn the rotary knob to choose Flowtube Resolution, and press the knob to highlight the choices. Turn the rotary knob to select either Change State or High Res, then press the knob to enter that selection.

High resolution provides a 3 digit resolution to the fresh gas digital flow displays - Low resolution provides 2 digit resolution.




Figure 29 Flowtube Resolution Screen

To return to the Flowmeter selection screen, scroll to the return arrow and press the knob.

7.2.2 O2 Whistle

When O2 Whistle is chosen, the screen will appear as shown in [Figure 30](#). Turn the rotary knob to select Change State, and press the knob to highlight the choices. Then turn the rotary knob to select either ENABLE or DISABLE, and press the knob to enter that selection.

By selecting Enable, the unit provides an electronic whistle which activates when Lo O2 supply is detected.




Figure 30 O2 Whistle Selection

To return to the Configure category screen, scroll to the return arrow and press the knob.

7.2.3 Alarms

When Alarms is selected, the screen shown in [Figure 31](#) will appear. Turn the knob to scroll down to the desired selection, and press the knob to enter that function.




Figure 31 Alarms Selection Screen

To return to the Configure category screen, scroll to the return arrow and press the knob.

7.2.3.1 No Fresh Gas

When No Fresh Gas is chosen, the screen will appear as shown in [Figure 32](#). Turn the rotary knob to select Change State, and press the knob to highlight the choices. Then turn the rotary knob to select either ENABLE or DISABLE, and press the knob to enter that selection.

Enabling No Fresh Gas provides a message to the operator when the fresh gas sensor detects no flow while in an active ventilator mode.




Figure 32 No Fresh Gas Selection Screen

To return to the Alarms category screen, scroll to the return arrow and press the knob.

7.2.3.2 Fresh Gas Low Alarm

When Fresh Gas Alarm is chosen, the screen will appear as shown in [Figure 33](#). Turn the rotary knob to select Change State, and press the knob to highlight the choices. Then turn the rotary knob to select either ENABLE or DISABLE, and press the knob to enter that selection.

An Enabled Fresh Gas Lo alarm provides a message to the operator when the available fresh gas is set such that -8 mbar auxiliary air valve opens to maintain the present Vt delivery.




Figure 33 Fresh Gas Alarm Selection Screen

To return to the Alarms category screen, scroll to the return arrow and press the knob.

7.2.3.3 Threshold Low Alarm

When Threshold Low Alarm is chosen, the screen will appear as shown in [Figure 34](#). Turn the rotary knob to select Change State, and press the knob to highlight the choices. Then turn the rotary knob to select either ENABLE or DISABLE, and press the knob to enter that selection.

A Threshold Low alarm that is enabled provides a corresponding alarm when the current threshold is set too low relative to the peak pressure.




Figure 34 Threshold Low Alarm Selection Screen

To return to the Alarms category screen, scroll to the return arrow and press the knob.

7.2.4 Pressure

When Flowmeter is selected, the screen shown in [Figure 35](#) will appear. Turn the knob to scroll down to the desired selection, and press the knob to enter that function.




Figure 35 Pressure Screen

To return to the Configure category screen, scroll to the return arrow and press the knob.

7.2.4.1 Ambient Pressure

When Ambient Pressure is chosen, the screen will appear as shown in [Figure 36](#). Turn the rotary knob to select Adjust, and press the knob. Then turn the knob to set the correct ambient pressure in mbar for the location the unit will reside, and press the knob to enter the value.

Using the Site Configurations table found in the PMS section, set the ambient pressure (mbar) in accordance with the local elevation. Refer to the Site Configurations portion of the PMS section for additional information.



NOTE

It is imperative the ambient pressure setting is properly set. Otherwise, erroneous Vt settings could result.

Changes in barometric pressure resulting from local weather conditions have minimal to no affect on Vt readings. Therefore, it is not necessary to change the ambient pressure setting with changing local weather conditions.




Figure 36 Ambient Pressure Screen

To return to the Pressure category screen, scroll to the return arrow and press the knob.

7.2.4.2 Plateau-Mean Display Screen

The breathing pressure monitor window can display Mean airway pressure or Plateau airway pressure.

When Plateau is chosen, the screen will appear as shown in [Figure 37](#). Turn the rotary knob to make the desired selection, and press the knob.




Figure 37 Plateau-Mean Display Screen

To return to the Pressure category screen, scroll to the return arrow and press the knob.

7.3 Secure Options

When the Secure Options category is selected, the screen shown in [Figure 38](#) will appear.

Turn the rotary knob to scroll down to select the desired function, then press the knob to confirm.

The Secure Options screen provides the DrägerService or Authorized Service Organization representative to enable or disable certain options.




Figure 38 Secure Options - Enabled

Inquiries regarding Secure Options or Ventilator Options shall be directed to: DrägerService - Technical Support, Phone: 1-800-4-Dräger, Phone: 215-721-5402, or e-mail to techsupport@draegermed.com.

8 Serial Port

When the Serial Port category is selected, the screen shown in [Figure 39](#) will appear.

Turn the rotary knob to select Parameters, then press the knob to confirm. The serial port communication parameters can then be set using the screen as shown in [Figure 40](#).




Figure 39 Serial Port Category Selected

8.1 Serial Port Parameters

With Parameters selected, the screen will appear as shown in [Figure 40](#).

Turn the rotary knob to select the parameter you want to set, and press the knob to highlight that parameter. Then turn the knob to select the desired value for that parameter. Press the knob again to enter the value. Scroll to the next parameter and repeat the process to enter a value.

The [Table 11](#) lists the values available for each parameter, along with the factory default values.




Figure 40 Serial Port Parameters Screen

Table 11 Values for Parameters

| Item | Selections | Factory Default |
|-----------|--------------------------------------|-----------------|
| BAUD RATE | 1200, 2400, 4800, 9600, 19200, 38400 | 9600 |
| PARITY | NONE, ODD, EVEN | NONE |
| STOP BITS | 1 or 2 | 1 |
| DATA BITS | 7, 8 | 8 |
| PROTOCOL | VITALINK or MEDIBUS | VITALINK |



NOTE

It is important to ensure that communication protocols selected on each host and external device are correct. Vitalink and Medibus protocols are similar and if not set identically on each device, inaccurate data may be displayed on the remote device.

To return to the main service screen, scroll to the return arrow and press the knob.

9 Pump On/Off

The pump on/off softkey provides the ability to toggle the pump ON or OFF as part of the troubleshooting and diagnostics process.

1. Access the 'Real Time Values' screen. Refer to [Real Time Values Screen Appearance with PCB 4118079](#).
2. Press the 'PUMP ON/OFF' softkey to toggle the pump between the On and Off state.



NOTE

The pump label located at the bottom left corner of the 'Real Time Values' screen will be highlighted in reverse when the pump has been activated. When first entering the screen, the initial state of the pump is On.

REPLACEMENT PROCEDURES

Replacement Procedures

This section outlines removal and replacement procedures for the field-replaceable assemblies of the Fabius Tiro Anesthesia System.

These procedures are to be performed only by a qualified Technical Service Representative (TSR).

The following are the only procedures authorized by DrägerService to be performed in the field. All other service procedures shall be referred to Draeger Medical, Inc., DrägerService.



NOTE

The PMS PROCEDURE must be performed after any replacement, removal, calibration or adjustment procedure.

1 Core Module Inversion

The following Core Module Inversion instructions is referenced throughout many of the replacement procedures found within the documentation for the Fabius Tiro. Its purpose is to provide general instructions to safely invert the unit upside down to perform certain service.

1.1 Core Module Inversion Disassembly

1. Disconnect all pipeline hoses and set the System Power switch to 'ON'.
2. Close all cylinder valves except the O2 valve.
3. Set the oxygen flow to 5 l/min.
4. Open the other gas flow control valves to deplete the pressure from the system.
5. Close the O2 cylinder valve and close the flow control valves. Press the O2 Flush valve to deplete oxygen pressure from the system.
6. Set the System Power switch to 'OFF' and disconnect AC power from the machine.
7. Detach the cylinder supply from the yoke.
8. Remove any third party monitoring devices and accessories from the unit.
9. Remove the vaporizer (refer to section [Vaporizers](#) for details).
10. Disconnect the necessary hoses and cables from the breathing system and remove the breathing system from the machine, including the GCX mount.
11. Loosen the set screws located in the rear GCX rails that secure the unit to its base. Refer to [Figure 1](#).
12. Wall Mount (If Applicable)
 - a) Remove the two swivel mount screws that secure the front of core module to the swivel mount. Refer to [Figure 2](#).



WARNING

Ensure all safety provisions have been made prior to performing the next step. Failure to do so may result in unit damage or personal injury.

13. Make necessary provisions to prevent opening of the ventilator (door) during step [14](#).
14. Using at least two people, invert the unit upside down on a table.
15. Perform the necessary service.

1.2 Core Module Inversion Reassembly

1. Ensure all necessary service is complete and unit is prepared for inversion to its upright position.
2. Using at least two people, invert unit to its upright position back onto its base, ensuring the unit sits properly on the rear mounting posts.
3. Tighten the GCX rail set screws previously loosened. Refer to [Figure 1](#).
4. Wall Unit (If Applicable)
 - a) Reinstall two swivel mount set screws previously removed.
5. Reinstall the breathing system and reconnect the previously removed hoses and cables from the breathing system to the unit.

6. Reinstall the vaporizer (refer to the [Vaporizers](#) section for details).
7. Reinstall any previously removed third party monitoring devices and accessories to the unit.
8. Reattach the cylinder supply to the yoke.
9. Reconnect all pipeline hoses and restore AC power to the machine.
10. Set the System Power switch to 'ON'.

SV88063




Figure 1 Core Module Inversion GCX Rail Location

Replacement Procedure

Fabius Tiro




Figure 2 Wall Mount Service Access

**2 Wall Mount Service
Access Procedure (If
Applicable)**

1. Remove the Vaporizer (Refer to the [Vaporizers](#) Section).
2. Make the necessary provisions to all monitoring devices, breathing system, and accessories on the Core Unit to ensure the unit can pivot up to 180°.
3. Remove the unit pivot knob and washer located under the wall mount bracket, then using an allen wrench, remove the service access screw. Refer to [Figure 2](#).
4. Remove the cylinder wrench from the holder, then carefully swivel the core unit to gain access to the rear of the unit.
5. Perform necessary service.
6. Pivot core module back into its original position and reinstall the service access screw, pivot knob, and washer.
7. Reinstall all monitoring devices, breathing system vaporizer or accessories to the core unit that were previously removed, then reinstall the cylinder wrench in its holder.

3 Cylinder Yokes (Fixed) and Regulators

Replacement of a yoke, check valve, or cylinder regulator requires that the yoke be removed from the anesthesia machine. [Figure 3](#) shows the cylinder yokes mounting arrangement. Access to the yoke mounting screws and tubing connections is at the underside of the table top at the back of the machine.

1. Invert unit by performing steps given in section [Core Module Inversion Disassembly](#).



WARNING

Ensure cylinder is stored in a safe place and is laid on its side.

2. Remove the two yoke mounting screws.
3. Hold the check valve assembly with a wrench, and carefully turn the yoke counter-clockwise to remove it.
4. Where applicable, remove the regulator by disconnecting the copper tubing at the HP side of the regulator, and the plastic tubing at the LP side of the regulator. The regulator assembly, with check valve attached, can then be removed from the machine.



NOTE

If a yoke is being replaced, verify that the pin indexing arrangement and the label are in agreement with the gas designation stamped on the mounting surface of the yoke. Refer to the parts list.



NOTE

Where a regulator is replaced, assemble the fittings and check valve to the regulator with Loctite #271 (red).

5. Install the regulator assembly in the machine. Hold the check valve assembly with a wrench and carefully tighten the yoke to the correct horizontal position.
6. Reattach the yoke to the machine with the hardware previously removed.
7. If a new cylinder is being installed, remove the old sealing washer from the gas inlet of the yoke and install a new washer.
8. Install the correct cylinder on the yoke, making sure that the index pins are properly engaged before tightening the handle bolt. The cylinder should hang vertically after the handle is tight.
9. Perform the following leak test on the yoke assembly:



NOTE

The cylinder used for this test must contain the following minimum pressure:

- a) Open the cylinder valve and check for a pressure indication on the corresponding cylinder pressure gauge.
 - O2: 1000 psi
 - N2O : 700 psi
 - b) Close the cylinder valve and remove the cylinder from the yoke.
 - c) For any gas, the pressure should not drop more than 50 psi in two minutes.
10. Invert unit into its upright position by performing steps given in the [Core Module Inversion Reassembly](#) Section.
 11. Perform the PMS Procedure.

Replacement Procedure

Fabius Tiro




Figure 3 Cylinder Yoke and Regulator Assemblies

R6020002T03_Replacement_Cylinder_Yokes.fm 17.06.04
All rights reserved. Copyright reserved.

4 Cylinder Pressure Gauges

Each pressure gauge is secured by two 6-32 kep nuts, as shown in [Figure 4](#). Access to gauge mounting nuts and tubing connections is at the underside of the machine.

1. Invert unit by performing steps given in section [Core Module Inversion Disassembly](#).



WARNING

Ensure cylinder is stored in a safe and is laid on its side.


2. Disconnect the tubing at the back of the gauge.
3. Remove the two nuts that secure the gauge mounting plate to the Core Unit, then remove bracket and gauge(s).
4. Remove the two kep nuts securing the gauge to its mounting plate.
5. Install the replacement gauge(s) in their proper location. Reinstall the bracket and hardware previously removed, ensuring the gauges are straight.
6. Reconnect the tubing.
7. Reinstall the cylinders.
8. Perform the following leak test:



NOTE

The cylinder used for this test must contain the following minimum pressure:

- a) Open the cylinder valve and check for a pressure indication on the corresponding gauge.
 - O₂: 1000 Psi
 - N₂O : 700 Psi
 - b) Close the cylinder valve and remove the cylinder from the yoke.
 - c) For any gas, the pressure should not drop more than 50 psi in two minutes.
9. Invert unit into its upright position by performing steps given in section [Core Module Inversion Reassembly](#).
 10. Perform the PMS Procedure.



SV88074

Figure 4 Cylinder Pressure Gauges

5 Auxiliary Oxygen Flow Meter

The auxiliary oxygen flowmeter is attached to the outside of the ventilator door on the left side of the machine, by a stud and nut arrangement - accessible by opening the door. A flexible O₂ supply tube from the flowmeter connects to a T-fitting below the floor of the ventilator housing. [Figure 5](#) shows the mounting and tubing arrangement.

1. Disconnect all pipeline hoses and close all cylinder valves.
2. Press the O₂ Flush button to drain oxygen pressure from the system.
3. Set the System Power switch to OFF and disconnect AC power from the machine.
4. Open the ventilator door.
5. Cut the tie strap on the flexible tube at the T-fitting, and disconnect the tube.
6. Cut the tie strap on the flexible tube just below the floor of the ventilator housing.
7. Remove the screw nuts securing the auxiliary O₂ flowmeter to the door, and remove the flowmeter.
8. Position the replacement flowmeter on the ventilator door (feed the flex tubing through the clearance hole) and secure the auxiliary O₂ flowmeter with the two nuts that were previously removed.
9. Connect the flex tubing to the T-fitting and secure it with a tie strap.
10. Install a tie strap on the flex tube just below the floor of the ventilator housing.
11. Close and secure the ventilator door.
12. Connect the pipeline hoses and restore AC power to the machine.
13. Perform the PMS Procedure.

SV00475




Figure 5 Auxiliary O2 Flowmeter

6 Vaporizers

The vaporizer mounting is shown in [Figure 6](#). Before removing a vaporizer from the machine for shipment, it must be completely drained and dried in accordance with the procedure given below. Be sure to have a suitable packing or storage container available in which to place the vaporizer.



CAUTION

The following steps must be performed in the sequence given.



WARNING

Do not inhale anesthetic vapors as this could result in personal injury.



WARNING

This procedure must be performed in a well ventilated area and without any other personnel present.



WARNING

For Vapor 19.x series, do not tilt a vaporizer that contains anesthetic agent more than 45 degrees. Failure to observe this precaution will render the handwheel calibration invalid.

1. Set the Main Power switch to ON.
2. Set all vaporizer handwheels to their Zero or OFF position.
3. Remove the filler and drain plugs, and drain the vaporizer into a suitable container. Dispose of the residual agent in an approved manner.
4. Make the necessary machine adjustments to direct gas flow through the breathing system to the scavenger.
5. Verify AGS scavenger is connected to vacuum manifold and AGS is active.
6. Turn the vaporizer handwheel to the maximum concentration setting.
7. Set the oxygen flow to 10 L/min. for at least 20 minutes.
8. Turn the vaporizer handwheel to 0 (zero), and replace the filler and drain plugs.
9. Turn the oxygen flow off, and set the Main Power switch to OFF.



NOTE

Should a 19.x series vaporizer containing anesthetic agent be accidentally tilted more than 45 degrees, it must be drained and flushed in accordance with instructions given in the manual supplied with the vaporizer.

10. Un-lock the vaporizer and carefully lift it from the manifold.
11. Place the vaporizer in a suitable container for transport or storage.
12. Set the handwheel on the replacement vaporizer to its Zero position.

13. Inspect the O-rings on the manifold and replace them if needed.
14. Install the replacement vaporizer on the manifold and turn the locking lever to secure the vaporizer to the manifold.
15. Perform the PMS Procedure.




Figure 6 Vaporizer Installation

7 O2 Flush Valve

The O2 flush valve is accessible from the underside of the table top. Valve removal requires disconnecting the plastic tubing and removal of the fittings. [Figure 7](#) shows mounting and tubing arrangement for the valve.

1. Invert unit by performing steps given in the [Core Module Inversion Disassembly](#) Section.
2. Disconnect the tubing from the fittings on the O2 flush valve.
3. Remove the fittings from the outlet of the O2 flush valve.
4. Un-screw the plastic retaining nut from the valve body, and withdraw the valve from the front of the unit.
5. Install the replacement valve in the unit, and secure it with the plastic retaining nut.
6. Install fittings from the removed O2 flush valve onto the replacement O2 flush valve.
7. Reconnect the tubing to the fittings as shown in [Figure 7](#).
8. Install the appropriate label on the replacement valve.
9. Invert unit into its upright position by performing steps given in the [Core Module Inversion Reassembly](#) Section.
10. Perform the PMS Procedure.

SV88006




Figure 7 O2 Flush Valve

R6020002T07_Replacement_O2_Flush_Valve.fm 17.06.04
All rights reserved. Copyright reserved.

8 Caster

The casters are mounted into the frame rails on the trolley via a nut and stem arrangement. A typical arrangement is shown in [Figure 8](#). Caster replacement requires that the machine be tilted to provide enough clearance for the caster stem to be withdrawn from the frame rail.




WARNING

Do not tilt the machine more than 10 degrees or raise the casters more than 3½ inches from the floor. Failure to observe this precaution may result in a tip-over, causing personal injury. Vaporizers containing anesthetic agent may also be damaged.

1. Obtain a brace capable of supporting one side of the machine with its casters two to three inches from the floor.
2. Remove all unsecured equipment and accessories from the machine.
3. Lock the front casters.
4. Using at least two people, tilt the machine until the casters on one side are raised two to three inches from the floor, and position the support brace under the trolley rail between the front and back casters.
5. Remove the appropriate end cap from the frame end to expose the casters' nut and stem.
6. Using the appropriate box wrench to secure the nut located inside the frame, unscrew the caster from the frame using a caster wrench, P/N S010055.
7. Apply a small amount of Loctite (#271 - Red) to the threaded stem of the replacement caster.
8. Install the replacement caster onto the frame and reinstall the hardware previously removed, then tighten.
9. Using at least two people, tilt the machine, remove the support brace and carefully lower the machine to the floor.
10. Check for proper operation of the caster and ensure that the front casters lock properly.
11. Perform the PMS procedure, including a vaporizer calibration verification.
12. Reinstall any unsecured equipment and accessories that were previously removed.

Fabius Tiro

Replacement Procedure



*R = LOCTITE #271 RED

Figure 8 Caster Replacement

R6020002108_Replacement_Caster.fm 17.06.04
All rights reserved. Copyright reserved.

9 Ventilator

The ventilator and door assembly are retained in the monitor housing by a hinge rod. [Figure 9](#) and [Figure 10](#) show the hinge arrangement and the disconnects needed to remove the ventilator from the machine. [Figure 14](#) shows the lower cylinder and motor arrangement. While it is possible to service some components without removing the ventilator assembly from the machine, the complete disassembly sequence is given below.

Ventilator Removal



NOTE

If servicing a wall mount unit, refer to [Wall Mount Service Access Procedure](#) prior to performing this procedure.

1. Disconnect all pipeline hoses and set the System Power switch to 'ON'.
2. Close all cylinder valves except the O2 valve.
3. Set the oxygen flow to 5 l/min.
4. Open the other gas flow control valves to deplete the pressure from the system.
5. Close the O2 cylinder valve and close the flow control valves. Press the O2 Flush valve to deplete oxygen pressure from the system.
6. Set the System Power switch to 'OFF' and disconnect AC power from the machine.
7. Detach the cylinder supply from the yoke.
8. Disconnect the ventilator hose at the side of the machine.
9. Remove the back panel to expose the controller PCB assembly. Provide support for the panel to prevent damage to the power supply cable connections.



NOTE



Make note of the positioning and routing of all cables and pneumatic lines for reassembly.

10. Swing open the ventilator door. Disconnect the ventilator cables from J11 and J12 on the controller PCB, and separate the flying connector on the cable attached to J10 on the controller PCB. Feed these three cables forward into the ventilator compartment.
11. Disconnect the ventilator chamber pressure sensor hose at patient assembly.
12. Disconnect the ventilator tubing from the T-fitting (see [Figure 9](#)) and feed this tube into the ventilator compartment.
13. If applicable, disconnect the O2 supply tube from the auxiliary O2 flow meter.
14. Disconnect the ventilator ground wires from the stud on the ventilator door.
15. Remove the E-ring from the hinge rod. While supporting the ventilator assembly by the top handle, slide the hinge rod down and out through the underside of the table:

16. Remove the ventilator from the door by removing Screw **A** (see [Figure 10](#)) and sliding out the upper and lower hinge bushings.

Replacement Procedure

Fabius Tiro




SV00523

Figure 9 Ventilator Replacement

R6020002T09_Replacement_Ventilator.fm 17.06.04
All rights reserved. Copyright reserved.

Fabius Tiro

Replacement Procedure



*P = LOCTITE # 222 (PURPLE)

Figure 10 Ventilator Assembly Details

R6020002T109_Replacement_Ventilator.fm 17.06.04
 All rights reserved. Copyright reserved.

Replacement Procedure

Fabius Tiro

Ventilator Servicing: Upper Diaphragm (patient side)

17. Release the three clamps holding the patient assembly to the upper cylinder; lift the patient assembly up and out.
18. Lift the diaphragm out of the upper cylinder as shown in [Figure 11](#) in the direction of the arrow.





Figure 11 Upper Diaphragm Removal

Lower Diaphragm: (motor & case assembly)

19. Lower the replacement diaphragm into the cylinder. Ensure that the lip of the diaphragm fits properly into the top groove on the cylinder.
20. Carefully reinstall the patient assembly, making sure the diaphragm seats properly. Secure the three clamps holding the patient assembly to the upper cylinder.
21. Noting the position of the clamp for reassembly, remove the bolt securing the band clamp holding the upper and lower cylinders together. Open the band clamp and move it down onto the lower cylinder. Separate the upper cylinder from the lower cylinder.
22. Using the proper 'Pozi' driver, remove the three screws from the top of the piston. See [Figure 14](#).
23. Lift the piston and diaphragm up and out of the lower cylinder.
24. Remove the O-ring off the piston.
25. Place the piston upside down on a table. Remove the retaining ring (cir-clip) securing the diaphragm at the inside of the piston, and pull the diaphragm up off the piston. See [Figure 12](#).



SV00545

Figure 12 Lower Diaphragm Removal

26. Install the replacement diaphragm down over the piston; fold the lower lip of the diaphragm (now facing up) over the piston and reinstall the retaining ring in the inner groove of the diaphragm inside the piston. See [Figure 12](#) and [Figure 14](#).



NOTE

The gap in the retaining ring (circlip) shall be aligned opposite the mark on top of the piston as shown in [Figure 14](#). This note does not apply to later units where the mark on the piston does not exist.




Figure 13 Lower Diaphragm Installation

27. Reinstall the O-ring in the groove on the outside of the diaphragm.
28. Slightly roll the diaphragm down the piston until the diaphragm (overlaps) covers the O-ring evenly around the bottom.
29. Install the piston over the motor and reinstall the three screws that were previously removed.



NOTE

The lip on the lower diaphragm is captured between the upper and lower cylinders, and must be correctly installed.

30. Place the upper cylinder on the lip of the diaphragm, and carefully fit the lip into the groove in the upper cylinder.
31. Carefully press the upper cylinder down until it is resting on the lower cylinder. Note that the cylinders are keyed to fit together in only one position. Look down through the upper cylinder and verify that the lower diaphragm shows a uniform roll around its circumference.




Figure 14 Lower Cylinder and Motor Assembly Details

32. While holding the cylinders together, move the band clamp into position and reinstall its bolt to join the cylinders. Make sure the band clamp is positioned as shown in [Figure 10](#) or the ventilator door may not close properly.

Replacement Procedure

Fabius Tiro

Motor/Ballscrew Assembly



NOTE

Upper and lower cylinders must be separated, and piston must be removed.

33. At the bottom of the lower cylinder, remove the ground plate, and remove the center cap.
34. Disconnect and remove the tach cable and positioning sensor cables. Note the arrangement of the cables so they can be reinstalled in the same manner. (The motor cable remains with the motor/ballscrew assembly.)
35. Remove the four screws securing the motor/ballscrew assembly to the lower cylinder. Note the orientation of the unit so the replacement assembly can be installed in the same manner.
36. Remove existing standoffs and shocks from removed motor/ballscrew assembly then reinstall onto the replacement motor/ballscrew assembly.
37. Install the replacement motor/ballscrew assembly with the hardware that was previously removed.
38. Reinstall the tach and sensor cables. Route all cables in the original manner.
39. Reinstall the ground plate and the center cap.
40. Reinstall the piston.



NOTE

The lip on the lower diaphragm is captured between the upper and lower cylinders, and must be correctly installed.

41. Place the upper cylinder on the lip of the diaphragm, and carefully fit the lip into the groove in the upper cylinder.
42. Carefully press the upper cylinder down until it is resting on the lower cylinder. Note that the cylinders are keyed to fit together in only one position. Look down through the upper cylinder and verify that the lower diaphragm shows a uniform roll around its circumference.
43. While holding the cylinders together, move the band clamp into position and reinstall its bolt to join the cylinders. Make sure the band clamp is positioned as shown in the illustration, or the ventilator door may not close properly.

Ventilator Reassembly

44. Position the ventilator against the door; slide the upper and lower hinge bushings in place, and reinstall Screw **A** (Figure 10) to attach the ventilator to the door.
45. Locate the hinge rod that was previously removed. Orient it with the E-ring groove toward its upper end. Hold the ventilator assembly by its top handle and position its lower hinge in the monitor housing (make sure the spacers are in place). Sight down through the hinge bushings to align the lower hinge bushing with the hole in the monitor housing. Slide the hinge

rod up through the table and up through the lower hinge bushing. Tilt the ventilator assembly into position and continue sliding the rod up and into the upper hole in the monitor housing until the E-ring groove on the rod is visible just above the upper hinge bushing. Reinstall the E-ring on the hinge rod.

46. Feed the three ventilator cables through the opening in the back wall of the ventilator compartment. Reconnect these to J11 and J12 on the controller PCB, and reconnect the remaining one to the flying connector on the cable attached to J10 on the controller PCB.
47. Feed the tubing from the ventilator patient assembly through the opening in the back wall of the ventilator compartment and reconnect it to the T-fitting.
48. Reconnect the vent pressure chamber pressure hose.
49. Re-attach the ground wires to the ventilator door.
50. If applicable, reconnect the O2 supply tube to the auxiliary O2 flow meter.
51. Reinstall the back panel of the controller PCB assembly.
52. Close the ventilator door and restore all external hose connections.
53. Restore AC power to the machine.
54. Perform the PMS procedure.

Patient Assembly

In case of repair, replace the high pressure valve, the negative pressure relief valve, or the complete patient system. See [Figure 10](#).

55. Disconnect the ventilator hose from the side of the machine.
56. Swing open the ventilator door.
57. Disconnect the ventilator chamber sensor hose.
58. Release the three clamps holding the patient assembly to the upper cylinder, and lift the patient assembly up and out.

High Pressure Valve

59. Unscrew the high pressure valve using plastic-jaw pliers.
60. Remove all residual thread sealant from the high pressure valve mount on the patient system.
61. Dispose of the old high pressure valve.
62. Apply an ample amount of thread sealant (DI-Wacker Silicon Rubber P/N 1202537) to threads of the new high pressure valve, and screw the high pressure valve into the patient system using plastic-jaw pliers.

Negative Pressure Relief Valve

63. Unscrew the negative pressure relief valve using plastic-jaw pliers.
64. Remove all residual thread sealant from the negative pressure relief valve mount on the patient system.
65. Dispose of the old negative pressure relief valve.
66. Apply an ample amount of thread sealant (DI-Wacker Silicone Rubber P/N 1202537) to threads of the new negative pressure relief valve, and screw the negative pressure relief valve into the patient system using plastic-jaw pliers.
67. Carefully reinstall the patient system, making sure the diaphragm seats properly.

Replacement Procedure

Fabius Tiro

68. Secure the three clamps holding the patient assembly to the upper cylinder.
69. Reconnect the ventilator chamber sensor hose.
70. Reconnect the ventilator hose to the patient system.
71. Close the ventilator door.
72. Reconnect the pipeline and cylinder supplies and restore power to the machine.
73. Perform the PMS procedure.

10 Gas Inlet Assembly (including O2 supply pressure switch)

The gas inlet assembly is attached to the side of the gas box assembly as shown in [Figure 15](#). This assembly includes regulators, check valves, filters, and the O2 supply pressure alarm switch. Refer to the parts list for details on internal parts locations for replacement.

1. Invert unit by performing steps given in the [Core Module Inversion Disassembly](#) Section.
2. Disconnect the wire harness from the O2 supply pressure alarm switch.



NOTE

Make a note of the connections and their destinations in the following step for later reassembly.

3. Disconnect the following tubing:
 - P1 (N2O pipeline pressure gauge)
 - P2 (Air pipeline pressure gauge)
 - P3 (O2 pipeline pressure gauge)
 - P4 (N2O cylinder regulator output check valve)
 - P5 (N2O flow control valve)
 - P6 (Air flow control valve)
 - P7 (O2 flow control valve)
 - P8 (O2 T-fitting)
4. Remove the screws securing the gas inlet assembly to the side panel and carefully lift the gas inlet assembly out.
5. Install the replacement assembly on the side panel in the same manner as the original, ensure there is no pinching or kinking of the attached connections.
6. Reconnect the tubing at P1 thru P8.
7. Reconnect the wire harness to the O2 supply pressure alarm switch.
8. Invert unit into its upright position by performing steps given in the [Core Module Inversion Reassembly](#) Section.
9. Perform the PMS procedure.

SV88009




Figure 15 Gas Inlet Assembly

11 Battery

The two 12 V rechargeable batteries are located in the vent & monitor assembly behind the front bezel. The batteries and their connections are accessible through the ventilator compartment as shown in [Figure 16](#).

1. Set the main power switch to OFF and disconnect AC power from the machine.
2. Remove the battery fuse from its holder (located on the back of the machine next to the main power switch).
3. Swing open the ventilator door.
4. Loosen the retainer screw securing the battery retainer bracket. Remove the bracket, and slide the batteries down to the lower level of the ventilator compartment.
5. Note the battery wire harness colors to the battery terminal markings for reconnection. Disconnect the battery wire harness from the battery terminals and remove the batteries.
6. Install the replacement batteries into the lower level of the ventilator compartment with their terminals oriented as shown in the illustration.
7. Reconnect the battery wire harness to the battery terminals. Make sure the connections are as shown in the illustration.
8. Move the batteries up to their original upper position.
9. Reinstall the battery retainer bracket and secure it with the retainer screw.
10. Close the ventilator door.
11. Reinstall the battery fuse and restore power to the machine.
12. Perform the PMS procedure.




Figure 16 Battery Replacement

12 Power Supply

The power supply is attached to the inside of the controller assembly back panel. [Figure 17](#) shows the mounting and connection arrangement for the power supply.



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount service Access Procedure](#) prior to performing the following step.

1. Disconnect all pipeline hoses and set the System Power switch to 'ON'.
2. Close all cylinder valves except the O2 valve.
3. Set the oxygen flow to 5 l/min.
4. Open the other gas flow control valves to deplete the pressure from the system.
5. Close the O2 cylinder valve and close the flow control valves. Press the O2 Flush valve to deplete oxygen pressure from the system.
6. Set the System Power switch to 'OFF' and disconnect AC power from the machine.
7. Detach the cylinder supply from the yoke.
8. Remove the screws securing the back panel on the controller assembly; while supporting the panel, carefully open it as shown in the illustration.



CAUTION

The power supply contains static sensitive devices. Use ESD protection while handling this assembly.

9. Disconnect the input cable from J1 on the power supply, and disconnect the output cable from J3 on the power supply.
10. Disconnect the ground wire from the back panel. Remove the wires from the cable clips, and separate the panel from the machine.
11. Remove the two power supply mounting screws and lock washers holding the power supply to the panel.
12. Install the replacement power supply on the panel - oriented as shown in the illustration.
13. While supporting the panel, reconnect the ground wire, and reconnect the input and output cables to the power supply at J1 and J3.
14. Place the wires back in their cable clips.
15. Reinstall the back panel on the controller assembly.
16. Reconnect the pipeline and cylinder supplies and restore power to the machine.
17. Set the main power switch to ON and verify that the power-up sequence is successfully completed.
18. Perform the PMS procedure.

SV00464




Figure 17 Power Supply

Replacement Procedure

Fabius Tiro

13 Control PCB Assembly

The control PCB assembly is accessible by removing the back cover of the controller housing. [Figure 18](#) shows the mounting arrangement and electrical connections to this assembly.



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount Service Access Procedure](#) prior to performing the following step.



NOTE

Verify the original software level of the PCB prior to removal. Replacement control PCBs may contain a different software level; a software update may be necessary. If applicable, perform the [Software Update Procedure](#) with the appropriate level software according to the applicable Software Upgrade Service Procedure after the replacement PCB is installed.

1. If possible, record the market settings on the unit for re-configuration of the replacement PCB. Refer to the Diagnostics Section to access the appropriate service screen. If these settings can not be displayed, refer to the Table 1 for the settings that correspond to your region.
2. If applicable, record the Device ID number, Machine Serial number, and Release Code number located in the Secure Options screen in the Main Service menu for all ventilator feature option(s) that are enabled, for reconfiguration of the replacement processor.
3. Set the System Power switch to 'OFF' and disconnect AC power from the machine.
4. Remove the back panel from the controller housing. Provide support for the panel to prevent damage to the power supply cable connections.



CAUTION

The control PCB contains static sensitive devices. Use ESD protection when handling this assembly.

5. Disconnect all of the cables from the board connectors indicated in the illustration. Also disconnect the tubing from the pressure sensors on the PCB.
6. Remove the screws holding the PCB assembly to the controller housing, and carefully remove the PCB assembly.
7. Install the replacement PCB assembly with the hardware that was previously removed.



NOTE

CONTROL PCB P/N 4116632:

Verify jumper at JP3 is installed on both pins. This activates the 3V lithium battery (BT1) for the non-volatile RAM.



NOTE:

CONTROL PCB P/N 4118079

Verify jumper at JP2 is installed on both pins. This activates the 3V lithium battery (BT1) for the non-volatile RAM.

8. Reconnect all of the cables and hoses that were previously disconnected.



NOTE

Replacement PCB's come with a new hose attached to the LS2 pressure sensor. Remove and discard the original 2-hose assembly from the breathing pressure interface panel connector inside the unit and attach the new hose to the breathing pressure interface connector.

9. Reinstall the back panel on the controller housing.
10. Restore power to the machine and observe the power-up diagnostic display (see the [Diagnostics](#) Section) to verify that the replacement PCB is working properly.
11. Access the appropriate service screen and update the machine serial number, clear service log and reset hours run. The machine serial number is located on the rear of the machine.
12. Perform the market settings per the original machine configuration or customer requirements.



NOTE

The default "Market Selection" for the replacement PCB is "U.S." To change the language and pressure measurement units, the "Market Selection" must be changed to "Non-U.S." (see the [Diagnostics](#) Section). Refer to [Table 1](#) for the market settings that correspond to your region.

13. If Applicable, perform [Pump Calibration](#) procedure in accordance with [Adjustment and Calibration](#) Section.
14. If applicable, enter any previously recorded release codes to restore any ventilator feature option(s), or contact DrägerService - Technical Support: Phone 1-800-4-Dräger, Phone: 215-721-5402, or e-mail to techsupport@draegermed.com and provide the following information to receive the necessary release codes:
 - Machine Type (Fabius GS or Fabius Tiro)
 - Feature Description and Part Number
 - Machine Serial Number
 - Device ID Number
15. Perform the PMS procedure, including all calibrations and site configurations.

Table 1 Market Settings

| | O2 Whistle | Language | Date Format | Time Format | O2 Position | Alarm Sounds | Pres. Units |
|-----|------------|----------|-------------|-------------|-------------|--------------|-------------|
| USA | Disabled | English | M/D/Y | 12 Hour | Right | Dräger | cmH2O |

Replacement Procedure

Fabius Tiro

| | O2 Whistle | Language | Date Format | Time Format | O2 Position | Alarm Sounds | Pres. Units |
|--------------------|-------------------|-----------------|--------------------|--------------------|--------------------|---------------------|--------------------|
| Germany | Enabled | German | D/M/Y | 24 Hour | Left | EN740 | mbar |
| UK/Ireland | Enabled | UK English | D/M/Y | 24 Hour | Left | EN740 | mbar |
| France | Enabled | French | D/M/Y | 24 Hour | Left | EN740 | mbar |
| Austria | Enabled | German | D/M/Y | 24 Hour | Left | EN740 | mbar |
| Belgium/Luxembourg | Enabled | French | D/M/Y | 24 Hour | Left | EN740 | mbar |
| Switzerland | Enabled | German | D/M/Y | 24 Hour | Left | EN740 | mbar |
| Canada | Enabled | English | M/D/Y | 12 Hour | Right | Drager | cmH2O |
| Spain | Enabled | Spanish | D/M/Y | 24 Hour | Left | EN740 | mbar |
| Netherlands | Enabled | Dutch | D/M/Y | 24 Hour | Left | EN740 | mbar |




Figure 18 Control PCB Assembly, P/N 4116632

R602002T13_Replacement_Control_PCB.fm 17.06.04
 All rights reserved. Copyright reserved.

Replacement Procedure

Fabius Tiro




Figure 19 Control PCB Assembly, P/N 4118079

14 Pneumatic (PEEP Control) Assembly

The Pneumatic assembly is accessible by removing the back cover of the controller housing. [Figure 20](#) shows the mounting, tubing and electrical connections for this assembly.



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount Service Access Procedure](#) prior to performing the following step.

1. Disconnect all pipeline hoses and set the System Power switch to 'ON'.
2. Close all cylinder valves except the O2 valve.
3. Set the oxygen flow to 5 l/min.
4. Open the other gas flow control valves to deplete the pressure from the system.
5. Close the O2 cylinder valve and close the flow control valves. Press the O2 flush valve to deplete oxygen pressure from the system.
6. Set the System Power switch to 'OFF' and disconnect AC power from the machine.
7. Detach the cylinder supply from the yoke.
8. Remove the back panel from the controller housing. Provide support for the panel to prevent damage to the power supply cable connections.
9. Disconnect the following items:
 - Pump cable at J8 on controller PCB
 - PEEP cable at flying connector (from J7)
 - Small dia tubing at LS1 sensor on controller PCB
 - PEEP hose at hose pneumatic assembly PEEP valve
 - Vacuum hose at pneumatic assembly reservoir
10. Loosen the captive mounting screw on the pneumatic assembly.
11. Remove the rear mounting screw and carefully slide the assembly out of the controller housing. (Make sure the back panel with the power supply attached is supported to prevent damage to its connections.)
12. Slide the replacement pneumatic assembly into the controller housing and secure it with the captive mounting screw. reinstall the rear mounting screw and cable clamp.
13. Reconnect the items that were previously disconnected:
 - Pump cable at J8 on controller PCB
 - PEEP cable at flying connector
 - Small dia tubing at LS1 sensor on controller PCB
 - PEEP hose
 - Vacuum hose
14. Reinstall the back panel on the controller assembly.
15. Reconnect the pipeline and cylinder supplies and restore power to the machine.

16. If Software version 2.N is installed, perform the [Pump Calibration](#) procedure in accordance with the Adjustment and Calibration Section.
17. Perform the PMS procedure.

SV00459




Figure 20 PEEP Control Assembly

R602002T14_Replacement_Pneumatic.fm 17.06.04
All rights reserved. Copyright reserved.

15 Flow Meter Bezel Assembly

Following is the procedure for removing and replacing the flow meter bezel assembly. See [Figure 21](#). This assembly comprises the fresh gas flow meter, flow control valves, flow sensors and filter, pipeline pressure gauges, SORC, and the fresh gas display PCB. Replacement procedures for these individual sub-assemblies are given in subsequent paragraphs.



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount Service Access Procedure](#) prior to performing the following step.

1. Disconnect all pipeline hoses and set the System Power switch to 'ON'.
2. Close all cylinder valves except the O2 valve.
3. Set the oxygen flow to 5 l/min.
4. Open the other gas flow control valves to deplete the pressure from the system.
5. Close the O2 cylinder valve and close the flow control valves. Press the O2 Flush valve to deplete oxygen pressure from the system.
6. Set the System Power switch to 'OFF' and disconnect AC power from the machine.
7. Detach the cylinder supply from the yoke.
8. Disconnect the pipeline and cylinder supplies, and deplete all pressures.
9. Remove the two screws that secure the Service Access Panel at the rear of the Core Unit.
10. Loosen the two captive screws inside the access panel securing the flow meter bezel assembly.
11. At the front of the machine, carefully pull the assembly forward about 2.5 inches. Disconnect the ribbon cable and the O2 alarm switch wire harness from the fresh gas display PCB, and disconnect each flow sensor cable. Continue to pull the assembly forward.

The individual sub-assemblies can now be serviced as described in subsequent paragraphs.

12. When reinstalling the flow meter bezel assembly, push the assembly into its compartment; reattach the ribbon cable and O2 Alarm Switch Wire Harness to the fresh gas display PCB, and reconnect each flow sensor cable. Continue to push the assembly into its compartment until it is fully seated.
13. Tighten the two captive screws that secure the Flow Meter Bezel Assembly.
14. Reinstall the two screws and Service Access Panel located at the back of the machine.
15. Reconnect the pipeline and cylinder supplies and restore power to the machine.
16. Perform the PMS procedure.




Figure 21 Flowmeter Bezel Assembly

16 Fresh Gas Flow Meter

Access to the fresh gas flow meter requires removal of the [Flow Meter Bezel Assembly](#). [Figure 22](#) shows the tubing and mounting arrangement for the fresh gas flow tube. For clarity, other sub-assemblies are not shown in the illustration.



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount Service Access Procedure](#) prior to performing the following step.

1. Disconnect all pipeline hoses and set the System Power switch to 'ON'.
2. Close all cylinder valves except the O2 valve.
3. Set the oxygen flow to 5 l/min.
4. Open the other gas flow control valves to deplete the pressure from the system.
5. Close the O2 cylinder valve and close the flow control valves. Press the O2 Flush valve to deplete oxygen pressure from the system.
6. Set the System Power switch to 'OFF' and disconnect AC power from the machine.
7. Detach the cylinder supply from the yoke.
8. Remove the screws securing the bezel assembly (refer to the [Flow Meter Bezel Assembly](#) instructions) and bring the assembly forward for access to the fresh gas flow meter.
9. Disconnect the tubing from the flow tube ports and remove the Luer fittings.
10. Remove the flow tube mounting nuts to release the flow tube.
11. Install the replacement flow tube with the mounting brackets and nuts supplied. Ensure that the float moves freely without sticking, and that the flow tube is oriented properly.
12. Reinstall the Luer fittings and reconnect the tubing
13. Reinstall the flow meter bezel assembly in the machine (refer to the [Flow Meter Bezel Assembly](#) instructions).
14. Reconnect the pipeline and cylinder supplies and restore power to the machine.
15. Perform the PMS procedure.

Fabius Tiro

Replacement Procedure




Figure 22 Fresh Gas Flow Meter

17 Flow Control Valves

Flow control valves can be replaced from the front panel while the flow meter bezel assembly is in place. [Figure 23](#) shows the hardware arrangement for a typical flow control valve. Should the valve manifold need to be replaced, its mounting arrangement is shown in the illustration. For clarity, other sub-assemblies are not shown in the illustration.


1. Set the main power switch to OFF.
2. Disconnect the pipeline and cylinder supplies.
3. Close all cylinder valves except the O2 valve.
4. Set the oxygen flow to 5 l/min.
5. Open the other gas flow control valves to deplete the pressures from the system.
6. Close the O2 cylinder valve and close the flow control valves. Press the O2 Flush valve to deplete oxygen pressure from the system.
7. Remove the cover from the knob of the valve to be replaced.
8. Loosen the nut securing the knob to the valve, and pull the knob from the valve.
9. Un-screw the flow control valve from the valve manifold.
10. Install the replacement valve in the manifold.



NOTE

Use plastic jaw pliers (P/N 7910296) to hold the knob while tightening the nut in the next step.

11. Reinstall the knob and tighten the nut to secure the knob to the valve.
12. Reinstall the cover on the knob.
13. Reconnect the pipeline and cylinder supplies and restore power to the machine.
14. Set the main power switch to ON and ensure that the gas flow can be controlled over its entire range.
15. Perform the PMS procedure.



SV00485

Figure 23 Flow Control Valves

18 Fresh Gas Flow Sensors and Filter Assembly

Access to the flow sensors and filter assembly requires access to the [Flow Meter Bezel Assembly](#). [Figure 24](#) shows the arrangement for this sub-assembly. For clarity, other sub-assemblies are not shown in the illustration.



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount Service Access Procedure](#) prior to performing the following step.



NOTE

The safety valve, filter, and flow sensors can be serviced without completely removing the sub-assembly from the flow meter bezel assembly. The illustration shows an exploded view to indicate connections that would not otherwise be visible.

1. Disconnect all pipeline hoses and set the System Power switch to 'ON'.
2. Close all cylinder valves except the O2 valve.
3. Set the oxygen flow to 5 l/min.
4. Open the other gas flow control valves to deplete the pressure from the system.
5. Close the O2 cylinder valve and close the flow control valves. Press the O2 Flush valve to deplete oxygen pressure from the system.
6. Set the System Power switch to 'OFF' and disconnect AC power from the machine.
7. Detach the cylinder supply from the yoke.
8. Loosen the screws securing the bezel assembly (refer to the [Flow Meter Bezel Assembly](#) instructions) and bring the assembly forward for access to the flow sensor sub-assembly.
9. Remove and retain the two top screws and two underside screws attaching the Flow Sensor Housing Cover to the Flow Sensor Housing. Carefully remove the housing cover to gain access to the safety valve, filter, and flow sensors.
10. The filter is accessible by un-screwing the safety valve from the manifold. The flow sensors and manifold are held in place by a spring clip arrangement. The SORC bypass valve is secured to the flow sensor bracket by two screws and kep nuts as shown in the illustration.
11. Apply Loctite #425 (Blue) to the adapter fitting on the end of the flowsensor, then install the replacement angle connector on the flowsensor hand tight.



CAUTION

While securing the plastic hex located on the flowsensor, tighten the angle connector to the flowsensor and align per original location.

12. Install the O-ring on the white manifold squarely by twisting the O-ring side to side until it is fully seated against the shoulder. Do not allow the O-ring to roll, it will not allow proper sealing once the flow valve is installed.

13. Place the flow meter opening squarely over the O-ring and twist the flow meter side to side until it is seated onto the manifold. Do not force the flow meter over the O-ring or it will not seat properly.
14. Reinstall the two top screws and two underside screws and carefully reattach the Flow Sensor Housing Cover to the Flow Sensor Housing.
15. Reinstall the flow meter bezel assembly in the machine (refer to the [Flow Meter Bezel Assembly](#) instructions).
16. Reconnect the pipeline and cylinder supplies and restore power to the machine.
17. Perform the PMS procedure including calibrations.




Figure 24 Fresh Gas Flow Sensors and Filter Assembly

19 Pipeline Pressure Gauges

Access to the pipeline pressure gauges requires removal of the [Flow Meter Bezel Assembly](#). [Figure 25](#) shows the tubing and mounting arrangement for the gauges. For clarity, other sub-assemblies are not shown in the illustration.



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount Service Access Procedure](#) prior to performing the following step

1. Disconnect all pipeline hoses and set the System Power switch to 'ON'.
2. Close all cylinder valves except the O2 valve.
3. Set the oxygen flow to 5 l/min.
4. Open the other gas flow control valves to deplete the pressure from the system.
5. Close the O2 cylinder valve and close the flow control valves. Press the O2 Flush valve to deplete oxygen pressure from the system.
6. Set the System Power switch to 'OFF' and disconnect AC power from the machine.
7. Detach the cylinder supply from the yoke.
8. Loosen the screws securing the bezel assembly (refer to the [Flow Meter Bezel Assembly](#) instructions) and bring the assembly forward for access to the pipeline pressure gauges.
9. Disconnect the tubing from the fitting on the gauge to be replaced.
10. Remove the retainer plate screws, and separate the retainer plate with the gauges from the bezel.
11. Un-thread the straight fitting from the back of the gauge to be replaced.
12. Remove the gauge mounting nut, and remove the gauge from the retainer plate.
13. Install the replacement gauge in the retainer plate - make sure it is oriented properly when tightening the gauge mounting nut.
14. Reinstall the straight fitting in the back of the gauge.
15. Reinstall the gauge retainer plate in the bezel.
16. Reconnect the tubing that was previously removed.
17. Reinstall the flow meter bezel assembly in the machine (refer to the [Flow Meter Bezel Assembly](#) instructions).
18. Reconnect the pipeline and cylinder supplies and restore power to the machine.
19. Perform the PMS procedure.

SV88012




Figure 25 Pipeline Pressure Gauges

R6020002T19_Replacement_Pipeline_Pressure_Gauges.fm 17.06.04
All rights reserved. Copyright reserved.

20 SORC (Sensitive Oxygen Ratio Controller)

Access to the SORC requires removal of the [Flow Meter Bezel Assembly](#). [Figure 26](#) shows the tubing and mounting arrangement for the SORC. For clarity, other sub-assemblies are not shown in the illustration.



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount Service Access Procedure](#) prior to performing the following step.

1. Disconnect all pipeline hoses and set the System Power switch to 'ON'.
2. Close all cylinder valves except the O2 valve.
3. Set the oxygen flow to 5 l/min.
4. Open the other gas flow control valves to deplete the pressure from the system.
5. Close the O2 cylinder valve and close the flow control valves. Press the O2 Flush valve to deplete oxygen pressure from the system.
6. Set the System Power switch to 'OFF' and disconnect AC power from the machine.
7. Detach the cylinder supply from the yoke.
8. Loosen the screws securing the bezel assembly (refer to the [Flow Meter Bezel Assembly](#) instructions) and bring the assembly forward for access to the SORC.
9. Disconnect the tubing from the SORC where shown in [Figure 26](#). Mark these tubes so they can be reattached to the replacement SORC in the same manner.
10. Remove the four SORC mounting nuts and remove the SORC from the flow meter bezel assembly.
11. Install the replacement SORC with the four mounting nuts that were previously removed.
12. Re-connect the O2 and N2O tubing to the SORC in the same manner as the original.
13. Reinstall the flow meter bezel assembly in the machine (refer to the [Flow Meter Bezel Assembly](#) instructions).
14. Reconnect the pipeline and cylinder supplies and restore power to the machine.
15. Perform the PMS procedure.

Fabius Tiro

Replacement Procedure




Figure 26 SORC

R6020002120_Replacement_SORC.fm 17.06.04
All rights reserved. Copyright reserved.

21 Fresh Gas Display PCB

Access to the fresh gas display PCB requires removal of the [Flow Meter Bezel Assembly](#). [Figure 27](#) shows the cable connections and PCB mounting details. For clarity, other sub-assemblies are not shown in the illustration.



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount Service Access Procedure](#) prior to performing the following step.

1. Disconnect all pipeline hoses and set the System Power switch to 'ON'.
2. Close all cylinder valves except the O2 valve.
3. Set the oxygen flow to 5 l/min.
4. Open the other gas flow control valves to deplete the pressure from the system.
5. Close the O2 cylinder valve and close the flow control valves. Press the O2 Flush valve to deplete oxygen pressure from the system.
6. Set the System Power switch to 'OFF' and disconnect AC power from the machine.
7. Detach the cylinder supply from the yoke.
8. Loosen the screws securing the bezel assembly (refer to the [Flow Meter Bezel Assembly](#) instructions) and bring the assembly forward.
9. Remove the two screws securing the flow sensor sub-assembly to the bezel assembly, and move the sub-assembly aside for access to the fresh gas display PCB. (Refer to the [Fresh Gas Flow Sensors and Filter Assembly](#) instructions.)
10. Disconnect the cables from the fresh gas display PCB where shown in [Figure 27](#).
11. Remove the PCB mounting screws and lock washers, and remove the PCB.
12. Install the replacement fresh gas PCB, with the insulator in place using the hardware that was previously removed.
13. Reconnect the cables to the PCB.
14. Reinstall the flow sensor sub-assembly.
15. Reinstall the flow meter bezel assembly in the machine (refer to the [Flow Meter Bezel Assembly](#) instructions).
16. Reconnect the pipeline and cylinder supplies and restore power to the machine.
17. Perform the PMS procedure.

SV88014




Figure 27 Fresh Gas Display PCB

R6020002T21_Replacement_Fresh_Gas_Display.fm 17.06.04
All rights reserved. Copyright reserved.

22 Monitor Bezel Assembly

The monitor bezel assembly comprises the keypad, flat panel display, rotary switch, and LED lamp PCB assembly. Access to the bezel assembly mounting screws is through the ventilator compartment as shown in [Figure 28](#).



NOTE

Replacement keypads do not include keypad inserts. Therefore, removal of the original keypad insert for later reuse is necessary. Refer to [Figure 28](#) for keypad insert location.

1. Set the main power switch to OFF and disconnect AC power from the machine.
2. Swing open the ventilator door.
3. While supporting the bezel assembly, remove the two screws securing the left side of the bezel.



CAUTION

Do not pull the bezel outward as this will break the locking tabs at the right side of the bezel.

4. Slide the bezel to the right to release its locking tabs from the machine. See [Figure 28](#).



NOTE

Verify that bezel locking tabs are clear of the machine before pulling the bezel away from the machine.

5. Disconnect the large ribbon cable from the keypad connector, and the smaller ribbon cable from the display panel.
6. Disconnect the ground wire.

Refer to the parts list for identification of individual components.

7. Reconnect the ground wire and the two ribbon cables to the replacement assembly.
8. Carefully fold the cables behind the bezel assembly as you bring it into position on the machine. Hold the bezel against the machine and slide it to the left to engage the locking tabs in the slots in the machine housing.
9. While supporting the bezel assembly, reinstall the two screws at the back of the bezel on the left side.
10. Close the ventilator door.
11. Restore power to the machine and verify all of the keypad and display functions.
12. Perform the PMS procedure.

SV00463




Figure 28 Monitor Bezel Assembly

23 Spirolog Sensor and Cable

1. Set machine power switch to OFF and disconnect AC power from the machine.
2. Remove the expiratory hose terminal from the breathing system, and slide the Spirolog sensor out of the breathing system.
3. Disconnect the cable from the sensor by squeezing the connector at the points opposite the latch while pulling the connector away from the sensor.



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount Service Access Procedure](#) prior to performing the following step.

4. Remove the back panel from the controller housing, and disconnect the Spirolog cable from J18 on the controller PCB.
5. Loosen the nut on the strain relief bushing at the interface panel. Remove the strain relief bushing from the panel while supporting the rear nut.
6. Install the replacement Spirolog strain relief bushing into the interface panel, and tighten the nut.
7. Connect the replacement Spirolog cable to J18 on the controller PCB.
8. Route the cable from the machine to the breathing system in the same manner as the original.
9. Connect the outboard end of the Spirolog cable to the flow sensor. Note that the sensor and connector are keyed.
10. Reinstall the back panel on the controller housing, and restore AC power to the machine.
11. Perform the PMS procedure.

ADJUSTMENT AND CALIBRATION PROCEDURES

Adjustment and Calibration Procedures

Fabius Tiro

1 Cylinder Pressure Regulator Adjustment

1. Set the System Power switch to OFF.



NOTE

Minimum cylinder pressures for this adjustment shall be N2O: 700 psi; O2 or Air: 1000 psi.



NOTE

[Figure 1](#) shows test connections for the O2 regulator adjustment. If you are adjusting the N2O regulator, connect the test gauge in series with the N2O pipeline hose.

2. Connect test pressure gauge (P/N 4114807) between machine's pipeline inlet connector and the pipeline supply hose.
3. Open the cylinder valves and set the System Power switch to ON.
4. Set the oxygen flow to 4 liters per min. (If you are adjusting the N2O regulator, also set the nitrous oxide flow to 4 liters per minute.)
5. Depress the push button on the test device.
6. Release the push button. After the pressure decay stabilizes, the gauge should indicate 32 - 40 psi.



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount Service Access Procedure](#) prior to performing the following step.

7. Remove the Service Access Panel at the rear of the unit. Refer to [Figure 1](#).
8. Remove the acorn nut from the regulator to expose the adjusting screw. For N2O, turn the screw until the test gauge indicates 35 psi. (50 psi for CSA machines). For O2 and Air, use the compensated regulator output setting based on the cylinder pressure given in the following table.



CAUTION

Based on information supplied by the cylinder regulator manufacturer, when the regulator is used for O2 or Air, its output pressure will decrease 0.5 psi for every 100 psi increase in cylinder pressure above 1000 psi. Currently, these pressure regulators are calibrated at 35 psi with a cylinder supply of 1000 psi. If a 2000 psi cylinder is then installed, the regulator output will be 42 psi. This change in output must be compensated for to provide accurate performance throughout the cylinder's working range.

SV88015




Figure 1 Cylinder Pressure Regulator Adjustment: O2 Connection Illustrated



NOTE
Cylinder pressure compensation for the N2O regulator is not required.

| Cylinder Pressure (psi) | Compensated Regulator Output Tolerances (-4/+2) |
|-------------------------|---|
| 2000 | 27 - 33 (*30 - 36) |
| 1800 | 28 - 34 (*31 - 37) |
| 1600 | 29 - 35 (*32 - 38) |
| 1400 | 30 - 36 (*33 - 39) |
| 1200 | 31 - 37 (*34 - 40) |
| 1000 | 32 - 38 (*35 - 41) |

* Canada settings

9. Verify the adjusted regulator maintains its compensated output tolerance.
10. If the O₂ cylinder regulator output was adjusted, perform the following test:
 11. Set the ventilator to Volume Mode ventilation.
 12. Set the Auxiliary O₂ and Fresh Gas O₂ flows to 10 L/min.
 13. Press the O₂ Flush and verify the Lo O₂ Supply alarm is not active. If the alarm is active refer to Section 3 [Oxygen Supply Pressure Alarm Switch Adjustment](#).
 14. Reinstall the acorn nut on the regulator.
 15. Close the cylinder valves and allow pressure to drain from the system.
 16. Close all of the flow control valves and set the System Power switch to OFF.
 17. Disconnect the test gauge.
 18. Connect the pipeline hoses.
 19. Perform the PMS Procedure.

2 Gas Inlet Regulator Output Adjustment

O₂ Inlet Regulator

1. Bleed all cylinder and pipeline pressures. Disconnect all cylinders and pipeline hoses.



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount Service Access Procedure](#) prior to performing the following step.

2. Wall Mount Access (If Applicable)
 - a) Regulator access can be obtained from the underside when the unit is rotated 180°.
3. Trolley Mount Access (If Applicable)

- a) Remove the top drawer by releasing the locking tabs located in the drawer slides (x2).
- b) Fully extend the writing tray to its out position.
4. Remove the pneumatic hose from the O₂ outlet of the gas inlet block and interconnect a digital manometer as shown in [Figure 2](#).
5. Reconnect the O₂ pipeline supply hose to the O₂ pipeline inlet connector and pressurize the O₂ supply.
6. Set the O₂ flow to 4 L/min.
7. After the digital manometer stabilizes, adjust the O₂ regulator output pressure to 30 psi.
8. Deplete O₂ pressure from the pipeline supply.
9. Close the flow control valve and disconnect the O₂ pipeline hose from the inlet.
10. Remove test equipment and reconnect the O₂ pneumatic hose leading from the O₂ flow control valve to the O₂ connector on the inlet block.
11. Perform the PMS procedure.
- N₂O Inlet Regulator
12. Disconnect the N₂O pneumatic hose from the outlet connector on the gas inlet block that connects to the N₂O flow control valve.
13. Interconnect a digital manometer between the N₂O outlet connector and the hose removed in the previous step. See [Figure 2](#).
14. Reconnect the N₂O pipeline hose to the inlet connector and activate the N₂O supply.
15. Set the O₂ and N₂O flows to 4 L/min. After the digital manometer stabilizes, adjust the N₂O regulator output pressure to 30 psi.
16. Deplete the O₂ and N₂O pipeline pressures.
17. Close the O₂ and N₂O flow control valves, and disconnect both pipeline hoses from the inlets.
18. Remove the test equipment and reconnect the N₂O hose from the flow control valve to the N₂O outlet on the inlet block.
19. Perform the PMS procedure.

Adjustment and Calibration Procedures

Fabius Tiro




Figure 2 Setting Inlet Regulators

3 Oxygen Supply Pressure Alarm Switch Adjustment



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount Service Access Procedure](#) prior to performing the following step.

1. Wall Mount Access (If Applicable)
 - a) Regulator access can be obtained from the underside when the unit is rotated 180°.
2. Trolley Mount Access (If Applicable)
 - a) Remove the top drawer by releasing the locking tabs located in the drawer slides (x2).
 - b) Fully extend the writing tray to its out position.
3. Invert unit by performing steps given in the [Core Mount Inversion Disassembly](#) Section.

4. Connect test pressure gauge (P/N 4114807) between machine's oxygen pipeline inlet connector and the oxygen pipeline supply hose.
5. Open the O2 cylinder valve and set the System Power switch to ON.
6. Set the oxygen flow to 200 mL per min.
7. Close the oxygen cylinder valve.
8. As the pressure drops, the O2 SUPPLY alarm should activate when the pressure is between 16 and 24 psi as shown on the test gauge.
9. If the alarm activates when the pressure is below 16 psi or above 24 psi, turn the adjustment set screw (see [Figure 3](#)). Repeat the test and adjust as necessary to bring the set point into the correct range.



NOTE

Turn the screw counter-clockwise to decrease the alarm activation.

10. Set the System Power switch to OFF.
11. Disconnect the test gauge.
12. Connect the pipeline hoses.
13. Perform the PMS Procedure.




Figure 3 Oxygen Supply Pressure Alarm Switch Adjustment

4 Sensitive Oxygen Ratio Controller (SORC) Adjustment



NOTE

If servicing a wall mount unit, refer to the [Wall Mount Service Access Procedure](#) prior to performing this procedure.

1. Set the System Power switch to OFF.
2. Remove the two screws that secure the Service access panel near the rear of the unit.
3. Loosen the two captive screws in the access panel securing the flow meter bezel assembly. Carefully pull the assembly out of the machine for access to the SORC.
4. Connect a test pressure gauge to the O2 port on the SORC shown in [Figure 4](#).
5. Connect an oxygen analyzer to the fresh gas outlet.
6. Connect the pipeline supplies.
7. Set the System Power switch to ON.
8. Set the oxygen flow to 200 mL/min.
9. The test pressure gauge should read 25 ± 1 cmH₂O. If needed, turn the O2 control pressure adjustment screw until the test gauge reads 25 ± 1 mbar. Apply a small amount of thread sealant to the screw head.

Setting O2 Control Pressure

Adjusting the Operating Point

10. Fully open the N2O flow control valve.
11. Set the oxygen flow to 200 mL/min.
12. The N2O flow should be 50 ± 30 mL/min. Make the following adjustment if needed:
 - a) Loosen the lock nut on top of the SORC.
 - b) Turn the operating point adjustment counter-clockwise to decrease the N2O flow, or clockwise to increase the N2O flow until it is at 50 ± 30 mL/min.
13. Tighten the lock nut and verify the N2O flow is 50 ± 30 mL/min. Apply a small amount of thread sealant where the lock nut meets the threads of the adjustment point.

Adjusting the O2/N2O Ratio

14. With the N2O flow control valve fully open, set the oxygen flow to 1 L/min.
15. The N2O flow should be 3 ± 0.3 mL/min, resulting in an O2 concentration of 25 ± 3 vol.% O2 at the fresh gas outlet.

If needed, turn the N2O flow adjustment screw until the flow is 3 ± 0.3 mL/min, resulting in an O2 concentration of 25 ± 3 vol.% O2 at the fresh gas outlet. Apply a small amount of thread sealant on the adjustment screw.

O2 Concentration Test:

16. Verify O2 concentration at the fresh gas outlet, with the O2 flow settings in the following table

| O2 Flow | N2O Flow Control Valve | O2 Concentration |
|-----------------------|------------------------|------------------|
| 0.5 ± 0.05 L/min. | Fully Open | 23 - 33 |
| 1.0 ± 0.1 L/min. | | 23 - 28 |
| 3.0 ± 0.2 L/min. | | 23 - 33 |


Flow Test

17. Fully open the N2O flow control valve.

Adjustment and Calibration Procedures

Fabius Tiro

18. Slowly open the O₂ flow control valve.
19. As the O₂ flow increases, the N₂O flow should also increase, to over 10 L/min. at the maximum.
20. Slowly close the O₂ flow control valve. The N₂O flow should decrease proportionally.
21. Close the flow control valves and remove all test equipment.
22. Set the System Power switch to OFF, and reinstall the flow meter bezel assembly in the machine.
23. Perform the PMS procedure.



SV88017

Figure 4 Sensitive Oxygen Ratio Controller (SORC) Adjustment

5 Oxygen Sensor Calibration

1. Set the main power switch to ON.
2. Access the main service screen (ref. [Main Service Screen](#) Section).
3. Scroll to Calibration and select O2 to bring up the O2 Zero calibration screen. See [Figure 5](#).

Adjustment and Calibration Procedures

Fabius Tiro

Zero Calibration

4. Remove the oxygen sensor capsule from its housing and allow several minutes for the displayed offset readings to stabilize.



NOTE

The difference between the Cell A and Cell B readings should be no greater than 8.

21% Calibration


5. Press the Store Zero soft key to store the current values as the new zero calibration.
6. Reinstall the sensor capsule in its housing.
7. Expose the sensor to ambient air only (away from any open part of the breathing system) and allow it to stabilize for several minutes.
8. Press the Standby key to return to the standby screen.
9. Press the Calibrate O2 Sensor softkey. An instruction window replaces the softkey labels.
10. Press the rotary knob to start the 21% O2 calibration. During calibration the 'Calibrate O2 Sensor' softkey shall illuminate. Following calibration the screen returns to normal.



NOTE

If an "O2 Sensor Calibration Failed" message appears, refer to the Oxygen Monitoring section of the Fabius Tiro Operator's Instruction Manual for further information.

11. Reinstall the O2 sensor assembly in the inspiratory valve dome.



SV88018

Figure 5 Oxygen Sensor




Figure 6 Oxygen Monitor Service Screen

6 Pressure Calibration

1. Set the main power switch to ON.
2. Access the main service screen (ref. [Main Service Screen](#) Section).
3. Scroll to Calibration and select Pressure to bring up the pressure calibration screen. See [Figure 7](#).

Zero Calibration

- a) Remove the pressure sample line and ventilator hose from the breathing system, and expose the sample line to air.
 - b) Let the current pressure value stabilize.
 - c) Touch the 'Store Zero' softkey to enter the current value as the calibration zero.
4. Re-connect the pressure sample line and the ventilator hose to the breathing system.
 5. Press the Standby key to return to the standby screen.




Figure 7 Pressure Calibration Screen

7 Fresh Gas Flow Calibration

1. Set the main power switch to ON.
2. Access the main service screen (ref. [Main Service Screen](#) Section).
3. Scroll to Calibration and select Fresh Gas Flow to bring up the flow calibration screen. See [Figure 8](#).
4. Close all flow control valves.
5. Verify no fresh gas flow is present.
6. Touch the "Store Zero" softkey.
7. At completion of calibration "Stored Zero" shall indicate Pass or Fail.
8. Press the Standby key to return to the standby screen.




Figure 8 Fresh Gas Flow Calibration

8 PEEP Valve Calibration



NOTE

Before performing the PEEP valve calibration, a valid pressure calibration and system leak test must be performed first. Otherwise an 'Inv pres. cal' message will appear.

1. Set the main power switch to ON.
2. Access the main service screen (ref. [Main Service Screen](#) Section).
3. Scroll to Calibration and select PEEP to bring up the PEEP calibration screen.
4. Attach test hose (P/N 9995112, 22 mm x 12 in long) from the expiratory valve terminal to the inspiratory valve terminal and remove breathing bag.
5. Touch the Calibrate softkey. (The calibration process can take several minutes.)
6. At completion of calibration a "Pass" or "Fail" message will appear.
7. Press the Standby key to return to the standby screen.
8. Remove the short hose from the inspiratory and expiratory valves.
9. Re-connect the breathing system hoses and breathing bag.

Adjustment and Calibration Procedures

Fabius Tiro




Figure 9 PEEP Valve




Figure 10 PEEP Valve Calibration Screen

9 Vacuum Adjustment

The vacuum adjustment is located on the pneumatic assembly and is accessible by removing the back cover of the controller housing.



NOTE

If the Core Unit is wall mounted, refer to the [Wall Mount Service Access Procedure](#) prior to performing the following step.

1. Remove the back panel from the controller housing. Provide support for the back panel to prevent damage to the power supply cable connections.
2. Set the Main Power switch to ON and set vent control to Volume mode.
3. Remove the APL bypass hose from the breathing system and attach it to a digital manometer.
4. Adjust the variable restrictor on the pneumatic assembly (see [Figure 11](#)) until the vacuum is between -150 and -240 cmH₂O (mbar) as shown on the digital manometer.
5. Disconnect the APL bypass hose from the digital manometer and reconnect it to the APL port on the breathing system.
6. Reinstall the back panel on the controller assembly.
7. Perform the PMS procedure.




Figure 11 Vacuum Adjustment

10 Vacuum Pressure Pump Assembly (Software Version 2.N)

10.1 Pump Calibration - If Applicable



NOTE

The ability to calibrate the pump is dependant upon the current PCB installed in the Core Unit.

1. Access the 'System Service Screen' screen and note the installed 'Hardware Part Number'. Refer to [System Service Screen](#) Section.
2. If the 'Hardware Part Number' listed is 4118079, proceed with the next step. If the 'Hardware Part Number' is 4116632, the pump calibration procedure is not available.
3. Access the 'Real Time Values' screen, Refer to [Real Time Values](#) Section.



NOTE

Upon initial entry into this screen, the calibration operation is not active.



NOTE

The pump calibration operation can be performed with the pump in either the ON or OFF state. However, after the calibration operation is complete, the pump will remain in the ON state unless manually turned off. Refer to [Pump On/Off](#) Section.

4. Press the 'PUMP CAL' softkey to activate the pump calibration operation. The 'PUMP CAL' label will be highlighted in reverse upon activation and will remain in this condition until the calibration operation is either complete or aborted.

Periodic Manufacturer's Certification

Refer to Appendix A

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

Spare and Replacement Parts List

Appendix B

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

SOFTWARE UPDATE PROCEDURE

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

Software Update Procedure

This section outlines the equipment needed, Service Procedures for specific software versions, and the BootStrap Download Procedure.



NOTE

The screen illustrations contained in this document are for reference only and therefore may or may not reflect the current equipment or software version being installed or exact Fabius Product.

1 Requirements

Equipment required:

- Cable Asm, Vitalink, Part No. 4110328 (9-pin to 9-pin). Alternate cable combination: 7901808 and 7901762
- IBM® PC or IBM PC Compatible configured with:
 - Windows 98, Windows NT, Windows 2000, Windows NT 4.0, or Windows XP
 - RS-232C Serial port connected to COM 1 or COM 2
 - Parallel port (LPT1 or LPT2)
 - CD-ROM and Floppy Drive
 - Standard sound card
 - 2 PCMCIA slots (on laptops)
 - Modem (or external modem) V34+ 28800 Baud

The following Service Procedure(s) is available via Lotus Notes and provides detailed instructions to install/upgrade specific software versions. Refer to the following Service Procedure or contact DrägerService - Technical Support: Phone 1-800-4-Dräger, Phone: 215-721-5402, or e-mail to techsupport@draegermed.com.

- SP00297 - Software Version 2.0N to 2.N Upgrade

2 Boot Strap Download Procedure

This procedure uses Version 2.N of the SEND.IMG program. V2 is a new Win32 program that adds the Fabius Tiro boot Strap load function to the original SENDIMG.EXE program that is currently used in existing products.

When a Fabius Tiro does not have a valid bootable image installed, it can not perform an Automatic Image Reload Operation. The same PC and files are used (including SENDIMG.EXE V2).

1. Record the market setting on the unit for reconfiguration of the controller PCB after download. Refer to the [General](#) portion of the Diagnostics Section to access the appropriate service screen.
2. Set the machine power switch to the OFF position.
3. Remove the back panel and change the JP2 (on Control PCB P/N 4116632) or JP3 (on Control PCB P/N 4118079) connector on the control PCB from the Boot Flash position ([Figure 1](#)) to the Boot Strap position ([Figure 2](#)). Refer to or for proper connector locations.




Figure 1 Boot Flash Position




Figure 2 Boot Strap Position

4. Connect the 3-pin plug end of the boot strap download cable (P/N 4117459) to J13 on the control PCB. See [Figure 6](#).
5. Connect the DB9 end of the cable to the Com1 port on the PC.
6. Power up the PC and double-click the “send_img” icon on desk top. The following screen will appear:




Figure 3 Communications Settings Screen

7. Choose “Boot Strap and Image transfer”, then click OK or press Enter.

8. Set the machine power switch to ON. The bootstrap load of the FLASH programmer load program then proceed according to the status messages on the PC. The following screen will be displayed:




Figure 4 Image File Transfer Execution Screen

9. Follow the same sequence as a normal download procedure to access the software. Refer to the appropriate Service Procedure. Once the load program is executing on the machine, the image load will proceed according to messages on the machine and the PC until a message on the machine indicates completion.
10. Close the PC Loader screen and power down the PC.
11. Disconnect the boot strap download cable.
12. Return the boot selector jumper at JP2 (on Control PCB P/N 4116632) or JP3 (on Control PCB P/N 4118079) to the Boot Flash position.




Figure 5 Boot Flash Position

13. Cycle machine power and verify all self-diagnostics indicate a PASS condition. NOTE: NVRAM settings will be set when step 16 is performed.
14. Reinstall the back panel of the control unit.
15. Access the appropriate service screens (see the [General](#) portion of the Diagnostics Section) and update the machine serial number, clear service log and reset hours run. The machine serial number is located on the rear of the machine.
16. Perform PMS test Steps Calibrations, Site Configurations, Oxygen Monitor, Oxygen Concentrations, Pressure Monitor and Ventilator. Verify the correct market kit settings are entered per original machine configuration or customer requirements.




Figure 6 Boot Strap Jumper and Download Connections on Control PCB P/N 4116632




Figure 7 Boot Strap Jumper and Download Connections on Control PCB P/N 4118079

PMS PROCEDURE

8.0 PMS PROCEDURE, FABIUS Tiro

The procedures in this section shall be performed in their entirety each time a component is removed, replaced, calibrated, adjusted and during all scheduled Periodic Manufacturer's Service (PMS) visits. A PMS Checklist form, available from DrägerService, shall be completed by the Technical Service Representative (TSR) each time a PMS is performed. Steps in the procedure marked with (✓) require a response at the corresponding line on the checklist form. Space is also provided on the PMS checklist form to record the results of a vapor concentration test. Contact DrägerService for vapor concentration verification procedures.

NOTE: Verify the dates on test equipment calibration labels. DO NOT USE any test equipment having an expired calibration date. Notify your supervisor immediately if any equipment is found to be out of calibration.

In the space provided at the bottom of the PMS checklist form, record the Model and ID number of all calibrated test equipment used. Also record the calibration due dates. Examples are: multimeter, digital pressure meter, Riken gas analyzer, safety analyzer, volumeter, trace gas analyzer, simulators.

| Description | Page |
|--|------|
| Test Equipment Required: | 8-2 |
| Fabius Tiro Parts Replacement Schedule: | 8-3 |
| Periodic Manufacturer's Certification General Instructions | 8-11 |
| 8.1 Electrical Safety | 8-16 |
| 8.2 System Diagnostics | 8-18 |
| 8.3 Battery Circuit | 8-19 |
| 8.4 Configuration | 8-19 |
| 8.5 Service Data | 8-19 |
| 8.6 Calibrations | 8-20 |
| 8.7 Site Configurations | 8-21 |
| 8.8 Scavenger | 8-24 |
| 8.9 Breathing System | 8-26 |
| 8.10 Yokes & Gauges | 8-37 |
| 8.11 Gas Inlet Regulator Output | 8-38 |
| 8.12 Cylinder Regulator & Pipeline Gauges | 8-41 |
| 8.13 High Pressure Leak | 8-45 |
| 8.14 Oxygen Supply Failure Protection | 8-46 |
| 8.15 Flowmeters | 8-48 |
| 8.16 Oxygen Monitor | 8-49 |
| 8.17 Oxygen Concentrations | 8-50 |
| 8.18 SORC | 8-51 |
| 8.19 Pressure Monitor | 8-52 |
| 8.20 Ventilator | 8-54 |
| 8.21 Volume Alarms | 8-58 |
| 8.22 Audio Silence | 8-59 |
| 8.23 Oxygen Flush Valve | 8-60 |
| 8.24 Final Tests | 8-61 |

Test Equipment Required:

- Multi-Meter (Fluke or equivalent)
- Electrical Safety Analyzer (Biotek 501 Pro or equivalent)
- Test Pressure Gauge, P/N 4114807 or equivalent
- Fresh Gas Test connector, P/N 4117361 or equivalent
- Fresh Gas Leak Test Device, P/N 4113119 or equivalent
- Test Cap, P/N M33972 or equivalent: two are required
- Adapter Assembly, Test Terminal, P/N 4104389 or equivalent: two are required
- Flowmeter Test Stand (0-250 cc), P/N S000081 or equivalent
- Breathing System Leak Test Device, P/N S010159 or equivalent
- Baromed Pressure Test Fixture or equivalent
- Test Minute Volume Meter, P/N 2212300 or equivalent
- Digital Pressure Manometer (SenSym PDM 200CD or equivalent)
- Riken Gas Indicator, Model 18H, or 1802D or equivalent
- Hose Asm, vapor testing (for use w/Riken gas indicator), 4117905
- Stop Watch
- Siemens Test Lung, P/N 8401892
- Rubber Plug, P/N 7901297
- Tube, Corrugated, 22 mm x 12 in. long, P/N 9995112

Accessories:

- T-connector, P/N M32803
- Washer, M4 white (3x) M31602 (for T-connector above)
- Fitting, Str 1/4 tube x 1/8 NPT, P/N 4109318
- Fitting, Str 0.130 ID hose x 1/8 MPT, P/N 4102963

Materials Required:

- Breathing Bag, 3 liter, P/N 9995330 or equivalent
- DI-Wacker Silicon Rubber, P/N 1202537
- Loctite, Purple 222, P/N 4118558-001
- Loctite, Red 271, P/N 4118558-003
- Loctite, Blue 425, P/N 4118558-008

Repair Tools:

- Plastic Jaw Pliers, P/N 7910296
- Nut Driver (modified), P/N 4117530
- Nut Driver (two-hole spanner). P/N 7910305
- Cable Assembly, flash load RAM, P/N 4117459
- Nut Driver, 3/8 in.
- Cable Assembly, Vitalink, P/N 4110328
- Hose cutter, P/N 7900894
- Pozi Drive Set
- Wrench, Caster, P/N S010055

Test equipment illustrations are shown on following pages.


| | | |
|--------------------|-------------------|----------------------------------|
| FABIUS TIRO | APPENDIX A | PMS PROCEDURE (continued) |
|--------------------|-------------------|----------------------------------|


Fabius Tiro Parts Replacement Schedule:

| Quantity | Description | Part Number | Alternate Part Number |
|--|----------------------------------|----------------------|-----------------------|
| As Required | | | |
| 1 | Flow Sensor (5 pack) | 8403735 | |
| 1 | O2 Sensor | 6850645 | |
| 2 | O-ring, Vent Hose | 2M08777 | |
| 2 | Valve Disc | 2123249 or M23225 | |
| 1 | Hose asm, PEEP/Pmax - APL Bypass | 4117027 | |
| 1 | Hose, Pressure | 1190520 | |
| 1 | Filter, AGS | M33294 | |
| Annually: Fabius GS/Tiro 1 year kit | | 4117360-001 | 4199912 |
| 2 | Filter, Pressure & Pneumatic asm | 8402868 | |
| 1 | Diaphragm, Patient | 2600650 | |
| 4 | O-ring, Vapor | 4115864 | |
| 1 | O-ring, for Diaphragm (piston) | 8604831 | |
| Every 3 Years: Fabius GS/Tiro 3 year kit | | 4117360-002 | 4199911 |
| 1 | Diaphragm, Piston | 2600651 | |
| 2 | Battery, 12V Rechargeable | 4114229 | |
| 1 | Canister Assembly | M29320 | |
| 1 | Lip Seal | M30455 | |
| 1 | Packing Ring | M30456 | |
| 1 | PM Kit; 1 Year | 4117360-001 | 4199912 |

Disposal of used batteries and O2 Sensors:


- Batteries must be disposed of in conformity with local waste disposal regulations.
- Expired O2 sensors can be returned to:
 - Dräger Medical AG & Co. KGaA
 - Moislinger Allee 53-55
 - Reparaturannahme
 - 23542 Lübeck
 - Germany
- Do not open forcibly: danger of chemical burns.
- Do not incinerate: danger of explosion.






8401892
SIEMENS TEST LUNG

SV00025




S010159
BREATHING SYSTEM LEAK TEST DEVICE



4104389
TEST TERMINAL
ADAPTER


SV00515




SV00559

S000081
FLOW METER
TEST STAND
0-250cc/MIN

3V00552




2212300
MINUTE
VOLUMETER




4115087 33mm x 22mm ADAPTER

M33972
TEST CAP




SV00516


4114807 PRESSURE TEST ASSEMBLY , WITH ADAPTERS




SV00376




4114830-002




4114830-001




4114830-004



4114830-003



4114830-006



4114830-005

Periodic Manufacturer's Certification General Instructions

The purpose of these procedures is to provide detailed instructions for performing a Periodic Manufacturer's Certification (PMC) inspection on the Fabius Tiro anesthesia machine.

A PMC consists of a complete Periodic Manufacturer's Service procedure and a certification level inspection based on Draeger Medical, Inc. Recommendations and equipment performance. Additional inspections are also performed to ensure proper product labeling.

Several additional documents have been created to assist the technician through the process. Following is a brief description of the purpose of each document.

Field Service Procedure:

Periodic Manufacturer's Certification Forms - Part Number SP00175.

This procedure illustrates the sample checklists with typical periodic maintenance items filled in, including vapor concentration verification tests, parts replaced, general comments and certification levels. Also included are sample PMC labels marked to show several levels of certifications. An excerpt from Draeger Medical, Inc.'s *Anesthesia System Risk Analysis and Risk Reduction* is included, and also a sample of an Executive Summary to be furnished to the hospital's Risk Manager or Chief of Anesthesia.

Field Service Procedure:

NAD Recommendation Guidelines Index Anesthesia Systems - Part Number S010250.

This Guideline was created to provide an assessment of each machine's certification. It contains various comprehensive overviews of possible equipment conditions and their associated certification levels.

The first list in the Recommendation Guidelines is a reference chart for machine certification based on equipment status. The second is an abbreviated summary of all NAD Recommendations and Failure Codes including the Condition Number, Equipment Condition, Recommended Corrections, Certification Code, and Tests Affected when applicable.

There is also a matrix classified as "Failure Codes" which identifies the correct manner in which to document equipment tests that fail, or were unable to be performed due to circumstances beyond the control of the service technician performing the inspection. (Ex: Air cylinder supply is unavailable to perform Air High Pressure Leak test.) The Failure Codes section also indicates suggested resolution of the situation. Failure Code numbers begin at 34 and use the same certification levels strategy, and carry the same weight as NAD Recommendation equipment condition codes.

The next section of the guideline lists all NAD Recommendations identified at a machine's major assembly level. This section is divided into subsections titled: "Anesthesia System", "Vaporizers", "Absorber System", "Ventilator" and "Scavenger System". The final matrix is the most comprehensive index sorted by machine model and includes Equipment Condition, Certification Code, and NAD Recommendations. It also specifies any suggested upgrade path including ordering information that should be taken such as installing a Bellows with Pressure Limit Control 4109664-S01 Kit, after market modification kit to a machine not equipped with pressure limit control.

The letters A, B, C, D and the Roman Numerals I, II are used as codes in the individual matrix for each model of anesthesia machine. The letters A, B, C, and D are used in descending order to indicate the certification level of the equipment. They are as follows:

- A = Certified
- B = Certified with Recommendations
- C = Conditionally Certified
- D = No Certification

Roman Numerals I and II do not affect the certification level but rather are provided to give further instructions to the end user as follows:

- I = The system in its present configuration shall only be used with a CO2 monitor incorporating an apnea warning. The operator of the system is advised to frequently scan the CO2 readings and alarm thresholds.
- II = The present configuration of equipment requires that the unit operate at all times with an oxygen analyzer that includes a low oxygen warning. The operator of the system is advised to frequently scan the oxygen readings and alarm limits.

Following is an explanation of machine certification levels:

Certified- No recommendations apply to machine being inspected. (Only item number 33 - "No Recommendations" shall apply for this certification level.)

Certified with Recommendations- A numbered recommendation with a code of B applies to the machine being examined.

Conditionally Certified- A numbered recommendation with a code of BC, BCI or BCII applies to the machine being examined.

No Certification- A numbered recommendation with a code of D applies to the machine being examined.

When multiple recommendations apply, use the above list in descending order from bottom to top. For example, "No Certification" would take precedence over "Conditionally Certified" and "Certified with Recommendations". "Conditionally Certified" would take precedence over "Certified with Recommendations".

For example:

A **Narkomed Standard** or **Compact** could have recommendation numbers 10, 11 and 15 apply.
10 - Breathing pressure monitor is not interfaced with anesthesia machine system power switch or ventilator. Code is B.

11 - Ventilator is equipped with descending bellows and lacks CO2 monitoring. Code is D1.

15 - Vaporizer mounting using a vaporizer selection valve. Code is B.

The correct certification for this machine is D, which means "NO CERTIFICATION", and additional recommendation I applies.

A **Narkomed AMIII** or **Narkomed II** could have recommendations 6 apply, and failure code 46.2.
6 - Oxygen ratio monitor instead of oxygen ratio monitor controller. Code D II.

46.2 - O2 analyzer malfunction. Code D.

Correct certification for this machine is D, which means "NO CERTIFICATION".

A **Narkomed 2A** could have recommendation numbers 12, 14 and 27 apply.

12 - No integrated exhaled CO2 monitor with user adjustable alarm limits. Code D1.

14 - CO2/Agent monitor exhaust port is not properly connected to waste gas disposal system. Code B.

27 - Ventilator bellows assembly has a PEEP valve assembly. Code B.

Correct certification for this machine is D1, which means "NO CERTIFICATION", and additional recommendation I applies.

A **Narkomed 3** could have recommendation number 21 and failure code 61.1 apply.

21 - No ventilator pressure limit control. Code is B.

61.1 - Enflurane agent is unavailable to test. Code is BC.

Correct certification for this machine is BC, which means CONDITIONALLY CERTIFIED WITH RECOMMENDATIONS.

A **Narkomed 4** could have recommendation numbers 14 and 21 apply.

14 - CO2/Agent monitor exhaust port is not properly connected to the waste gas scavenger. Code B.

21 - No ventilator pressure limit control. Code B.

The correct certification for this machine is B, which means "CERTIFIED WITH RECOMMENDATIONS".

A **Narkomed 2B, 2C** or **GS** could have recommendation 30 apply.

30 - Anesthesia machine is equipped with inhalation anesthesia vaporizers without an agent analyzer in the breathing system. Code B.

The correct certification for this machine is B, which means "CERTIFIED WITH RECOMMENDATIONS".

A **Narkomed 6000** could have no NAD recommendations or failure codes apply. The correct certification level for this machine is Code A, "CERTIFIED".

Code D, which means "NO CERTIFICATION", also means the machine shall not receive a Periodic Manufacturer's Certification label. The machine shall also receive a "WARNING - This System is Not Certified" label, P/N 4114857. This label shall be placed at a prominent location on the right side of the machine after all other previous PM and "Vigilance Audit® Validation" labels have been removed.

PM Certification Procedure for Fabius Tiro Anesthesia System

1. Use the PM Certification form P/N 4118460 for the Fabius Tiro Anesthesia System.
2. Completely fill in the header information.
3. Perform the vapor concentration test on all Dräger vapor vaporizers every six months in accordance with SP00073 at a six month maximum interval. Perform the vaporizer concentration test on all Desflurane vaporizers in accordance with SP00091 for fixed mount vaporizers and SP00189 for user removable D-tec vaporizers at a six month maximum interval. For every vaporizer tested, fill out a "VAPOR VAPORIZER CALIBRATION CHECK" label (part # S010016). Information on this label shall include your signature, type of agent, date tested, test results @ 1%, 2.5%, 4% for H, E, I, or S vaporizers, or @ 4%, 6%, 10%, 12%, 16% for Desflurane vaporizers, and a PASS or FAIL indication. This label shall be attached to the upper right side of the vaporizer. If vaporizer fails the concentration verification, internal leak, or exclusion system tests, check "NO" in the "RECOMMENDED FOR USE" section on the PM Certification form.

Place a "CAUTION DO NOT USE" label (part # 4114327) on the vaporizer, and issue a departmental alert. The TSR shall also seek permission from the equipment operator to remove the failed vaporizer from the machine and install a replacement vaporizer or an adapter block onto the mount. All nonfunctional Dräger vaporizers must be removed from service for machine to receive certification.

4. Proceed with PM Certification in accordance with Section 8.0. If any tests fail refer to the "Failure Codes" listing in NAD Recommendations Guidelines Index, P/N S010250, to determine correct certification level starting point. Failure codes shall be documented on the "RECOMMENDATIONS / GENERAL COMMENTS" section of the PM Certification form and on the Executive Summary if applicable. If a test fails that has not been identified by the "Failure Codes" list, consult with Draeger Medical, Inc. to assess the proper certification level.

PM Certification Procedure for Fabius Tiro Anesthesia System

5. Based on the "EQUIPMENT CONDITION" inspect the machine for any "NAD RECOMMENDATIONS" that would apply. Use the Fabius Tiro section of the "NAD RECOMMENDATION GUIDELINES INDEX", P/N S010250. Note all applicable NAD recommendations on the Executive Summary.

NOTE: If using a carbon form, indicate the Equipment Condition number and to see reverse side under "RECOMMENDATIONS / GENERAL COMMENTS" section of the form.

6. Determine the correct certification level of the machine based on the combined lowest common denominator of "Equipment Conditions" and "Failure Codes". If the machine is at least conditionally certified, fill out the "PM CERTIFICATION" label. Check the box(s) on the validation label where appropriate. Write the month and year, (six months from date of PM Certification) next to "NEXT VISIT DUE:." If certification level is "D", machine shall not receive a "PM CERTIFICATION" label. Any machine not receiving a PM Certification label shall receive a "WARNING NOT CERTIFIED" label. This label shall be placed at a prominent location on the left side of the machine after all other previous PMS and Vigilance Audit Validation labels have been removed.
7. In the "CERTIFICATION LEVEL" section of the PM Certification form, record the last visit certification level, the current certification level and the next visit due month and year, (six months from date of PM Certification) in the spaces provided.
8. If applicable, remove the previous PM CERTIFICATION VALIDATION label and attach the new label in a prominent location on the rear of the anesthesia machine.
9. Check the appropriate boxes on the "PM CERTIFICATION NOTICE" label, (part # S010011). If the machine is not certified, the last box of this notice label shall be checked. Attach this notice to the flow shield of the anesthesia machine.
10. Have the customer sign each PM Certification form or the Executive Summary, and review the equipment conditions and recommendations with the customer.
11. Return the top copy to Draeger Medical, Inc. Service Department, keep middle copy for service organization records, give bottom copy to customer.

NOTE: The following procedure will require making adjustments and calibrations with the possibility of removal of accessories and/or additional monitoring equipment. Therefore, take note of these devices for reconfiguration and/or reinstallation after the PMS procedure is complete.

8.1 Electrical Safety

(✓) 8.1.1 Protective Ground Continuity

8.1.1.1 Turn the System Power switch to STANDBY.

Plug the unit into the safety analyzer, and plug the power cord of the safety analyzer into an AC receptacle. Turn analyzer to ON position.

NOTE: The BIOTECH 501 PRO will automatically test the source outlet for open ground (or ground resistance of 31 Ω or higher), reverse polarity, open neutral and open line. (The latter two conditions will prevent the analyzer from powering up.)

8.1.1.2 Set the safety analyzer function switch to the GROUND WIRE RESISTANCE position. Attach the test lead to the SINGLE LEAD connector of the analyzer. Connect the other end of the red test lead to the ground socket of the front panel outlet on the safety analyzer. Verify a displayed resistance of 0.000 Ω or, if necessary, press the CALIBRATE key on the front panel of the analyzer to zero the device.

8.1.1.3 Set the safety analyzer GROUND switch to NORMAL. Set the POLARITY switch to OFF.


8.1.1.4 The safety analyzer shall indicate 0.1 Ω or less with its test lead applied to the following points:

- Cylinder yoke (if applicable)
- Serial port cable retaining nut on processor
- Vaporizers
- Push/pull handle
- GCX rails

(✓) 8.1.2 Circuit Isolation

NOTE: This test is not applicable on machines with serial numbers ≤ 10240 .

8.1.2.1 With a multimeter set to its highest resistance range, check for continuity between the serial port cable retaining nut on processor and the O2 sensor connector pin shown in the illustration. There shall be no continuity between these points.



O2 SENSOR

(✓) 8.1.3 Auxiliary Outlet Strip (if applicable)

NOTE: This test will check the auxiliary strip outlets for fault conditions such as open ground ($>31 \Omega$), reverse polarity, open line and open neutral. This is done each time the BIOTECH 501 PRO is powered up and allowed to cycle through its self test.

8.1.3.1 Shut off and unplug the safety analyzer. Remove the anesthesia machine plug from the analyzer, and plug it into the same outlet that was being used by the analyzer.

8.1.3.2 Plug the safety analyzer into the first outlet to be tested, and turn the analyzer power switch ON. Allow the analyzer to cycle through its Auto Test sequence. If no wiring fault is indicated, shut off the analyzer and move its plug to the next outlet. Test this outlet in the same manner, and continue until all convenience outlets on the auxiliary strip outlets are tested.

8.1.4 Chassis Leakage Current

8.1.4.1 Turn the anesthesia machine System Power switch to ON and set the safety analyzer to the CHASSIS LEAKAGE CURRENT position.

8.1.4.2 Attach the safety analyzer test lead to a rear GCX rail.


- (✓) 8.1.4.3 Record the total leakage current with the Polarity and Ground switches set as follows:

| | |
|--------|----------|
| Ground | Polarity |
| Open | Normal |
| Normal | Normal |
| Open | Reversed |
| Normal | Reversed |

Verify that the leakage current is 300 microamps or less in each of the switch positions (500 microamps or less for the 220/240 volt power supply option).

(✓) **8.2 System Diagnostics**

- 8.2.1 Connect the pipeline supply or open the cylinders.
- 8.2.2 Push the System Power switch to ON (⊙).
- 8.2.3 Verify that the following is displayed and system diagnostics indicate a "Pass" condition:

| SYSTEM DIAGNOSTICS | System Status |
|--|---------------|
| Watch Dog | Pass |
| System RAM | Pass |
| Program Memory | Pass |
| Video Test | Pass |
| Interrupts | Pass |
| A/D Converter | Pass |
| NV RAM | Pass |
| Serial Port | Pass |
| Clock | Pass |
| Speaker | Pass |
| Main Power | Pass |
| Battery | Pass |
|  Fabius Tiro SW XX | |

- (✓) 8.2.4 Record the machine software version on the header of the checklist form.

(✓) **8.3 Battery Circuit**

- 8.3.1 Press the Man Spont key, and press the rotary dial to confirm.
- 8.3.2 Unplug the AC power cord and verify the "Power Fail" message appears on the display within one minute.
- 8.3.3 Restore AC power to the machine and verify the "Power Fail" message disappears within one minute.

(✓) **8.4 Configuration**

- 8.4.1 Press the Standby key, then press the rotary key to confirm.
- 8.4.2 Press the Setup key, and press the rotary dial to confirm.

Default Settings:

- 8.4.3 Verify all the parameters listed in the 'Default Settings' screen are adjustable.
- 8.4.4 Adjust the Alarm Volume to its highest setting and press the rotary dial to confirm.
- 8.4.5 Press the rotary dial once to return to the Setup Screen.

Configuration:

- 8.4.6 Using the rotary dial, scroll to Configuration and press the rotary dial to confirm.
- 8.4.7 Verify all parameters listed in the 'Configuration Settings' screen are adjustable and correct per local configurations. Press to confirm two times to exit Setup mode.

8.5 Service Data

- 8.5.1 Press and hold the Home and Standby keys, and press the rotary dial. Verify the System Service screen appears.
- (✓) 8.5.2 Verify on-screen serial number with machine serial number located on rear of unit, and record on header of PMS form.
- (✓) 8.5.3 Record the Last Service Date, Hours Run Since Last Service, Total Hours Run, and Total Ventilator Hours Run on the PMS form.
- 8.5.4 Press the rotary knob to go to the Main Service screen.
- 8.5.5 Scroll to Preventative Maintenance, then to Activate.

- 8.5.6 Change the Preventative Maintenance date to reflect next PM due date. Press the rotary dial once to return to the main service screen.
- 8.5.7 Inspect the condition of the ventilator (piston diagram) O-ring. Is the O-ring in good condition? (Y)
- 8.5.8 If applicable, install the appropriate PM kit(s) and record in the Parts Replacement section of the PM form. See [“Fabius Tiro Parts Replacement Schedule:” on page 3.](#)

8.6 Calibrations

Fresh Gas Flow:

- 8.6.1 Scroll to Calibration and press to confirm.
- 8.6.2 Scroll to Fresh Gas Flow and confirm.
- 8.6.3 Follow on-screen instructions and perform Fresh Gas Flow calibration.
- (✓) 8.6.4 Verify Stored Zero indicates Pass (Y)
- 8.6.5 Press the Exit key.

Pressure:

- 8.6.6 Scroll to Pressure and confirm.
- (✓) 8.6.7 Follow the on-screen instructions and perform the Zero calibration. Reconnect pressure line and ventilator hose to breathing system.
- 8.6.8 Press the Exit key, then scroll to O2 and confirm.

O2 Offset:

- (✓) 8.6.9 Follow the on-screen instructions and perform the O2 zero calibration. Reassemble O2 housing and install the inspiratory dome plug in inspiratory dome.

PEEP:

NOTE: Before performing the PEEP valve calibration, a valid pressure calibration must be performed. Otherwise, an 'Inv pres. cal' message will appear.

- 8.6.10 Press the Exit key; scroll to PEEP and confirm.
- (✓) 8.6.11 Follow the on-screen instructions to calibrate PEEP. Verify PEEP calibration indicates PASS (Y).
- 8.6.12 Press Exit, then confirm (x2).

8.7 Site Configurations

8.7.1 Scroll to Configure and confirm.

O2 Position (virtual flowtubes):

- (✓) 8.7.2 Scroll to Flowmeter and confirm. Scroll to O2 Position and confirm. Select position of virtual flowmeter (left or right) in accordance with local requirements or customer demand. Press to confirm. record position on PMS form.

Gas Selection:

- (✓) 8.7.3 Press to confirm, scroll to Gas Selection and confirm. Verify Gas Selection type (2 or 3) is in accordance with machine configuration. Press to confirm (x2). Record selection type on PMS form.

Flowtube Resolution:

8.7.4 Scroll to Flowtube Res. and press to confirm. Scroll to change state and confirm. Scroll to change state in accordance with customer demand.

O2 Whistle:

- (✓) 8.7.5 Press to confirm, scroll to O2 whistle. Verify O2 whistle selection (enabled or disabled) is in accordance with local requirements or customer demand. Press to confirm (x2). Record position on PMS form.

Alarms:

- (✓) 8.7.6 Scroll to No Fresh Gas and confirm. Verify state of No Fresh Gas alarm is enabled. Press to confirm (x2). Record position on PMS form.
- (✓) 8.7.7 Scroll to Fresh Gas Low Alarm. Verify Fresh Gas Low Alarm position (enabled or disabled) is set to enabled position. Press to confirm (x2). Record position on PMS form.
- (✓) 8.7.8 Scroll to Threshold Low and confirm. Verify state of Threshold Low Alarm is enabled. Press to confirm (x2). Record position on PMS form.

Pressure:

- (✓) 8.7.9 Scroll to Ambient Pressure and confirm. Turn rotary dial to adjust. Press to confirm (x2). Record pressure setting on PMS form.

NOTE: Using the table found on the following page, set the ambient pressure (mbar) in accordance with the local elevation. To verify the local elevation, contact the nearest regional airport, or contact:

DMI - Technical Support - 1-800-543-5047
DMT - Tech Line - 49-451-882-4222

- 8.7.10 Scroll to Plateau/Mean Dis and verify setting is in accordance with customer demand. Press to confirm (x3).

- 8.7.11 Scroll to the 'System Settings' screen and confirm. Verify 'Model Type' listed is consistent with actual unit.

NOTE: Machines from the factory are configured with the appropriate model type. Other than performing a software download or PCB replacement, the model type should not be changed. Changing the model type will disable some ventilator options and will require configuration.

Serial Ports:

- (✓) 8.7.12 Scroll to Serial Ports and confirm. Verify parameters are adjustable to the following protocols. If necessary, set the protocols for any third party monitoring device connected to the machine. Refer to the Fabius Tiro Operator's Manual and third party Operator's Manual.

| | |
|------------|--------------------------------------|
| Baud Rate: | 1200, 2400, 4800, 9600, 19200, 38400 |
| Parity: | NONE, ODD, EVEN |
| Stop Bits: | 1, 2 |
| Data Bits: | 7, 8 |
| Protocol: | Vitalink - Medibus |

NOTE: It is important to ensure that communication protocols selected on each host and external device are correct. Vitalink and Medibus protocols are similar and if not set identically on each device, inaccurate data may be displayed on the remote device.

- 8.7.13 Press Standby key to exit.

| | | |
|--------------------|-------------------|----------------------------------|
| FABIUS TIRO | APPENDIX A | PMS PROCEDURE (continued) |
|--------------------|-------------------|----------------------------------|

| Elevation Range in Feet | Elevation Range in Meters | Barometric Pressure Setting |
|-------------------------|---------------------------|-----------------------------|
| -200 - 199 | -60 - 60 | 1013 |
| 200 - 599 | 61 - 182 | 1000 |
| 600 - 999 | 183 - 304 | 985 |
| 1000 - 1399 | 305 - 426 | 970 |
| 1400 - 1799 | 427 - 548 | 960 |
| 1800 - 2199 | 549 - 670 | 945 |
| 2200 - 2599 | 671 - 792 | 930 |
| 2600 - 2999 | 793 - 914 | 920 |
| 3000 - 3399 | 915 - 1036 | 905 |
| 3400 - 4199 | 1037 - 1158 | 890 |
| 3800 - 4199 | 1159 - 1280 | 880 |
| 4200 - 4599 | 1281 - 1401 | 865 |
| 4600 - 4999 | 1402 - 1523 | 855 |
| 5000 - 5399 | 1524 - 1645 | 840 |
| 5400 - 5799 | 1646 - 1767 | 830 |
| 5800 - 6199 | 1768 - 1889 | 820 |
| 6200 - 6599 | 1890 - 2011 | 805 |
| 6600 - 6999 | 2012 - 2133 | 795 |
| 7000 - 7399 | 2134 - 2255 | 785 |
| 7400 - 7799 | 2256 - 2377 | 770 |
| 7800 - 8199 | 2378 - 2499 | 760 |
| 8200 - 8999 | 2500 - 2620 | 750 |
| 8600 - 8999 | 2621 - 2742 | 740 |
| 9000 - 9399 | 2743 - 2864 | 730 |

(✓) **8.8 Scavenger**

8.8.1 AGS Scavenger

AGS Scavenger Cleaning

- 8.8.1.1 Remove all scavenger hoses one at a time and drain all accumulated moisture. Inspect all scavenger hoses for deterioration and replace any worn hoses.
- 8.8.1.2 Disconnect the hospital vacuum source from the scavenger.
- 8.8.1.3 Remove the reservoir canister from the scavenger by partial turn counter-clockwise of the canister.
- 8.8.1.4 Remove and inspect the silencer; replace if needed.
- 8.8.1.5 Remove the flowmeter from the scavenger by turning the mounting nut counter-clockwise. Inspect the tube and clean with compressed air if needed.
- 8.8.1.6 Reassemble the scavenger assembly, and reactivate the vacuum source.

AGS Scavenger Functional Test

- 8.8.1.7 Verify vacuum hose is connected to one of the two gas supply connections.
- 8.8.1.8 Verify the additional gas surplus connection is plugged.
- 8.8.1.9 Verify all hose connections to the scavenger are complete at all destinations.
- 8.8.1.10 Verify vacuum waste disposal system at manifold is active.
- 8.8.1.11 Reinstall canister and adjust scavenger flow control valve; verify float moves freely between the upper and lower limits.

8.8.2 Passive Scavenger

Passive Scavenger Cleaning

- 8.8.2.1 Inspect all scavenger hoses for signs of wear and deterioration. Replace any worn hoses.
- 8.8.2.2 Remove the anti-occlusion cage from the passive scavenger by unscrewing it.
- 8.8.2.3 Inspect the filter.

- 8.8.2.4 If necessary, remove the filter for cleaning. Brush any accumulated lint or dust off the filter. The filter can be further cleaned with a low flow of clean air or oxygen.
- 8.8.2.5 Place the filter back, making sure that it lays flat on the valve port orifice.
- 8.8.2.6 Reinstall the anti-occlusion cage on the scavenger body, making sure the filter is properly seated.

Passive Scavenger Functional Test

- 8.8.2.7 Check for moisture accumulation in the breathing and scavenger hoses. Remove and moisture found.
- 8.8.2.8 Short-circuit the COSY breathing system inspiratory and expiratory valves with a 2 mm breathing hose.
- 8.8.2.9 Install a breathing bag on the COSY breathing system.
- 8.8.2.10 Set the COSY breathing system APL valve to SPONT.
- 8.8.2.11 Open the oxygen flow control valve to a flow of 10 L/min and occlude the passive scavenger exhaust port connection.
- 8.8.2.12 After the breathing bag inflates, the absorber systems' breathing pressure gauge must indicate a pressure of less than 5 cm H₂O.
- 8.8.2.13 Systems that are not equipped with a pressure gauge, observe the generated pressure waveform, set the threshold to 5 cm H₂O. The pressure waveform shall not rise and remain above the dotted threshold trace during this performance test.

8.9 Breathing System

- (✓) 8.9.1 Breathing System Inspection
 - (✓) 8.9.1.1 Record the serial number (located on the side of the COSY) on the PMS form.
 - 8.9.1.2 Remove the inspiratory and the expiratory valve domes.
 - 8.9.1.3 Is there a broken or bent pin on the valve assembly?
Inspiratory ___ (N) Expiratory ___ (N)
 - 8.9.1.4 Is there a broken pin on the valve domes?
Inspiratory ___ (N) Expiratory ___ (N)
 - 8.9.1.5 Is the valve disc in good condition?
Inspiratory ___ (Y) Expiratory ___ (Y)
 - 8.9.1.6 Is there excessive wear on the valve craters?
Inspiratory ___ (N) Expiratory ___ (N)
 - 8.9.1.7 Are the valve dome washers in good condition? ___ (Y)
 - 8.9.1.8 Reinstall the inspiratory and expiratory valve domes.
 - 8.9.1.9 Un-screw the 22 mm expiratory port. Is the O-ring inside the port and 22mm taper in good condition? ___(Y)
 - 8.9.1.10 Remove the flow sensor from the breathing system housing. Is the O-ring inside housing in good condition? ___(Y)
 - 8.9.1.11 Reinstall the flow sensor and expiratory port.
 - 8.9.1.12 Remove the canister and inspect the canister, gaskets, and condition of soda lime. Are these components in good condition? ___(Y)
 - 8.9.1.13 Remove the breathing pressure gauge from its mount and inspect the O-rings in the mount and on the 90° hose connection - if applicable.
 - 8.9.1.14 Is the breathing pressure gauge in good condition and set at zero (0)? ___ (Y) (if applicable)
 - 8.9.1.15 Inspect the O-rings on the inspiratory dome plug assembly and O₂ sensor housing. Are the O-rings in good condition? ___(Y)

- 8.9.1.16 Examine all pneumatic hoses connecting from the interface panel to the breathing system. Are the hoses kink free and in good condition? ___(Y)
- 8.9.1.17 Inspect the breathing circuit and manual breathing bag. Is the breathing circuit and manual breathing bag in good condition? ___(Y)
- 8.9.1.18 Inspect the ventilator hose and associated O-rings connected between breathing system and machine. Are the hose and O-rings in good condition? ___(Y)
- 8.9.1.19 Inspect the fresh gas hose connected to the breathing system. Is the fresh gas connector and washer in good condition? ___(Y)
- 8.9.1.20 Inspect APL valve labeling for legibility. Are all markings on APL valve easy to see and legible? ___(Y)
- 8.9.1.21 Inspect the scavenger gas connection on the breathing system. Is the scavenger connector in good condition? ___(Y)

8.9.2 Fresh Gas Leak

- 8.9.2.1 Push the System Power switch OFF ($\overset{\bullet}{\circ}$).
- 8.9.2.2 Remove the fresh gas connector from the breathing system.
- 8.9.2.3 Connect the fresh gas hose from the breathing system to the fresh gas leak test fixture (4113119 modified) or equivalent, via the fresh gas test connector 4117361.
- 8.9.2.4 Connect a digital pressure manometer to the fresh gas leak test device.
- 8.9.2.5 Apply 50 cm H₂O of pressure to the system.
- (✓) 8.9.2.6 After thirty (30) seconds, what is the pressure on the manometer? ___ (>40 cm H₂O)
- 8.9.2.7 Turn on the vaporizer to the first graduated marking.
- 8.9.2.8 Apply 50 cm H₂O of pressure to the system.
- (✓) 8.9.2.9 After thirty (30) seconds, what is the pressure on the manometer? ___ (>40 cm H₂O)
- 8.9.2.10 Turn off the vaporizer.
- 8.9.2.11 Remove the test equipment from the fresh gas hose.

- 8.9.2.12 Push the System Power switch ON (⊙).
- 8.9.2.13 Open the O₂ flow control valve to 5 L/min., purge the system for 5 seconds, then close the O₂ flow control valve.
- 8.9.2.14 Push the System Power switch OFF (⊙).
- 8.9.2.15 Reconnect the fresh gas connector from the machine to the freshgas outlet connector on the breathing system.


(✓) 8.9.3 Breathing System Leak

- 8.9.3.1 Verify the System Power switch is On.
- 8.9.3.2 Press the Leak/Compl Test softkey.
- 8.9.3.3 Close all flow control valves.
- 8.9.3.4 Perform the on-screen instructions.
- 8.9.3.5 Verify all applicable tests indicate a Pass condition.

8.9.4 APL Valve Verification

- 8.9.4.1 Connect a digital manometer between the inspiratory and expiratory ports of the breathing system (using two test terminals and modified 4113119 fresh gas leak test device). See following illustration.
- 8.9.4.2 Attach a breathing bag to the bag terminal connector on the breathing system and install a test cap on the ventilator port.

SV88055



FRESH GAS LEAK TEST DEVICE
4113119 (MODIFIED)

| | Lever Style APL Valve | Rotary Knob Style APL Valve |
|-----|--|--|
| | 8.9.4.3 | Verify the fresh gas hose is connected to the breathing system. |
| | 8.9.4.4 | Verify the APL valve is in the MAN position and set to 10 cmH ₂ O. |
| | 8.9.4.5 | Activate the flush to inflate the reservoir bag, then release. |
| (✓) | 8.9.4.6 | Adjust the O ₂ and Air flows to 10 L/min. (20 L/min. total). Verify pressure remains between 8 and 12 cmH ₂ O. |
| (✓) | 8.9.4.7 | Repeat Steps 8.9.4.4 thru 8.9.4.6 with the following settings, and verify APL valve accuracy: |
| | APL Valve Setting | APL Valve Setting |
| | Spec (cm H ₂ O) | Spec (cm H ₂ O) |
| | Low High | Low High |
| | 30 27 33 | 40 34 46 |
| | 50 45 55 | |

8.9.4.8 Close the O₂ and Air flow control valves.

8.9.4.9 Remove the digital manometer and associated test equipment from the inspiratory & expiratory valve ports and bag terminal.


8.9.5 Inhalation and Exhalation Valves

Inhalation:

8.9.5.1 Install test terminal on the 22 mm inspiratory connector.


8.9.5.2 Connect 0 - 250 cc flowmeter (S000081) with a test terminal as shown, to the ventilator port of the breathing system.

SV02002




- 8.9.5.3 Remove the freshgas hose from the breathing system and connect it to the freshgas leak connector (4117163).
- 8.9.5.4 Connect the fresh gas leak connector to the fresh gas leak test device (4113119 modified) as shown.
- 8.9.5.5 Increase O₂ flow to maintain 30 cmH₂O on the digital manometer.
- (✓) 8.9.5.6 What is leakage as shown on 0 - 250 cc flowmeter? (<60 cc/min.)
- 8.9.5.7 Remove all test equipment from breathing system and disconnect digital manometer from fresh gas hose.
- Exhalation:
- 8.9.5.8 Occlude ventilator and inspiratory port, each with test cap M33972.
- 8.9.5.9 Interconnect digital manometer and 0-250 cc test flowmeter to bag terminal as shown.
- 8.9.5.10 Connect Fresh gas hose from machine to fresh gas port of breathing system.
- 8.9.5.11 Slowly increase O₂ flow to build up and maintain 30 cmH₂O/mbar on digital manometer.
- NOTE:** Verify pressure on digital manometer. Pressure decay is not visible on system cmH₂O pressure gauge.
- (✓) 8.9.5.12 Verify leakage is <60 cc/min. as indicated on test flowmeter.
- 8.9.5.13 Remove all test equipment except for test cap on inspiratory port.

SV02003




8.9.6 Reverse Flow Test (Fresh Gas Decoupling Valve)

- 8.9.6.1 Verify the breathing bag port, fresh gas connection and expiration port connections are open and not occluded.
- 8.9.6.2 Connect a digital pressure meter between the COSY pressure port and the pressure hose from the machine as shown.
- 8.9.6.3 Verify occlusion of the inspiration port with test cap M33972.
- 8.9.6.4 Connect the fresh gas hose to the ventilator port on the COSY with the test connector as shown.



- 8.9.6.5 Verify the APL valve is in the MAN position and set to 70 cmH₂O.
- 8.9.6.6 Slowly increase O₂ flow to maintain 40 cmH₂O on the digital pressure meter.
- (✓) 8.9.6.7 Verify the flow on the 0-250 cc test flowmeter indicates less than 10 mL/min. (0.01 L/min.).
- 8.9.6.8 Close the O₂ flow control valve.
- 8.9.6.9 Remove the occlusion plugs and test equipment from the breathing system.
- 8.9.6.10 Reconnect fresh gas hose to breathing system; connect ventilator hose to vent port and breathing system.
- 8.9.7 Leakage Control Port
 - 8.9.7.1 Disconnect the APL and PEEP control hoses from the interface panel at rear of machine.
 - 8.9.7.2 Connect the PEEP control hose to the input of the 0-250 cc test flowmeter.

SV02005



- 8.9.7.3 Short circuit the inspiratory and expiratory ports, attach a test terminal to the breathing bag port and attach a digital manometer to the test terminal.
- 8.9.7.4 Increase O₂ flow to maintain 40 cmH₂O on the digital manometer.
- (✓) 8.9.7.5 Verify the test flowmeter indicates < 10 cc/min.
- 8.9.7.6 Remove the 0-250 cc test flowmeter, digital manometer and short circuit hose from machine.
- 8.9.7.7 Return all hose connections to their normal positions.

8.10 Yokes & Gauges

- (✓) 8.10.1 Yokes & Check Valves (if applicable)
 - 8.10.1.1 Turn the System Power switch Off.
 - 8.10.1.2 Disconnect all pipeline hoses and close all cylinder valves.
 - 8.10.1.3 Remove cylinder or yoke plug from each yoke assembly.
 - 8.10.1.4 Do all the yoke handles adjust smoothly? ___ (Y)
 - 8.10.1.5 Are the two (2) yoke pins installed securely in each yoke? ___ (Y)
 - 8.10.1.6 Is there only one (1) cylinder washer on each yoke assembly? ___ (Y)
 - 8.10.1.7 Is the proper gas I.D. label affixed to each yoke assembly? ___ (Y)
 - 8.10.1.8 Attach a cylinder to each yoke assembly, open the cylinder valve, let the pressure stabilize, close the cylinder valve, and remove the cylinder from the yoke assembly.
 - 8.10.1.9 Does the yoke check valve assembly prevent the escape of excessive pressure? ___ (Y)
 - 8.10.1.10 Attach the cylinders to the yokes.
- (✓) 8.10.2 Cylinder Gauges (if applicable)
 - 8.10.2.1 Are the pressure gauges correct for the gases indicated by the flowmeters? ___ (Y)
 - 8.10.2.2 Bleed all pressure from the cylinder circuits using the flow control valves.

8.10.2.3 Are the cylinder gauges at zero (0) PSI? ___ (Y)

8.10.2.4 Open the cylinder valves.

8.10.2.5 Do the cylinder pressure gauges respond properly? ___ (Y)

8.11 Gas Inlet Regulator Output

8.11.1 O2 Inlet Regulator


8.11.1.1 Trolley or Wall Mount - Service Access

8.11.1.1.1 Fully extend the writing tray to its out position.

8.11.1.1.2 Remove the top drawer by releasing the locking tabs located in the drawer slides (x2). (Trolley Mount Only)


8.11.1.2 Bleed all cylinder and pipeline pressures. Disconnect all cylinders and pipeline hoses.

8.11.1.3 Remove the pneumatic hose from the O2 outlet of the gas inlet block and interconnect a digital manometer as shown.



- 8.11.1.4 Reconnect the O2 pipeline hose to the O2 pipeline inlet connector and pressurize the O2 supply.
- 8.11.1.5 Set the O2 flow to 4 L/min.
- (✓) 8.11.1.6 After the digital manometer stabilizes, what is the regulator output pressure? ___psi (28 - 32 psi)
- 8.11.1.7 Deplete O2 pressure from the pipeline supply.
- 8.11.1.8 Close the flow control valve and disconnect the O2 pipeline hose from the inlet.
- 8.11.1.9 Remove test equipment and reconnect the O2 pneumatic hose leading from the O2 flow control valve to the O2 connector on the inlet block.
- 8.11.2 N2O Inlet Regulator
 - 8.11.2.1 Disconnect the N2O pneumatic hose from the outlet connector on the gas inlet block that connects to the N2O flow control valve.
 - 8.11.2.2 Interconnect a digital manometer between the N2O outlet connector and the hose removed in the previous step. See illustration on previous page.
 - 8.11.2.3 Reconnect the O2 and N2O pipeline hoses to the inlet block and activate the pipeline supplies.
 - (✓) 8.11.2.4 Set the O2 and N2O flows to 4 L/min. After the digital manometer stabilizes, what is the N2O regulator output pressure? ___psi (28 - 32 psi)
 - 8.11.2.5 Deplete the O2 and N2O pipeline pressures.
 - 8.11.2.6 Close the O2 and N2O flow control valves and disconnect both pipeline hoses from the inlet block.
 - 8.11.2.7 Remove test equipment and reconnect the N2O hose from the flow control valve to the N2O outlet on the inlet block.
- 8.11.3 Pipeline Check Valves
 - O2 Pipeline Check Valve:
 - 8.11.3.1 Attach a Swagelock fitting (P/N 4109318) to the inlet of the S000081 test flowmeter.

- 8.11.3.2 Detach the hose from the Pressure Test adapter (P/N 4114807) gauge assembly and attach it to the inlet of the test flowmeter as shown.
- 8.11.3.3 Attach the appropriate O2 adapter to the O2 connector on the inlet block, and connect the other end of the hose with the quick disconnect fitting to the adapter as shown.



- (✓) 8.11.3.4 Open an O2 (reserve) cylinder valve. What is the flow as indicated on the test flowmeter? ___cc ≤ 5 cc/min.

N2O Pipeline Check Valve (if applicable)

- 8.11.3.5 Attach the appropriate N2O adapter to the N2O connector on the inlet block, and transfer the end of the hose with the quick disconnect fitting to the N2O adapter.

- (✓) 8.11.3.6 Open the N₂O (reserve) cylinder. What is the flow as indicated on the test flowmeter? ___cc ≤ 5 cc/min.


8.12 Cylinder Regulator & Pipeline Gauges

Minimum cylinder pressure requirements for this test are:

N₂O: 600 psi;
O₂, Air: 1000 psi.

- (✓) 8.12.1 N₂O Cylinder Regulator (if applicable)
- 8.12.1.1 Configure test gauge 4114807 using the appropriate adapters between the pipeline inlet of machine and central supply hose. A typical test connection is shown on following illustration (O₂ illustrated).
 - 8.12.1.2 Connect the test fixture hose to the machine's nitrous oxide pipeline inlet.
 - 8.12.1.3 Does the side panel correctly identify the nitrous oxide inlet? ___(Y)
 - 8.12.1.4 Connect the nitrous oxide pipeline supply hose to the test fixture.
 - 8.12.1.5 Open the nitrous oxide and the oxygen cylinder valves.
 - 8.12.1.6 Set the oxygen and nitrous oxide flows to 4 L/min.
 - 8.12.1.7 Depress the push button on the test device.
- (✓) 8.12.1.8 Release the push button. After the pressure decay stabilizes, what is the regulator output pressure? ___psi (32 - 40)
- NOTE:** If a pressure decrease does not occur, either the hospital's supply pressure is too low or the regulator pressure is set too high.
- (✓) 8.12.2 N₂O Pipeline Gauge Accuracy
- 8.12.2.1 Close the nitrous oxide cylinder valve and drain all nitrous oxide pressure.
 - 8.12.2.2 Depress and hold the push button on the test device.
 - 8.12.2.3 Does the nitrous oxide gauge respond properly? ___(Y)
- (✓) 8.12.2.4 After the pressure stabilizes, are the indicated pressures on the test gauge and the nitrous oxide pipeline pressure gauge within the normal operating range? 41 - 87 psi (non-USA), 50 - 55 psi (USA).

SV88025



8.12.3 N2O Pipeline Leak

8.12.3.1 Close the nitrous oxide flow control valve.

8.12.3.2 Release the test device push button.

(✓)

8.12.3.3 After 30 seconds, what is the pressure loss? ___psi (<5)

8.12.3.4 Remove the test equipment and reconnect the nitrous oxide pipeline hose.

(✓) 8.12.4 Air Cylinder Regulator

8.12.4.1 Configure the test gauge using the appropriate adapters between the pipeline inlet of machine and central supply hose.

8.12.4.2 Connect the test fixture hose to the machine's air pipeline inlet.

8.12.4.3 Does the side panel correctly identify the air inlet? ___(Y)

8.12.4.4 Connect the air pipeline supply hose to the test fixture.

8.12.4.5 Set the air flow to 4 L/min.

8.12.4.6 Depress the push button on the test device.

(✓) 8.12.4.7 Release the push button. After the pressure decay stabilizes, what is the regulator output pressure? ___psi (based upon cylinder pressure as given in the following table)

| Cylinder Pressure (psi) | Compensated Regulator output tolerances (-4/+2) |
|-------------------------|---|
| 2000 | 27 - 33 (*30 - 36) |
| 1800 | 28 - 34 (*31 - 37) |
| 1600 | 29 - 35 (*32 - 38) |
| 1400 | 30 - 36 (*33 - 39) |
| 1200 | 31 - 37 (*34 - 40) |
| 1000 | 32 - 38 (*35 - 41) |

* Canada settings

NOTE: If a pressure decrease does not occur, either the hospital's supply pressure is too low or the regulator pressure is set too high.

(✓) 8.12.5 Air Pipeline Gauge Accuracy

8.12.5.1 Close the air cylinder valve and drain all air pressure.

8.12.5.2 Depress the push button on the test device.

8.12.5.3 Does the air gauge respond properly? ___(Y)

8.12.5.4 After the pressure stabilizes, are the indicated pressures on the test gauge and the air pressure gauge within the normal operating range? 41 - 87 psi (non-USA), 50 - 55 psi (USA).

8.12.6 Air Pipeline Leak

8.12.6.1 Close the air flow control valve

8.12.6.2 Release the test device push button

(✓) 8.12.6.3 After thirty (30) seconds, what is the pressure loss? ___psi (<5)

8.12.6.4 Remove the test equipment and reconnect the air pipeline hose.

(✓) 8.12.7 O2 Cylinder Regulator

8.12.7.1 Configure a test gauge 4114807 using the appropriate adapters between the pipeline inlet of machine and central supply hose.

8.12.7.2 Connect the test fixture hose to the machine's oxygen pipeline inlet.

8.12.7.3 Connect the oxygen pipeline supply hose to the test fixture.

8.12.7.4 Does the back panel correctly identify the oxygen inlet? ___(Y)

8.12.7.5 Set the oxygen flow to 4 L/min.

8.12.7.6 Depress the push button on the test device.

(✓) 8.12.7.7 Release the push button. After the pressure decay stabilizes, what is the regulator output pressure? ___psi (based upon cylinder pressure as given in the following table)

| Cylinder Pressure (psi) | Compensated Regulator output tolerances (-4/+2) |
|-------------------------|---|
| 2000 | 27 - 33 (*30 - 36) |
| 1800 | 28 - 34 (*31 - 37) |
| 1600 | 29 - 35 (*32 - 38) |
| 1400 | 30 - 36 (*33 - 39) |
| 1200 | 31 - 37 (*34 - 40) |
| 1000 | 32 - 38 (*35 - 41) |

* Canada settings

NOTE: If a pressure decrease does not occur, either the hospital's supply pressure is too low or the regulator pressure is set too high.

- (✓) 8.12.8 O2 Pipeline Gauge Accuracy
- 8.12.8.1 Close the oxygen cylinder valve and drain all oxygen pressure.
 - 8.12.8.2 Depress the push button on the test device.
 - 8.12.8.3 Does the oxygen gauge respond properly? ___(Y)
 - 8.12.8.4 After the pressure stabilizes, are the indicated pressures in the test gauge and the oxygen pressure gauge within the normal operating range? 41 - 87 psi (non-USA), 50 - 55 psi (USA).
- 8.12.9 O2 Pipeline Leak
- 8.12.9.1 Release the test device push button
- (✓) 8.12.9.2 After thirty (30) seconds, what is the pressure loss? ___psi (<5)

8.13 High Pressure Leak

- 8.13.1 Oxygen High Pressure Leak
- 8.13.1.1 Verify the System Power switch is OFF.
 - 8.13.1.2 Open the oxygen cylinder valve.
 - 8.13.1.3 Let the pressure stabilize.
 - 8.13.1.4 Close the oxygen cylinder valve and remove the cylinder.
 - 8.13.1.5 Observe the oxygen cylinder pressure gauge.
- (✓) 8.13.1.6 After two (2) minutes, what is the pressure loss? ___ PSI (<50)
- 8.13.1.7 Attach the cylinder.
- 8.13.2 Nitrous Oxide or Air High Pressure Leak (if applicable)
- 8.13.2.1 Turn the System Power switch to ON.
 - 8.13.2.2 Open one (1) oxygen cylinder valve and one (1) nitrous oxide or air cylinder valve.
 - 8.13.2.3 Adjust the oxygen flow to 4 L/min.
 - 8.13.2.4 Let the pressure stabilize.


- 8.13.2.5 Close the nitrous oxide or air cylinder valve and remove the cylinder.
- 8.13.2.6 Observe the nitrous oxide or air cylinder pressure gauge.
- (✓) 8.13.2.7 After two (2) minutes, what is the pressure loss? ___ PSI (<50)
- 8.13.2.8 Attach the cylinder.
- 8.13.2.9 Close the oxygen flow control valve.

8.14 Oxygen Supply Failure Protection

- (✓) 8.14.1 Nitrous Oxide (if applicable)
 - 8.14.1.1 Open and close the oxygen cylinder valve.
 - 8.14.1.2 Open the nitrous oxide cylinder valve.
 - 8.14.1.3 Set the O₂ and N₂O flows to 4 L/min.
 - 8.14.1.4 Does the flow of nitrous oxide cease when the oxygen flow is depleted? ___ (Y)
 - 8.14.1.5 Connect the O₂ pipeline supply.
 - 8.14.1.6 Close the nitrous oxide cylinder valve and bleed the pressure from the circuit.
 - 8.14.1.7 Connect the N₂O pipeline supply.
 - 8.14.1.8 Disconnect the O₂ pipeline supply.
 - 8.14.1.9 Does the flow of nitrous oxide cease when the oxygen pressure is depleted? ___ (Y)
 - 8.14.1.10 Close the nitrous oxide flow control valve.
 - 8.14.1.11 Disconnect the N₂O pipeline supply.
- 8.14.2 Oxygen Supply Pressure Alarm
 - 8.14.2.1 If not already connected, connect the Pressure Test Assembly (P/N 4114807) with O₂ adapters between the central supply hose and the machine.
 - 8.14.2.2 Open and close an oxygen cylinder.
 - 8.14.2.3 Set the oxygen flow to 0.5 L/min.

- 8.14.2.4 Depress and release the test device push button.
- (✓) 8.14.2.5 What is the pressure on the test gauge when the "O₂ SUPPLY LOW!!!" message and associated red indicator appear? ___ PSI (16 - 24)
- (✓) 8.14.2.6 Does the "O₂ SUPPLY LOW!!!" message appear on the display? ___ (Y)
- 8.14.2.7 Bleed the remaining O₂ pressure from the system, then close the flow control valve.
- 8.14.2.8 Remove the test gauge from the machine and reconnect the pipeline supply hose.

SV88056




8.15 Flowmeters

- (✓) 8.15.1 Oxygen Flowmeter
 - 8.15.1.1 Open the O₂ cylinder valve.
 - 8.15.1.2 Is it possible to adjust the flow of oxygen to 10 L/min.? ____ (Y)
 - 8.15.1.3 Close the O₂ cylinder valve and bleed the pressure.
 - 8.15.1.4 Connect the O₂ pipeline supply, and verify the operation of the oxygen flowmeter.
 - 8.15.1.5 Is the correct flow control knob and label attached to the oxygen flow control valve? ____ (Y)
 - 8.15.1.6 Close the oxygen flow control valve.

- (✓) 8.15.2 Nitrous Oxide Flowmeter (if applicable)
 - 8.15.2.1 Set the oxygen flow to 4 L/min.
 - 8.15.2.2 Open the nitrous oxide cylinder valve.
 - 8.15.2.3 Is it possible to adjust the flow of nitrous oxide to 10 L/min.? ____ (Y)
 - 8.15.2.4 Close the nitrous oxide cylinder valve and bleed the pressure.
 - 8.15.2.5 Connect the N₂O pipeline supply, and verify the proper operation of the N₂O flow.
 - 8.15.2.6 Is the correct flow control knob attached to the N₂O flow control valve? ____ (Y)
 - 8.15.2.7 Close the oxygen and nitrous oxide flow control valves.

- (✓) 8.15.3 Air Flowmeter
 - 8.15.3.1 Connect the Air pipeline supply (if applicable).
 - 8.15.3.2 Is it possible to adjust the flow of air to 10 L/min.? ____ (Y)
 - 8.15.3.3 Close the air flow control valve and disconnect the Air pipeline supply.
 - 8.15.3.4 Is the correct flow control knob attached to the air flow control valve? ____ (Y)

- (✓) 8.15.4 Auxiliary Oxygen Flowmeter (if applicable)
- 8.15.4.1 Verify the auxiliary oxygen flowmeter flow control valve is closed.
 - 8.15.4.2 Connect a pressure manometer to the Aux O₂ outlet.
 - 8.15.4.3 Is there an increase in pressure? ____ (N)
 - 8.15.4.4 Slowly open the Aux O₂ flow control valve to apply a pressure of 50 cm H₂O to the manometer, then close control valve and start timer.
- (✓) 8.15.4.5 After 30 seconds, what is the pressure on the manometer? ____ (40 - 50 cm H₂O)
- 8.15.4.6 Remove the manometer test fixture from the outlet.
 - 8.15.4.7 Is it possible to adjust the flow over the full range of the auxiliary oxygen flowmeter? ____ (Y)
 - 8.15.4.8 Set the flow rate to 5 L/min.
 - 8.15.4.9 Hold the calibrated oxygen sensor at the auxiliary oxygen flowmeter outlet.
- (✓) 8.15.4.10 After the value stabilizes, what is the oxygen concentration? ____% (97-100)
- 8.15.4.11 Remove the oxygen sensor from the auxiliary oxygen flowmeter, and insert it into the inspiratory valve dome adapter.
 - 8.15.4.12 Close the auxiliary oxygen flow control valve.
- (✓) **8.16 Oxygen Monitor**
- 8.16.1 Press the Man Spont key and press the rotary knob to confirm.
 - 8.16.2 Verify the ventilator hose is connected to the breathing system.
 - 8.16.3 Disconnect the oxygen sensor cable from the Oxygen Sensor interface.
 - 8.16.4 The following message shall appear on the display: O₂ Sensor Fail!.
 - 8.16.5 Reconnect the oxygen sensor.
- NOTE:** Make sure that the sensor has stabilized in ambient air for several minutes.
- 8.16.6 Press the Setup key, then perform O₂ calibration per on-screen instructions.

- 8.16.7 After calibration is completed, verify the "O2 Sensor Fail" message disappears.
- 8.16.8 What is the oxygen concentration? ___% (21)
- 8.16.9 What is the low oxygen alarm default? ___% (20)
- 8.16.10 Press the  key and select the OXYGEN LOW alarm limit. Does the low alarm limit illuminate? ___ (Y)
- 8.16.11 Verify that the low alarm limit has a range from 18 to 99%. Adjust the alarm limit above current O2 monitor value.
- (✓) 8.16.12 INSP O2 LOW!!! shall appear on the display and the heading shall be flashing with a corresponding audible alarm.
- 8.16.13 Place the oxygen sensor into the inspiratory valve dome adapter and set the APL valve to 70 and MAN position. Attach a 12-inch hose to the inspiratory port and occlude the bag mount.
- 8.16.14 Set the oxygen flow to 4 L/min.
- 8.16.15 Verify that the INSP O2 LOW!!! message has cleared.
- 8.16.16 Select the OXYGEN HIGH alarm limit. Does the high alarm limit illuminate? ___ (Y)
- 8.16.17 What is the high oxygen alarm default? ___% (100)
- 8.16.18 Verify that the high alarm limit has a range from 100 to 19%.
- 8.16.19 Set the high alarm limit to below the correct O2 reading.
- (✓) 8.16.20 Does the message INSP O2 HIGH!! appear on the display? ___ (Y)
- 8.16.21 Return the high alarm limit to 100 and confirm.
- 8.16.22 The INSP O2 HIGH message shall disappear.
- (✓) 8.16.23 Within 3 minutes, what is the oxygen concentration? ___% (97-100)

8.17 Oxygen Concentrations

- (✓) 8.17.1 Oxygen + Nitrous Oxide Concentration (if applicable)
- 8.17.1.1 Verify the oxygen flow is at 4 L/min.
- 8.17.1.2 Depress the O₂ Flush button for 5 seconds.

- 8.17.1.3 Does the oxygen monitor read 97-100% after the value stabilizes? __ (Y)
- 8.17.1.4 Set the nitrous oxide flow to 2 L/min.
- (✓) 8.17.1.5 After the value stabilizes on the oxygen sensor, what is the O₂ concentration? ___% (62-72)
- 8.17.1.6 If measured value is within the range, close the N₂O flow control valve and proceed to Section [8.18.2](#).
- (✓) 8.17.1.7 If measured value is not within range, close the nitrous oxide flow control valve and adjust the oxygen flow until total flow meter reads 4 L/min. What is the corresponding display flow rate? ___ L/min (3.9 - 4.7)
- (✓) 8.17.1.8 Set nitrous oxide flow to 2.5 L/min and oxygen flow to 2.5 L/min on the displays. What is the corresponding flow visually approximated on the total flow meter? ___ L/min (4.5 - 5.5)
- 8.17.1.9 Close the nitrous oxide flow control valve.
- 8.17.2 Oxygen + Air Concentration
 - 8.17.2.1 Depress the O₂FLUSH button for 5 seconds.
 - 8.17.2.2 Does the oxygen monitor read 97-100% after the value stabilizes? __ (Y)
 - 8.17.2.3 Set the air flow to 2 L/min.
 - (✓) 8.17.2.4 After the value stabilizes, what is the O₂ concentration? ___% (71-77)
 - 8.17.2.5 Close the air and O₂ flow control valves.

8.18 SORC

- 8.18.1 Depress the O₂FLUSH for 5 seconds.
- 8.18.2 Set the O₂ and N₂O flow to 10 L/min.
- 8.18.3 Set the O₂ flow control valve to 0.8 L/min.
- (✓) 8.18.4 What is the oxygen concentration after the value stabilizes? ___% (≥23%)
- 8.18.5 Adjust the oxygen flow to 1.5 L/min.
- (✓) 8.18.6 What is the oxygen concentration after the value stabilizes? ___% (≥23%)

- 8.18.7 Adjust the oxygen flow to 2 L/min.
- (✓) 8.18.8 What is the oxygen concentration after the value stabilizes? ___% ($\geq 23\%$)
- 8.18.9 Adjust the oxygen flow to achieve 10 L/min. N₂O.
- (✓) 8.18.10 What is the oxygen concentration after the value stabilizes? ___% ($\geq 23\%$)
- (✓) 8.18.11 Reduce the O₂ flow to 0.5 L/min. Verify that the N₂O flow is greater than or equal to 0.6 L/min. ___ L/min.
- 8.18.12 Reduce O₂ flow to zero, then slowly increase O₂ flow until N₂O flow is present.
- (✓) 8.18.13 What is O₂ flow when N₂O starts to flow (0.1 - 0.2 L/min.)
- 8.18.14 Close the O₂ and N₂O flow control valves.
- (✓) **8.19 Pressure Monitor**
- 8.19.1 Verify unit is in Man Spont mode.
- 8.19.2 Disconnect the breathing pressure sensor line from the breathing system.
- 8.19.3 Connect a test pressure gauge and syringe to the breathing pressure sensor line.
- 8.19.4 Press the "Alarms" key.
- 8.19.5 What is the pressure apnea threshold default? ___ cm H₂O (8)
- 8.19.6 Verify that the pressure apnea threshold limit has a range from 5 to 30 cm H₂O.
- 8.19.7 Adjust the threshold limit to 10 cm H₂O.
- 8.19.8 Select the HIGH alarm limit, and confirm.
- 8.19.9 What is the high alarm limit default? ___ cm H₂O (40)
- 8.19.10 Verify that the high alarm limit has a range from 10 to 70 cm H₂O.
- 8.19.11 Set the high alarm limit to 40 cm H₂O, and confirm. Set the ventilator to Volume Control mode and confirm.
- 8.19.12 Increase the pressure to 25 cm H₂O, then decrease the pressure to 0 cm H₂O and start a stop watch.

- (✓) 8.19.13 What is the time when APNEA PRESSURE!! appears as a medium alarm tone with associated flashing yellow LED? ___ sec (26-34)
- (✓) 8.19.14 What is the time when the APNEA PRESSURE!!! appears as a high alarm with flashing red LED and corresponding audible alarms? ___ sec (56-64)
- NOTE:** With the threshold low alarm enabled a corresponding threshold low alarm shall appear.
- 8.19.15 After the APNEA PRESSURE!!! alarm is displayed as a high, slowly increase the test pressure.
- (✓) 8.19.16 At what pressure does the APNEA PRESSURE alarm deactivate? ___ cm H₂O (7-13)
- 8.19.17 Adjust the threshold to 18 cm H₂O.
- 8.19.18 Increase the pressure to 20 cm H₂O, maintain the pressure, and start a stopwatch.
- (✓) 8.19.19 What is the time when CONTINUOUS PRESSURE appears as a high alarm? ___ sec (12-18)
- (✓) 8.19.20 Decreasing the pressure slowly, what is the pressure at which the CONTINUOUS PRESSURE alarm deactivates? ___ cm H₂O (15-21)
- 8.19.21 Increase the test pressure.
- (✓) 8.19.22 At what pressure does the Airway Pressure High!!! alarm activate? ___ cm H₂O (38-42)
- 8.19.23 Bleed the pressure.
- 8.19.24 Using a syringe in place of the squeeze bulb, slowly create a sub-atmospheric pressure.
- (✓) 8.19.25 At what pressure does the Pressure Negative!!! alarm activate? ___ cm H₂O (-2 to -8)
- 8.19.26 Return the pressure to zero.
- 8.19.27 Does the PRESSURE NEGATIVE alarm deactivate? ___(Y)
- 8.19.28 Using the syringe, increase the pressure to 20 cmH₂O; then decrease the pressure to zero.
- 8.19.29 Press the Autoset soft key.
- 8.19.30 Does the threshold limit adjust to within four (4) cmH₂O of peak pressure? ___ (Y)

8.19.31 Press the Man Spont key and confirm.

8.19.32 Reconnect the breathing pressure line to the breathing system.

8.20 Ventilator

(✓)

8.20.1 Manual Ventilation

8.20.1.1 Cycle machine power-using On/Off switch and allow completion of the self-test.

8.20.1.2 Verify the ventilator switches to standby mode after completion of self-test.

8.20.1.3 Connect a test lung (8401892) to the breathing circuit Y-piece, and attach to breathing system.

8.20.1.4 Attach a breathing bag to the breathing bag connector of the breathing system.

8.20.1.5 Select Man Spont.

8.20.1.6 Set the fresh gas flow to 3 L/min.

8.20.1.7 Set the APL valve to MAN, 30 cmH₂O.

8.20.1.8 Verify that manual breathing is possible by manually squeezing the breathing bag. ____ (Y)

(✓)

8.20.2 Spontaneous Breathing

8.20.2.1 Set the APL valve to SPONT.

8.20.2.2 Verify that spontaneous breathing is possible with the test lung. ____ (Y)

8.20.2.3 Close O₂ flow control valve.

(✓)

8.20.3 Flow Sensor Zeroing

8.20.3.1 Select Standby mode and confirm.

8.20.3.2 Disconnect the expiratory hose from the breathing system.

8.20.3.3 Press the Calibrate Flow Sensor key and perform the on-screen instructions. At completion of calibration verify "Flow Calibration Completed" message appears.

- 8.20.3.4 Reconnect the breathing circuit to the inspiratory and expiratory ports of the breathing system.
- (✓) 8.20.4 Ventilator Delivery
- 8.20.4.1 Set APL to Man position and set to 30 cmH₂O.
- 8.20.4.2 Switch to Volume Control mode and set the fresh gas flow to 3 L/min.
- 8.20.4.3 Interconnect a test volumeter (P/N 2212300 w/4115087 connectors) between the expiratory limb of breathing system and breathing circuit.
- 8.20.4.4 Is Volume Mode displayed and flashing? ___ (Y)
- 8.20.4.5 Depress Flush momentarily to inflate bag.
- 8.20.4.6 Set: pmax = 70 mbar (cmH₂O)
- Vt = 380 mL
Freq = 12 BPM
TI:TE = 1:1
TIP:TI = 10%
PEEP = 0 cmH₂O
- 8.20.4.7 Press the rotary knob.
- 8.20.4.8 Is Volume Control displayed? ___ (Y)
- 8.20.4.9 Verify that ventilation starts. ___ (Y)
- 8.20.4.10 Does drive run quietly and smoothly ___ (Y)?
- (✓) 8.20.4.11 Verify display Vte = 300 - 450 mL (___ mL/min.)
- (✓) 8.20.4.12 Verify Minute Volume on test volumeter and MV on display are within 20% of each other.
- 8.20.5 PEEP Pressure Accuracy
- 8.20.5.1 After 10 breaths, press PEEP parameter button.

(✓) 8.20.5.2 Verify displayed values for the following PEEP settings:

| Setting | Displayed Value |
|---------|--|
| 0 | 0 + 2 cmH ₂ O |
| 10 | 10 ± 2 cmH ₂ O or 20% of setting, which ever is greater |
| 15 | 15 ± 3 cmH ₂ O or 20% of setting, which ever is greater |

8.20.5.3 Return PEEP setting to Zero and confirm.

8.20.6 Pmax Accuracy

8.20.6.1 Connect a test pressure gauge in line with the pressure port connector of breathing system.

8.20.6.2 Remove test lung and seal off Y-piece.

(✓) 8.20.6.3 Verify test gauge display and breathing gauge values for the following Pmax settings:

| Setting | Test Gauge Value |
|---------|---------------------------|
| 25 | 25 ± 5 cmH ₂ O |
| 40 | 40 ± 5 cmH ₂ O |
| 60 | 60 ± 5 cmH ₂ O |
| 70 | 70 ± 5 cmH ₂ O |

8.20.7 APL Valve Man/Spont

8.20.7.1 Set ventilator to Man Spont.

| | Lever APL Valve | Rotary Knob APL Valve |
|-----|--|---|
| | 8.20.7.2 Set O ₂ and Air Fresh Gas flows to 10L/min ea. (20 L/min. total) | |
| | 8.20.7.3 Set APL valve to MAN, 30 cmH ₂ O. | Set APL Valve to 40 cm H ₂ O |
| (✓) | 8.20.7.4 Verify pressure is 27 - 33 mbar (cmH ₂ O) as shown on digital manometer. | Verify pressure is 34 - 46 mbar (cmH ₂ O) as shown on digital manometer. |
| | 8.20.7.5 Set APL valve to SPONT. | |
| (✓) | 8.20.7.6 Verify pressure display is 0 - 3 mbar. | |

8.20.7.7 Set fresh gas flow to zero.

8.20.7.8 Reconnect test lung to breathing circuit and remove the test pressure gauge from the breathing system.

8.20.8 Pressure Limiting Valve Test

8.20.8.1 Set ventilator to Volume Control mode and confirm.


8.20.8.2 Set: $p_{max} = 40$ mbar (cmH₂O)

$V_t = 600$ mL
Freq = 6 BPM
TI:TE = 1:1
TIP:TI = 10%
PEEP = 0 mbar (cmH₂O)

8.20.9 Remove ventilator hose from ventilator and connect a digital manometer to the ventilator outlet port.

(✓) 8.20.10 Verify relief valve opens between 70 and 80 mbar (cmH₂O) when ventilator is in inspiratory phase.

SV88022



8.20.11 Auxiliary Air Valve Test

- (✓) 8.20.11.1 Close all flow control valves. Verify vacuum valve opens at -7.0 to -8.0 mbar (cmH₂O) pressure when ventilator is in expiratory phase.

8.20.12 Vacuum Pressure

8.20.12.1 Set Fabius Tiro to Volume mode.

8.20.12.2 Remove APL bypass hose from breathing system and attach to digital manometer.

- (✓) 8.20.12.3 Verify vacuum is between -150 and -240 cmH₂O (mbar) (-110 and -175 mmHg) as shown on digital manometer. ___ cmH₂O

8.20.12.4 Disconnect APL bypass hose from digital manometer and reconnect to APL port on breathing system.

8.21 Volume Alarms

8.21.1 Set ventilator to Volume Control mode and confirm.

8.21.2 Set: pmax = 40 mbar (cmH₂O)
Vt = 380 mL
Freq = 12 BPM
TI:TE = 1:1
TIP:TI = 10%
PEEP = 0 mbar (cmH₂O)

8.21.3 Press the Alarms key, scroll to MV High Alarm and confirm. Does a box appear around the Minute Volume High Alarm Limit? ___(Y)

8.21.4 What is the High Minute Volume Alarm default? ___ (12) Press to confirm.


8.21.5 Scroll to Minute Volume Low Alarm and confirm. What is the minute Volume Low default? (3.0)

- (✓) 8.21.6 Adjust the low minute volume alarm to 5.0 liters and confirm. Does the MINUTE VOLUME LOW!! alarm appear as a mid alarm with a corresponding yellow LED? ___ (Y).

8.21.7 Adjust minute volume alarm to 2.0 and confirm. Does the Minute Volume Low alarm disappear? ___ (Y).

8.21.8 Disconnect breathing hose from expiratory connection and start a stop watch.

- (✓) 8.21.9 What is the time when APNEA FLOW!! appears as a mid alarm? ___ sec (13 - 17)

- (✓) 8.21.10 What is the time when APNEA FLOW!!! appears as a high alarm? ___ sec (26 - 34)
- 8.21.11 Disconnect the respiratory volume sensor cord from the Flow sensor.
- (✓) 8.21.12 Does the "FLOW SENSOR FAIL!" message appear as an Advisory with a single audible tone? ___ (Y)
- 8.21.13 Does the APNEA FLOW!!! alarm disappear? (Y)
- 8.21.14 Do the numerical values associated with volume disappear? ___ (Y)
- 8.21.15 Reconnect the respiratory volume sensor cord to the flow sensor, and expiratory limb of breathing circuit to breathing system.
- 8.21.16 Press Standby, confirm, and recalibrate the flow sensor.
- 8.21.17 Fresh Gas Low
- 8.21.17.1 Select Volume Control mode and confirm.
- 8.21.17.2 Disconnect the hose attached to the breathing system's inspiratory port.
- (✓) 8.21.17.3 After 30 seconds verify the APNEA PRESSURE!!!, APNEA FLOW!!!, MINUTE VOLUME LOW!!!, FRESH GAS LOW!!, and NO FRESH GAS!!! alarm messages appear on the display with a corresponding audible alarm.
- 8.21.17.4 Reconnect the patient Y connector to the patient inspiratory port and activate the flush.
- 8.21.17.5 Verify the FRESH GAS LOW!! alarm disappears after two (2) seconds following the next ventilator cycle. ___ (Y)
- (✓) **8.22 Audio Silence**
- 8.22.1 Press the  key.
- 8.22.2 Does the LED on the Silence Alarms key light? Is the audio alarm silenced for 120 sec.?
- 8.22.3 Press Standby key and confirm.
- 8.22.4 Cycle the system power and wait for the Standby screen to appear.
- 8.22.5 Press Volume Control and confirm.
- 8.22.6 Verify the 120-sec. delay starts at power-up.

8.22.7 Disconnect all test equipment.

(✓) **8.23 Oxygen Flush Valve**

8.23.1 Press and release the O₂FLUSH button.

8.23.2 Does the flow of oxygen stop immediately? ___ (Y)

8.23.3 Connect a 12-inch hose to the inspiratory valve.

8.23.4 Occlude the bag mount with test plug.

8.23.5 Insert the sensor from a calibrated O₂ Med into the valve dome adapter on the inspiratory valve.

8.23.6 Close all flow control valves.

8.23.7 Press the O₂FLUSH button.

(✓) 8.23.8 What is the O₂ concentration after the value stabilizes? ___%O₂ (97-100)

8.23.9 Remove the O₂ sensor from the inspiratory dome and install the plug.

8.23.10 Disconnect the fresh gas hose from the breathing system. Using the fresh gas test connector (4117361), test terminal (4104389), and appropriate adapter (4115087), connect to fresh gas hose and attach test terminal to in port of test volumeter.

8.23.11 Press and hold the O₂ FLUSH button for 15 seconds; multiply the value by 4.

(✓) 8.23.12 What is the oxygen flush flow rate? ___ L/min.

Minimum: 35 L/min. @ 50 psi (domestic)

Minimum: 25 L/min. @ 40 psi (2.8 bar)

Maximum: 75 L/min. @ 87 psi (6.0 bar)

8.23.13 Remove the test minute volumeter and test fixture, and reconnect the fresh gas hose.

8.23.14 Turn the System Power switch to ON.

8.24 Final Tests

- (✓) 8.24.1 Operator's Instruction Manual
 - 8.24.1.1 Verify that the availability/location of the machine's Operator's Instruction Manual is in close proximity of the machine.
- (✓) 8.24.2 Lamp Test
 - 8.24.2.1 Verify that the table lamp is working properly.
- (✓) 8.24.3 Final Check
 - 8.24.3.1 Verify that all cylinder pressure gauges indicate zero (if applicable).
 - 8.24.3.2 Verify that the pipeline hoses are connected to the hospital pipeline.
 - 8.24.3.3 Verify that the APL valve knob is fully open.
 - 8.24.3.4 Verify that the O₂ sensor is removed from the valve dome and the plug is inserted in the inspiratory valve dome.
 - 8.24.3.5 Verify that the machine is plugged into a live outlet.
 - 8.24.3.6 Verify all test equipment is removed from machine.
 - 8.24.3.7 Return all machine controls and settings to their original state.

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

Spare and Replacement Parts List

Appendix B

| | | |
|-------------|------------|-----------------------------|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS |
|-------------|------------|-----------------------------|


SPARE AND REPLACEMENT PARTS

Part numbers for field-replaceable items on the Fabius Tiro anesthesia system are listed on the following pages, along with part numbers for related hardware and cables.

The item numbers are keyed to the accompanying illustrations to aid in identifying the item and its location.

NOTE: Certain items in the following parts listings include an “alternate part number” corresponding to those original part numbers that exceed seven (7) digits. Use the original part number whenever possible when ordering spare or replacement parts.

| ASSEMBLY/PART | PAGE |
|---|------|
| Figure 1. Trolley Mount and Wall Mount, Fabius Tiro | 2 |
| Figure 2. Spirolog Flow Sensor, O2 Sensor, and Fresh Gas Hose | 4 |
| Figure 3. Breathing System Assembly w/ APL Bypass | 6 |
| Figure 4. Breathing Bag Extensions | 8 |
| Figure 5. Absorber Assembly, Top, COSY | 10 |
| Figure 6. Canister Assembly | 12 |
| Figure 7. Breathing System Mount Assembly, Patient Airway Pressure Gauge Assembly | 14 |
| Figure 8. Flowmeter Tube Assembly, Flow Control Valves, and Valve Manifold (part of Flowmeter Bezel Assembly 4118287 or 4118288) | 16 |
| Figure 9. Fresh Gas Flow Sensors and Filter Assembly (part of Flowmeter Bezel Assembly 4118287 or 4118288). | 18 |
| Figure 10. Pipeline Pressure Gauges and Fresh Gas Display PCB (part of Flowmeter Bezel Assembly 4118287 or 4118288) | 20 |
| Figure 11. SORC | 22 |
| Figure 12. Cylinder Yokes, Regulators, and Gauges | 24 |
| Figure 13. Gas Inlet Assembly | 26 |
| Figure 14. Controller PCB and Related Parts | 28 |
| Figure 15. Pneumatic (PEEP Control) Assembly | 30 |
| Figure 16. Bezel Assembly, Monitor | 32 |
| Figure 17. Ventilator | 34 |
| Figure 18. Auxiliary O2 Flowmeter | 36 |
| Figure 19. Vapor Mounting System | 38 |
| Figure 20. Scavenger, AGS | 40 |
| Figure 21. Scavenger, AGS, w/ Adjustment Valve | 42 |
| Figure 22. O2 Flush Valve | 44 |



TROLLEY MOUNT


WALL MOUNT

SV88026

Figure 1. Trolley Mount and Wall Mount, Fabius Tiro

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|----------------------|---|-------------|
| Figure 1. | | |
| | Trolley Assembly | 4118056 |
| 1. | Caster, locking (2x) | 8411945 |
| 2. | Caster (2x) | 8411946 |
| 3. | Drawer (3x) | 4118258 |
| 4. | Drawer Slide (6x) | 4118397 |
| | Pop-Rivet (2X/side) 5/31 x 1/8 grip (not shown) | 4118470 |
| 5. | Writing tray | 4118259 |
| 6. | Drawer Slide (2x) | 4118359 |
| 7. | Rail (2x) | 4118388 |
| 8. | Rail Ring (4x) | 2600239 |
| 9. | Cap, end (4x) | G60455 |
| 10. | Screw, socket head (2x) | 1329952 |
| 11. | Cap, end (2x) | 4118396 |
| | | |
| | Wall Mount Assembly | 4118141 |
| 12. | Frame, wall mount | 4118232 |
| 13. | Pin (x2) | 4118230 |
| 14. | Screw, cap socket head (x2) | HW01137 |
| 15. | Screw, socket head (x4) | 1329952 |
| 16. | Nut, kep (x4) | HW55011 |
| 17. | Rail Ring (x4) | 2600239 |
| 18. | Rail (x2) | 4118388 |
| 19. | Plate, dovetail | 4118439 |
| 20. | Cap, end (x4) | G60455 |
| 21. | Washer, lock (x2) | HW65005 |
| 22. | Screw, set | HW04017 |
| 23. | Spacer | 4110792-106 |
| 24. | Spacer | 4110792-092 |
| 25. | Plate, swivel | 4118333 |
| 26. | Screw, shoulder | 4118334 |
| 27. | Knob | 4115758 |
| 28. | Screw, flat head (x2) | HW03058 |
| 29. | Washer, flat (x2) | HW66005 |
| 30. | Washer, lock | HW65006 |
| 31. | Nut, hex | HW50010 |
| 32. | Screw, cap socket head | 1338196 |
| 33. | Tape, 3/4" UHMW Polyethylene | 4118387 |



SV88040

Figure 2. Spirolog Flow Sensor, O2 Sensor, and Fresh Gas Hose

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 2.

Spirolog Flow Sensor

- | | | |
|----|---|-------------|
| 1. | Spirolog flow sensor, 5 pack | 8403735 |
| 2. | Spirolog cable assembly, COSY | 4117081-001 |
| | Alternate part number | 4199951 |

O2 Sensor


- | | | |
|----|--|---------|
| 3. | Sensor capsule | 6850645 |
| 4. | Housing assembly, less capsule | 4117454 |
| 5. | Adapter (cover) | 4106385 |
| 6. | O-ring, #025 (Viton) | 4105791 |
| 7. | O-ring, #016 (Viton) (2x) | 4106388 |

Fresh Gas Hose

- | | | |
|----|--|-------------|
| 8. | Fresh gas hose assembly | 4108577-002 |
| | Alternate part number | 4199941 |
| | Fresh gas hose assembly, Right hand COSY | 4108577-003 |
| | Alternate part number | 4199910 |

Miscellaneous Items (Not Shown)

- | | | |
|--|---|-------------|
| | Hose Asm - Silicon - Ventilator | 4117590 |
| | Hose Asm, Silicon - Ventilator, Right hand COSY | 4117152-001 |
| | Hose Asm - Peep/Max - APL Bypass | 4117027 |



SV02006

NOTE A - DEPENDANT UPON THE CONFIGURATION OF THE COSY CURRENTLY INSTALLED, ITEM NUMBER 19 MAY BE LOCATED IN EITHER POSITION SHOWN ON ILLUSTRATION.

Figure 3. Breathing System Assembly w/ APL Bypass

| | | |
|-------------|------------|---|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|-------------|------------|---|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 3.

| | | |
|-----|--|--|
| 1. | Dome, insp valve | 4107960 |
| 2. | Plug assembly, O2 sensor | 4115265 |
| 3. | Plunger. | 4117399 |
| 4. | Breathing System Assembly w/APL bypass, Service Exchange (new style). | SE4118378 |
| | Alternate part number | 4199939 |
| | Breathing System Assembly w/APL bypass, Service Exchange | SE4116398 |
| | Breathing System Assembly - Non-US, (new style). | 4118379 |
| | Breathing System Assembly - Non-US | 4117529 |
| | NOTE: Verify the style of the original breathing system for its equivalent replacement. See the following figure for details. | |
| 5. | Scavenger port assembly. | 4116715 |
| 6. | APL valve assembly (rotary knob style). | MK00625 |
| 6a. | APL valve assembly (level style) | 4117149 |
| 6b. | Retaining nut | M30448 |
| 7. | Crater asm - Kit (includes O-rings and shoulder bolt) | 4118168 |
| 8. | Shoulder Bolt | 4118166 |
| 9. | Union nut (2x). | M30453 |
| 10. | Packing ring (2x) | M09231 |
| 11. | Valve disk (2x) | M23225 |
| 12. | Exp. hose terminal (incl. O-ring). | 8603275 |
| 13. | Label - Exp. | 4118054-002 |
| 14. | Bag mount port. | 8603276 |
| 15. | Label - Breathing Bag | 4118054-003 |
| | O-ring | M20622 |
| 16. | Insp. hose terminal. | 8603249 |
| 17. | Label - Insp. | 4118054-001 |
| | O-ring | M20622 |
| 18. | Absorber top & canister | Refer to Figure 5. and Figure 6. |
| | Fresh gas decoupling valve: | |
| | Cap. | 2600734 |
| | O-ring, cap | 8603336 |
| | Spring Cross. | M17450 |
| | Diaphragm, Viton | 8604105 |
| | Crater | 2600737 |
| | PEEP/Pmax valve: | |
| | Cap & hose barb assembly | 4116714 |
| | O-ring | 8603336 |
| | Diaphragm asm | 4117885 |
| | Crater | 2600731 |
| | Sealing ring | 2600479 |
| | APL bypass valve: | |
| | Cap & hose barb assembly | 4116929 |
| | Spring. | 4116925 |
| | Spacer, spring holder | 4116926 |
| | Crater assembly. | 4116933 |
| 19. | Sample line gas return fitting (non-U.S. breathing systems) | 2600751 |




Figure 4. Breathing Bag Extensions

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 4.

| | | |
|----|---|-----------|
| 1. | Pole Assembly, Breathing Bag Extension (for use with item #3) | 4117753 |
| 2. | Pole Assembly, Breathing Bag Extension (for use with item #4) | 4118380 |
| 3. | Breathing System Configuration (for reference only) | N/A |
| 4. | Breathing System Configuration | SE4118378 |




Figure 5. Absorber Assembly, Top, COSY

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 5.

| | | |
|----|--------------------|---------|
| 1. | Absorber top | M33013 |
| 2. | Lip seal | M30455 |
| 3. | O-ring | R50313 |
| 4. | O-ring | R18352 |
| 5. | O-ring (3x) | 4302469 |
| 6. | Washer (9x) | 2600459 |
| 7. | Screw (3x) | 1266225 |
| 8. | Washer (3x) | 1329472 |
| 9. | Packing ring | M30456 |

SV00419




Figure 6. Canister Assembly

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 6.

| | | |
|----|--|---------|
| 1. | Canister assembly | M29320 |
| 2. | Screen insert | M29999 |
| 3. | Canister | M29994 |
| 4. | Packing ring | M30563 |
| 5. | Bead chain | M30562 |
| 6. | Screw, oval hd AM4 x 6 DIN 85 - A4 | 1315838 |
| 7. | Acorn nut, M4 DIN 1587 - M A4/051 | 1334751 |
| 8. | Screw, M4 x 8 DIN 921 - A2 | 1336029 |




Figure 7. Breathing System Mount Assembly, Patient Airway Pressure Gauge Assembly

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 7.

| | | |
|----|--|-------------|
| 1. | Breathing system mount assembly | .4117400 |
| | Breathing system mount assembly, Right hand COSY | 4117400-002 |
| | Alternate part number | 4199772 |

Gauge assembly, patient airway pressure

| | | |
|-----|--|-------------|
| 2. | Pressure gauge assembly | .4114290 |
| 3. | O-ring, #010 (neoprene)..... | 8604831 |
| 4. | Gauge mount adapter | .4117072 |
| 5. | Hose barb coupling, 90 deg | .4117070 |
| 6. | Hose, 4 x 1.5 SI NF clear | .1190520 |
| 7. | Screw, 8-32 x 1/2 in. btn hd skt | HW09002 |
| 8. | Lock Washer, #8 int-t | HW67000 |
| 9. | O-ring, #109 (EPDM) | 4112628-001 |
| 10. | Filter | 8402868 |
| 11. | Restrictor | .4118440 |




Figure 8. Flowmeter Tube Assembly, Flow Control Valves, and Valve Manifold (part of Flowmeter Bezel Assembly 4118287 or 4118288)

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| | | |
|------|-------------|-------------|
| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|

Figure 8.

1. Flowmeter tube assembly (incl. brackets and mounting nuts)4116950

Flow control valves

[part of flow meter bezel assembly 4118287 (2-gas) or 4118288 (3-gas)]

2. Flow control valve (3x)4118483
 3. Spacer (3x) 4110792-054
 Alternate part number 4199752
 4. Collet (3x)4116969
 5. Chuck (3x) 2600308
 6. Knob, N2O & Air (2x) M24452
 7. Knob, O2 M24900
 8. Washer, DIN 1751 (3x) M19190
 9. Spring washer (3x) M24963
 10. Nut (3x) 2600305

Valve Manifold Assembly

[part of flow meter bezel assembly 4118287 (2-gas) or 4118288 (3-gas)]

11. Manifold4116988
 12. Washer, 4 mm blue M30937
 13. Angle connector M30953
 14. Straight connector, M5 M30952
 15. Ferrule, pilot to Divan M30963
 16. Washer, M4 white M31602

Knob caps part numbers:

| Country | O2 | | Air | | N2O | |
|------------------------------|-------|---------|--------|--------|-------|--------|
| | Color | P/N | Color | P/N | Color | P/N |
| Germany, Austria | Black | M34305 | Black | M34307 | Black | M34306 |
| USA | Green | M25147 | Yellow | M25797 | Blue | M24901 |
| France, Spain, UK/Ireland | White | 4115657 | Blk/wh | M26205 | Blue | M24901 |




Figure 9. Fresh Gas Flow Sensors and Filter Assembly (part of Flowmeter Bezel Assembly 4118287 or 4118288)

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|


| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 9.

Fresh gas flow sensors and filter assembly

[part of flow meter bezel assembly 4118287 (2-gas) or 4118288 (3-gas)]

| | | |
|-----|---|-------------|
| 1. | Flow sensor housing | 4118279 |
| 2. | Flow sensor assembly (3x) (w/o angle connector (11) and O-ring (8)) | 4116938 |
| 3. | Cable assembly, gas flow sensor | 4118281 |
| 4. | Safety valve | 4115950 |
| 5. | O-ring, #113 (silicone) | 4107961 |
| 6. | Filter assembly | M16201 |
| 7. | Flow sensor manifold, 3-gas | 4116978 |
| | Flow sensor manifold, 2-gas | 4116978-001 |
| | Alternate part number | 4199937 |
| 8. | O-ring, #106 (neoprene) (3x) | 4109322 |
| 9. | Collar, 6 mm white | M31603 |
| 10. | Angle connector, M5 | M30961 |
| 11. | Angle adapter, 1/4 MPT to M5 (3x) | 4116986 |
| 12. | Straight connector, M5 | M30952 |
| 13. | Washer, M4 white | M31602 |
| 14. | Y-plug | M30962 |
| 15. | Ferrule, Pilot to Divan | M30963 |
| 16. | Washer, 4 mm blue | M30937 |
| 17. | Input valve, complete | 4117293 |
| 18. | Washer, lock (6x) | HW06700 |
| 19. | Screw, button head (6x) | HW09110 |
| 20. | Flow sensor housing cover | 4118280 |
| 21. | Clip, power cable | 4107047 |



SV88030

Figure 10. Pipeline Pressure Gauges and Fresh Gas Display PCB (part of Flowmeter Bezel Assembly 4118287 or 4118288)

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 10.


Pipeline pressure gauges and fresh gas display PCB
[part of flow meter bezel assembly 4118287 (2-gas) or 4118288 (3-gas)]

Pipeline pressure gauges

| | | |
|----|---|---------|
| 1. | Pipeline pressure gauge (3x) | 4117144 |
| 2. | Screw, M4 x 12 mm self tap (4x) | 4117021 |
| 3. | Washer, 4 mm blue | M30937 |
| 4. | Straight connector, M5 | M30952 |
| 5. | Ferrule, pilot to Divan | M30963 |
| 6. | Washer, M4 white | M31602 |

Fresh gas display PCB

| | | |
|-----|---|-------------|
| 7. | Insulator, flowmeter display | 4116965 |
| 8. | PCB assembly, fresh gas display | 4116623 |
| 9. | Cable assembly, display | 4112149-010 |
| | Alternate part number | 4199753 |
| 10. | Cable assembly, O2 pressure | 4117054 |



SV88031

Figure 11. SORC

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

[Figure 11.](#)

SORC

[part of flow meter bezel assembly 4118287 (2-gas) or 4118288 (3-gas)]

| | | |
|----|-------------------------------------|--------|
| 1. | SORC (proportional valve) | M32940 |
| 2. | Resistor, O2 red | M27543 |
| 3. | O-ring | E20274 |
| 4. | O-ring, DIN 3771 | M30308 |
| 5. | Resistor, N2O blue | M27539 |
| 6. | Filter (inside unit) | M32922 |




Figure 12. Cylinder Yokes, Regulators, and Gauges

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 12.

Cylinder Yokes, Regulators, and Gauges:

| | | |
|-----|---|-------------|
| 1. | O2 yoke (2x) | 1101620 |
| | N2O yoke | 1101621 |
| | AIR yoke | 1101625 |
| 2. | Index pin - screw, sltd 0.157 OD 6-32 x 0.718 L (2x per yoke) | 4105929 |
| 3. | Label, 3/8 DOT O2 | 4111266-007 |
| | Alternate part number | 4199935 |
| | Label, 3/8 DOT N2O | 4111266-008 |
| | Alternate part number | 4199934 |
| | Label, 3/8 DOT AIR | 4111237-004 |
| | Alternate part number | 4199930 |
| 4. | Yoke handle | 4113536 |
| 5. | Check valve assembly (2x) | 4113932-002 |
| | Alternate part number | 4199943 |
| 6. | O2 cylinder pressure regulator, 36 psi | 4117244 |
| 7. | N2O cylinder pressure regulator, 36 psi | 4117245 |
| 8. | Straight fitting, 1/4 MPT x 1/8 FPT (2x) | 4102906 |
| 9. | Check valve (2x) | 4117686 |
| 10. | Washer, 4 mm white | M31602 |
| 11. | Washer, 4 mm blue | M30937 |
| 12. | Connector, 90 degree | M30935 |
| 13. | L-fitting, 3/16 tube x 1/4 MPT | 4109409-001 |
| | Alternate part number | 4199928 |
| 14. | Cylinder pressure gauge, 3000 psi, NMM (2x) | 4114247-002 |
| | Alternate part number | 4199944 |
| 15. | Straight fitting, 3/16 tube x 1/8 FT (2x) | 4109402 |




Figure 13. Gas Inlet Assembly

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|


| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 13.

| | | |
|-----|---|---------|
| 1. | Connector, NIST N2O male | 4111382 |
| | Connector, DISS N2O | 4111384 |
| 2. | Connector, NIST Air male | 4110386 |
| | Connector, DISS Air | 4102886 |
| 3. | Connector, NIST O2 male | 4110388 |
| | Connector, DISS O2 | 4102563 |
| 4. | Gas Inlet assembly, 3-gas/3 connector | 4118228 |
| 5. | Inlet regulator (2x) | 4116946 |
| 6. | 90 deg connector G 1/8 (3x) | M30935 |
| 7. | Fitting, pilot to Divan (4x) | M28816 |
| 8. | Contact 04 mm (2x) | M31348 |
| 9. | Ferrule, Pilot to Divan (2x) | M30963 |
| 10. | Pressure switch, O2 supply | 4117143 |
| 11. | Washer, M4 white (4x) | M31602 |
| 12. | Washer, M4 blue (2x) | M30937 |
| 13. | Filter holder (4x) | 4116942 |
| 14. | Filter (4x) | 4102532 |
| 15. | O-ring, #008 (EPDM) (4x) | 4115960 |
| 16. | Check valve (4x) | 4117759 |
| 17. | O-ring, #022 (Viton) (3x) | 4101882 |
| 18. | Lock washer, #10 split (4x) | HW65003 |
| 19. | Screw, socket head (4x) | HW01128 |
| 20. | Connector, straight (4x) | M30952 |

Related parts (not shown):

| | |
|---|---------|
| Hose, 2.7 x 0.65 PAE, black/white | 1210181 |
| Hose, 2.7 x 0.65 PAE, white | 1210157 |
| Hose, 2.7 x 0.65 PAE, blue | 1210211 |



NOTE A: CONFIGURATION SHOWN IN ILLUSTRATION IS BASED UPON CONTROL PCB P/N 4116632. THE CONFIGURATION OF CONTROL PCB P/N 4118079 WILL DIFFER.

Figure 14. Controller PCB and Related Parts

| | | |
|-------------|------------|---|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|-------------|------------|---|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 14.

- | | | |
|----|--|-----------|
| 1. | PCB assembly, Service Exchange | SE4116632 |
| | Alternate part number | 4199927 |

NOTE: The SE4116632 Control PCB is not interchangeable in a unit originally configured with a 4118079 Control PCB.

- | | | |
|--|--|---------|
| | PCB assembly, Service Exchange (for use only with Software Version \geq 2.X) | 4118744 |
|--|--|---------|

NOTE: Replacement PCBs may contain a previous software version. Therefore, a software update may be necessary after installation. If applicable, perform the Software Update Procedure outlined in Section 8 with the appropriate level software. Refer to applicable Technical Service Bulletins for details.

- | | | |
|--|---|-------------|
| | Fuse, 5 x 20 mm time lag, 1.6 A | EC01200-001 |
| | Alternate part number | 4199955 |

NOTE: The following fuses are used only on PBC assembly SE4116632 and 4199927:

- | | | |
|--|---|-------------|
| | Fuse, 5 x 20 mm time lag, 2.5 A | EC01200-002 |
| | Alternate part number | 4199954 |
| | Fuse, 5 x 20 mm time lag, 4.0 A | EC01200-004 |
| | Alternate part number | 4199952 |

- | | | |
|----|---|-------------|
| 2. | Fuse Holder | 1832662 |
| 3. | Fuse, 5 x 20 mm time lag, 3.15 A | EC01200-003 |
| | Alternate part number | 4199953 |
| 4. | Cable assembly, power switch | 4117053 |
| 5. | Power inlet & mains fuse holder | 4117013-001 |
| 6. | Fuse, 5 x 20 mm time lag, 2.5 A (2x) | EC01200-002 |
| | Alternate part number | 4199954 |
| 7. | Power supply | 4117060 |
| 8. | Cable assembly, power supply AC input | 4117050-002 |

NOTE: (new style) Replacement of this cable may also require replacement of the AC Power Inlet, Item #5.

- | | | |
|----|--|-------------|
| 9. | Cable assembly, power supply output (incl. capacitors & inductor) | 4117051 |
| | Cable assembly, power supply output (w/o capacitors & inductor) (new style) | 4117051-001 |

NOTE: Not to be used in conjunction with the Control PCB 4116632 Series.

- | | | |
|-----|--|-------------|
| 10. | Tie strap (2x) | 1101732 |
| 11. | Spacer, R | 4110972-106 |
| 12. | Cable assembly, charge capacitors (incl. capacitors) | 4117089 |

Miscellaneous Items (Not Shown)

Power cord assembly:

- | | | |
|--|--|-------------|
| | USA, 15 ft. | 4117266 |
| | UK/Ireland | 4115377-001 |
| | Alternate part number | 4199926 |
| | Spain/France/Germany/Austria | 4115367-001 |
| | Alternate part number | 4199925 |
| | Battery, 12V Rechargeable (x2) | 4114229 |




Figure 15. Pneumatic (PEEP Control) Assembly

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 15.

| | | |
|-----|--|-------------|
| | Pneumatic (PEEP Control) Assembly | 4117004 |
| 1. | Top Gasket | 4116156 |
| 2. | Muffler | 4116158 |
| 3. | PEEP valve assembly | 2600573 |
| 4. | Filter | 8402868 |
| 5. | Top plate | 4116155 |
| 6. | Straight fitting, 1/16 ID x 10-32 w/seal | 4114732-001 |
| 7. | Straight fitting, 3/32 ID hose x 10-32M (4x) | 4114732 |
| 8. | Adapter fitting for 1/16 ID hose x 1/8 MPT | 4111446 |
| 9. | Pump assembly | 4118160 |
| 10. | L-fitting, 3/32 ID hose x 10-32M | 4117971 |
| 11. | Variable restrictor | 4107995 |
| 12. | Hose, 2.5 x 2 SI (0.025 M) | 4115747 |
| 13. | Tie strap, 0.09 x 4 1/8 | 4106068 |
| 14. | Hose, 2.5 x 2 SI (0.127 M) | 1190520 |
| 15. | Hose, natural, 0.075 ID x 0.030 W (0.178 M) | ML08003 |
| 16. | Hose, 4 x 1.5 SI NF clear (0.012 M) | 1190520 |
| 17. | Hose, 2.5 x 2 SI (0.040 M) | 4115747 |
| 18. | Hose, 2.5 x 2 SI (0.127 M) | 4115747 |
| 19. | Hose, 4 x 1.5 SI NF clear (0.178 M) | 1190520 |




Figure 16. Bezel Assembly, Monitor

| | | |
|-------------|------------|---|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|-------------|------------|---|


| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 16.

| | | |
|-----|--|-------------|
| 1. | Bezel Assembly, Monitor | 4117019 |
| 2. | Keypad (control panel) 4 key | 4118119-001 |
| | Alternate part number | 4199863 |
| | Keypad (control panel) | 4118119-002 |
| | Alternate part number | 4199862 |
| | NOTE: A Switch Panel Insert may need to be ordered with replacement of the Keypad or Bezel Assembly. Refer to item 15 for various inserts based on the original Keypad configuration. | |
| 3. | Bezel | 4116919 |
| 4. | Flat panel display | 4113755 |
| 5. | Knob | M29655 |
| 6. | LED lamp PCB assembly | 4117017 |
| 7. | Rotary switch (incl. mtg. hardware) | 4117012 |
| 8. | Flat washer, 0.168 ID x 0.624 OD | 4117020 |
| 9. | Screw, self tap, M4 x 12 mm | 4117021 |
| 10. | Spacer, 0.162 ID x 0.205 L (4x) | 4113200-001 |
| | Alternate part number | 4199923 |
| 11. | Kep nut, M3 x 0.5 (4x) | HW55008 |
| 12. | Kep nut, M4 x 0.5 (5x) | HW55010 |
| 13. | Lock washer, #8 ext-t. | HW68001 |
| 14. | Gasket, EMI, 0.275 x 0.080, 3.125 long. | 4113906-002 |

Switch panel inserts:

| | | |
|-----|---|-------------|
| 15. | Non-US IPPV, PCV | 4118120-002 |
| | Alternate part number | 4199860 |
| | IPPV | 4118120-004 |
| | Alternate part number | 4199858 |
| | US Volume Control, Pressure Control | 4118120-006 |
| | Alternate part number | 4199856 |
| | Volume Control | 4118120-008 |
| | Alternate part number | 4199854 |



SV88052


Figure 17. Ventilator

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 17.

| | | |
|---|--|---------------|
| Ventilator | | 4117080-001 |
| Ventilator, Service Exchange | | SE4117080-001 |
| Alternate part number | | 4199791 |
| 1. Motor/ball screw assembly (includes encoder) | | 8603150 |
| 2. Diaphragm, patient side | | 2600650 |
| 3. Diaphragm, piston | | 2600651 |
| 4. O-ring | | 8604831 |
| 5. Band clasp | | 2600785 |
| 6. Screw, M6 x 80 hex hd. | | HW08025 |
| 7. Patient assembly | | 2600775-003 |
| Alternate part number | | 4199890 |
| 8. High Pressure valve | | 2600680 |
| 9. Negative Pressure Relief valve (-8 mbar) | | 8604217 |
| 10. Cable assembly, light barrier (7-pin) | | 4117056 |
| 11. Cable assembly, incrementer (6-pin) | | 4117055 |
| 12. Dosing Cage | | 2600676 |
| 13. Gear Case | | 2600677 |
| 14. Rapid Fastener (x3) | | 2600704 |
| 15. Screw, Pan Head, M3 x 8 - DIN 562 (x6) | | 1338544 |
| 16. Washer, Flat, M3-DIN 9201 (x6) | | 1335359 |
| 17. O-Ring, Vent Hose (x2) | | 2M08777 |
| 18. Window - Vent | | 4116991 |



SV00453


Figure 18. Auxiliary O₂ Flowmeter

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 18.

| | | |
|-----|---|-------------|
| 1. | Auxiliary O2 Flow Meter, USA, complete | 4109310-003 |
| | Auxiliary O2 Flow Meter, ISO | 4109310-004 |
| | Alternate part number | 4199906 |
| 1A | Auxiliary O2 Flow Meter, Global (flowmeter and valve asm) | .4117687 |
| 2. | Set screw, cup pt 6-32 x ¼ in. | HW04003 |
| 3. | Knob, auxiliary oxygen | .4111442 |
| 4. | Screw, btn hd skt 6-32 x ¼ in. (2x) | HW09004 |
| 5. | Label, 5/8 DOT - green/wht rings (USA) | 4109373 |
| | Label, 5/8 DOT - white (ISO) | 4105981 |
| 6. | Set screw, cup pt 10-32 x 7/8 in. (2x) | .HW04011 |
| 7. | Label "AUXILIARY OXYGEN" | 4109381 |
| 8. | Tie Strap 0.09 W x 4 1/8 L (3x) | 4106068 |
| 9. | Hose, 0.171 ID x 0.040 W 20 in. | ML08006 |
| 10. | Hose, 0.075 ID x 0.030 W 24 in. | ML08003 |
| 11. | Coupling, 1/8 NPT (F) | 4103668 |
| 12. | T-fitting, 1/8 BSPP x 4 mm OD | .4117057 |
| 13. | Adapter, 1/16 ID hose x 1/8 MPT | .4111446 |
| 14. | Label, "O2" | 4109871 |
| 15. | Kep nut, 10-32 (2x) | HW55002 |
| 16. | Hose clamp, press-on | 4104161 |
| 17. | L-fitting, NY 1/16 ID hose x 10-32 M | .4110173 |



SV88034

SELECTATEC VAPOR MOUNTING SYSTEM
(VIEW LOOKING DOWN)

*P = LOCTITIE #222 (PURPLE)

Figure 19. Vapor Mounting System

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 19.

| | | |
|----|--|---------|
| 1. | Vapor mounting system | 4118227 |
| 2. | O-rings, vaporizer mount (2x) | U04314 |
| 3. | Holder, lifter | M28935 |
| 4. | Screw (2x) | 1334514 |
| | | |
| 5. | Vapor mounting system (Selectatec) | 4118262 |
| 6. | O-ring, vapor mount (2x) | U04314 |
| 7. | Clip | M35122 |
| 8. | Screw, slotted, M2.5x6L9 (4x) | 1330535 |

Miscellaneous Items

| | |
|--|-------------|
| Loctite, Purple, #222 | 4118558-001 |
| Loctite, Blue, #425 | 4118558-008 |
| Touch-up paint: Euro white | 7900380 |
| Touch-up paint: Euro blue | 7901261 |
| Spray paint: Euro white | S000094-002 |
| Spray paint: Euro blue | S000094-003 |
| Silicone Adhesive (Walker-Elastosil) | 1202527 |




Figure 20. Scavenger, AGS

| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 20.

| | | |
|----|----------------------|--------|
| | Scavenger, AGS | M33300 |
| 1. | Container | M33292 |
| 2. | Filter | M33294 |
| 3. | Flow tube | M33293 |
| 4. | Screw plug | M33291 |
| 5. | Coupling | M33149 |

SV88076




Figure 21. Scavenger, AGS, w/ Adjustment Valve


| | | |
|--------------------|-------------------|--|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|--------------------|-------------------|--|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 21.

| | | |
|----|--|---------|
| | AGS Scavenger w/adj valve | 4117494 |
| 1. | Kit, Scavenger mount. | 4117515 |
| | Rail clamp, AGS. | M32967 |
| 2. | L-fitting, Q-disc adapter | M33151 |
| 3. | Metering valve | AF00539 |
| 4. | Adapter, 30 mm | M29430 |
| 5. | Hose, scavenger, 1.0 M | M33295 |
| 6. | Hose, 4 x 1.5 SI NF clear, 1.5 M | 1190520 |

SV88060



UNDERSIDE VIEW OF MACHINE

Figure 22. O2 Flush Valve

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

[Figure 22.](#)

O2 Flush Valve

| | | |
|----|---------------------------------|---------|
| 1. | Valve assembly, fresh gas | 4117226 |
| 2. | Connector, straight (x2) | M30952 |
| 3. | Washer, M4, white (x2) | M31602 |
| 4. | Label, O2 flush: | |
| | Germany | 2600494 |
| | USA | M28934 |




Figure 23. Scavenger, Passive

| | | |
|-------------|------------|---|
| FABIUS TIRO | APPENDIX B | SPARE AND REPLACEMENT PARTS (continued) |
|-------------|------------|---|

| ITEM | DESCRIPTION | PART NUMBER |
|------|-------------|-------------|
|------|-------------|-------------|

Figure 23.

| | | |
|----|--|-------------------------|
| 1. | Passive Scavenger | 4118605 |
| 2. | Mounting Bracket | supplied with scavenger |
| 3. | Mounting Screws | supplied with scavenger |
| 4. | Adapter, 19 mm | 1101312 |
| 5. | Hose, Corrugated | 9995230 |
| 6. | Cage, Anti-Occlusion | 4118551 |
| 7. | Filter, Lint (10 pcs.) | 1199 |
| | Container, Lint Filter (Not Shown) | 14014 |

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

RETURN TO THIS MANUAL'S TABLE OF CONTENTS
RETURN TO CD-ROM TABLE OF CONTENTS

Drägermedical

A Dräger and Siemens Company

DrägerService is a division of
Dräger Medical, Inc.
3122 Commerce Drive
Telford, PA 18969
Tel: (215) 721-5402
(800) 543-5047
Fax: (215) 721-5784
Web: www.draegermedical.com
Printed in the U.S.A.

Fabius Tiro Service Manual

Rev. A summary of changes

| Page | Description |
|--|--|
| Cover | Updated revision level and date |
| TOC | Updated as required |
| 6 | Added Fabius Tiro to Trademark Notices |
| 177 | Changed first NOTE, Changed last NOTE |
| 178 | Added first NOTE, changed step 14, Updated Table 1 |
| 181 | Updated Figure 19 |
| 191 | Added new step 9, renumbered subsequent steps |
| 192 | Added new step 14, renumbered subsequent steps, updated Figure 24 |
| 229 | Software Update Procedure: Removed the following sections: Requirements Installing Send Image (utility) file on laptop Updating Send Image (utility) file All information contained in the removed sections can be found in existing Service Procedures. The existing Service Procedures are now listed in the Requirements Section of the Software Update Procedure |
| PMS Section, Appendix A | |
| 24 | Changed title of Section 8.8, Added Passive Scavenger cleaning and functional test |
| 45 | Added AIR to steps 8.13.2 and 8.13.2.2 |
| 46 | Added AIR to steps 8.13.2.5 and 8.13.2.6 |
| Spare and Replacement Parts Section, Appendix B | |
| 7 | Changed information under item #4, removed P/N from item #11 |
| 17 | Changed P/N of item #2 |
| 25 | Added AIR yoke to item #1, Added AIR Label to item #3 |
| 29 | Changed P/N of SE PCB Assembly (w/ SW Version 2.X) |
| 35 | Changed P/N of item #1, deleted alternate P/N |
| 46, 47 | Added Passive Scavenger Figure #23 |

Rev. A summary of changes

NOTE: As you will see, the manual has changed dramatically since its prior release. This was due to the request for development of structured documentation that is compatible with an XML format.

During the restructuring of this document, it was required that we no longer maintain our existing format style. This has resulted in the relocation of the PMS and Spare and Replacement Parts Sections of the manual.

The PMS and the Spare and Replacement Parts Sections of the manual are now contained as Appendix A and Appendix B, respectively. In order to maintain integrity between the PMS and its associated PM form, Appendix A paragraph numbering will remain unchanged within the appendix. However, the paragraph number no longer will relate to the location within the manual.

Below are the technical changes that occurred in this manual since its last distribution on the Anesthesia Service CD, 1Q04:

| Page | Description |
|------------|--|
| Cover..... | Updated revision level and date |
| 14..... | Updated Figure 6 caption |
| 15..... | Added new Figure 7 |
| 23..... | Updated Figure 14 caption |
| 24..... | Added new Figure 15 |
| 25..... | Updated Figure 16 caption |
| 26..... | Added new Figure 17, changed 1st paragraph |
| 27..... | Changed Figure 18 caption, added 3 paragraphs |
| 28..... | Added new Figure 19 |
| 29..... | Changed Figure 20 caption |
| 30..... | Added new Figure 21, changed first paragraph |
| 31..... | Changed Figure 22 caption, added last 3 paragraphs |
| 32..... | Added first 2 paragraphs, added new Figure 23 |
| 34..... | Changed Figure 24 caption |
| 35..... | Added new Figure 25 |
| 36..... | Changed Figure 26 caption |
| 37..... | Added new Figure 27 |
| 43..... | Changed Figure 32 (removed PCB ASM P/N) |
| 46..... | Changed Figure 35 (added 4th softkey and renumbered) |
| 47..... | Added new #4 in table, resequenced subsequent |
| 64..... | Changed Figure 1 (removed PCB ASM P/N) |

| | |
|----------|---------------------------|
| 66..... | Added new Figure 3 |
| 141..... | Changed Figure 2 |
| 164..... | Updated NOTE |
| 166..... | Updated Figure 14 |
| 177..... | Added step 2 |
| 182..... | Changed Figure 20 |
| 184..... | Updated Figure 21 |
| 197..... | Added 1st NOTE |
| 198..... | Changed Figure 28 |
| 217..... | Changed Figure 10 caption |
| 218..... | Changed Figure 11 |

APPENDIX A - PMS

| | |
|---------|---------------------------------|
| 2..... | Changed P/N of Test Lung |
| 5..... | Changed P/N of Test Lung |
| 17..... | Added 8.1.2 - Circuit Isolation |
| 53..... | Changed step 8.20.1.3 |

APPENDIX B - Spare and Replacement Parts Section

| | |
|---------|---|
| 6..... | Changed Figure 3 |
| 7..... | Added NOTE |
| 8..... | Added new Figure 4 |
| 9..... | Added parts list for new Figure 4 |
| 30..... | Changed Figure 15 |
| 31..... | Changed P/N's of items #6, #10, and #14 |
| 32..... | Added item #10 |
| 33..... | Added NOTE, changed item #15 |

Rev. - summary of changes (4118302-002)

| Page | Description |
|------|-------------|
|------|-------------|

The part number for this manual has changed due to the restructuring of its format. There were no technical changes to the manual for this dash revision.

Rev. C summary of changes (4118302-001)

| Page | Description |
|-------|--|
| Cover | Updated revision level and date |
| 5-32 | Changed first NOTE, first and second paragraphs, added 'Machine Type' to listing |
| 5-46 | Updated third paragraph, added fourth paragraph |
| 6-38 | Changed second NOTE, changed first paragraph, replaced step 6.13.2 |
| 6-39 | Changed steps 6.13.13 and 6.13.14, changed first paragraph |
| 6-40 | Changed figure caption |
| 6-41 | Added new figure, renumbered subsequent |
| 6-43 | Changed step 6/14/6 |
| 7-19 | Changed title of steps 7.10 and 7.10.1 |
| 8-2 | Changed Loctite P/N's, removed 5 repair tools |
| 8-3 | Changed P/N of O-Ring, Diaphragm |
| 8-17 | Changed step 8.1.2 |
| 8-22 | Changed second NOTE, updated steps 8.7.10 and 8.7.11 |
| 8-50 | Changed step 8.17.2.4 |
| 8-52 | Changed steps 8.19.13 and 8.19.14 |
| 9-1 | Changed first paragraph. changed first NOTE |
| 9-3 | Changed step 9.2.8 |
| 9-4 | Changed step 9.3.3, deleted NOTE |
| 9-5 | Changed step 9.4.2 |
| 9-10 | Changed NOTE and steps 9.4.11, 9.4.12, and 9.4.13 |
| 9-12 | Deleted NOTE, changed step 9.5.3 |
| 9-13 | Changed step 9.5.12 |
| 9-14 | Changed figure caption on figure 9-5, added new figure 9-6 |
| 10-3 | Added new parts to item #4 |
| 10-5 | Added new parts to item #8, added miscellaneous items |
| 10-13 | Added new parts to item #1, changed P/N of item #3 |
| 10-15 | Added alternate P/N to item #3 |
| 10-17 | Switched description and P/N of item #18 and #19 |
| 10-19 | Added alternate P/N to item #9 |
| 10-26 | Added NOTE A to illustration |
| 10-27 | Added first, fourth and fifth NOTEs, changed second NOTE, changed P/N for item #5 and #8, added P/N to item #9 |
| 10-30 | Updated illustration |
| 10-31 | Renumbered all items |
| 10-33 | Changed P/N for item #4 |
| 10-37 | Changed Loctite P/Ns |

Rev. B summary of changes

| Page | Description |
|-----------------------|--|
| Cover..... | Updated revision level |
| 4-10..... | Corrected illustration |
| 5-2..... | Added Secure Options line |
| 5-3..... | Relocated Figure 5-2 to page 5-3 |
| 5-4..... | Relocated Figure 5-2 onto page |
| 5-32..... | Updated NOTE and added second statement to NOTE |
| 5-37..... | Updated second paragraph of step 4.7.2.1.3 |
| 5-38, 5-39, 5-40..... | Moved to pages 5-43, 5-44, and 5-45 |
| 5-40..... | Changes last statement on page |
| 5-41..... | Page change |
| 5-42..... | Page change |
| 5-43..... | Page change |
| 5-44..... | Page change |
| 5-45..... | Updated both the first and second paras in step 4.7.2.4.2 |
| 5-48..... | Changed Serial Port Parameters Screen |
| 6-26..... | Changed wording in step 6.9.28 |
| 6-27..... | Updated illustration |
| 6-44..... | Changed para 6.15.11 |
| 8-2..... | Added repair tools and P/N |
| 8-3..... | Added (piston) to O-ring, for Diaphragm |
| 8-20..... | Added new steps 8.5.7 and 8.5.8 |
| 8-29..... | Replaced steps 8.9.4.3 through 8.9.4.6 |
| 8-51..... | Changed step 8.19.2 |
| 8-52..... | Changed step 8.19.26 |
| 8-53..... | Added step 8.19.32 |
| 8-55..... | Replaced steps 8.20.7.2 through 8.20.7.6 |
| 9-5..... | Corrected step 9.4.2 |
| 9-10..... | Corrected step 9.4.13 |
| 10-6..... | Updated illustration |
| 10-7..... | Added item numbers 6a and 6b |
| 10-27..... | Added additional PCB board numbers, distinguished between fuse usage, updated first NOTE, added second NOTE. |
| 10-32..... | Added item 18, removed Motor and Case Asm P/N |
| 10-33..... | Added item 18 |
| 10-37..... | Added miscellaneous parts |

Rev. A summary of changes

| Page | Description |
|------------|--|
| Cover..... | Updated revision level |
| 1-2..... | Changed illustration |
| 1-3..... | Changed illustration |
| 2-1..... | Changed illustration |
| 2-2..... | Changed illustration |
| 2-3..... | Changed illustration |
| 2-4..... | Changed illustration |
| 2-5..... | Changed illustration |
| 2-6..... | Changed illustration |
| 2-8..... | Changed illustration |
| 2-9..... | Changed illustration, removed 'if' from second paragraph |
| 2-10..... | Changed illustration |
| 2-11..... | Removed last sentence from first paragraph, deleted 'inlet valve vent' and parenthesis from fourth paragraph, changed illustration |
| 2-13..... | Changed 'three' to 'various' in first paragraph |
| 2-14..... | Changed item 3 in Key to state 'MAN/SPONT APL Valve' |
| 2-15..... | Changed illustration |
| 2-16..... | Changed illustration |
| 2-17..... | Changed illustration |
| 2-18..... | Changed illustration |
| 2-19..... | Added '/Pressure' to heading, made last sentence of last paragraph its own paragraph and added 'During ventilation' |
| 2-20..... | Added '/Pressure' to heading, inserted '/' (slashes) between PEEP and Pmax, changed illustration |
| 2-21..... | Added '/Pressure' to heading, changed illustration |
| 2-26..... | Changed illustration |
| 2-27..... | Added '/' (slash) to header between PEEP and Pmax |
| 2-28..... | Changed illustration |
| 2-30..... | Changed illustration |
| 4-2..... | Changed illustration |
| 4-3..... | Added P/N to figure caption |
| 4-4..... | Changed illustration |
| 4-7..... | Changed illustration |
| 4-8..... | Changed illustration |
| 4-9..... | Changed illustration |
| 4-10..... | Changed illustration |
| 4-11..... | Changed illustration |

| | |
|------------------|--|
| 4-13 | Changed illustration |
| 4-14 | Changed illustration |
| 4-15 | Changed illustration |
| 4-16 | Changed illustration |
| 5-2 | Added quick reference tree |
| 5-3 through 5-31 | Page numbers changed due to added page |
| 5-3 | Changed Figure 5-2 caption |
| 5-4 | Deleted figure |
| 5-5 | Figure number changed |
| 5-6 | Figure number changed |
| 5-7 | Figure number changed |
| 5-8 | Figure number changed and changed dates in figure |
| 5-13 | Figure number changed |
| 5-14 | Figure number changed |
| 5-17 | Figure number changed |
| 5-18 | Figure number changed |
| 5-19 | Figure number changed |
| 5-20 | Figure number changed |
| 5-21 | Figure number changed |
| 5-22 | Figure number changed |
| 5-23 | Figure number changed |
| 5-24 | Figure number changed |
| 5-25 | Figure number changed |
| 5-26 | Figure number changed |
| 5-27 | Figure number changed |
| 5-28 | Figure number changed |
| 5-29 | Figure number changed and corrected spacing within figure |
| 5-30 | Figure number changed |
| 5-31 | Figure number changed and deleted figure |
| 5-32 | Deleted existing data on page 5-32, moved prior data to current page, Figure number changed |
| 5-33 | Figure number changed |
| 5-34 | Figure number changed |
| 5-35 | Figure number changed |
| 5-36 | Figure number changed |
| 5-37 | Figure number changed |
| 5-38 | Figure number changed |
| 5-39 | Figure number changed and added text to address changing weather conditions |
| 5-40 | Figure number changed |
| 5-41 | Figure number changed |

5-42 Figure number changed
5-43 Figure number changed
5-44 Figure number changed
5-45 ... Figure number changed and changed enabled to lower case in second paragraph
5-46 Figure number changed
5-47 Figure number changed
5-48 ... Figure number changed and changed default PARITY value to NONE
5-49 Added text to support operation of Pump On/Off
6-2 Changed 'Core Mount Conversion' to 'Core Module Conversion' - all places, changed 'proved' to 'provide' in first paragraph,
6-3 .. Added new step 6.1.1.13, re sequenced subsequent steps, changed 'Core Mount Inversion Reassembly' to "Core Module Inversion Reassembly"
6-5 Changed illustration
6-7 Changed first warning
6-8 Changed illustration
6-10 Reversed existing steps 6.4.3 and 6.4.4, moved existing step 6.4.8 to after step 6.4.9
6-11 Changed illustration
6-18 Added new step 6.8.6
6-19 Changed illustration
6-23 Changed illustration
6-24 Changed 'Roll' to 'Remove' in step 6.9.24
6-26 Deleted reference to 'mark' in step 6.9.29
6-27 Changed illustration
6-28 Added new step 6.9.36
6-29 ... Deleted existing steps 6.9.43 through 6.9.48, re sequenced subsequent steps
6-29 through 6-end. All page numbers have changed resulting from the deletion of steps removed as referenced above
6-38 ... Added new NOTE after 6.13.7, changed cross reference in the second NOTE, added third NOTE, added new step 6.13.2
6-39 Text from previous page rolled, added 'If applicable' to paragraph 6.13.12, added new step 6.13.14, deleted step 6.13.16
6-40 ... Text from previous page rolled. All text on pages hereafter have been rolled.
6-41 Added 'If applicable' to the beginning of step 6.14.16
6-45 Changed illustration
6-50 Changed last NOTE to a step, added a CAUTION
6-51 ... Added new steps 6.18.11 and 6.18.12, re sequenced subsequent steps
6-54 Moved last step from previous page to this page

6-55 Moved step on page to previous page
7-2 Changed illustration
7-3 Changed reference in step 7.1.13
7-10 Relocated steps previously found on this page
7-11 Added step previously relocated from page 7-10
8-3 Made 'Use of Batteries' lower case
8-11 Changed 'machines' to 'machine' in first paragraph
8-14 Added PMC form part number
8-17 Removed reference to 'Convenience' outlet strip, moved ending step and
table portion to next page
8-18 Added step and table portion from previous page
8-19 Deleted step 8.4.3 and renumbered steps, changed text in steps 8.4.6,
8.5.2, and 8.5.3, re sequenced
8-20 Changed steps 8.6.7 and 8.6.9
8-22 Changed second NOTE, removed - (dash) and placed comma between 1
and 2 Stop Bits in table
8-24 Added checkmark by 8.8
8-26 Made 8.9.2 Freshgas two words, deleted checkmark from 8.9.2.11
8-27 Corrected 8.9.2.12 symbol
8-29 Deleted steps 8.9.4.4, 8.9.4.6, and 8.9.4.7 and renumbered, changed text
in steps 8.9.4.8 and 8.9.4.9, added row of data to table, deleted
checkmark from 8.9.4.8
8-41 Updated illustration
8-49 Deleted checkmark from step 8.16.8, changed text in step 8.16.11
8-50 Deleted checkmark from step 8.17.1.9
8-54 Added (cmH₂O) to step 8.20.4.6
8-55 Removed space between cm H₂O in step 8.20.7.3, added (cmH₂O) to end of
step 8.20.7.6
8-56 Added (cmH₂O) to first and sixth line in step 8.20.8.2
8-57 Added (cmH₂O) to step 8.20.11.1, added (cmH₂O) to first and sixth line of
step 8.21.2
8-58 Changed text in step 8.21.17.3, deleted step 8.21.17.6
9-1 Deleted second paragraph
9-4 Changed NOTE
9-5 Replaced step 9.4.2
9-6 Changed illustration
9-10 Replaced step 9.4.13
9-11 through end Text and illustrations on pages rolled onto following pages
10-2 Changed illustration
10-3 Change Qty for item number 4, added new item number 6, re sequenced
subsequent item numbers

10-5 Added alternate part number for item 8
10-7 Removed 'SE' from items 1,2,and 3, changed part number for item 4,
reversed order of alternate part number and fitting under item 4
10-12 Changed illustration
10-17 Added detailed text to item number 2
10-22 Changed illustration
10-23 Added alternate part number for item 13
10-26 Changed illustration
10-31 Added Non-US keypad P/N and alternate P/N
10-33 Added alternate part number for item 7
10-40 Changed illustration
10-41 Added Draeger Mount, under item number 1