Fisher
Bioblock Scientific

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Operating instructions

DÜPERTHAL – Under-Bench Cabinets in accordance with EN14470-1

Model: **Type 90 - UTS line:**

29-060664-xyz, 29-060660-xyz(L), 29-061160-xyz, 29-061163-xyz, 29-061165-xyz,
29-061164-xyz, 29-061463-xyz, 29-061465-xyz

Type 90 - UTS-line with Special Depth of 500 mm:

29-060654-xyz, 29-060650-xyz(L), 29-061150-xyz, 29-061153-xyz, 29-061155-xyz,
29-061154-xyz, 29-061453-xyz, 29-061455-xyz

Type 30 - UTS BASIC line: 23-061160-xyz

Please enter the following details:

(This information will be important if there are queries.)

Cabinet Model: _____ Manufacturing No. : _____

Key No.: _____ Location: _____

Accountable Person: _____

TABLE OF CONTENTS for Operating Instructions

Section			Page	Section			Page
1	Manufacturer, Technical Data for Under-Bench Cabinet and Safety Instructions		15	8	Opening the Under-Bench Cabinet after a Fire		21
2	Installing the Under-Bench Cabinet		17	9	Inspection/ Faults		21
3	Transporting the Under-Bench Cabinet		17	10	Disposal		22
4	Commissioning the Under-Bench Cabinet		17	11	Spare Parts		22
5	Interior of the Under-Bench Cabinet		18	12	EG Declaration of Conformity		55
6	Ventilation		19	13	GS Certificate		56, 57
7	Operating the Under-Bench Cabinet		20				

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Fax +49 (0) 61 88 / 7 81-1 54

1.2 Description / Application

Under-Bench Cabinet in accordance with EN 14470-1 - Type 90 / Type 30 with 90 minute / 30 minute fire resistance (FWF 90 / FWF 30), for the storage of flammable liquids in workplaces.

1.3 Construction

1.3.1 UTS line

Outer body: Powder-coated steel sheet

Insulation: Insulating plates with sandwich construction.

Interior: Melamine resin coated decorative panels.

1.3.2 UTS BASIC line

Body: Impervious B1 material with impervious material

Key x

0 = Light grey body RAL 7035, zinc-yellow door(s) RAL 1018

4 = Light grey body RAL 7035, light grey door(s) RAL 7035

8 = Body and door(s) in special colour

9 = Light grey body RAL 7035, door(s) special colour

The manufacturing number (serial number), e.g. 40001, is printed on the labelling sticker.

1.4 Fittings

	Model-size								
	XLD(-5)	XLK(-5)	L(-5)	LD(-5)	LK(-5)	LT(-5)	ST(-5)	S(-5)	BLT
UTS line With Wing doors Special depth 500 mm						29-061160-x00 29-061160-x01	29-060660-x00(L) 29-060660-x01(L)		
UTS line with Drawers Special depth 500 mm	29-061465-x00 29-061465-x01		29-061164-x00 29-061164-x01	29-061165-x00 29-061165-x01				29-060664-x00 29-060664-x01	
UTS line Comb. with Wing Doors / Drawers Special depth 500 mm	29-061455-x00 29-061455-x01		29-061154-x00 29-061154-x01	29-061155-x00 29-061155-x01				29-060654-x00 29-060654-x01	
UTS line Comb. with Wing Doors / Drawers Special depth 500 mm		29-061463-x00 29-061465-x01			29-061163-x00 29-061163-x01				
BASIC-line With Wing doors									23-061160-x00 23-061160-x01

(Key to variables: x = Colouring, yz = Definition of interior Fittings, L = optional door opening to left)

1.5 General Technical Data

1.5.1 UTS line and BASIC UTS line - Depth 595 mm

Model	UTS	UTS	UTS	UTS	UTS	UTS	UTS	UTS	UTS BASIC
Model-size	S	ST	LT	LK	LD	L	XLK	XLD	BLT
Type	90	90	90	90	90	90	90	90	30
<u>External dimensions:</u>									
Length app. mm	595	595	1100	1100	1100	1100	1400	1400	1100
Width app. mm	595	595	595	595	745	595	595	595	595
Height (without Pedestal) app. mm	600	600	600	600	600	600	600	600	--
Height (with Pedestal) app. mm	635	635	635	635	635	635	635	635	635
Empty weight with interior fittings app. kg	163	158	206	216	221	221	251	276	97
Max. Total weight with full load app. kg	193	233	281	321	281	251	356	336	172
Volume/storage space app. m ³	0,11	0,11	0,22	0,22	0,22	0,22	0,28	0,28	0,17

1.5.2 UTS line - Depth 500 mm

Model	UTS	UTS	UTS	UTS	UTS	UTS	UTS	UTS
Model-size	S-5	ST-5	LT-5	LK-5	LD-5	L-5	XLK-5	XLD-5
Typ	90	90	90	90	90	90	90	90
<u>External dimensions:</u>								
Length app. mm	595	595	1100	1100	1100	1100	1400	1400
Width app. mm	500	500	500	500	700	500	500	500
Height (without Pedestal) app. mm	600	600	600	600	600	600	600	600
Height (with Pedestal) app. mm	635	635	635	635	635	635	635	635
Empty weight with interior fittings app. kg	155	150	195	205	210	210	239	262
Max. Total weight with full load app. kg	185	225	270	310	270	240	344	322
Volume/storage space app. m ³	0,09	0,09	0,18	0,17	0,17	0,18	0,23	0,23

1.6 Guarantee

Our liability for delivery faults is defined in our Conditions of Delivery. We accept no liability for damage arising from non-compliance with the Operating and Maintenance Instructions and Conditions of Use.

Only products in perfect condition with the features we guarantee leave our factory. Proper operation is guaranteed if these Operating and Maintenance Instructions are observed.

1.7 Safety

These Operating and Maintenance Instructions contain basic instructions to be observed during operation and maintenance/commissioning. It is therefore essential that the operator's specialist staff read these Operating and Maintenance Instructions before commissioning, that the operator draw up operating instructions based on these ones, but taking into account local and company-specific conditions and that appropriate instructions are given to the accountable, specialist staff before they start any activity.

These Operating Instructions contain safety precautions, non-observance of which could be hazardous to humans and the environment and which are specifically identified with the general warning sign in accordance with the BGV A 8 accident prevention regulation.

This sign is specially used to warn of electrical voltage.

For safety precautions, were non-observance could damage technical equipment or cause a malfunction, the word

WARNING is added.

Instructions attached directly to the Under-Bench Cabinet must be observed to the letter and maintained in a fully legible condition.

1.8 Dangers of Non-observance of Safety Precautions

Non-observance of safety precautions may be hazardous to humans and the environment as well as to technical equipment. Such non-observance may also result in any guarantees and compensation claims being forfeited.

The safety precautions listed in these Operating and Maintenance Instructions, current legal regulations, accident prevention regulations and the operator's own operating and safety regulations must be complied with.

Safety Precautions for the Operator

- The ban on storing hazardous material and preparations together must be observed.
- The operator must prohibit unauthorized access to hazardous materials stored in the under-bench cabinet.
- The operator must make rules about smoking and handling of naked flames, as well as working with naked flames and on spark producing jobs in and around the under bench cabinet.
- Safety equipment must be maintained in working order at all times.
- The incoming and exhaust-air vents of the technical ventilation system must be kept in unimpaired working order.

- Hazardous materials must be stored in such a way that the catch pan can be checked visually every [working] day for contaminants.
- In accordance with accident prevention regulations, any spills must be contained immediately with absorbers, wiped up and disposed of in a regulation-compliant manner

Safety Precautions for Maintenance, Inspection and Assembly Procedures

The operator must ensure that authorized and qualified specialist staff carry out all maintenance and assembly procedures.

Work on the electrical system should always be undertaken only by qualified electricians at zero current. For more information, please refer to the applicable accident prevention regulations, VDE regulations and the regulations of your local power supplier.

Improper operating methods

The safety equipment such as mechanical ventilation system, catch pans etc built into the DÜPERTHAL under-bench cabinet must be maintained in perfect working order at all times.

Improper operation also means allowing unqualified or insufficiently qualified personnel and any other unauthorized persons free access to hazardous materials stored in the under-bench cabinet.

Operation is unlawful, if the under-bench cabinet operator fails to comply with legal requirements such as:

- Not issuing the required operating instructions and/or
- Failing to carry out a risk assessment or write up an explosion prevention document and/or
- Failing to perform required maintenance and inspection work or carry out periodical tests and/or
- Failing to provide staff with documented training, which must be done at least once a year, on the proper storage of hazardous materials and/or materials and preparations that are hazardous to water.

1.9 Applications

The DÜPERTHAL under-bench cabinet is approved for the storage of flammable liquids in workplaces in accordance with the Technical Regulations for flammable Liquids (TRbF 20, Appendix L).

2. Transport

2.1 Safety measures for In-House Transport

2.1.1 Under-Bench Cabinet UTS-line / UTS BASIC-line

The unloaded under-bench cabinet may be lifted up only from the underside, using suitable means of transport and must be transported horizontally. Pick-up points are on the side or the back, never lift holding the drawers / or door panels. The receiving surface of the transporting vehicle e.g. lift truck, must be at least as large as the surface area of the under-bench cabinet. Damage to the underside of the cabinet from lift trucks or such is to be avoided as is tilting while the cabinet is being lifted for transportation. Tilting may damage the cabinet and impair or destroy the functioning of the built-in safety features.

3. Installation

3.1 Installation site requirements

The ground/floor must be level and able to bear the weight of the under-bench cabinet when it is fully loaded. Further, the installation must be chosen so as to ensure that the under-bench cabinet:

- cannot not be damaged by vehicles;
- is not be exposed to any direct or indirect source of heat;
- Is protected from humidity;
- is not be operated at temperatures below -5°C or above +45°C.

3.2 Installing the Under-Bench Cabinet

The under-bench cabinet must be installed in a level position. Alignment of the under-bench cabinet is the only way to ensure that it will not overturn and likewise only under this condition is the perfect functioning of the self-closing door mechanism guaranteed.

4. Commissioning

4.1 Drawer / Wing Door

4.1.1 Locking ram for the UTS line

The drawer/door can be locked by means of the integrated locking cylinder. The key number, for example A007, is stamped on the locking ram and on the keys supplied with it.

4.1.2 Locking ram for the UTS line

The door can be locked by means of the integrated cylinder lock. The key number, for example 0001, is embossed on the locking ram and on the keys supplied with it.

4.2 Door fixing device / Closing Mechanism

4.2.1 Door fixing device for the UTS line Models

The drawer or/wing door is opened by pulling the handle and remains open in any position. This ensures safe storage or removal. At a temperature of about 50 °C, the open drawer or/door is closed by the safety mechanism.

4.2.2 Door Closing Mechanism (Door fixing device optional) for the UTS BASIC line models

The wing door is open by pulling the handle. The door closes automatically when the handle is released. If the cabinet has a fixing device (optional) fitted on it, pull the door as far as it can go; it stops automatically in the factory-installed door fixing device. The door fixing device holds the door open for as long as is necessary to ensure safe storage or removal. At a temperature of about 50 °C, the door held open by the fixing device is closed by the safety mechanism.

WARNING

The drawer / door of the under-bench cabinet must not be held open by using wooden wedges or such objects. This will block the safety mechanism!

4.3 Noise emission

Noise emitted by the door closure is < 70dB(A).

4.4 Earthing options

The interior is conductively connected to an equipotential bonding terminal on the outer body of the cabinet. Earthing is therefore possible in accordance with BGR 132 (Prevention of Ignition Risks). Please refer to sections 4.4.1 or 4.4.2 according to model. If an isolated contact is required on each wall inside the cabinet, it can be connected directly to the connecting terminal in the cabinet roof (inside). If an additional isolated contact is required, the appropriate DÜPERTHAL 2-000700 potential equalizer bar must be installed.

4.4.1 Earthing connection for the UTS line model

Earthing can be done with an equipotential bonding terminal on the rear wall of the cabinet on right and left.

4.4.2 Earthing connection for the UTS BASIC line model

Earthing can be done with an equipotential bonding terminal on the cabinet rear wall.

5. Interior fittings UTS line / UTS BASIC line

5.1 Inserting shelves (optional)

Position the brackets supplied in the delivery to the desired height. Place the shelves on the particular support brackets and screw them tightly 2x each on to the side walls of the cabinet. Thread the earthing bolt to the equipotential bonding terminal through the appropriate drill hole on the shelves. When adjusting the height loosen the earthing screw and the retaining screw and remove the brackets. Reposition the brackets at the desired height, secure the shelves with the retaining screws and thread in the earthing bolt to the equipotential bonding terminal.

5.2 Positioning the tray shelves

The shelves are permanently mounted during production. Reconfiguration is not possible.

5.3 Technical Specifications for Fittings in UTS line / UTS BASIC line

see section 5.3.1 and 5.3.2

5.4 Base bottom tray

5.4.1 Drawer Side

Shelves are fitted on to the drawers in the under-bench cabinet, so no additional spill tray is required under the shelf on the bottom drawer. It is not an additional storage surface for containers. You can however make full use of the surface by using a perforated sheet, see DÜPERTHAL catalogue Section 2. The perforated sheet can be dispensed with, if you store only one container with the max. permitted volume.

5.3.1 Depth 595 mm

Model		UTS	UTS	UTS	UTS	UTS	UTS	UTS	UTS	UTS BASIC
Model-size		S	ST	LT	LK	LD	L	XLK	XLD	BLT
Load capacity per shelf with uniform weight distribution (Max. Loadability / storage) in kg		--	30	30	30	--	--	30	--	30
Load capacity per tray shelf with uniform weight distribution (Max. Load / storage) in kg		30	--	--	30	30	30	30	30	--
Max. Volume of the largest individual container which may be stored in the cabinet in l	left side	10	10	21	10	10	20	16	15	17
	right side				10	10		10	10	
Bottom tray Volume in l	left side	11	11	24	11	11	22	18	17	19
	right side				11	11		11	11	
Max. Total weight with full load in approx. kg		193	233	281	321	281	251	356	336	172

5.3.2 Tiefe 500 mm

Model		UTS	UTS	UTS	UTS	UTS	UTS	UTS	UTS
Model-size		S-5	ST-5	LT-5	LK-5	LD-5	L-5	XLK-5	XLD-5
Load capacity per shelf with uniform weight distribution (Max. Loadability / storage) in kg		--	30	30	30	--	--	30	--
Load capacity per tray shelf with uniform weight distribution (Max. Load / storage) in kg		30	--	--	30	30	30	30	30
Max. Volume of the largest individual container which may be stored in the cabinet in l	left side	7	7	17	8	7	15	12	12
	right side				7	7		7	7
Bottom tray Volume in l	left side	8	8	19	9	8	17	14	14
	right side				8	8		8	8
Max. Total weight with full load in approx. kg		185	225	270	310	270	240	344	322

5.4.2 Door Side

The spill tray in the floor area of the under-bench cabinet is for the collection of spillages inside the cabinet. It is not an additional storage surface for containers. However, this area can be utilized with the insertion of a perforated sheet, see DÜPERTHAL catalogue Section 2.

The bottom tray pan is also bolted to the equipotential bonding terminal. Please ensure that the bottom tray pan is connected back to the equipotential bonding terminal after dismantling, e.g. for cleaning.

5.5 Control of the bottom tray

WARNING

The leakages spilled into the spill tray are to be removed immediately and disposed of in a regulation-compliant manner. We therefore recommend keeping an emergency set handy, e.g. the DÜPERTHAL 11.30.000.

5.6 Clamping Device for Piping (optional), only UTS-line Type 90

The clamping device item 29 -70000-044 is available as an optional extra, for the DÜPERTHAL under-bench cabinet Type 90.

It is meant for securing piping leading from the cabinet. The fixture is suitable for max. 3 steel pipes with 6 to 10 mm Ø and is approved by the Material Testing Institute. Their report states that correct use of the clamping device will in no way adversely affect the fire resistance properties. Open pipe connections which are not in use must be blocked. Please contact the manufacturer for this.

Pipe feed through sockets varying in size e.g. larger than size 10 mm Ø, can be requested from DÜPERTHAL. A range of special sizes is also available with valid approvals.

6. Ventilation

6.1 General

Under-bench cabinets can be connected to a mechanical venting system leading to a non-hazardous location in the open, as per TRbF 20 Appendix L. In ventilated under bench cabinets there must be a minimum 10-fold air-exchange / hour, with the doors closed and the drop in pressure must not be greater than 1.5 mbar (150 Pa). If the under-bench cabinet is connected to existing technical ventilation, ensure that the ventilation is operating 24 hours a day.

NOTE

If the under-bench cabinet is connected to separate, technical ventilation, the fan should be mounted as close as possible to the exhaust-air line exiting into the open.

6.1.1 Pressure Differential Table Depth 595 mm

Model Type	Model-size	Type	Exhaust-Air Connection Back Wall of Cabinet	Volume Flow 10-fold air exchange/hr
UTS	S	90	left	approx. 1,5 m³/h
UTS	ST	90	left	approx. 1,5 m³/h
UTS	LT	90	left	approx. 2,5 m³/h
UTS	LK	90	left	approx. 2,5 m³/h
UTS	LD	90	left	approx. 2,5 m³/h
UTS	L	90	left	approx. 2,5 m³/h
UTS	XLK	90	left	approx. 3,0 m³/h
UTS	XLD	90	left	approx. 3,0 m³/h
UTS BASIC	BLT	30	right	approx. 2,0 m³/h

Drop in pressure in cabinet (all models): < 1 Pa.

6.1.2 Pressure Differential Table Depth 500 mm

Model Type	Model-size	Type	Exhaust-Air Connection Back Wall of Cabinet	Volume Flow 10-fold air exchange/hr
UTS	S-5	90	left	approx. 1,0 m³/h
UTS	ST-5	90	left	approx. 1,0 m³/h
UTS	LT-5	90	left	approx. 2,0 m³/h
UTS	LK-5	90	left	approx. 2,0 m³/h
UTS	LD-5	90	left	approx. 2,0 m³/h
UTS	L-5	90	left	approx. 2,0 m³/h
UTS	XLK-5	90	left	approx. 2,5 m³/h
UTS	XLD-5	90	left	approx. 2,5 m³/h

Drop in pressure in cabinet (all models): < 1 Pa.

6.2. Connection to an Exhaust-Air System

6.2.1 UTS-line

The openings for air supply and exhaust-air extraction are on the back of the under-bench cabinet. An incoming and exhaust-air connections NW Ø 75 are included in the standard delivery. This facilitates the connection to a ventilation system, e.g. DÜPERTHAL 2.00.315 for control of the exhaust-air with a fan/ventilator (ATEX compliant).


The exhaust-air pipe is connected to the exhaust-air connecting piece – **left connecting piece, seen from the front**. The pipeline must be connected with a sleeve or something similar to the exhaust air connecting piece. The air supply can be taken from the room and needs no additional external air supply is required under normal working conditions. After installation of the under-bench cabinet, check the connection to the extraction system, e.g. with smoke pipes.

The installation of technical ventilation or connection to an existing exhaust-air system is not part of the DÜPERTHAL supply package. Please contact a ventilation company near you.


WARNING

Fitting a ventilator entails a combination of devices. The fan/ventilator, as a Zone 2 operating device must comply with the requirements of RL 94/9/EG (ATEX) Equipment Group II, Category 3 G.


S, ST

LT, LK, LD, L

XLD, XLK

Rear view

6.2.2 UTS BASIC-line


The exhaust-air and air supply opening is situated at the rear of the under-bench cabinet. An exhaust and air supply connection NW Ø 100 is included in the standard delivery. A reducing piece of NW 75 Ø is available from DÜPERTHAL as an optional extra and this makes it possible to connect to a ventilation system, say the DÜPERTHAL 2.00.315 for exhaust-air control with fan/ventilator (ATEX compliant).

The exhaust-air line is to be connected to the exhaust-air connecting piece – **right connecting piece, seen from the front**. The pipeline must be connected with a sleeve or something similar to the exhaust-air connecting piece. The air supply can be taken from the room and no additional external air supply is required under normal working conditions. After installation of the under-bench cabinet, check for proper connection to the extraction system e.g. with smoke pipes.

The installation of technical ventilation or connection to an existing exhaust system is not part of the DÜPERTHAL supply package. Please contact a ventilation company near you.

WARNING

Fitting a ventilator entails a combination of devices. The fan/ventilator, as Zone 2 operating device must comply with the requirements of RL 94/9/EG (ATEX) Equipment Group II, Category 3 G.

6.3 Air Circulation inside the Cabinet

The ventilation is effective directly above the bottom tray.

6.4 Fan/Ventilator Selection

Only those fans/ventilators may be used which meet the technical specifications defined in the VDMA Unit Sheet 24 169 Part 1.

In non-explosive areas, this means, according to VDMA: Interior: Zone 2 – Outside: Non-explosive area

WARNING

Fitting a ventilator entails a combination of devices. The fan/ventilator, as a Zone 2 operating device must comply with the requirements of RL 94/9/EG (ATEX) Equipment Group II, Category 3 G.

Also refer to DÜPERTHAL Catalogue Range section PA 2- Fans/Ventilators.

6.5 Operating without mechanical ventilation

In accordance TRbF 20, Appendix L, the DÜPERTHAL under-bench cabinet can be installed and operated without connection to a technical ventilation system.

If the under-bench cabinet is operated as an operating device without connection to an exhaust-air system, a notice to this effect must be displayed on the under-bench cabinet.

Under these operating conditions, the interior of the under-bench cabinet becomes a Zone 1 explosive area. When operating the under-bench cabinet without connection to a technical extraction system, the following must be observed:

- In particular the relevant regulations of the explosion-proofing guidelines.
- Determine the Ex-Zone as per the hazardous materials ordinance for all possible operating conditions and mark this permanently and visibly.
Handling of naked flames and smoking is prohibited in the explosion endangered area. Similarly, no spark-producing tools may be used in this area. It must further be ensured that no electrostatic charges build up in these zones and no operating resources are used with surface temperatures of higher than the igniting temperatures of the stored flammable liquids. Electrical equipment may be operated only if it conforms to the requirements of the VDE 0165.
- Under-bench cabinets without technical ventilation may be operated in working spaces only under the following conditions:
 - A radius of at least 2.5 m around the under-bench cabinet is a Zone 2 explosive area to a height of at least 0.5 m above the floor.
 - In technically ventilated workplaces with a minimum 5-fold air-exchange/ hour, the Zone 2 explosive area can be reduced to 1 m from the front of the under-bench cabinet and 0.5 m from the sides of the cabinet, to a height of 0.3 m above the floor.

WARNING

Mark the Zone 2 explosive area as per the Accident Prevention Regulations (BGV A 8) in a clearly visible manner with the warning symbol shown below:

The dimensions of all symbols and signs should be in proportion to the size of the under-bench cabinet.

7. Operation

7.1 Basic Principles

In accordance with TRbF 20 Appendix L No. 4, the DÜPERTHAL under-bench cabinet may be used for the storage of flammable liquids in work spaces and can be operated with – as well as without connection to a technical ventilation system.

The following information must be displayed visibly on the front of the under-bench cabinet:

- The instruction: Close Door.
- W 1 warning sign (warning of fire-hazardous materials) in accordance with DIN 4844-1.
- P02 prohibiting sign (fire, naked flame and smoking prohibited) in accordance with DIN 4844-1.

Even in cases where the operational situation would seem to require it, unauthorized persons must be prohibited from accessing the hazardous materials stored in the under-bench cabinet. This ban must be highlighted by means of a prohibiting sign or a sign with the inscription „Operation by Unauthorized Persons Prohibited“ on the front of the front of the cabinet as per the accident prevention regulations BGV A 8, P06.

WARNING

Observe the mark for maximum loading height on the inside of the drawer / door. Stacking higher than this mark poses the risk of the stored containers being hit by the top inside edge of the cabinet when the drawer is closed. The door of the under-bench cabinet must not be held open with wooden wedges or the like. This will block the safety mechanism!

As per TRbF 20 No. 3.1.3, the following quantities of flammable liquids may be stored in the DÜPERTHAL under-bench cabinets with minimum **90-minute** fire resistance

Highly and easily flammable liquids with a flash point less than 21 °C

- in breakable containers, e.g. glass, porcelain etc. maximum of up to 60 l
- in other containers e.g. metal etc. maximum of up to 450 l

Flammable liquids with a flash point higher than 21 °C and equal to/less than 55 °C

- in breakable containers, e.g. glass, porcelain etc. maximum of up to 200 l
- in other kinds of containers e.g. metal etc. maximum of up to 3,000 l

Flammable liquids with flash point higher than 55 °C

maximum of up to 5.000 l

For an under-bench cabinet with 90-minute fire resistance, the conversion conditions set out in TRbF 20 No. 3.1.1 Section 2 apply for calculation of the total storage quantity when storing together flammable liquids with a flash point higher than 55 °C and flammable liquids with lower flash points. The lower quantities specified in No. 3.1.3 must not be exceeded

The DÜPERTHAL under-bench cabinet can be installed and operated directly next to other 90-minute fire resistant under-bench cabinets without any need to leave some distance between them.

In accordance with TRbF 20 Appendix L No 3.3, Inflammable liquids to a total quantity of 300 litres can be stored in under-bench cabinets with fire resistance (FWF) of **minimum 20 but less than 90 minutes**, in a work space of up to 100 m².

Highly and easily inflammable liquids with a flash point less than 21 °C

- in breakable containers, e.g. glass, porcelain etc. maximum of up to 60 l and
- in other containers e.g. metal etc. maximum of up to 200 l

The total storage quantity of highly and easily flammable liquids must not exceed 200 litres.

In areas larger than 100 m², the storage quantity in such cabinets can be increased, taking into consideration the maximum quantities in proportion to the room size, if the free space around the installation site of the under-bench cabinets is at least 10 m.

The storage quantities permitted in under-bench cabinets with fire resistance (FWR) of at least 20 but less than 90 minutes may be doubled taking into consideration the above-mentioned allowable maximum quantities, if

- an automatic fire alarm system and a team of certified industrial fire fighters is available with a maximum call out time of 5 minutes after the alarm has been raised.
- there is an automatic fire extinguishing system in the workplace

7.2 Storage

In the DÜPERTHAL under-bench cabinets the following hazardous materials can also be stored together with flammable liquids:

- Noxious materials and preparations
- Irritant hazardous materials
- Class 2 and 3 fire-accelérant materials to a maximum of 200 kg – see technical regulations for hazardous materials 515 (TRGS 515) No. 1.1 (2) and No. 3.3.3.1 outline.
- Poisonous materials and preparations up to a maximum of 200 kg and extremely poisonous materials to a maximum of 50 kg, as long as the poisonous and the extremely poisonous materials are both flammable – see TRGS 514 No. 1.2 (4) and No. 3.2.2 (2) and/ or No. 3.2.6.

The following must **not** be stored in the under-bench cabinet:

- Materials and preparations with an ignition temperature below 100 °C - e.g. carbon disulphide – unless the materials and preparations are stored in ventilated under-bench cabinets in packages which prevents them from igniting.
- Self-igniting materials, including all materials and preparations, which are classified under Class 4.2 of the regulations for road transportation of hazardous materials.
- Materials and preparations with unstable properties which can lead to the outbreak of fires and explosions.
- Unsealed containers and receptacles
- Aggressive materials, products and preparations.

9.4 Faults

Fault	Possible Cause	Possible solution
Drawer / door(s) not closing	Under-bench cabinet is not in plumb	Align the under-bench cabinet
No extraction	Fire protection flaps are closed, because seal mechanism is no longer working	Replacement of seal mechanism by technical customer service
Drawers/ door(s) heavy, not running smoothly	Air exchange inside the cabinet is too high Corrosion on moving parts (hinges etc.)	Remove rust, lubricate parts, remove aggressive materials from cabinet, inform technical customer service

8. Opening the Under-Bench Cabinet after a Fire

After a fire, the under-bench cabinet cannot be opened until at least 24 hours have elapsed.

WARNING

Depending on the duration and intensity of the fire, an explosive vapour/air mix may have built up inside the under-bench cabinet. Before opening the under-bench cabinet, all ignition sources, naked flames lights etc must be removed from within a 10 m radius of the cabinet. If the under-bench cabinet has to be opened with tools, only non-spark producing tools must be used.

9. Inspection / Faults

The under-bench cabinet must be inspected for visible external faults or defects after installation, and before commissioning, thereafter at the subsequently stipulated intervals and after alterations or repairs.

9.1 Daily Operating Checks

In accordance with water balance regulations, the base [collection] pan must be checked every day of operation. Liquid spillages must be mopped up straightaway and disposed of in the regulation-compliant manner.

9.2 Monthly inspection

The monthly inspection encompasses the following points:

- Drawer / Door closure – Open the door and check the closure.
- Ventilation – check the effectiveness of the ventilation with a wool thread, smoke pipe etc in the inside of the cabinet in front of the left exhaust-air duct.
- Sealing strips – check that the sealing strips are in place in the frame of the frame of the cabinet body and the front edge of the door(s). If any damage is apparent, the sealing strips must be replaced straightaway.

9.3 Annual inspection

The annual inspection of the under-bench cabinet must undertaken by a qualified technician.

In the vent of faults, you can assist technical customer service (for address and telephone number see cover page) by stating the cabinet model number, manufacturing and key number as well as providing a description of the fault.

10. Disposal

The under-bench cabinet can be dismantled completely. The individual components e.g. metal, insulation panels, door closing mechanisms etc can be recycled separately. The national and local disposal regulations must be observed. To protect natural resources, parts of the under-bench cabinet and/or the entire cabinet should not be put in the inorganic waste or domestic refuse collections.

11. Spare Parts and Accessories

Product:

- Sealing Strips
- Fire protection Flaps
- Tray Bottoms/Drip Trays
- Base Pan
- PP-Insert
- Anti-slip mat (rubber)
- Cylinder lock
- Handle
- Holed Steel Insert

General

These Operating instructions have been drawn up, with reference to the requirements of the applicable legal regulations as of May 2004.

WARNING

We will not update these Operating Instructions to reflect changes to applicable laws, regulations etc, unless specifically instructed to do so.

These Operating Instructions do not replace the user guide for the **Under-Bench Cabinet for Inflammable Liquids** which is required by law to be drawn up by the operator.

Technical dimensions given are approximate.
Subject to technical alterations. Printing errors excepted.

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All rights reserved. Reprinting, even in parts, prohibited. No part of this work may be reproduced or edited using electronic systems, duplicated or disseminated without written permission from M/s DÜPERTHAL in any form whatsoever (photocopy, microfilm or any other process), even for purposes of training.

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EG-DECLARATION OF CONFORMITY

in accordance with Machinery Directive 98/37/EG, Appendix II A

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hereby declare, that the design, construction and commercialized execution of the below mentioned machine complies with the essential health and safety requirements of the EG Machinery Directive. The validity of this declaration will cease in case of any modification or supplement not being agreed with us previously. Furthermore, validity of this declaration will cease in case that the machine will not be operated correctly and in accordance to the operating instructions and/or not to be inspected regularly.

Machine description:

safety cabinet

Machine type:

29-201260-xyz / 29-200960-xyz / 29-200660-xyz(L) /
29-200670-xyz(L) / 29-130660-xyz(L)

29-200970-xyz

23-201260-xyz / 23-200660-xyz(L) / 23-061160-xyz

29-060664-xyz / 29-060660-xyz(L) / 29-061160-xyz /
29-061163-xyz / 29-061165-xyz / 29-061164-xyz /
29-061463-xyz / 29-061465-xyz

29-060654-xyz / 29-060650-xyz(L) / 29-061150-xyz /
29-061153-xyz / 29-061155-xyz / 29-061154-xyz /
29-061453-xyz / 29-061455-xyz

variables:

x = colour

yz = definition of the interior

L = optional door hang left

Relevant
EG Directives:

Machinery Directive (98/37/EG)

Transposed national
standards and technical
specifications in
particular:

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Certification Mark:

Product: Safety cabinets
Under-bench Cabinet

The product was tested on a voluntary basis and complies with the essential requirements. It meets the requirements of January 06, 2004 of the German Equipment and Product Safety Act. The certification marks shown above can be affixed on the product. See also notes overleaf.

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„x“ colour description, stands for number 0, 4, 5, 8, 9

„yz“

description of interior fitting, stands for number:

00 without interior fittings

01 with interior fittings made of sheet steel

02 with interior fittings made of stainless steel

03 with additional interior fittings

made of polyethylen or polypropylen

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© Copyright 2004