

SOLVENT DELIVERY MODULE LC-10AT
SHIMADZU HIGH PERFORMANCE
LIQUID CHROMATOGRAPH
INSTRUCTION MANUAL

SHIMADZU CORPORATION
CHROMATOGRAPHIC INSTRUMENTS DIVISION
KYOTO, JAPAN

Precautions

High performance liquid chromatography, which often requires large volumes of organic solvent, requires proper care in handling and an appropriate installation site.

Be sure to pay due attention to the precautions noted throughout this manual, in addition to those below.

1. Ventilation

Solvents used in high performance liquid chromatography are often inflammable and toxic, therefore the room in which they are used must be well ventilated.

2. Fire

Never use an open flame in the vicinity of a high performance liquid chromatograph. Do not install other equipment that may generate sparks in the same room.

Install a fire extinguisher at a conspicuous location nearby.

3. Static Discharge

Static discharge also may cause a fire. Appropriate measure should be taken. See "Section 12.1 Precautions on Static Electricity".

4. Protective Goggles

Wear protective goggles when using solvents.

5. Other Facilities Required

A sink should be located as close as possible to the unit so that toxic solvents splashed in the eyes or on the skin may be immediately rinsed off.

6. Power Source

This instrument employs various voltages for each part, in addition to that of the commercial power supply.

As maintenance may in general be carried out with the main cover closed, do not open the cover. If you need to open the cover, be sure to contact the field representative or service person.

7. Other Considerations

Avoid installation of the instrument in direct sunlight or where variations of ambient temperature are extreme.

In earthquake-prone areas, use the supplied tie plates when stacking the units, to prevent them from toppling.

CONTENTS

Precautions

Chapter 1 General

1.1 Outline	1-2
1.2 Features	1-3

Chapter 2 Check Upon Delivery

2.1 Lists of Supplied Parts, etc.	2-2
--	-----

Chapter 3 Construction and Functions

3.1 Front panel	3-2
3.2 Rear Panel	3-4
3.3 The Right Side and the Bottom	3-5

Chapter 4 Installation

4.1 Conditions for Installation Site	4-2
4.2 Examples of System Configurations	4-3
4.3 Mounting of Multiple Units	4-6
4.4 Electrical Connections	4-7
4.5 Piping of the Unit	4-8
4.6 Connection of a Degasser	4-12
4.7 Mounting the Mixer	4-16
4.8 Connection of the Injector, Column and Detector	4-20
4.9 Installation of a High Pressure Gradient Elution System	4-23
4.10 Installation of a Low Pressure Gradient Elution System	4-26

Chapter 5 Operation

5.1 Basic Operation	5-2
5.2 Constant Flow Solvent Delivery	5-8
5.3 Constant Pressure Solvent Delivery	5-11
5.4 Time Program Creation and Execution	5-12
5.5 Auxiliary Functions (AUX.FUNC)	5-20
5.6 High Pressure Gradient Elution Mode	5-29
5.7 Low Pressure Gradient Elution Mode	5-32
5.8 Plunger Washing Method	5-33
5.9 Change of Mobile Phase	5-35

Chapter 6 Performance Checking

6.1	Operation Checking for Simple Isochratic System	6-2
6.2	Checking Concentrations in High Pressure Gradient Elution	6-5
6.3	Checking Concentrations in Low Pressure Gradient Elution	6-7

Chapter 7 Control from External Equipment

7.1	Connections of External Equipment to the Input/Output Terminals	7-2
7.2	Control from the SCL-10A	7-4

Chapter 8 Maintenance

8.1	Replacement of Plunger Seals	8-2
8.2	Replacement of Plungers	8-5
8.3	Replacement of Washing Seal	8-7
8.4	Cleaning and Replacement of Check Valves	8-9
8.5	Replacement of Line Filter	8-10
8.6	Replacement of Fuse	8-11
8.7	Cleaning	8-12

Chapter 9 Troubleshooting

9.1	Symptoms and Countermeasures for Troubles	9-2
9.2	Error Messages	9-3

Chapter 10 Specifications

10.1	LC-10AT Unit Specifications	10-2
10.2	High Pressure Gradient Elution Specifications	10-3
10.3	Low Pressure Gradient Elution Specifications	10-4

Chapter 11 Spare Parts and Optional Units

11.1	Consumable Parts List	11-2
11.2	List of Other Necessary Parts	11-3
11.3	Optional Units List	11-4

Chapter 12 Reference Material

12.1	Precautions on Static Electricity	12-2
12.2	Mobile Phase Characteristics	12-6

Chapter 1 General

Contents

1.1 Outline	1-2
1.2 Features	1-3

The LC-10AT is a solvent delivery unit with a dual-plunger, tandem-flow (i.e., 2-stage) pump which has been developed for improving accuracy and sensitivity of analysis in high-performance liquid chromatography.

In addition to the solvent delivery unit(s), a high performance liquid chromatography system requires an automated or manual injector, column, injector holder, detector and so on, which may be separately ordered through your sales representative.

This instruction manual covers operations for the LC-10AT and relevant accessories. For use of the other modules and special accessories, please refer to the instruction manuals for each.

- 1. Stable Solvent Delivery with Little Flow Pulsation** The LC-IOAT is a solvent delivery unit with a dual-plunger, tandem -flow pump for use in high-performance liquid chromatography. It features stable solvent delivery with little flow pulsation.
- 2. A Variety of Functions for Gradient Elution** Two gradient elution modes, low pressure and high pressure, are available. The high pressure gradient elution mode is highly accurate and allows a minimum time lag, while the low pressure gradient elution mode can handle solvent delivery of up to four liquids with a single pump. In addition, two types of gradient elution control are available depending upon the system setup. One is to use the SCL-IOA System controller, and the other allows control of gradient elution from the LC-IOAT alone.
- 3. Long Life Plunger Seal** The plunger seal is made of ultra-high-molecular-weight Polyethylene (UHMW-PE) which offers a low degree of wear.
- 4. Seals and Plungers Protected Against Buffer Salts** The plunger can be washed at the rear side of the plunger seal. When a buffer solution is used as the mobile phase, the mechanism effectively allows a longer service life of the plunger seal by preventing the formation of salt deposits.
- 5. Easy Maintenance** The LC-IOAT has a simple construction with a small number of components for easy maintenance. In addition, the replacement of plunger seals, plungers and check valves can be carried out from the front side of the LC-IOAT.

Chapter 2 Check Upon Delivery

Contents

2.1 Lists of Supplied Parts, etc.	2-2
--	-----

The LC-10AT is composed of the following parts.

- 1. LC-10AT Main Body
- 2. Standard Parts and Accessories

Part Name	Part No.	Qty
100V power cord or 200V power cord	071-60814-01 071-60814-06	1
Suction filter	228-18740-91	1
Instruction manual (Japanese version) or Instruction manual (English version)	228-30055 228-30056	1
Accessory kit (See the following page.)	228-32133-91	1
Drain tube kit	228-28161-91	1

3. Accessory Kit

Category	Part Name	Part No.	Qty
Tools	Wrench, 8 x 10mm	086-03006	2
	Wrench, 13 x 17mm	086-03017	1
	Allen wrench, 3mm	086-03804	1
	Seal remover	228-25142	1
	File (for cutting SUS pipes)	670-18928-02	1
Parts	Tie plate	228-18751	1
	Optical cable	070-92025-51	1
	Clip (for fixing tube)	046-00994-03	1
	Male nut, 1.6MN	228-16001	2
	Ferrules, 1.6F	228-16000	2
	Drain tube	228-25495-91	1
	SUS pipe, 1.6 x 0.3	670-10006-02	2m
	Syringe needle (for disposable syringe)	228-18216-91	1
	Syringe, 20ml	046-00038-01	1
	Pipe clamp	670-11610-01	1
	Spiral wrap	018-26002	0.2m
	Lid	228-17644	1
	Bottle cap	228-18887	1
	Bushing 1.6MN PEEK	228-18565	2
	PTFE Tube 1 x 0.5	016-37502	2m
	Polyvinyl Tube	016-31401	0.12m
	Remote cable	228-28253-91	1
Consumables	Plunger seal	228-21975	1

The above are packed together as the accessory kit (P/N 228-32133-91).

4. Drain Tube Kit

Category	Part Name	Part No.	Qty
Parts	Drain tube for solvent leakage	228-25162-03	1
	Drain tube, elbow	228-28094	1
	T joint	228-28162	1
	Straight joint	228-28163	1
	Pipe clamp	670-11610-01	1
	Instruction manual	228-10913	1

Chapter 3 Construction and Functions


Contents

3.1	Front Panel	3-2
3.2	Rear Panel	3-4
3.3	The Right Side and the Bottom	3-5

3. 1

Front Panel

1. Front Panel



d

Fig. 3.1

No.	Part Name	Description
①	Operator's panel	For operation, see item 5.1.
②	Front cover	Covers the pump head and flow lines, etc.
③	Drain valve knob	Turn this knob to open/close the flow line drain valve.
④	Front cover latch button	For opening the front cover.
⑤	Power switch	Push this button to turn on/off the instrument

2. Inside the Front Cover




Fig. 3.2

No.	Part Name	Description
①	Head holder	For mounting pump head. Includes a washing flow line for cleaning the plungers.
②	Pump head	The plunger reciprocates in the pump head to deliver solvent.
③	Drain pipe	Outlet of the washing flow line for the plungers
④	Check valve IN	Inlet check valve
⑤	Pump inlet	The reservoir suction filter (or selector valve outlet) is connected to this inlet.
⑥	Pump outlet	Usually for connecting the tubing to the injector.
⑦	Line filter	Traps dirt, etc. to prevent column clogging
⑧	Drain valve	The valve has a built-in pressure sensor, and is used to replace mobile phase or to release air from the flow line.
⑨	Drain port	Drain tube is connected to this port.
⑩	Check valve OUT	Outlet check valve
⑪	Tray	If washing solution leaks from the washing seal, this tray catches the washing solution.




Fig. 3.3

No.	Part Name	Description
①	External input/output terminal block	Used to make connection with external equipment.
②	Earth terminal	Used for grounding the instrument
③	Power cord connector	Used to connect a power cord.
④	Fuse holder	Two fuses are in the holder.
⑤	SOL.V connector	Used to connect a solvent selector valve unit FCV-I0AL/FCV-I1AL.
⑥	REMOTE connector	Used to connect with SCL-IOA or an additional LC-IOAT. See Section 4.9 "Installation of the High Pressure Gradient Elution System."




Fig. 3.4

①	Drain port of the pan in the front panel	For connecting the drain tube for solvent leakage. See Section 4.5.5. "Connection of the Drain Tube for Solvent Leakage."
②	Mounting screws for the gradient mixer	A mixer or a manual injector can be mounted here.
③	Fan vent	For cooling the interior. Do not block this opening.

Chapter 4 Installation

Contents

4.1	Conditions for Installation Site	4-2
4.2	Examples of System Configurations	4-3
4.3	Mounting of Multiple Units	4-6
4.4	Electrical Connections	4-7
4.5	Tubing of the Unit	4-8
4.6	Connection of the Degasser	4-12
4.7	Mounting the Mixer	4-16
4.8	Connection of the Injector, Column and Detector	4-20
4.9	Installation of a High Pressure Gradient Elution System	4-23
4.10	Installation of a Low Pressure Gradient Elution System	4-26

For correct and safe use of the instrument, proper care about the installation site should be taken.

1. **Environmental Conditions** To assure long service life and good performance of the instrument, avoid installation in a place exposed to corrosive gases or dust.
2. **Precautions on Ventilation and Fire**
 - Provide adequate ventilation when using flammable or toxic solvents as mobile phase.
 - Never use an open flame in the room particularly when flammable solvents are used.
3. **Electromagnetic Noise** Avoid installation in the vicinity of such equipment that generates a strong magnetic field. Use an additional noise filter if the power line has much noise.
4. **Installation Space Requirements** The LC-10AT is designed to be used on a table or stand, preferably a **solid** and flat surface with a depth of 60cm or more. See Section 4.2 “Examples of System Configurations” for typical configurations of systems and installation spaces.

<<Cautions>>

In selecting the installation site, due care should be taken with regard to the following items in order to assure the optimum performance of the instrument.

- (1) Ambient temperature should be within 5 to 35°C, and without extreme fluctuations.
- (2) Do not expose the instrument to the direct output from a heater or a cooler.
- (3) Do not expose the instrument to direct sunlight.
- (4) The installation site should be vibration-free.
- (5) Relative humidity should be within 45% to 85%.

5. Use in the Patient Environment

LC-10AT is not designed in consideration of the safety regulation for medical equipments.

Therefore LC-10AT can not be used in the patient environment.”

- * Patient environment: the area within a distance of 2.5m from patients.

The following shows examples of HPLC system configurations incorporating the LC-10AT.

1. Simple System

An example of comparatively simple system configuration. It includes a column oven, detector, etc. that make a minimum set of modules required for stable isocratic analyses.




Fig. 4.1

2. Low Pressure Gradient Elution System 1

An example configuration of a low pressure gradient elution system with manual sample injection. The system is highly economical and allows handling solvent delivery of up to four liquids for quaternary gradient elution.




Fig. 4.2

3. Low Pressure Gradient Elution System 2

A low pressure gradient elution system including an automatic injector controlled by the SCL-10A System Controller. Centralized control of each module is available through the use of the system controller allowing great ease of operation. The system can also readily be used for automatic analyses.




Fig. 4.3

4. High Pressure Gradient Elution System 1

An example configuration of a high pressure gradient elution system with manual sample injection. The LC-10AT high pressure gradient elution system allows highly accurate gradient analysis with a small time lag.




Fig. 4.4

5. High Pressure Gradient Elution System 2

A high pressure gradient elution system including an automatic injector controlled by the SCL-10A System Controller. In contrast with High Pressure Gradient Elution System 1, functions such as automatic sample injection, sample cooling, and selecting of mobile phase are available with this system, allowing flexible use for various types of analyses.




Fig. 4.5

Mounting of Multiple Units

The LC-10AT may be stacked one on another for use. The SPD-10A detector, CTO-IOA column oven, SIL-10A automatic injector, etc. may also be stacked on an LC-IOAT. When stacking the units, it is possible to **fix** units together to prevent them from falling over during an earthquake or the like. The units should be fixed using the accessory tie plate.

- (1) Unscrew the screws that are used to fixed the unit cases.
- (2) Use the same screws to fix the tie plate.




Fig. 4.6

Before connecting the power source, confirm that the following conditions are met:

- Voltage and capacity of the power outlet

Part number	228-31900-91, 92
90 – 130V~	100VA 50/60Hz
Part number	228-31900-93
200 – 250V~	100VA 50/60Hz

The instrument will not exhibit satisfactory performance if operated on unstable line voltage or insufficient power capacity. In addition, the power capacity required from the overall system should be considered when preparing the power source.

- Make sure that the power switch on the main unit is turned off.

1. Connection to the power outlet

Plug in the female end of the accessory power cord to the power connection on the back of the unit. Connect the male end to the power outlet.

<<Caution>>

- (1) The LC-IOA employs a three-pronged power cord including a grounding wire. Be sure to connect the power cord to a three-pronged power outlet including a protective conductor terminal so as to ensure proper grounding.
- (2) For prevention of electric shock and to ensure stable operation, be sure to ground the instrument.

Piping of the Unit

1. **Preparation of Reservoir** Prepare a reservoir of a capacity of 500ml or more.
2. **Removing the Front Cover**

Before connecting the flow lines of the unit, remove the front cover.

 - (1) Press the latch button to open the cover.
 - (2) The cover may be removed by pulling it while pressing the side of the cover as shown in the illustration below.




Fig. 4.7

3. **Connection of the Suction Filter**
 - (1) Unscrew the transportation bushing dosing the pump inlet by hand.

Replace this bushing when the unit is not to be used for a long period of time in order to prevent the entrance of dust into the flow line.
 - (2) Remove the 3D bushing from the suction filter tube.
 - (3) Cut the tube to an appropriate length according to the distance between the reservoir and the pump inlet.
 - (4) Put the filter in the reservoir as shown in the following illustration and run the tubing through the lid with four holes followed by the cap. Replace the 3D bushing back on the tube.
 - (5) Connect the 3D bushing of the suction filter line to the pump inlet.
 - (6) Fix the tube using the tubing clip as shown in the next illustration.





Fig. 4.8

<<Cautions>>

- Keep the filter element clean to prevent clogging.
- For constantly stable analysis, the mobile phase in the reservoir should be degassed. See Section 4.6 "Connection of the Degasser."

4. Connection of the Drain Tube

- (1) Remove the transportation cap mounted in the drain port.
- (2) Screw the 1.6MN bushing (PEEK) of the accessory drain tubing into the drain port.
- (3) Prepare a waste container and put the other end of the drain tube in it.
- (4) Fix the drain tube using the tubing clip as shown in the next illustration.




Fig. 4.9

5. Connection of the Drain Tube for Solvent Leakage

If liquid leaks in the unit, it accumulates in the pan on the front panel. Connect the drain tube for solvent leakage to lead the solvent to a waste container.

- (1) Fit the L-shaped drain tube to the drain port for solvent leakage at the right side of the unit.
- (2) Connect the drain tube for solvent leakage and the L-shaped drain tube using the straight joint.
- (3) Put the other end of the drain tube for solvent leakage in a waste container.

<<Caution>>

If the flow line of the drain tube for solvent leakage is located above the drain port for solvent leakage, the leaked solvent will not be discharged. Be sure to put the waste container lower than the unit, and direct the drain tube downward. Pour some water in the pan to check that waste will not back up but flow freely to the waste container,




Fig. 4.10

6. Tubing for Solvent Leakage in the LC-10A System

Each component of the LC-IOA series has a drain port for solvent leakage at the front right side of the unit, which discharges liquid if leakage occurs in the unit. Tubing for solvent leakage should be conducted when installing the system.

- (1) Fit the accessory L-shaped drain tube to the drain port for solvent leakage on each component.
- (2) Connect the L-shaped drain tube fitted on the lowest component with a drain tube using a straight joint, and lead it to a waste container.
- (3) For the next higher component, use a T joint to connect on L-shaped drain tube and the tubing from upper component, and lead the other to the waste container. For the third and higher components, make interconnection with the other components using a T joint and a drain tube cut in an appropriate length. Be sure to incline the L-shaped drain tubes downward and fix them on the sides of the units using the accessory pipe clips, if necessary.

Place the waste container lower than the lowest unit.




Fig. 4.11

Degassing of mobile phase is required for constantly stable analysis for the following reasons.

- (1) It prevents troubles that would occur from generation of bubbles.
- (2) It prevents unstable measurements that would occur due to variations in the concentration of dissolved gases.

Helium degassing and vacuum degassing through a resin membrane are popular in HPLC. Choose from these according to your needs. See Section 11.3 "Optional Units List."

As examples, the connection of DGU-2A and DGU-3A to LC-10AT will be described in this chapter. The former is for helium degassing, and the latter is for vacuum degassing using a resin membrane.

1. **Connection of the DGU-2A (helium degassing)** Use helium gas of high purity (99.995% or higher) for degassing.
 - (1) Set a pressure regulator (optional) to the helium gas cylinder and connect the pressure regulator and the DGU-2A with a carrier gas pipe (optional).
 - (2) For the reservoir, prepare a type LSI glass container. (A commonly used glass reagent container of a capacity of 500 or 1000ml).




Fig. 4.12

<<Caution>>

Never use a container with a crack or other flaw for a reservoir.

(3) Assemble the reservoir cap as shown in the illustration.




Fig. 4.13

Category	Part Name	Description
①	Lid	Outer lid of the reservoir (B3)
②	Bottle cap	Inner lid of the reservoir
③	Teflon packing	Packing for the inner lid
④	Teflon tubing	3.0 x 2.0 x 230
* ⑤	Teflon tubing	3.0 x 2.0, Connect to the pump inlet.
: ⑥	Filter element of the suction filter	Used to filter the mobile phase.
⑦	Helium gas bubbler frit	Helium gas bubbler frit
⑧	Teflon tubing	3 x 2 x 600 To be connected to an "OUT" (A, B, C, D) port of the DGU-2A.
⑨	Teflon tubing	3 x 2 x 600 To be connected to an "IN" (A, B, C, D) port of the DGU-2A.

: Parts ⑤ and ⑥ above are accessories supplied with the LC-10AT. Parts with other numbers are accessories supplied with the DCU-2A.

- (4) Put the mobile phase in the reservoir. Also put the filter element assembled in (3) into the reservoir, then firmly close the cap.




Fig. 4.14

- (5) Remove the caps set to DGU-2A ports to which you will make connection. Connect Teflon tubing led from the reservoir to OUT (A, B, C, or D) and IN (A, B, C, and D) ports of the DGU-2A, as shown in the illustration. Keep the unused ports covered with the caps.

Save the removed caps to use again to protect the degasser from entry of dust into flow lines when it is not to be used for a long time.




Fig. 4.15

2. **Ventilation of Helium Gas**
- (1) Prepare a glass container. Fill it with a solvent that is miscible with the mobile phase.
 - (2) Connect one end of the accessory Teflon tubing to the **VENT** port of the **DGU-2A**. Put the other end in the container and loosely cap the container with aluminum foil.
 - (3) Place the container near a fume hood, ventilating fan or window so as to release the gas from the **DGU-2A** to outside the room.




Fig. 4.16

3. **Plumbing of the DGU-3A**
- (1) Carry out plumbing from the reservoir to the **DGU-3A** and from the **DGU-3A** to the **LC-10AT** as shown in the illustration below. First, remove the caps set to **DGU-3A** ports which you want to use. Keep unused ports covered with caps. Save the removed caps to use again to protect the degasser from entry of dust into flow lines when it is not to be used for a long period.




Fig. 4.17

For gradient elution using the LC-IOAT, a special mixer (P/N 228-28000-91, optional) that is excellent in mixing solvents is available. Follow the instructions below for mounting the mixer and selecting mixer capacity.

1. Mounting Position

The mixer may be mounted in any of the following positions.

- (1) The right side panel of the LC-10AT
- (2) The inside left of the CTO-10A
- (3) The inside right of the CTO-10A

2. Mixer Capacity

The following three capacities can be selected. Select the right mixer capacity according to the analysis requirements.

- (1) 2.6ml
- (2) 1.7ml
- (3) 0.5ml

3. Selecting Mixer Capacity

Mixer capacity can be changed by altering the plumbing. When it is shipped, it is set to 2.6ml. To change the setting, follow the instructions below:

- (1) Remove the mixer cover.




Fig. 4.18

- (2) According to the desired capacity, change the plumbing as shown in the following illustration:

(A) Plumbing for 0.5m/

(B) Plumbing for 1.7m/

(C) Plumbing for 2.6m/
(Factory default)




Fig. 4.19

4. Mounting the Mixer on the LC-10AT

- (I) Move the pre-mixer section forward and fasten it by screws as shown in the following illustration:




Fig. 4.20

- (2) Unscrew the screws on the right side of the pump, and fix the mixer with the screws as shown below:





Fig. 4.21

- (3) Carry out piping between the mixer and injector referring to Section 4.9 "Installation of the High Pressure Gradient Elution System."
(4) Remount the mixer cover.



The pre-mixer joint sticks out of the cover.

Fig. 4.22

5. Mounting the Mixer in the CTO-10A

- (1) Mount the mixer on the inner wall of the CTO-IOA, using the accessory screws as shown in Fig. 4.24. The mixer is **usually** mounted at the position indicated in (A). Note, however, that when mounting the 7725 manual injector with position sensor on the CTO-10A, the mixer should **be** mounted at the position shown in (B).




Fig. 4.23

(A) Mounting on the inside left wall

(B) Mounting on the inside right wall (In contrast with mounting it on the left side, the mixer is mounted upside down.)





Fig. 4.24

The following example of connections apply to the basic system shown in the illustration below. Connections should be made according to each system, referring to the examples.




Fig. 4.25

For the method to mount the model 7725 manual injector and a column in the CTO-IOA column oven, refer to the instruction manual for the column oven. If not using the column oven, an injector holder and column holder (optional) are available for mounting the 7725 manual injector and column on the LC-10AT.

1. Connection of the Manual Injector

As an example, the connection between the LC-10AT and model 7725 manual injector will be described.

- (1) Cut the accessory SUS pipe 1.6×0.3 to the length required for piping from the pump outlet to the injector port 2.
- (2) Fit a 1.6MN bushing and 1.6F ferrule (that are supplied with the unit) at one end of the SUS pipe and connect it to the pump outlet. Also fit a bushing and ferrule (that are supplied with the manual injector) at the other end of the SUS pipe and connect that end to the model 7725 manual injector.

<<Cautions>>

- The discharging pipes are connected to ports 5 and 6 of the 7725 manual injector. The other ends of the discharging pipes should be located at the same height as the injector.
- For pipe connections to No. 1 to 6 of the manual injector model 7725, use the bushings and ferrules that are supplied with the model 7725.




Fig. 4.26

2. Connection between the Injector and Column

The following is an example of a typical connection between the manual injector model 7725 and column.

- (1) Cut the accessory SUS pipe 1.6 x 0.3 to the length required for tubing from the injector to the column.
- (2) Fit bushings and ferrules at both ends of the SUS pipe.
- (3) Connect these ends of SUS pipe to the injector and column.

<<Cautions>>

- In order to minimize sample broadening, make the tubing between the injector and column as short as possible.
- In order to eliminate dead volumes from the flow line, cut the pipe in such a manner that the cut face is perpendicular to the pipe axis.




Fig. 4.27

3. Connection between Column and Detector

As an example, the connection between column and the SPD-10A detector will be described.

- (1) Cut the ETFE tube (1.6 x 0.3) that is supplied with the SPD-10A to the length required to make connection from the column outlet to the cell inlet of the SPD-10A.
- (2) Fit bushings (1.6MN PEEK) at both ends of the ETFE tube.
- (3) Connect the ends of the ETFE tube to the column outlet and cell inlet (with a blue mark) of the detector as shown in the following illustration.

See cautions for "Section 4.8.2 Connection between the Injector and Column."




Fig. 4.28

High pressure gradient elution systems using LC-IOAT are possible in two types of configuration characterized by the control method.

- (1) A high pressure gradient elution system of two or three solvents controlled from the SCL-IOA system controller.
- (2) A two-solvent high pressure gradient elution system controlled by the LC-IOAT.

Installation for both types of configuration will be described in the following paragraphs.

1. Wiring for a System Controlled by the SCL-10A

- (1) Connect the REMOTE connector on the rear of the unit and the REMOTE connector on the rear of the SCL-IOA using the optical cable included in accessories for the unit. Channel 1 and 2 of the REMOTE connector are reserved for the SIL-IOA and FRC-10A. Channels 3 to 8 can be used for the pump(s).
- (2) The ADRS parameter of each LC-IOAT must be set to the channel number with which the unit is connected to the SCL-10A. (See Section 5.5 "Auxiliary Functions (AUX. FUNC).")




Fig. 4.29

2. Wiring for a System Controlled by the LC-10AT

- (1) Connect two LC-10ATs by the REMOTE connectors on the rear panels, using the accessory optical cable.




Fig. 4.30

- (2) Set 85 to the ADRS parameter of each LC-10AT.
(See Section 5.5 "Auxiliary Functions (AUX.FUNC).")

3. Piping of the High Pressure Gradient Elution System

As an example, piping for a high pressure binary gradient elution system with two LC-10ATs, a degasser and a mixer will be described.




Fig. 4.31

Carry out piping referring to the following connection method.

- (1) Prepare a reservoir. Degassing of mobile phases is necessary for stable gradient elution of good reproducibility. Carry out piping up to the pump inlet referring to Section 4.6 "Connecting the Degasser."

- (2) Mount the mixer referring to Section 4.7 "Mounting the Mixer."
- (3) After mounting a mixer on the lower LC-10AT, make connection of the LC-10ATs and mixer inlets as shown in the following illustration.

First remove the plugs on inlets A and B of the pre-mixer section, and carry out piping from the outlets of each LC-10AT to inlet A and B using the accessory SUS pipe (1.6 x 0.3).

For a ternary gradient system, remove the stop joint from inlet C and carry out piping from the inlet C to the pump outlet of the third LC-10AT.




Fig. 4.32

- (4) Connect the injector and the mixer outlet with the SUS pipe (1.6 x 0.3.)
- (5) Fix on the cover of the mixer referring to Section 4.7 "Mounting the Mixer."

This pump allows configuration of low pressure gradient elution systems that can handle up to four liquids, using the FCV-IOAL low pressure value module. The installation of the system is described in the following paragraphs.

1. Configuration of pressure gradient elution system

The simplest low pressure gradient elution system consists of the units and parts listed below:

- ① LC-10AT 1 ea.
- ② Mixer 1 ea.
- ③ FCV-IOAL Low pressure value module 1 ea.

2. Wiring to the FCV-10AL

- (1) Make sure that the power switch of the module is turned Off.
- (2) Using the cable included in the accessories for the FCV-IOAL, connect the FCV-IOAL and the SOL.V connector on the rear of the pump module.
- (3) Fasten the fixing screws on the connector of the module with a screw driver. The FCV-IOAL end of the cable has a flat cable connector; As such, fixing by screws is not done for the FCV-10AL connector.




Fig. 4.33

<<Cautions>>

- In order to reduce dead volume in the flow line and allow a minimum time lag for gradient elution, install the FCV-10AL as close to the LC-IOAT as possible.

3. Piping of the Low Pressure Gradient Elution System

Carry out piping referring to the following connecting method.

- (1) Prepare reservoir containers. Because degassing of the mobile phases is needed, carry out piping to the inlet joints of the FCV-10AL referring to Section 4.6 "Connecting the Degasser."

- (2) Connect the outlet of the 5-way branching block and the pump inlet of the module with the accessory tube (40cm length).

The FCV-10AL has four inlet ports (from A to D). If not every port is used, disconnect the pipe corresponding to the unused flow line from that inlet of the 5-way branching block. Then, be sure to fit the unused inlet with a plug (included in the accessories for the FCV-10AL).




Fig. 4.34

- (3) Remove the normal drain value line filter with the supplied wrench. Mount the low pressure line filter included with the FCV-10AL.




Fig. 4.35

- (4) Install the mixer, referring to the instruction in 4.7 "Mounting the mixer".
- (5) Connect the LC-IOAT and the inlet of the mixer as shown in the figure below.

Remove the plugs of inlets A and B of the pre-mixer, and connect the pump outlet of the LC-IOAT with inlet B by using the 1.6 x 0.3 SUS pipe provided with the mixer.

Cover inlet A with the stop joint included with the mixer.




Fig. 4.36

- (6) Connect the injector and the mixer outlet with the 1.6 x 0.3 SUS pipe.
- (7) Mount the cover, referring to the instruction in 4.7 "Mounting the mixer".

Chapter 5 Operation


Contents

5.1	Basic Operation	5-2
5.2	Constant Flow Solvent Delivery	5-8
5.3	Constant Pressure Solvent Delivery..	5-11
5.4	Creation of Time Programs and Execution	5-12
5.5	Auxiliary Functions (AUX.FUNC)	5-20
5.6	High Pressure Gradient Elution Mode	5-29
5.7	Low Pressure Gradient Elution Mode	5-32
5.8	Plunger Washing Method	5-33
5.9	Change of Mobile Phase	5-35

The operation of the unit is performed through the keys on the operator's panel. The status of operation may be checked on the display at any time. When using the SCL-IOA for control, see Section 7.2 "Control from the SCL-IOA."

1. Start-up

- (1) Press the power switch located at the lower left of the front panel to turn on the unit. The unit is turned off by repressing the switch.



Power switch in OFF state Power switch in ON state

Fig. 5.1

- (2) When the power is switched on, a memory check is carried out automatically. If there is nothing irregular in the memory, the ROM version number is displayed for several seconds and then the display shown in the following illustration will appear to enable operation. (Values displayed vary according to setting.) This is the initial state.




Fig. 5.2

<<Caution>>

If an error message is displayed after turning on the unit, see Section 9.2 "Error Messages" for appropriate measures.

2. **Description of the Display** The display has a screen and LED lamps. The function of each part is as follows:




Fig. 5.3

No.	Indication or Name	Function
①	flow/press	Indicates the set flow rate (ml/min) when the module is in the constant flow solvent delivery mode, or the set pressure value ($\times 10^5$ Pa) in the constant pressure solvent delivery mode.
②	pressure	The pressure value measured by the pressure sensor ($\times 10^5$ Pa)
③	p.max	Indicates the set upper limit pressure ($\times 10^5$ Pa)
④	p.min	Indicates the set lower limit pressure ($\times 10^5$ Pa)
⑤	pump	Pump operation indicator lamp The lamp is lit when the pump is in operation.
⑥	program	Programmed operation indicator lamp The lamp is lit when a program is running.
⑦	remote	REMOTE mode indicator lamp Flashes when the pump is controlled by the SCL-10A.
⑧	c.flow	Constant flow mode indication lamp The lamp is lit when the module is in the constant flow solvent delivery mode.
⑨	c.press	Constant pressure mode indication lamp The lamp is lit when the module is in the constant pressure solvent delivery mode.

Note. 10^5 Pa = 1 bar = 1kgf/cm²

3. Description of Keys

The module has 20 keys on the front panel which are used for operation and setting. Functions of each key are as follows.

pump	run	7	8	9
purge	CE	4	5	6
del	back	1	2	3
edit	func	0	•	Enter

pump Pump key

Key for start/stop of the pump. Press this key to start the pump. By pressing this key during operation, the pump stops. If it is pressed when purging, purging stops.

purge Purge key

Key for stadstop of purging. Press this key to start purging. Purging stops automatically after a period which can be set with the P.TIMER parameter. Purging can be **also** stopped manually by pressing this key during operation.

(del) Delete key

Deletes one line from the program during time program editing.

(edit) Edit key

Shifts to the time program edit mode.

run Run key

Key for stadstop of a time program. (If no time program has been set, this key has no effect.)

CE Clear entry key

- Sets the display screen to the initial state.
- If it is pressed during entry of numeric values, values entered so far are cleared.
- If it is pressed when an error message is on the display, the screen is cleared and the alarm is canceled.

back Back key

- Press this key to select the previous parameter during the editing of a time program.
- If this key is pressed during the setting of the parameters for basic operation, entry position for setting of the parameter moves to the left, i.e., to the next left parameter.
- If it is pressed during AUX.FUNC setting, the setting display shows the previous parameter in the list.

[func] Function key

- If this key is pressed during the setting of a parameter for basic operation, the entry position for setting of a parameter moves to the right, i.e., to the next right parameter.
- Press this key to select the next parameter for entry during the editing of a time program.
- If it is pressed during AUX.FUNC setting, the setting display shows the next parameter in the list.

[Enter] Enter key

Pressing this will register the value keyed in from the numeric keys.

0 - 9 Numeric keys

Used to key in numeric values for each settable parameter.


4. Settable Parameter List Settable parameters, their ranges, minimum values, initial values, and flow modes for basic operation are as follows:

Settable parameter	Settable range	Settable step	Default or Initial value	Mode
flow	0 - 9.999ml/min	0.001ml/min	0ml/min	Constant now solvent delivery
press	10 - 400 [x 10 ⁵ Pa]	1 [x 10 ⁵ Pa]	10 [x 10 ⁵ Pa]	Constant pressure solvent delivery
p.max	10 - 440 [x 10 ⁵ Pa] ^{Note}	1 [x 10 ⁵ Pa]	100 [x 10 ⁵ Pa]	
p.min.	0 - 440 [x 10 ⁵ Pa] ^{Note}	1 [x 10 ⁵ Pa]	0 [x 10 ⁵ Pa]	


Note: The maximum pressure value is limited to 220 [x 10⁵ Pa] if the flow is set to 5.001ml/min or higher.

5. Basic Setting Procedures The basic procedures used for entry of the values in the setting display are as follows:

- (1) Set the unit to the initial state (that appears when power is turned on). If the unit is not in the initial state, press [CE] to set it to the initial state. Then numerics corresponding to FLOW/PRESS in the display will flicker, indicating that the entry of numeric values to this parameter is enabled.



- (2) Key in a value with the numeric keys then press [Enter] to register the value.
- (3) If you want to set another parameter, press [func] to move the flashing field to the right.



- (4) To cancel a value just keyed in, press **[CE]** before pressing **[Enter]**.
- (5) To return the display to the initial state after setting, press **[CE]**.

6. Setting of the Maximum Pressure Value

To protect columns and other components in the flow line, set an upper pressure limit. If the pressure value measured by the pressure sensor exceeds the upper pressure limit, the limiter will be activated to stop solvent delivery automatically.

When the limiter is activated, an alarm beeps and the following message will be displayed.

flow/press	pressure	p.max	p.min
E R R O R P - M A X			

The setting procedures are as follows:

<<Setting Example>>

Example: To set p.max (upper pressure limit) to 150 [$\times 10^5$ Pa]:

- (1) Press **[func]** to enable the entry of p.max (the entry field flashes).

flow/press	pressure	p.max	p.min
0 . 0 0 0	0	1 0 0	0

- (2) Press **[1]** , **[5]** , **[0]** , and **[enter]**.

The following display appears, and the setting is complete.

flow/press	pressure	p.max	p.min
0 . 0 0 0	0	1 5 0	0

7. Setting the Lower Limit Pressure

To prevent drawing air into the flow line when the mobile phase in a reservoir has been exhausted, or for safety measures in the case solvent is leaking from the flow line, a lower pressure limit should be set. The limiter will be activated if the pressure is still lower than the minimum pressure even one minute after starting solvent delivery. When the limiter is activated, the alarm beeps, the solvent delivery stops automatically and the following message will be displayed.

flow/press	Pressure	p.max	p.min
------------	----------	-------	-------

ERROR P - MIN			
---------------	--	--	--

The setting procedures are as follows:

<<Setting Example>>

Example: To set p.min (minimum pressure) to 20 [$\times 10^5$ Pa]:

- (1) Press **func** to enable the entry of p.min (the entry field flashes).

flow/press	Pressure	p.max	p.min
0 . 0 0 0	0	1 5 0	0

- (2) Press **2** , **0** , and **Enter** .

The following display will appear, and the setting is complete.

flow/press	Pressure	p.max	p.min
0 . 0 0 0	0	1 5 0	2 0


<<Caution>>

A value of 0 disables the low pressure limiter function.

1. An Operation Example of Constant Flow Solvent Delivery

The following is an example of operation for constant flow solvent delivery:

- (1) Prepare a mobile phase and pour it into a reservoir followed by the suction filter.
- (2) Carryout piping of flow line referring to Section 4.8 “Connection of the Injector, Column and Detector.”
- (3) Turn the drain valve knob counterclockwise by 180° to open the drain valve. Make sure that the drain tube is put into a waste container.



Drain valve knob

Fig. 5.4

- (4) Press the power switch.
- (5) Check that the pressure value indicated on the display is in the range from -3 to 3 [$\times 10^5$ Pa]. If not, carry out zero adjustment for the pressure sensor. (See Section 5.5 “AUX.FUNC Functions.”)
- (6) Check that the upper pressure limiter and lower pressure limiter are appropriate. (See Section 5.1 “Basic Operation.”)
- (7) Press **[purge]**, then the pump will operate with a flow rate of about 9.9ml/min.

<<Note>>

In normal operation, solvent delivery at the preset flow rate starts immediately after the **[pump]** or **[purge]** key is pressed.

In the two cases (a. and b.) below, however, solvent is delivered at a low rate for a few seconds until the home position of the pump is detected. After the home position is detected, the solvent delivery rate increases to the preset value.

- a. The **[pump]** or **[purge]** key is pressed for the first time after power-on.
- b. The **[pump]** or **[purge]** key is pressed for the first time after the high pressure limiter is activated with a limit value of 220 [$\times 10^5$ Pa] or higher.

- (8) Observe the condition of solvent discharged from the drain tube. If the liquid is being discharged smoothly without bubbles, the unit is operating normally.




Fig. 5.5

<<Caution>>

If solvent does not flow, use the disposable syringe and syringe needle to draw the mobile phase through to clear the flow path as shown in the illustration below.




Fig. 5.6

- (9) Press **purge** or **(pump)** to stop the pump.
- (10) Set the desired flow rate.
Example: When you want to set 1ml/min, press **1** and **Enter**.

flow/press	pressure	p.max	p.min
1.000	0	100	0

- (11) Turn the drain valve knob clockwise to the stop to close the drain valve.
- (12) Press **pump** to operate the pump.
- (13) Check that the pressure at the pump outlet increases and becomes stable.

<<Caution>>

If the drain valve is left open, mobile phase may siphon out from the drain port.

Keep the drain valve closed except for the time when draining or purging the system.

1. **An Example of Operation for Constant Pressure Solvent Delivery** The following is an example of the operation for constant pressure solvent delivery.

- (1) Carry out the operating procedures from (1) to (9) of Section 5.2 "Constant Flow Solvent delivery."
- (2) Change the solvent delivery mode from the constant flow solvent delivery mode to the constant pressure solvent delivery mode. First, press **func** until the following display appears on the screen.

flow/press	pressure	p.max	p.min
0.000	MODE	CHANGE	

- (3) Each time **Enter** is pressed, the solvent delivery mode changes. When the unit is in the constant flow solvent delivery mode, the **c.flow** lamp is lit, and in the constant pressure solvent delivery mode, the **c.press** lamp. Now set the unit so that the **c.press** lamp is lit.

pump	prog run	remote	cflow	cpress
☐	☐	☐	☐	☒

- (4) Press **CE** after shifting the mode to return to the initial setting display. Check that the flow/press entry field is ready for entry of values (flashing).

flow/press	pressure	p.max	p.min
10	0	100	0

- (5) Set a desired pressure with the numeric keys.
Example: When you want to set 20 [$\times 10^5$ Pa], press **2**, **0**, and **Enter**.

flow/press	pressure	p.max	p.min
20	0	100	0

- (6) Turn the drain valve knob clockwise to the stop to close the drain valve.
- (7) Press **pump** to operate the pump.
- (8) Check that the pressure increases and becomes constant at approx. 20 [$\times 10^5$ Pa].

Using the time program function, parameters such as flow rate may be changed automatically at user-defined time intervals.

The created time program may be stored as a file and it can be combined with other time program files for sequential running.

1. Command List

The following commands can be used in a time program.

Command	Description	Settable range	Remarks
FLOW	Flow rate (effective only in the constant flow solvent delivery mode)	0 – 9,999ml/min	The minimum step is 0.001ml/min
PRESS	Pressure (effective only in the constant pressure solvent delivery mode)	10 – 400 [$\times 10^5$ Pa]	The minimum step is 1 [$\times 10^5$ Pa]
SV	Open/close solenoid valve in the FCV-IOAL or FCV-11AL(S) (optional).	0, 1, 2, 3, 123	See Note 1.
EVENT	Event output ON/OFF	0, 1, 2, 12	Select one of the four numeric values.
LOOP	Repeats a program from the first step.	from 0 to 255 Value 0 repeats a program 256 times.	The minimum step is 1.
STOP	Ends a program.	—	
GOTO	Runs programs by linking files. (Up to 10 files)	0 – 9	
BCONC	Concentration of solvent B (effective only when SYS = 2 or 4)	0 – 100% Note 2	The minimum step is 0.1%.
CCONC	Concentration of solvent C (effective only when SYS=4)	0 – 100%	Same as above
DCONC	Concentration of solvent D (effective only when SYS=4)	0 – 100%	Same as above

Note 1. If the FCV-IOAL has been specified by setting the FCV TYPE parameter, the settable valve is 1, 2, 3, or 4. If the FCV-11AL(S) has been specified, the settable value is 0, 1, 2, 3 or combinations of those values.


Note 2. The concentration of solvent A is determined by subtracting BCONC, CCONC and DCONC from 100.

$$A\text{CONC} = 100 - (B\text{CONC} + C\text{CONC} + D\text{CONC})$$

2. Description of the Display

To create a time program, set the unit to the edit mode, then proceed with programming. The procedures are as follows:


(1) Press **[edit]** and the following display will appear.



- ① Number of steps already programmed
- ② Abbreviation of the word "used"
- ③ Number of free steps remaining
- ④ Abbreviation of the word "left"


The above example shows that 10 steps of the time program have already been written and that there are 310 remaining steps.

(2) Now press **[Enter]** and the following display will appear.




- ① Elapsed time (minutes and decimal fraction) from the start of program
- ② Command or function name
- ③ Set value or argument

(3) Press **[Enter]** once again; the first step of the program will appear.



For detailed creation of a program, see the following Section 5.4.3 "Creation of Time Program."


3. **Creation of Time Program** The sequence of creating a time program starting from the initial state is as follows. Arbitrary values are shown for example only.



4. An Example of Creating a Time Program

An example of creating a time program will be described taking the case of varying the flow rate. (Let's assume the initial flow rate is 1ml/min.)

TIME	FUNC	VALUE
15.0	FLOW	2.00
20.0	FLOW	5.00
30.0	STOP	



- (1) Press **CE** to return to the initial display.

flow/press	pressure	p.max	p.min
1.000	0	100	0

- (2) Press **edit**

flow/press	pressure	p.max	p.min
0	U	320	LE

- (3) Press **[Enter]**.

flow/press	pressure	p.max	p.min
TIME	FUNC	VALUE	

- (4) Enter the time (minutes, tenths and hundredths) for the initial step. Press **1**, **5** and **[Enter]**.

flow/press	pressure	p.max	p.min
15.00	FLOW		

- (5) Select an instruction. When an instruction other than FLOW is displayed, press **func** until FLOW is displayed and then press **(Enter)**.

flow/press	pressure	p.max	p.min
15. 00	FLOW	V	U
		A	E
		L	

- (6) Enter the flow rate. Press **2** and **(Enter)**

flow/press	pressure	p.max	p.min
15. 60	FLOW	2.	00

- (7) Like the procedure from (4) to (6) above, set the second step. Press **2**, **0**, **(Enter)**, **(Enter)**, **5** and **(Enter)**.

flow/press	pressure	p.max	p.min
20. 00	FLOW	5.	00

- (8) Enter a step for the STOP instruction. After entering a time value by the same procedure as (5), press **func** to display STOP, then press **(Enter)**.

flow/press	pressure	p.max	p.min
30. 00	STOP		

- (9) Press **CE** to end the editing of the time program. A time program having three steps has now been completed in this example.

flow/press	pressure	p.max	p.min
3	U	317	LE

- (10) By pressing **CE** once again, the initial display shown in (1) will return.

5. Deleting a Step

Display the step you want to delete and press **[del]**.

The following is an example of deleting the first step of the program created in Section 5.4.4 “An Example of Creating a Time Program.”

- (1) As in the creation of a program, display the program step you want to delete.

flow/press	pressure	p.max	p.min
1 5. 0 0	FLOW	2.	0 0

- (2) Press **[del]**.

flow/press	pressure	p.max	p.min
2 0. 0 0	FLOW	5.	0 0


Now the first step of the program is deleted and the second step is displayed. (This is now the first step of the program.)

6. LOOP Instruction

Use of the LOOP instruction permits repeating a program at a user-defined time interval for a specified number of times.

	TIME	FUNC	VALUE
①	15.00	FLOW	2.00
②	20.00	FLOW	5.00
	30.00	LOOP	3

With this setting the program will be executed four times; The initial execution plus 3 looped executions for a total time of 120 minutes as illustrated here:



<<Cautions>>

- All steps set (timewise) after the LOOP command are ignored, with the exception of a GOTO command step.

- After the LOOP is completed, the program stops automatically. If a GOTO instruction exists after LOOP instruction, the pump parameter conditions are maintained until the time of the GOTO instruction; at which time the program executes the GOTO instruction and stops. (See below)
- A number up to 255 may be set as VALUE for the LOOP instruction. Note, however, that setting 0 means LOOP 256 times. Delete the LOOP step to disable looping.
- When setting multiple steps, it is not necessary to enter steps in the correct time sequence. The unit sorts the steps automatically.
- Note that a STOP instruction at the last step is usually required except when you want to execute the time program endlessly, or link program files by the GOTO instruction.

7. Start/Stop

To start or stop a time program that has been completed, follow the procedures below.

(1) Starting the Program

Press **run**.

The **prog run** lamp lights up and the program starts.

(2) Stopping the Program

To stop the program, two methods are available. One is to forcibly stop it by pressing the STOP key, and the other is to stop by means of a STOP instruction step set in the program.

The case of a forced stop is described here.

Press **run**.

The **prog run** lamp is turned off and the program stops.

8. GOTO Instruction (File chaining via Time Program execution)

The GOTO instruction stops the program and switches the active file from the present file to another file whose number is specified in the GOTO instruction argument.

After the file is switched, the instrument parameters are set to the initial conditions specified in the new file.

The automatic start of the new file's time program (if extant) is possible by carrying out the following setting.

- (1) Make connection between the external Input/Output terminals as illustrated here;




Fig. 5.7

- (2) Set the EVENT parameter to 0 in the initial conditions of the temporary file. Or carry it out at sometime in the time program *before* the GOTO instruction is executed.
- (3) Set the EVENT parameter to 1 in the initial conditions of the new file.

As soon as the new file is installed and activated (**by** the GOTO command), the relay closure (EVENT 1) will trigger the start of that file's time program.

Thus, each of the 10 files available can be installed and activated in any desired sequence ("chaining").

The unit has a number of auxiliary functions (AUX.FUNC). Use of these functions permits a variety of operations to be performed by this device. It is also possible to control the unit from an external device via some of these functions

1. AUX.FUNC List

The following table lists the AUX.FUNC functions. Details of each are on following pages

Type	Name	Function	Factory default value
1	BCONC	Concentration of solvent B in the gradient elution mode.	0
1	CCONC	Concentration of solvent C in the gradient elution mode.	0
1	DCONC	Concentration of solvent D in the gradient elution mode.	0
1	FILE NUM	Selects file.	0
1	FILE COPY	Copies a file.	—
2	FILE DEL	Deletes a file.	—
1	SV	Sets the condition of solenoid valve unit FCV-10AL or FCV-1IAL.	1
1	EVENT	Sets the state of relay contact point output.	0
1	EXT-S	Lets the EVENT 1 output work as a start signal for a time program. It also lets the EVENT 2 output work as a stop signal for external equipment in occurrence of an error.	0
1	S-PROT	Reduces flow rate without stopping the pump when the pressure exceeds the set value of p.max.	0
2			
1	MON FLOW	Monitors solvent delivery flow rate in the constant pressure solvent delivery mode.	0
1	MON TIME	Monitors elapsed time when running a time program.	0
1	MON REV	Monitors accumulated number of pump revolutions.	0
1	MON ID	Monitors pump's ID.	0
1	P-SET	Used when replacing the plunger and plunger seal.	—
1	COMP	Used for fine adjustment of compensation of solvent compressibility.	0.45
2	ZERO ADJ	Carries out zero adjustment for monitoring pressure.	—
2	CLOSE KEY	Disables key entry.	—
1	RANGE	Sets the span for pressure signals in the recorder output.	10
1	SYS	Specifies isocratic or gradient system in use.	1
1	FCV TYPE	Sets the type of solvent selector valve connected to the SOL.V connector on the rear of the unit.	0
1	LOCAL	Selects whether control is made from the system controller (=1) or by the pump (=0).	0
1	P TIMER	Sets the period of purging	3
1	ADRS	Sets the address to which the unit is connected.	1

- BCONC is displayed only when SYS = 2 or 4.
- CCONC and DCONC are displayed only when SYS = 4.

The type column in the list shows the type of operating method.
Type 1: Requires a number argument to execute; key in a value with the numeric keys then press **Enter** .

Type 2: Press **Enter** to directly execute the function.

Type 3: Monitor; The current value (condition) is displayed.

2. Setting Procedures for

AUX.FUNC**BCONC, CCONC, DCONC****(Concentration setting)**

flow/press	pressure	p.max	p.min
0.000	B C O N C	0	0

The concentration (volumetric ratio) of solvent B, C and/or D is set for the gradient elution mode. Key in the concentration of solvent B, C or D with the numeric keys, then press **Enter**. (Unit: %)

FILE NUM**(File selection)**

The unit permits creation of up to 10 files of programs that can be stored in the memory. A file may be selected (i.e., installed) through this parameter.

flow/press	pressure	p.max	p.min
0.000	F I L E	M U M	1

Key in the desired file number with the numeric keys, then press **Enter** to install that file as the active file.

FILE COPY**(Copy of a file)**

flow/press	pressure	p.max	p.min
0.000	F I L E	C O P Y	1

The program contents of the currently selected file (source file) are copied to a (destination) file of selected number.

Key in the destination file number with the numeric keys, and press **Enter**.

FILE DEL**(Deleting a program)**

flow/press	pressure	p.max	p.min
0.000	F I L E	D E L	

The time program of the currently selected file is deleted.

Press **Enter**. A new program can now be created in this file.

SV**(Setting the state of solenoid valves)**

flow/press	pressure	p.max	p.min
0.000	S V	1	

Any one of mobile phases can be selectively delivered by using the optionally available FCV-10AL or FCV-11AL as a selector valve

at the pump inlet. Key in the value corresponding to the desired mobile phase with a numeric key, then press **Enter** to switch the value.

Solenoid Valve type	Argument	Selected mobile phase
FCV-10AL	1	Solvent A
	2	Solvent B
	3	Solvent C
	4	Solvent D

Solenoid valve type	Argument	Selected mobile phase
FCV-11AL	0	All channels are set to the A side.
	1	Channel 1 is set to the B side.
	2	Channel 2 is set to the B side.
	3	Channel 3 is set to the B side.
In use, combine the arguments to select multiple channels. Example: By setting SV 123, channels 1, 2 and 3 are all set to the B side.		

<<Note>>

Before setting SV, specify the kind of solenoid valve unit installed with the FCV TYPE parameter. See next.

FCV.TYPE

(Specification of type of solenoid valve module)

flow/press	pressure	p.max	p.min
0.000	FCV	TYPE	1

The optional solvent selector valve to be connected to the SOL.V connector on the rear of the unit is specified with this parameter. Key in the appropriate value with a numeric key, then press **Enter** to register the flow control value type.

Argument	Solenoid valve module
0	FCV-10AL
1	FCV-11AL

EVENT

(Setting EVENT relay states)

flow/press	pressure	p.max	p.min
0.000	EVENT	0	

ON (contacts closed) and OFF (contact open) of the EVENT relay contacts (on the rear of the module) are set.

Key in the desired value with the numeric keys then press **Enter** to activate/deactivate the relay(s).

Argument	EVENT 1 output	EVENT 2 output
0	Relay 1 OFF	Relay 2 OFF
1	Relay 1 ON	Relay 2 OFF
2	Relay 1 OFF	Relay 2 ON
12	Relay 1 ON	Relay 2 ON

EXT-S (External signals)
(Function setting for the
EVENT relay terminals)

flow/dress pressure p.max p.min

0.000 EXT - S 1

Sets control mode for the EVENT output (relays 1 and 2).
key in the argument with a numeric key, then press **[Enter]**.

Argument	Function
0	Relays are controlled by the value set to EVENT. (Normal default condition)
1	Relay 1 (EVENT 1) is used as a start output signal when the time program starts. (Event 2 operates normally)
2	Relay 2 (EVENT 2) is used as an error output signal. (Event 1 operates normally)
3	Combination of functions 1 and 2. (Normal EVENT 1/2 operations are disabled.)

<<Note>>

The use of the EXT-S function disables normal operation of the EVENT parameter relevant to the corresponding EVENT terminal. Use care.

S-PROT
(Setting the system
protection function)

1 s pressure p a p.min

0.000 S - PROT 0

S-PROT controls solvent delivery in such a manner that when the P.MAX limiter is activated, flow rate is reduced by one half without stopping the pump, until the flow rate becomes lower than the P.MAX value. (Normally, exceeding P.MAX causes the pump to stop.)

Key in the argument with a numeric key, then press **[Enter]**.

Argument	Function
0	Cancels the system protection.
1	Activates the system protection.

To cancel the P.MAX error alarm, press **ICE**).

MODE CHANGE

(Selecting the solvent delivery mode)

flow/press	pressure	p.max	p.min
0.000	MODE	CHANGE	

Switches between the constant flow solvent delivery mode and the constant pressure solvent delivery mode.

Press **(Enter)** to select the desired solvent delivery mode. The selected mode is indicated by lighting of either the C.PRESS or C.FLOW LED.

<<Note>>

The pump should not be in operation when switching the solvent delivery mode.

MON FLOW

(Displaying flow rate in the constant pressure solvent delivery mode)

flow/press	pressure	p.max	p.min
0.000	MON	FLOW	0

Displays the approximate flow rate when in the constant pressure solvent delivery mode.

Key in the argument with a numeric key, then press **(Enter)**.

Argument	Function
0	Cancels the flow rate display function.
1	Activates the flow rate display function.

MON TIME

(Monitoring the elapsed time of the time program)

flow/press	pressure	p.max	p.min
0.000	MON	TIME	0

Argument	Function
0	Cancels monitoring of the elapsed time of the time program.
1	Activates monitoring of the elapsed time of the time program.

MON REV

(Monitoring the accumulated number of pump revolutions)

flow/press	pressure	p.max	p.min
0.000	MON	REV	0

The accumulated number of pump revolutions is monitored. It counts up to 16,777,215, then resets to 0.

Key in the argument with a numeric key, then press **(Enter)**.

Argument	Function
0	Cancels monitoring the accumulated number of pump revolutions.
1	Activates monitoring the accumulated number of pump revolutions.

The accumulated number of pump revolutions serves as an indicator for seal replacement.

It is recommended to replace the seals every 2.5 million revs, i.e., approx. every 2000 hours (at 1m³/min, 100 x 10⁵ MPa using water).

MON ID
(ID indication for remote control)

flow/press	pressure	p.max	p.min
0.000	MON	ID	0

If the unit is connected to the SCL-IOA, the pump connection address (A, B or C on the SCL-IOA) is indicated in the initial display. Key in the argument with a numeric key, then press **Enter**.

Argument	Function
0	Disables the ID indication.
1	Enables the ID indication.

flow/press	pressure	p.max	p.min
0.000	PUMP	A	

If the SCL-IOA is not connected, the following will be displayed.

flow/press	pressure	p.max	p.min
0.000	NOT	LINKED	

P-SET
(Plunger set)

flow/press	pressure	p.max	p.min
0.000	P - SET	0	

Moves the plungers and stops them in the home position for seal replacement.

Key in the argument with a numeric key, then press **Enter**.

Argument	Function
1	The left and right plungers move until they reach the same position, and stop. Replace the plunger seals at this position.

COMP

(Fine adjustment of compensation for fluid compressibility)

flow/press	pressure	p.max	p.min
0.000	COMP	0.45	

The LC-10AT employs compressibility compensation in order to reduce the fluctuation of the pressure due to solvent compressibility. Fine adjustment of compensation is possible by setting the compressibility factor for the solvent used.

For low pressure solvent delivery, changing of the parameter is not usually required, but for delivering solvents of greater compressibility such as hexane or methanol under a pressure of 200 [$\times 10^5$ Pa] or higher, this parameter should be adjusted accordingly. Key in the compressibility factor (GPa)⁻¹ as the set value with the numeric keys, then press **Enter**. The default value is 0.45.

Mobile phase	Compressibility (GPa) ⁻¹
Water	0.45
Acetone	1.24
Methanol	1.25
Hexane	1.6

<<Note>>

Compressibility factors vary slightly with temperature; The values shown are for room temperature (~25°C). However, correction is usually not required for the usual range of temperatures encountered in HPLC.

ZERO ADJ

(Zero adjustment for pressure indication)

flow/press	pressure	p.max	p.min
0.000	ZERO ADJ		

Output of the pressure sensor is reset to zero.

Press **Enter**. The zero adjustment is carried out. Before using this function, open the drain valve and set the conditions so that no pressure is applied.

CLOSE KEY

(Disabling key entry)

flow/press	pressure	p.max	p.min
0.000	CLOSE	KEY	

Disables key entry.

Press **[Enter]**. The unit is set to the state in which key depression is ineffective.

To cancel this key lockout effect, press **[CE]** and **[func]** simultaneously.

RANGE

(Changes the output pressure signal FS range)

flow/press	pressure	p.max	p.min
0.000	RANGE	10	

Sets the fullscale value for the output signal (through the “pressure” terminal on the rear of the unit). Key in the factor with the numeric keys, then press **[Enter]**.

Full scale = (factor) × 50[× 10⁵ Pa]. Full scale voltage is 1mV.

Examples: When factor = 1, the full scale is 50[× 10⁵ Pa].

When factor = 10, the full scale is 500[× 10⁵ Pa].

<<Note>>

The factor must be an integer value from 1 to 10; Fractional values are not accepted.

SYS

(Setting the system parameter)

flow/press	pressure	p.max	p.min
0.000	SYS	2	

Sets a system parameter.

Enter the set value with a numeric key, then press **[Enter]**

Argument	Function
1	Set this value when using the unit individually or when control is made externally.
2	Set this value when using the unit as a control pump in a two-module high pressure gradient elution system.
4	Set this value when using the pump in a low pressure gradient elution system.

<<Caution>>

Values other than the above should not be set.

LOCAL
(Setting the control mode)

flow/press	pressure	p.max	p.min
0 . 0 0 0	L O C A L	0	

By setting this function when the unit is linked with SCL-IOA, the unit can be operated independently.
Key in the argument with a numeric key, then press **Enter**.

Argument	Function
0	Control is made by the SCL-IOA.
1	The module is operated independently (local mode).

P TIMER
(Setting the period of purging)

flow/press	pressure	p.max	p.min
0 . 0 0 0	P T I M E R	3 . 0	

Sets the period of purging (min).
Key in the required purge time with numeric keys, then press **Enter**.
A value from 0.1 to 9.9 (min) can be selected.

ADRS
(Setting the REMOTE address)

flow/press	pressure	p.max	p.min
0 . 0 0 0	A D R S	3	

If the unit is connected to the SCL-10A system controller or another LC-10AT, the connection address should be stated with this parameter.
Key in the address number with the numeric keys, then press **enter**.

Address	Function
3 - 16	In the case when the unit is connected to the SCL-IOA, ADRS should be set to the port number of the SCL-IOA to which the unit is connected.
85	In the case when the unit is connected to another LC-IOAT, ADRS should be set to 85.

Values other than the above should not be set.

High pressure gradient elution is possible by controlling the 2 or 3 LC-10ATs in the system from the SCL-10A system controller, or by connecting two LC-10ATs with each other in such a manner that one of them (the “master”) controls the other (the “slave”). The latter setup will be described here.

For the operating method for control by the system controller, see the System Controller Instruction Manual.

1. Preliminary Setting

- (1) Set system parameters for the master LC-10AT (called pump A) that controls the other (pump B) as follows:

SYS = 2

ADRS = 85

- (2) Set system parameters for the slave LC-10AT (called pump B) that is to be controlled as follows:

SYS = 1

ADRS = 85

For the procedures for setting the system parameter, see Section 5.5 “Auxiliary Functions (AUX.FUNC).”

- (3) Set the upper limit pressure (P.max) and lower limit pressure (P.min) for each pump. Values P.max and P.min set at each pump are effective.

From now on operation is carried out from the master pump A only.

It is advised to lock the keys of pump B to prevent misoperation.

2. Setting Initial Conditions

- (1) Set the total flow rate of the mobile phase (through the column, i.e., total flow rate of pump A and B) as the flow rate of pump A.

- (2) Refer to Section 5.5 “Auxiliary Functions (AUX.FUNC),” enter the concentration value for solvent B (unit: %) to parameter BCONC.

The concentration of solvent of pump A (ACONC) is as follows.

$$ACONC(\%) = 100(\%) - BCONC(\%)$$

For the BCONC parameter, settable range is from 0 to 100%, and the minimum step is 0.1%.

- (3) By pressing pump on pump A, pumps A and B start solvent delivery simultaneously.

3. High Pressure Gradient Delivery Program

By setting the concentration of solvent B in a time program, it is possible to change concentration of solvent A and B linearly (inverse percentage ratio) on a time basis.

4. **An Example of Program Setting**
- (1) Call the time program editor screen. (Press **edit** .)
 - (2) Enter a time.
 - (3) Press **func** until BCONC is displayed on the setting screen, and then press **Enter** .
 - (4) Enter concentration.


The following is an example of setting for changing the concentration through a time program.

(See Section 5.4 "Time Program Creation and Execution".)

Initial value for BCNC: 0.0

TIME	FUNC	VALUE
10.00	BCONC	25.0
20.00	BCONC	100.0
30.00	BCONC	100.0
30.10	BCONC	0.0
40.00	STOP	

The program above results in a binary gradient elution in which the concentration of the solvent B can be graphically illustrated, as follows.



Note that the concentration delivered at any given time is based on a linear interpretation of proportionality between each defined point. (Nonlinear curves are available by control with the SCL-10A.)

5. **Operation of Pump B**
- (1) Pump A controls flow rate of pump B, but all parameters of pump B except the flow rate (FLOW) are effective.
 - (2) The **pump** key of each pump is always effective, so that pressing the **pump** key starts or stops delivery even if it is controlled by a time program of pump A.
 - (3) When not delivering, pressing the **purge** key will allow for purging.

<<Caution>>

Do not press the [RUN] key on pump B.

If [RUN] key is pressed, the time program (if any) in pump B will also **start**.

This may result in faulty solvent delivery.

The module by itself allows running of a binary, ternary or quaternary low pressure gradient elution. (It is also possible to carry out the setting from the SCL-10A system controller.) The following is a description of setting a quaternary (A, B, C, and D) gradient elution using a single LC-1OAT equipped with a selector valve and mixer, as described in Section 4.10.

1. Preliminary Setting

- (1) Enter 4 at the system parameter **SYS** of **AUX.FUNC**.
- (2) Enter total flow rate of the mobile phase (sum of flow rates of four liquids A, B, C and D) for the set value of the flow parameter.

Set the flow rate to 5ml/min here as an example.

- (3) Next, set concentration parameters. As an initial operation is being done here, setting is made so that each flow rate is equal to fill **all** the flow lines with the liquids. Set each of **BCONC**, **CCONC** and **DCONC** to 25%. (See Section 5.5 “(AUX.FUNC)”.

Concentration of liquid A (**ACONC**) is as follows.

$$\text{ACONC}(\%) = 100(\%) - \text{BCONC}(\%) - \text{CCONC}(\%) - \text{DCONC}(\%)$$

For parameters **BCONC**, **CCONC**, and **DCONC**, the setting range is from 0 to 100%, and the minimum step is 0.1%.

<<Note>>

Key input is automatically rejected if an attempt is made to enter a concentration to **BCONC**, **CCONC** or **DCONC** which causes the total conc. to exceed 100%.

2. Initial Operation

- (1) Open the drain valve.
- (2) Press (pump) to start the pump.
- (3) Draw the mobile phase from the pump outlet to fill the flow line with the mobile phase.
- (4) Make sure that no bubbles are in the flow line.
- (5) Purge the flow line with the mobile phase for about 10 minutes.
- (6) Press (pump) to stop the pump, and close the drain valve.

3. Low Pressure Gradient Operation



- (1) Set the (total) flow rate and the concentration(s) to the prescribed values, then press (pump) to start solvent delivery.

4. Low Pressure Gradient Delivery Program

By setting concentration of liquids B, C, and D at certain times in a time program, it is possible to vary the concentrations of liquids A, B, C, and D in linear proportions on a time basis.

1. Connections

- (1) Connect the accessory Teflon tubes with the two Teflon tubes at the bottom of the head holder by using the accessory vinyl tubes.
- (2) Place the open end of one of the Teflon tubes into a waste container.

**Fig. 5.0****Fig. 5.9**

2. Operation

- (1) Fill the disposable syringe with washing solution and insert the needle end into the rinsing tube.
- (2) Force the washing solution through the washing flow line. It is possible to conduct the washing irrespective of whether the pump is running or not.
- (3) Typical periods for manual washing are as follows.
Phosphoric acid buffer solution (low concentration)
: A few times/day
Boric acid buffer solution (low concentration)
: Several times/day

<<Caution>>

Use of a buffer solution as the mobile phase will produce crystalline salts upon evaporation of the solution, and which can damage the plunger and plunger seal and shorten their service life. When using a buffer solution, accordingly, it is recommended to wash the flow line frequently. Continuous washing can be done using an inexpensive peristaltic pump or the like to pump water through the wash line at a few ml/min.

As a mobile phase may be or may not be miscible with another, or as the use of buffer solutions requires extra care in handling, change of a mobile phase should be conducted correctly according to the following description. If you also change the column to a different type, remove the first column and carry out the procedure before installing the next column.

1. Change between Miscible Mobile Phases

- (1) Put approx. 100ml of the fresh mobile phase into a container of 200ml capacity.
- (2) Take out the suction filter from the first reservoir and swirl it around in the 100ml of fresh mobile phase to remove as much of the old mobile phase from the accessible surfaces of the filter and its tube.
- (3) Put the suction filter into a newly prepared reservoir filled with fresh mobile phase, open the drain valve, and press **purge** to purge the previous mobile phase completely out of the flow line.




Fig. 5.10

- (4) Break open the connection at the outlet of the manual injector and pipe the outlet to a waste container. Set the flow rate to 2-3ml/min then close the drain valve. Switch the injector between the load and inject positions and hold each position long enough to allow the sample loop and inner passages to be flushed. Thus, the mobile phase in the manual injector will be replaced.




Fig. 5.11

- (5) Reconnect the outlet of the manual injector, then replace the mobile phase in components in the flow line that follow the manual injector. Flow rate should be so adjusted that the delivery pressure does not exceed the permissible pressure rating of the column.




Fig. 5.12

2. Change between Immiscible Mobile Phases

- (1) Prepare an intermediate washing solution (such as isopropyl alcohol) that is miscible with both previous and next mobile phases.
- (2) Replace the old mobile phase with the intermediate washing solution. See procedures (1) to (5) of Section 5.9.1 "Change between Miscible Mobile Phases."
- (3) Replace the intermediate washing solution with the fresh mobile phase. See procedures (1) to (5) of Section 5.9.1 "Change between Miscible Mobile Phases."

3. Replacement of Buffered Mobile Phase

Use of a buffered mobile phase may result in crystalline substances being deposited upon dehydration or evaporation of the solution. Accordingly, replace the buffer solution as follows.

Note that if an organic solvent such as isopropyl alcohol is delivered, strongly bound salt deposits may form, so care should be taken.

- (1) Prepare distilled or deionized water.
- (2) Feed 200ml or more of the water to wash out the buffer solution. See procedures (1) to (5) of Section 5.9.1 "Change between Miscible Mobile Phases."
- (3) Then replace the water with the fresh mobile phase. See procedures (1) to (5) of Section 5.9.1 Change between Miscible Mobile Phases."

4. Use of Washing Flow Line

Use of a buffered mobile phase may result in the deposit of crystalline substances upon dehydration or evaporation of the solution: These salts can damage the plungers and plunger seals and shorten their service life. When using a buffer, it is advisable to frequently or continuously wash the back side of the plunger seals and plunger surfaces by pumping water through the built-in wash line. For the washing method, see Section 5.8 "Plunger Washing Method."

<<Caution>>

If you are not going to use the unit for some time (hours +) after feeding a buffer solution, replace the contents of the flow line with distilled or deionized water to prevent salt formation and deposition.

Chapter 6 Performance Checking

Contents

6.1 Operation Checking for Simple Isochratic System	6-2
6.2 Checking Concentrations in High Pressure Gradient Elution	6-5
6.3 Checking Concentrations in Low Pressure Gradient Elution	6-7

Operation checking for a simple isochratic system, as shown in the following illustration, is described here.




Fig. 6.1

1. Preparation

- (1) Prepare as many reservoirs as required.
- (2) Prepare an appropriate mobile phase and sample.

2. Connection

Connect the input terminal of the Chromatopac and AUX terminal of SPD-10A by using the signal cable and terminal block included in the accessories for the Chromatopac and the signal cable included in the accessories for the SPD-IOA.

3. Operation

- (1) Switch power on for the whole system.
- (2) Purge the flow line (from the suction filter to the pump drain).
- (3) Close the drain valve and operate the pump at 1 ml/min.
- (4) Check that the pressure is stable and liquid is flowing from the outlet of the SPD-IOA.
- (5) Set the temperature of the CTO-IOA to 40°C.
- (6) On the SPD-IOA, set the wavelength to 254nm, response to 4, and AUX RNG to 3. (The AUX output range becomes 2AU/V.)

- (7) On the Chromatopac, set ATTEN to 4.
- (8) On the Chromatopac, press **ZERO**, **5**, **0** and **ENTER** to set the pen position at the center.
- (9) On the Chromatopac, press **PL** and **ENTER** to carry out plotting.
- (10) Wait until the baseline becomes stable.




Fig. 6.2

- (11) On the Chromatopac, set ATTEN to 8.
(The Chromatopac's full scale becomes equivalent to 0.5AU.)
- (12) On the Chromatopac, press **ZERO**, **2**, **0**, and **ENTER**.
- (13) On the Chromatopac, press **PL** and **ENTER** to stop plotting.
- (14) Inject the sample. (see next page)
- (15) At the same time, press **START** of the Chromatopac.

4. Example of Analysis
(conditions)


The following is an example of the performance check under the conditions below:

Mobile phase : H₂O/Acetonitrile = 15/85
 Column : Shim-pack CLC-ODS 6.0mmø x 15cm
 Flow rate : 1ml/min
 : Wavelength 254nm, response 4
 Sample: : Naphthalene 10mg
 : Biphenyl 5mg
 : Acenaphthene 30mg
 : Anthracene 1mg
 : Pyrene 10mg

in 100ml
 Methanol
 Inject 5μl

(results)

START 1



CHROMATOPRC C-R6A
 SAMPLE NO 0
 REPORT NO 51 5

FILE 0
 METHOD 41

PKNO	TIME	RRER	MK	IDNO	CONC	NAME
1	4.857	214319			8,9471	Naphthalene
2	5.44	593455	SU		24,7171	Biphenyl
3	6.358	241558	U		10.8683	Acenaphthene
4	6.882	751233	U		31,2885	Anthracene
5	8.463	599924			24.9865	Pyrene
TOTAL		2488987			100	

The concentration of each solvent in high pressure gradient elution can be tested without a column or sample in the following method.

1. Conditions

Solvent A : 0.3% acetone in water
 B : Distilled water
Flow rate : 1ml/min
Load pressure : 10 – 30 [$\times 10^5$ Pa]
Measurement wavelength : 254nm

2. Specification


At liquid A concentration ratios of 10, 50 and 90%, the concentration error should be within $\pm 1\%$.

3. Precautions for Measurement

- (1) Be sure to degas the mobile phase.
- (2) Open the drain valve and wash the flow line for 5 minutes, with the flow rate set to 9ml/min and BCONC parameter to 50%.
- (3) Then set the flow rate to 3ml/min. Close the drain valve and wash the flow line for 20 minutes.
- (4) Last, set the flow rate to 2ml/min and the BCONC parameter to 100%. and check that the baseline has become stable, then start measurement.

4. Measurement Example

Initial parameter setting
Flow rate : 1ml/min
BCONC : 100



Time program		
TIME	FUNC	VALUE
0.10	BCONC	90.0
10.00	BCONC	90.0
10.10	BCONC	50.0
20.00	BCONC	50.0
20.10	BCONC	10.0
30.00	BCONC	10.0
30.10	BCONC	0.0
40.00	BCONC	0.0
40.10	BCONC	100.0

Calculation

$$\frac{A_{10}}{A_{100}} = 10\% \text{ actual concentration}$$

$$\frac{A_{50}}{A_{100}} = 50\% \text{ actual concentration}$$

$$\frac{A_{90}}{A_{100}} = 90\% \text{ actual concentration}$$


It is possible to read values of A10, A50, A90 and A100 directly from the recorder chart, however, this method yields serious errors. Accordingly, it is recommended to run the following BASIC program on the Chromatopac, to print out the LEVEL.

BASIC	PROGRAM
10	PRINT LEVEL;
20	WAIT 60
30	GOTO 10

1. **Conditions**

Solvent A : 0.3% acetone water solution
 Solvent B, C, and D : Distilled water
 Flow rate : 1ml/min
 Load pressure : 10 - 30 [$\times 10^5$ Pa]
 Measurement wave length : 254nm
2. **Specification**

At Solvent A mixing ratios 10, 50, and 90%, concentration error $\pm 2\%$ maximum
3. **Precautions for Measurement**
 - (1) Be sure to degas the mobile phase.
 - (2) Open the drain valve and wash the flow line for 5 minutes, setting the flow rate to 5ml/min and BCONC, CCONC, and DCONC parameter to 25%. In this case, be sure to connect the unused inlets with the reservoir and conduct replacement with either solution used.
 - (3) Then set the flow rate to 2ml/min. Close the drain valve and wash the flow line for 20 minutes.
 - (4) Last, set the flow rate 1ml/min and the BCONC parameter to 100%, CCONC and DCONC parameters to 0% and check that the base line has become stable, then start measurement.
4. **Measurement Example** Low pressure gradient Elution using Solvent A and B



Initial parameter setting

Flow rate: 1ml/min

B.CONC : 100.0

C.CONC : 0.0

D.CONC : 0.0

Time program

TIME	FUNC	VALUE
0.10	BCONC	90.0
10.00	BCONC	90.0
10.10	BCONC	50.0
20.00	BCONC	50.0
20.10	BCONC	10.0
30.00	BCONC	10.0
30.10	BCONC	0.0
40.00	BCONC	0.0
40.10	BCONC	100.0

Calculation

$$\frac{A_{10}}{A_{100}} = 10\% \text{ actual concentration}$$

$$\frac{A_{50}}{A_{100}} = 50\% \text{ actual concentration}$$

$$\frac{A_{90}}{A_{100}} = 90\% \text{ actual concentration}$$

It is possible to read values of A10, A50, A90 and A100 directly from the recorder chart, however, this method yields serious errors. Accordingly, **it is** recommended to run the following BASIC program on the Chromatopac, to print out the LEVEL.

BASIC	PROGRAM
10	PRINT LEVEL;
20	WAIT 60
30	GOTO 10

Chapter 7 Control from External Equipment

Contents

7.1	Connections of External Equipment to the Input/Output Terminals	7-2
7.2	Control from the SCL-IOA	7-4

1. The Input/Output Terminals




Fig. 7.1

Signal name	Description	Remarks
EVENT 1	Reed relay contact point output. The contact points close and open according to EVENT parameter values set in a program or AUX.FUNC setting.	The rated values of the reed relay contact points are 30VDC and 0.1A.
EVENT 2	Reed relay contact point output. The contact points close and open according to EVENT parameter values set in a program or AUX.FUNC setting.	
ON (PUMP)	The pump in the unit can be started by the external contact point signal through this terminal.	Control is made by short-circuiting these signal terminals and "COM" terminal. The closure time (tc) should be as follows: 0.5sec < tc < 10sec, i.e., 500ms or more
OFF (PUMP)	The pump in the unit can be stopped by the external contact point signal through this terminal.	
PRG START	The time program for the unit can be started by the external contact point signal through this terminal.	
COMMON	Common terminals for ON, OFF , START, and STOP signals.	
PRESSURE MONITOR	Outputs pressure values. The output scale may be set in ten steps within an output range of (0 - 50 [$\times 10^5$ Pa]) to (0 - 500 [$\times 10^5$ Pa]). See Section 5.5 "Auxiliary Functions (AUX. FUNC)"	Full scale voltage = 1mV

2. Wiring

- (1) Strip the insulation from the end of the connecting wire to leave approx. 10mm of exposed conductor.
This treatment is not necessary for the accessory remote cable.
- (2) If the core wire is solid, just insert it into the hole of a terminal.
If it is stranded, twist the end well, then insert it into the hole while pressing the button on the right side using a small screwdriver or the like. To disconnect the cable, press the button and pull the cable.




Fig. 7.2

<<Caution>>

The accessory remote cable can be used for the wiring.
If additional wire is required, carry out wiring with the following wires.

Solid wire: Q0.4- ϕ 1.0 (AWG26 - 18)

Stranded wire: 0.3mm^2 - 0.75mm^2 (AWG22 - 20)

Element wire diameter: ϕ 0.18 at the minimum

To protect the wire from breaking, stranded wire is recommended.
Twisted-pair cable composed of stranded wire with the exposed ends well tinned with solder is most highly recommended.

For the basic operation of the SCL-10A system controller, see the SCL-10A System Controller Instruction Manual.

1. **Switching Power on and the CONFIG Display** When the SCL power switch is turned on, power is supplied and a few seconds later the CONFIG display will appear.

If the solvent delivery unit and the SCL is connected correctly, **LINKED** is shown at the right of each unit name **PUMP.A**, **PUMP.B** and **PUMP.C**.

(In the case that three units are connected to the SCL.)

In this display, connection of units should be defined to control the solvent delivery unit from the SCL. If this definition is not carried out for a unit, link is not checked for the unit before an analysis even if it is connected to the SCL.

This definition is required only once after installation. The SCL stores the information in itself, so reassignment is not required unless you change the system configuration.

Define connection of the solvent delivery unit as shown below.

SVSTEM CONFIGURATION	
PUMP.A	LINKED
PUMP.B	LINKED LOCAL
PUMP.C	LINKED
DET.A	LINKED
DET.B	LINKED LOCAL
A.INJ	LINKED
FRC	LINKED
OVEN	LINKED LOCAL
SUBC	LINKED
press MENU	
PUMP.A	PUMP.B PUMP.C DET.A DET.B
time	a. fil 0 smpl / inj fr 0
9999.9	READY

unit name

Function line

- (1) Call up display "a" in the following illustration at the function line.
 - If "b" is displayed at the function line, press **(func)** once, then display of the line changes.

a	PUMP.A PUMP.B PUMP.C DET.A DET.B
b	A.INJ FRC OVEN SUBC KEYLOK

- (2) With the function keys, select the solvent delivery unit connected to the SCL. The selected unit will be highlighted on the screen.
The highlighted units on the screen are recognized by the SCL-10A as connected units.

<<Note>>

- If a unit is marked with **LOCAL**, control from the SCL is not accepted even if the name of the unit is highlighted. Set its local mode to 0, according to Section 5.5 (AUX.FUNC).

2. The Main Menu Screen

From the CONFIG screen, press **[menu]** to enter into the main menu screen.

It is possible to return to the main menu display whenever you press **[menu]** from operation in any other display.

```

                                MENU
ANALYSIS FILE PARAMETERS      : 0
                                TIME PROGRAM : 1
                                PUMP CONTROL  : 2
                                DET CONTROL   : 3
FRACTION COLLECTOR            : 4
AUTO INJECTOR/ANAL SEQUENCE : 5
MONITOR                       : 6
SYSTEM                        : 7

time a. file 0  smpl / inj 0  frc 0  KEYLOCK status
9999.9  READY

```

3. System Screen

Solvent selector valve unit FCV 10AL or FCV-11AL can be connected either to the pump unit or to the sub-controller. In any case the valve unit can be controlled by SCL-10A.

If the valve unit is connected to the pump unit or the sub-controller, state the connection in System Screen.

Select SYSTEM from the Main Menu to call the System Screen. In the System Screen the connection of optional units is stated. All information is stored in memory.

The CONFIG screen can be called from this screen.

```

SYSYTEM
RELAY1 : start stop error event1
RELAY2 : start stop error event2
RELAY3 : start stop error event3
A. INJ : online offline
FRC : online offline
P.LINE : on off
P.EXEC : on off
P.FRAC : on off
S. PROT : on off
AUTOZ. A : on off
AUTOZ. B : on off
EXT. S : yes no
PAC. CH : ch1 ch2 ch1&2
ALARM : enable disable
PRESS : kgf/cm2 psi MPa
PDERR : enable disable
SV : SUBC PUMP. A PUMP. B PUMP. C


CONFIG time a. file 0  smpl / inj 0  frc 0  KEYLOCK status
9999.9  READY

```

- (1) Select SV with cursor keys (Δ , ∇).
- (2) Select with cursor key ($\triangleright$, $\triangleleft$) the unit name to which the valve unit is connected.

4. Setting the Initial Parameters

- (1) Press $\boxed{0}$ or move the flashing position to 0 using the arrow keys in the main menu screen, then press $\boxed{\text{Enter}}$.
The analysis file parameter control screen will be displayed.
In addition to those for the solvent delivery units, parameters for any module may be displayed and updated from the analysis file parameter control screen.



- The $\boxed{\text{PUMP}}$ column in the figure shows parameters for the solvent delivery units.

The solvent delivery mode in use is indicated inside the parentheses.

If an FCV-IOAL or FCV-11AL is connected to the pump unit, the parameter SV for the sub-controller works as the parameter SV of the pump unit.

■ Operation

- (1) Press $\boxed{\text{func}}$ to display the following functions at the function line.

ISO (f1) B.GE (f2) T.GE L.P.GE KEYLOK

- (2) Select a solvent delivery mode for the solvent delivery units using the corresponding function key.

$\boxed{\text{ISO}}$ ($\boxed{\text{f1}}$) : Flow rates for up to three solvent delivery modules may be controlled independently.

$\boxed{\text{B.GE}}$ ($\boxed{\text{f2}}$) : Conducts high pressure binary gradient elution using two solvent delivery units (A, B). Flow rate of solvent delivery unit C is controllable independently.

T.GE (**f3**) : Conducts high pressure ternary gradient elution using three solvent delivery units.

LP.GE (**f4**) : Low pressure quarternary gradient elution is carried out using solvent delivery unit **A**. Flow rates of B and C are controllable independently.

<<Note>>

The number and contents of the parameters displayed on the screen will vary according to the solvent delivery mode.

- (3) Move the cursor with the arrow keys (, , , and ) to the parameter you want to change. Key in a new value and press **Enter** .
- If the **activate** lamp is lit, set parameters are sent to each unit at once for immediate execution.
 - If the function display field is switched to the help information by pressing **func**, the settable range and minimum step for the parameter at which the cursor is will be displayed.

Examples:

```
A. FLOW 0 - 9.99 (ml/min) STEP 0.01
EVENT 0. 1. 2. 3. 12. 13. 23. 123
```

- Parameters displayed on the screen are as follows (for the solvent delivery module only):

Parameter	Description
A.FLOW	Flow rate for pump A (ml/min)
B.FLOW	Flow rate for pump B (ml/min)
C.FLOW	Flow rate for pump C (ml/min)
T.FLOW	Total flow rate (ml/min) in gradient elution mode
A.PRES	Pressure of pump A ($\times 10^5$ Pa)
B.PRES	Pressure of pump B ($\times 10^5$ Pa)
C.PRES	Pressure of pump C ($\times 10^5$ Pa)
B.CONC	Concentration of solvent B (%)
C.CONC	Concentration of solvent C (%)
D.CONC	Concentration of solvent D (%)
B.CURV	Gradient curve of solvent B
C.CURV	Gradient curve of solvent C
D.CURV	Gradient curve of solvent D
P.MAX	Upper pressure limit ($\times 10^5$ Pa)
P.MIN	Lower pressure limit ($\times 10^5$ Pa)

5. Time Program

When the setting of the initial parameters in the parameter control screen is complete, press **[menu]** to return to the main menu screen and select **[TIME PROGRAM]**.

From the time program screen, a time program, as described in Section 5.4 "Time Program Creation and Execution," may be set through the SCL.

In addition, program steps set in this time program editing screen can include parameters for other modules such as detectors.

Program entry lines

TIME PROGRAM 0 (B. GE)					300 STEPS LEFT	
FILE	#	TIME	FUNC	VALUE		
0	1	30.00	B. CONC	100		
1	2	50.55	WAVE	A 350		
2	3	6000.0	STOP			
3	4					
4	5					
5	6					
6	7					
7	8					
8	9					
9	1					
FUNCTION FOR (PUMP DET. A DET. B OTHER)						
0	A. FLOW	5	B. CONC	8	B. CURV	11 A. PRES
1	B. FLOW	6	C. CONC	9	C. CURV	12 B. PRES
2	C. FLOW	7	D. CONC	10	D. CURV	13 C. PRES
4	T. FLOW					
TIME 0.01 - 9999.90 (min) STEP 0.01						
9999.9	time	a. file	0	amp	inj	fr c status
						READY

Instruction table

■ Operation

- (1) If a time program does not exist, the cursor is in the **[TIME]** field at step No. 0.

First at, using the numeric keys, **key** in the desired time for the parameter you want to change. After keying in values, press **[Enter]**, then the cursor moves to the **[FUNC]** field.

PRINT CLEAR COPY TABLE TABLE

- (2) Press **[func]** to select the function display that follows:
- (3) Use the **[TABLE←]** (**[f4]**) or **[TABLE→]** (**[f5]**) key to move the highlighted box and select the module name in the FUNCTION FOR column.

FUNCTION FOR (PUMP DET. A DET. B OTHER)

- (4) The contents of the parameter table vary with the module name which is displayed in the FUNCTION FOR column. The following is the parameter table for **[PUMP]**

0	A. FLOW	5	B. CONC	8	B. CURV	11	A. PRES
1	B. FLOW	6	C. CONC	9	C. CURV	12	B. PRES
2	C. FLOW	7	D. CONC	10	O. CURV	13	C. PRES
4	T. FLOW						

As for the parameters for another modules such as detectors, refer to the instruction manual of the SCL-IOA.

- (5) If you enter a number shown in the parameter table using the numeric keys, the **func** field will show the number. By pressing **Enter**, the number disappears and the corresponding parameter will appear. The parameter table is displayed for ease of reference for parameter numbers. If a parameter has a number and it is not displayed on the screen, the entry of the value is possible anyway.
- (6) If the selected parameter requires an argument (the entry of values), the cursor will move to the **VALUE** field. Key in an appropriate value and press **Enter**. If you select the help information by pressing **func**, the settable range and the minimum step of the argument required for the selected parameter will be displayed.
If the parameter requires no argument, i.e., is a command, the cursor will move to the **TIME** field of the next step.
- (7) Repeat the procedures (2) to (6) above and complete the time program.

<<Notes on Entry of program steps>>

The time you enter in the **TIME** field must be equal to or later than that in the previous step.

If you enter an earlier value, the alarm beeps when you press **Enter**, and the entry is not accepted. (This is different than the time program editor of the LC-10AT; See the SCL-10A manual for instructions on how to insert steps into a time program, etc.)

- (8) The following operation may be carried out from the time program screen using the function keys.
By pressing **func**, the function field switches and displays each function.

PRINT CLEAR COPY TABLE←TABLE→

PRINT (**f1**)

The time program on the screen is output to the printer connected to the Chromatopac.

CLEAR (**f2**)

All the steps of the time program corresponding to the file number on the screen will be deleted.

COPY (**f3**)

The time program corresponding to the file number on the screen will be copied to a file of another number.

TABLE← (**f4**)

Moves the box in the module name field of the parameter table to the left.

TABLE→ (**f5**)

Moves the box in the module name field of the parameter table to the right.

KEYLOK


KEYLOK (**f5**)

All key entry is disabled except **KEYLOK** (releases the key lock).

6. Pump Control Screen

The Gradient profile or flow rate profile in a time program can be displayed graphically on the pump control screen so that it is possible to visualize and check the time program during its creation. To display the pump control screen, select **PUMPCONTROL** from the main menu.

Plot of change in solvent delivery parameters



- (1) The horizontal axis of this graph corresponds to time of the time program. The full scale for the horizontal axis is determined automatically so that the STOP instruction is located adjacent to the right end of the axis (if a STOP instruction is not included, the time of the last step).
- (2) The vertical axis corresponds to pump flow rates, concentration, or pressures. For flow rates and pressures, the maximum flow rate and pressure of each solvent delivery unit are shown. In the case of concentration, the full scale value is 100%.

- (3) The solvent delivery mode field at the right shows the currently selected mode in highlight. Below it, in addition, the functional parameter names for the solvent delivery modules in that mode are shown.

■ Operation


- (1) Press **func** to switch the function display as follows:

ISO. B.GE T.GE LP.GE PARAM

- (2) Select the desired solvent delivery mode with the function key.
- The selected solvent delivery mode is highlighted in the solvent delivery field at the upper right on the screen. Below it, parameter names corresponding to the mode are shown.
- (3) By pressing function key **PARAM** (**f5**), the highlighted box moves down one line in the parameter name field. Move it to the desired parameter name for setting.
- Below the **TIME** field, the cursor flashes.
 - It is the same as **func** of the time program.
- (4) Enter the desired time with the numeric keys and press **Enter**.
- The entered time is displayed, and the cursor moves to the position below the **VALUE** field.
 - It is the same as **TIME** in the time program.
- (5) Enter a parameter value with the numeric keys and press **Enter**.
- The entered value is displayed and the cursor will move back to the **TIME** field.
 - It is the same as the parameter value in the **VALUE** field of the time program.
 - If you set the solvent delivery mode to the high pressure gradient elution and selected **CONC** in (3) for parameter name, proceed to (8).
- (6) By the procedures from (3) to (5) the entry of 1 step in time programming is completed.
- Changes in flow rate (pressure) of the solvent delivery unit or concentration, which is caused by the entry of a new step, are displayed graphically on the screen.
- (7) Complete the time program by repeating the entry procedures from step (3).
- (8) By entering a parameter value for **CONC** with the numeric keys, the cursor moves further down by one step and flashes.

Now enter the value for the shape of the gradient curve that starts from the current time.

- This is the value for the **CURV** parameter in time programming.
- If the value for the gradient curve is 0 (linear), just press **Enter**.
- For the value and shape of the gradient curve, refer to the instruction manual for the SCL-10A.
- When the entry is complete, the screen returns to (3).



<<Notes on Entry>>

- The entry of steps does not need to be carried out according to the order of time. Each new step is inserted in the time program in the correct order by an automatic sort function.
- If you change solvent delivery mode from the pump control screen, the solvent delivery mode on the parameter control screen will also be changed.
If a parameter for some other solvent delivery mode is entered while programming, it will be disregarded in the actual analysis.
- It is impossible to correct a wrong entry from this screen. In such a case, press **(menu)** to return to the main menu and correct the wrong entry with the editing function from the **TIME PROGRAM** screen.

■ Switching the Screen Mode

The graphics display has two screen modes:

OVERLAY: Information about **all** parameters listed at the right of the screen are displayed on the screen at once.

SINGL: The plot line of one parameter (selected with the **PARAM** key) is displayed.

Each time function key **G.MODE** (**f1**) is pressed, the screen mode switches.

The selected screen mode is highlighted at the right of the top line on the screen.

The following operation is available using the function keys from the pump control screen.

ISO. B. GE T. GE LP. GE PARM

ISO (f1)

How rates for up to three solvent delivery units may be controlled independently.

B.GE (f2)

Conducts high pressure binary gradient elution using two solvent delivery units (A, B). How rate for solvent delivery unit C is controlled independently.

T.GE (f3)

Conducts high pressure ternary gradient elution using three solvent delivery units.

LP.GE (f4)

Low pressure quaternary gradient elution is carried out using solvent delivery A. Flow rates for unit B and C is controlled independently.

G.MODE G.MODE G.MODE G.MODE KEYLOK

PARM (f5)

Selects parameter name.

G.MODE (f1)

Switches the screen display mode.

KEYLOK (f5)

All key entry is disabled except KEYLOK (releases the key lock).

Chapter 8 Maintenance

Contents

8.1	Replacement of Plunger Seals	8-2
8.2	Replacement of Plungers	8-5
8.3	Replacement of Washing Seal	8-7
8.4	Cleaning and Replacement of Check Valves	8-9
8.5	Replacement of Line Filter	8-10
8.6	Replacement of Fuse	8-11
8.7	Cleaning	8-12

It is recommended that the plunger seal be replaced approx. every 2000 hours of solvent delivery. (Water delivery at 1ml/min, 100×10^5 Pa.) See MON REV in Section 5.5 "Auxiliary Functions (AUX.FUNC)". When the plunger seal has been worn, the pressure becomes unstable, or solvent leaks from the gap between pump head and head holder or from the washing solution outlet. For replacement of the plunger seal follow the procedures below:

- (1) Assign 1 to the P.SET parameter. As for the setting method of the P.SET parameter, see Section 5.5 (AUX.FUNC).
- (2) The PUMP lamp will light up. Wait until this lamp goes off. Now the left and right plungers are set at the same position.
- (3) Remove the SUS pipe and teflon tube which are provided at the upper right and lower left areas of the pump head, respectively.




Fig. 8.1

<<Caution>>

When removing the suction filter, keep the top end of the suction filter tube at a position higher than the level of the solution in the reservoir, or empty the tube.

- (4) Loosen the three screws which fix the pump head gradually and alternately with the accessory Allen wrench. Then pull out the pump head.




Fig. 8.2

- (5) Two seals are used in the pump head. Remove the seal which is causing the liquid leakage by using the accessory seal remover. Insert the rimmed end of the seal remover into the seal to be removed, and pull it out.
- (6) Two white spacers made of Teflon are provided at the end of the sealing area in the pump head. Remove these spacers by using tweezers, etc.. Be careful not to damage (scratch) the pump head.
- (7) Wet the new seal and the two spacers with alcohol.
- (8) Mount the large spacer to the original position in the pump head.
- (9) Fit the new seal onto the non-rimmed end of the seal remover. Then, fit the small spacer to the spring of the seal. (Being wetted with alcohol, both the seal and the spacers stick together easily when they contact.)
- (10) Insert the parts assembled in step (9) above into the pump head. Be careful that the assembly is not cocked. Insert it squarely into the head.




Fig. 8.3




Fig. 8.4

- (11) Mount the pump head as before and tighten the three screws alternately and uniformly. Finally, fit the short end of the hex wrench to each screw, grab the long handle and securely tighten each screw. Pay attention to the up/down orientation of the pump head. The left side has the check valve, while right side does not.
- (12) Reconnect the tubing at the top and bottom of the pump head as before.

<<Caution>>

It may take some time after replacing the plunger seal until solvent delivery become stable. Accordingly, carry out running for about **an** hour using the solvent of usual with a flow rate setting of 1 to 5ml/min and pressure of 30 to 100×10^5 Pa before start of analysis.

<<Note>>

If only one seal leaks after both seals have been in service for some time, it is advisable to replace both since the undamaged seal will probably fail next, and possibly soon.

If leakage persists after replacement of the plunger seal or the service life of the new plunger seal is very short, the plunger surface may be damaged. Defects may be caused by foreign matter contained in the pump head or crystalline substances of the buffer solution. In that case, the plungers need to be replaced.

For replacement of the plungers, follow the procedures below:

- (1) Remove the pump head in the same manner as procedures (1) to (4) of Section 8.1 "Replacement of Plunger Seals."
- (2) Remove the seal holder assembly by pulling it out.




Fig. 8.5

- (3) Release and remove the two screws on the left and right sides of the head holder with the accessory Allen wrench. Then, release the two screws at the center alternately and gradually, and remove the head holder and the springs gently so that the head holder does not become cocked and jamb.
- (4) Remove the plunger by turning it counterclockwise with the accessory wrench (13 x 17) as shown in Fig. 8.6.
- (5) Mount the new plunger and **fix** it with the wrench.




Fig. 8.6

- (6) Thoroughly wipe the plunger surface with a piece of soft cloth or tissue paper.
- (7) Replace the head holder and spring to the original positions and tighten the two screws at the center alternately and uniformly until secure.
Then securely tighten the two screws on the left and right.
- (8) Clean the engaging section between the seal holder and head holder with a piece of soft cloth or tissue paper. Mount the seal holder to the original position.
- (9) Mount the pump head as before and tighten the three screws alternately and uniformly. Fit the short end of the hex wrench to each screw and securely tighten them.
- (10) Reconnect the tubing at the top and bottom of the pump head as before.

When the washing seal gets worn and washing solution leaks, the washing solution drips into the tray and is lead to the panel drain. If leakage of the solution is noticed, replace the washing seal.




Fig. 8.7

- (1) Remove the pump head in the same manner as procedures (1) to (4) of Section 8.1 “Replacement of Plunger Seals.”
- (2) Remove the seal holder by pulling it out.




Fig. 8.8

- (3) Remove the washing seal from the seal holder in the same manner as the replacement of the plunger seal, using the seal remover.
- (4) Slip a new seal over the non-rimmed end of the seal remover and push it vertically into the seal holder.




Fig. 8.9

- (5) Thoroughly wipe the engaging section between the seal holder and head holder with a soft cloth or a piece of tissue paper, then remount the seal holder gently without cocking it.
- (6) Mount the pump head and tighten the three screws gradually **and** alternately with an equal torque.
- (7) Reconnect tubing at the top and the bottom of the pump head **as** it was before.

A defective check valve may result in poor reproducibility of retention time or unstable pressure in solvent delivery. In such a case clean or replace check valves in the following procedure.

1. Cleaning the Check Valves

Replace the mobile phase with 2-propanol. Connect a resistance tube (0.1 x 2m) in place of the column, then feed the solution at a flow rate of 2ml/min for one hour or more.

2. Replacement of check valve

- (1) Remove the SUS pipe and teflon tube at the top and the bottom of the pump head.




Fig. 8.10

- (2) Remove the check valves with the accessory wrench.
- (3) Mount the new check valves.
- (4) Reconnect the tubing at the top and the bottom of the pump head as before.

A

<<Cautions>>

- Do not disassemble the check valves. If it is done, we cannot guarantee the functions of the check valve.
- Be careful not to set the IN check valve upside down, otherwise the parts inside the valve may fall out.

If a back pressure exists with the tubing at the pump outlet disconnected, a possible cause is a clogged line filter.

To replace the line filter follow the procedures below:

- (1) Disconnect the tubing of the pump outlet.
- (2) Unscrew the line filter with the accessory wrench.
- (3) Take out the gasket that is in the drain valve main body.
- (4) Wipe off dirt in the connecting port for the line filter with a soft cloth or a piece of tissue paper.
- (5) Mount a new line filter and a gasket.
- (6) Connect the outlet tubing as before.




Fig. 8.11

When the fuse is broken, replace it as follows.

Be sure to use the same type and the same rating of fuses.

Type and rating of fuses;

250V 2AT (5 × 20mm) (for 90-130V~)

250V 1AT (5 × 20mm) (for 200-250V~)

- (1) Turn the power switch off.
- (2) Remove the power cord from the power cord connector.
- (3) With a slotted-head screwdriver, catch the cover on the fuse holder and slide it out.




Fig. 8.12

- (4) After replacing the fuse, push the fuse holder cover back in until it clicks.




Fig. 8.13

Cleaning

If the cover or the front panel of the unit is dirtied, wipe away the dirt or dust with a soft cloth or a piece of tissue paper. If necessary, use a synthetic detergent.

Chapter 9 Troubleshooting

Contents

9.1	Symptoms and Countermeasures for Troubles	9-2
9.2	Error Messages	9-3

Pump operation troubles result from various causes all of which cannot be covered in this manual. This section describes the troubleshooting measures on the basis of the symptoms. As for the countermeasures, detailed information is given in each relevant section.

Symptom	Cause	Remedy
Solvent delivery is unstable. Fluctuation in pressure is large.	1. Bubbles in the pump chamber.	1. Flush liquid through the drain to purge bubbles. 2. Suck up bubbles through the drain tubing with a syringe.
	2. Old solvent remains in the pump chamber.	1. Increase flow rate and flush solvent through the drain to purge the old mobile phase completely.
	3. Bubbles in the suction filter enter the pump.	1. Increase the flow rate to purge the old solvent through the drain tubing. 2. Shake the suction filter to release bubbles. 3. If the suction filter is clogging, carry out ultrasonic cleaning. If the ultrasonic cleaning is not effective, replace the suction filter. 4. Degas the mobile phase.
	4. Malfunction of the check valves.	1. Clean or replace the check valves.
	5. Leak from the plunger seal.	1. Replace the plunger seal. 2. Replace the plunger.
	6. Leak in the flow line.	1. Tighten the joints in the flow line. 2. Replacement of defective parts.
	7. Flow line clogged	1. Clean or replace the line filter. 2. Check the flow line.
The pump operates but solvent is not delivered.	1. Bubbles in the pump chamber.	1. Flush liquid through the drain to purge bubbles. 2. Suck up bubbles through the drain tubing with a syringe.
	2. Air enters from the joint between the suction filter and inlet pipe.	1. Firmly connect the joint.
Measured flow rate is lower than the set flow rate.	1. Malfunction of the check valves.	1. Clean or replace the check valves.
	2. Clogged suction filter.	1. Clean or replace the suction filter.
Solvent is not delivered. (Pump is not operating.)	1. An error message is displayed. (P.MAX, P.MIN, etc.)	1. Press [CE]. Remove the cause of the error.
Pressure does not rise.	1. The drain valve is open.	1. Close the drain valve.
	2. Leakage in flow line.	1. Tighten joints. 2. Replacement of defective parts.
Pressure rises excessively. (Confirm it by removing the column.)	1. Line filter clogged	1. Clean line filter. 2. Replace line filter.
	2. Flow line clogged	1. Identify the clogged component and replace it.
	3. The inside diameter of tubing is excessively small.	1. Replace the tubing with another having an appropriate inside diameter.

This unit has multiple self-diagnostic testing functions, which beep an alarm and show an error message if they detect any abnormalities. Except for the case (1) and (2), pressing **CE** will stop the alarm beep and return the display to the initial screen.

1. ROM Error

ROM FAILURE

This error message will be displayed when there is something irregular in the **ROM**.

If this error message is displayed, pressing **CE** will not release the system from the error.

Turn off the power source for the unit, and contact our office or agent.

2. RAM Error

RAM FAILURE

This error message will be displayed when there is something irregular in the **RAM**. If this error message is displayed, pressing **CE** will not release the system from the error.

Turn off the power source for the module, and contact our office or agent.

3. Lost RAM File Error

NOT PROTECTED

This error message will be displayed when the contents of the **RAM** have been lost or destroyed due to power source abnormalities. By pressing **CE**, the module may become operable. However, if this error message is displayed every time the power switch is turned on, the problem may be caused by dead or low power backup batteries. Contact our office or agent for information on exchanging the backup batteries.

4. Overstepping Error

OUT OF MAX STEPS

This error message is displayed when the total number of steps used in the time programs exceeds 320.

If this error message is displayed, the time program can not be prepared any further. Continue entering program data after deleting unnecessary files.

5. Maximum Pressure Error

ERROR P — MAX

This error message will be displayed when the upper limit pressure, p-max, has been exceeded by the actual delivery pressure. In this case solvent delivery stops. (However, if S.PROT function is activated, the unit only reduces the flow rate by one half without stopping the pump.) When this error message is displayed, check the flow line for clogging, and if there is no problem, set the **P-MAX** value to an appropriate value.

6. Minimum Pressure Error

ERROR P — MIN

This error message will be displayed if the actual delivery pressure is less than lower limit pressure P.MIN after 1 min. of operation. In this case the limiter stop the solvent delivery. If this error occurs, check the flow line for leakage, and if there is no problem, set the p.min value to an appropriate value. This error trapping will not function in the first one minute after the start of the solvent delivery in order to allow for pressure equilibration when just beginning to pump.

7. Home Position Error

ERROR HOME POS

This error message is displayed when the home position of the motor can not be detected. It is also displayed when the motor slips.

In this case, press the **ICE** key to cancel the error, then press **pump** to operate. If the same error appears again, turn off the power for the unit, then turn it on again after 10 or 15 seconds and try it again. If this doesn't work, contact our office or agent.

8. Purge Error

O P E N D R A I N V A L V E

This error message is displayed if the actual delivery pressure of the pump exceeds 5kgf/cm^2 during purge operations started with the **(purge)** key. The **CE** key will release the system from the error. Start the purge operation again after opening the drain valve. If the error message is displayed even with the drain valve open, execute the zero adjustment of the pressure sensor (see **ZERO ADJ** in 5.5 "AUX.FUNC").

Chapter 10 Specifications

Contents

10.1	LC-10AT Unit Specifications	10-2
10.2	High Pressure Gradient Elution Specifications	10-3
10.3	Low Pressure Gradient Elution Specifications	10-4

1. Pump type A 2-stage (tandem-flow) double-plunger reciprocating pump.
2. Modes of solvent delivery Constant flow or constant pressure solvent delivery
 - 2.1 Constant flow delivery
 - Flow rate range: 0.001 - 5ml/min (10 - 400 [$\times 10^5$ Pa])
5.001 - 9.99ml/min (10 - 200 [$\times 10^5$ Pa])
 - Flow rate accuracy: Within +2% or $\pm 2\mu$ /min of set value, whichever is higher (0.1 - 5ml/min, 10 - 400 [$\times 10^5$ Pa], with water and room temperature at 20°C)
 - Flow rate precision: Within $\pm 0.3\%$ (0.1 - 5ml/min, 10 - 400 [$\times 10^5$ Pa] with water and room temperature at 20°C)
 - 2.2 Constant pressure delivery
 - Pressure setting range: 10 - 400 [$\times 10^5$ Pa] (in 1 [$\times 10^5$ Pa] steps)
 - Pressure setting accuracy: Within $\pm 10\%$ or ± 10 [$\times 10^5$ Pa] of set value, whichever is larger
3. Pressure limits: User-settable High/Low pressure limits
4. Materials in contact with solvent: SUS316, UHMW-PE, Elgiloy[®], Kalrez[®], PEEK
ETFE, PCTFE, PTFE
5. Suction filter: 10 μ m mesh
6. Line filter: 2 μ m mesh, liquid volume approx. 200 μ l
7. Time program: Flow rate, pressure, EVENT, LOOP (program repeat), 10 files, total 320 steps
8. Pressure indication accuracy: Within +2% or ± 5 [$\times 10^5$ Pa], whichever is larger
9. Plunger cleaning flow line: Built-in Cleaning of plunger seals is possible with the connection of syringe or pump.
10. Size: 260W x 140H x 420D (except protruding components)
11. Weight: 11kg
12. Operating temperature range: 5 - 35°C
13. Power requirements: (Power requirements vary depending upon the line voltage.)
P/N 228-31900-91, 92: 90-130V~ 50/60Hz 100VA
P/N 228-31900-93. : 200-250V~ 50/60Hz 100VA

1. LC-10AT Control - High Pressure Gradient Elution Specifications

- 1) Number of solvents: 2
- 2) Gradient profile: Stepwise and linear, multiple processes are possible in a single program.
- 3) Number of program files: 10 files with a total of 320 program steps
- 4) Program mable time: 0.01 min - 999.99 min, 0.01 min steps
- 5) Settable Concentration range: 0 - 100%, 0.1% steps
- 6) Concentration setting accuracy: $\pm 1\%$ (0 - 100%, 0.3 - 3ml/min, 10 - 400 [$\times 10^5$ Pa], with aqueous acetone/water)

2. SCL-10A Control - High Pressure Gradient Elution Specifications

- 1) Number of solvents: 2 or 3
- 2) Gradient profiles: Stepwise, linear, and exponential functions, multiple processes are possible in a single program.
- 3) Number of program files: 20 files, with a total of 400 program steps
- 4) Programmable time: 0.01 min - 9999.9 min, 0.01 min steps
- 5) Settable Concentration range: 0 - 100%, 0.1% steps
- 6) Concentration setting accuracy: $\pm 1\%$ (0 - 100%, 0.3 - 3ml/min, 10 - 400 [$\times 10^5$ Pa], with aqueous acetone/water)

1. LC-10AT Control - Low Pressure Gradient Elution Specifications

- 1) Number of solvents
mixed: up to 4
- 2) Gradient profiles: Stepwise and Linear, multiple processes are possible in a single program.
- 3) Number of program
files: a total of 10 files and 320 steps
- 4) Programmable time: 0.01 - 999.99 min, 0.01 min steps
- 5) Settable Mixing ratio
range: 0 - 100%, 0.1% steps
- 6) Concentration accuracy $\pm 2\%$
- 7) Flow rate range 0.1 - 2ml/min
(In the case of a flow rate of more than 2 ml/min, the life time of the valves of the low pressure gradient unit become short due to high speed cycling.)

2. SCL-10A Control - Low Pressure Gradient Elution Specifications

- 1) Number of solvents
mixed: up to 4
- 2) Gradient profiles: Stepwise, linear and exponential functions; multiple processes are possible in a single program.
- 3) Number of program
files: a total of 20 files and 400 steps
- 4) Programmable time: 0.01 - 9999.9 min
- 5) Mixing ratio setting
range: 0 - 100%, 0.1% steps
- 6) Concentration accuracy $\pm 2\%$
- 7) Flow rate range 0.1 - 2ml/min
(In the case of a flow rate of more than 2 ml/min, the life time of the valves of the low pressure gradient unit become short due to high speed cycling.)

Chapter 11 Spare Parts and

Contents

11.1 Consumable Parts List.....	11-2
11.2 List of Other Necessary Parts	11-3
11.3 Optional Units List	11-4

11. 1

Consumable Parts List

	Part name	Part number	Remarks
I	Plunger seal (yellow)	228-21975	Seal for pump head

1. Mechanical Parts

	Part name	Part number	Remarks
1	Plunger assy	228-32654-91	Seal for washing solution
2	Head holder	228-31903	
3	Seal PE (white)	228-28499	
4	Pump head	228-31905	
5	Check valve OUT	228-32531-92	Including line filter and drain valve
6	Check valve IN	228-32166-91	
7	Pressure sensor assy	228-32252-92	Low volume type → 228-32693-91
8	Filter assy standard for Low Pressure Gradient	228-12642-93 (228-32467-91)	
9	Gasket	228-12564	For line filter replacement
10	Suction filter assy	228-18740-91	For detection of the home position of pump
11	Photosensor assy	228-28288-91	
12	Drain valve shaft 10AS ASSY	228-28312-92	For plunger seal
13	Spacer (large)	228-32325-01	
14	Spacer (small)	228-32401	

	Part name	Part number	Remarks
1	PB-I assy ROM included	228-24785-95	ROM included
2	PB-I assy ROM not included	228-23685-93	ROM not included (Photocoupler included)
3	ROM LC-IOAT	228-24783-95	
4	Fuse IAT 5 × 20	072-01652-16	For 200V
5	Fuse 2AT 5 × 20	072-01652-19	For 100V

The following is a list of the optional units that may be used in combination with this unit. For details on the available optional units or for information on those not listed here, contact our office or agency.

Optional unit identification	Part number	Description
DGU-2A	228-21996-91	A helium degasser for four solvents. Up to four solvents may be degassed.
DGU-3A	228-24200-91	A degasser for three solvent. It permits continuous degassing by pressure reduction through a resin membrane for up to 3 flow lines.
DGU-4A	228-24201-91	A degasser for two solvent. It permits continuous degassing by pressure reduction through a resin membrane for up to 2 flow lines.
DGU-IOB	228-32235-91	A full automatic helium degasser for four solvent. Controlled by the LC-IOAT (or the SCL-IOA through option box L).
FCV-10AL	228-24861-91	A solenoid valve unit for low pressure gradient elution analysis. - Also usable as a four-solvent selector solenoid valve unit. - Controlled by the LC-IOAT (or by the SCL-IOA through option box L).
FCV-11AL	228-24813-91	A solenoid valve unit for selecting solvents (for 3 flow lines). - Three units of three way solenoid valves are built in, permitting selection of solvents for three flow line. - Controlled by the LC-IOAT (or by the SCL-IOA through option box L).
FCV-11AL(S)	228-24813-92	A solenoid valve unit for selecting solvents (for 1 flow line). - A three way solenoid valve is built in, permitting selection of two solvents. - Controlled by the LC-IOAT (or by the SCL-IOA through option box L).
FCV-13AL	228-24914-91	A valve unit for selecting solvents. - Having 6 setting positions, this valve unit permits switching 6 solvents for one flow line. - Controlled by the SCL-IOA through option box L.
A set of 3 suction filters	228-18907-91	A set of three suction filters. This set is for the FCV-13AL module. (The suction filters are not included with the FCV-13AL.)
Option Box L	228-25025-91(100V) 228-25025-92(200V)	This unit has 3 functions. - For control of the FCV-11AL, FCV-11AL(S), and FCV-13AL. However, setting is made through the SCL-IOA. - Up to 6 of the following units can be contained: FCV-11AL, FCV-11AL(S), FCV-13AL, DGU-1A, DGU-2A, and DGU-IOB. - Supplies AC power for other modules included in a system. 100V x 8, or 200V x 6.
Option Box S	228-25306-91	Up to 2 of the following units can be contained: FCV-11AL, FCV-11AL(S), FCV-13AL, DGU-1A, DGU-2A, DGU-IOB, etc. This unit can be placed under the LC-IOAT so it requires very little space.
Reservoir Box	228-25038-91	This box provides housing space for the reservoirs, and permits installation of a manual injector.
Mixer	228-28000-91	A static mixer for the high and low pressure gradient elution analysis.

Manual Injector Model 7725	228-32210-91	A manual injector for general analyses. The standard sample loop of 20 μ l is included.
Manual Injector Model 7725 i	228-32210-93	Manual Injector 7725 with the position sensing switch. A signal synchronized with sample injection is sent to the system controller and Chromatopac.
Manual Injector Model 8125 for semi-micro liquid chromatography	228-23200-91	A manual injector for semi-micro liquid chromatography. The standard sample loop of 5 μ l is included. The unit has a built-in position sensing switch.
Tubing Parts Kit for semimicro liquid chromatography	228-23198-91	This unit consists of a SUS pipe with an inside diameter of 0.2mm and a joint for the manual injector.
Injector Holder	228-25468-91	This unit is used for mounting the manual injector on the right side of the LC-IOAT.
Column Holder	228-15418-93	This unit is used when fixing a column on the right side of the LC-IOAT.
Air Trap	228-23675-91	This is attached to the tube for the suction filter and used to let the air out of the flow line.

Chapter 12 Reference Material

Contents

12.1 Precautions on Static Electricity	12-2
12.2 Mobile Phase Characteristics	12-6

Chapter 12 Reference Material

Contents

12.1 Precautions on Static Electricity	12-2
12.2 Mobile Phase Characteristics	12-6

Liquid chromatography using flammable organic solvents as mobile phase requires proper care against fire, explosion, etc. Particularly, among various possible accidents, those caused by static electricity are difficult to anticipate, and tend to occur only with unexpected conditions which often make countermeasures insufficient.

At a site where preparative liquid chromatography is practiced, a large amount of flammable substances may be used. Therefore, once an accident happens, it could lead to tremendous damage.

The mechanism of accident caused by static electrical discharge and preventive measures are described below. Take due care in safety measures in handling of equipment.

1. Mechanism of Static Electrical Discharge Accident (Example)

Occurrence of Static Electricity

Accidents caused by static electricity take place through the following processes.

When liquid is fed at high speed through a small-diameter tube like the pipe of a liquid chromatograph, static electrical charge occurs by friction between solid and liquid as shown in Fig. 12.1.




Fig. 12.1 Occurrence of Static Electricity by Friction between Solid and Liquid

Charging and storage

When the charged liquid is collected in an insulated vessel, the

Energy release by discharge

If some other conductive object is brought near the vessel, electricity is discharged at a certain distance from the vessel releasing heat energy.

Ignition of combustible substances

If flammable gas of sufficient concentration exists nearby, ignition is caused by this energy.




Fig. 12.2 Conditions for Accidents

2. Preventive Measures against Accidents

The principal preventive measure is the prevention of “charging and storage of static electricity” among those items shown in “Occurrence mechanism of static electrical discharge accidents.” The preventive measures are shown below. It is recommended to exercise two or more measures simultaneously.

- * Particularly when a large quantity of flammable solvent is held in a large vessel, be sure to observe the preventive measures 1, 2, and 3.

Preventive measure 1.

Use metallic (conductive) waste liquid vessel which is well grounded. This releases the charge of the waste liquid and vessel to ground.

The following items are available.

- (1) Grounding wire with clip P/N 228-21353-91
- (2) Metallic 18 liter can P/N 038-00044
- (3) Metallic 4 liter can P/N 038-00043-01

- * Be sure to ground the vessel properly. Disconnecting of grounding wire or poor grounding defeat the purpose of using a metallic vessel.
- * There are some metallic cans which have no conductivity due to an oxidized coating or lacquer on their surface. Be sure to confirm the grounding of vessels by a tester before application.
- * When a liquid with almost no conductivity (of 10^{-10} S/m or less) is discharged into the vessel, it is necessary to mix it

with another liquid with some conductivity. (The other liquid can be placed in the vessel in advance.)

Preventive measure 2. Minimize the clearance of both inlet and outlet of vessel to prevent flame from entering the vessel.

- (1) Cap with three holes for 18 liter and 4 liter cans (P/N 228-21354-91) is available.




Fig. 12.3 Anti-Static Electricity Measures for Vessel

Preventive measure 3. Do not approach the vessel with charged objects including the human body.

Charging prevention measures for human body

- Prevention of charging of shoes and clothes
- Grounding of human body
- Make working floor conductive

Suitable products to be used for those measures a), b), and c) are available on the market.

- ✱ When persons who use no charge prevention measures approach dangerous sections, they have to be grounded beforehand. (For example, they should contact grounded metal by hand.)

Preventive measure 4. Use pipes with inner diameter of 2mm or more for waste liquid line for large flow rates.

- ✱ Inclusion of bubbles in the tube may increase the amount of charging by ten times. Check that there is no inclusion of air via tube joints.

Preventive measure 5. When it is impossible to use a conductive vessel, use caution in the following points.

- a) Set the vessel **so** that the pipe outlet will be placed below the liquid level in the vessel. Or, dip **a** grounded metal (ex. pipe connected to the main body of device) in the liquid.
- ✱ This method is not effective for liquid with small conductivity (10^{-10} s/m or less).
- b) Use a vessel of the smallest possible capacity to minimize the damage by fire if it should occur.
- c) Prevent the room from being dry. Humidity of 65% or more has charge prevention effects.

(1) Solvent	*7<.5cP,>45° **7<.5cP,<45°	(2) Source	(3) UV Cutoff	(4) R.I. 25°	Boiling Point (°C)	Viscosity (cP,25°C)	(5) p'	(6) e _a [°]	(7) Water Solubility %* in 20°C Solvent	(8) Dielectric Constant ²⁰	(9) P'+ 0.25 e
1.	FC-78 (*) FC-75 (Fluorescent solvent) FC-43	(Particular to LC)	210 nm 210 (opaque or under)	1.267 1.276 1.291	50 102 174	0.4 0.8 2.6	< -2 < -2 < -2	-.25 -.25 -.25		1.88 1.86 1.9	p' and Dielect. const. (function in proportion to intensity)
2.	Isooctane (*) (2,2,4- tri methylpentane)	LC	197	1.389	99	0.47	0.1	0.01	0.011	1.94	0.1
3.	n-Heptane (*)	LC	195	1.385	98	0.40	0.2	0.01	0.010	1.92	0.5
4.	n-Hexane (*)	LC	190	1.372	69	0.30	0.1	0.01	0.010	1.88	0.5
5.	n-Pentane (**)	LC	195	1.355	36	0.22	0.0	0.00	0.010	1.84	0.5
6.	Cyclohexane	LC	200	1.423	81	0.90	-0.2	0.04	0.012	2.02	0.5
7.	Cyclopentane (*)	LC	200	1.404	49	0.42	-0.2	0.05	0.014	1.97	0.6
8.	1-Chlorobutane (*)	LC	220	1.400	78	0.42	1.0	0.26		7.4	2.8
9.	Carbon disulfide	LC	380	1.624	46	0.34	0.3	0.15	0.005	2.64	1.7
10.	2-Chloropropane (**)	LC	230	1.375	36	0.30	1.2	0.29		9.82	3.7
11.	Carbon tetrachloride	LC	265	1.457	77	0.90	1.6	0.18	0.008	2.24	2.3
12.	n-Butyl ether		220	1.397	142	0.64	2.1	0.25	0.19	2.8	2.4
13.	Trichthylamine			1.398	89	0.36	1.9	0.54		2.4	2.4
14.	Bromoethane (*)			1.421	38	0.38	2.0	0.35		9.4	4.3
15.	i-Propyl ether (*)		220	1.365	68	0.38	2.4	0.28	0.62	3.9	3.2
16.	Toluene	LC	285	1.494	110	0.55	2.4	0.29	0.046	2.4	2.9
17.	p-Xylene		290	1.493	138	0.60	2.5	0.26		2.3	3.0
18.	Chlorobenzene			1.521	132	0.75	2.7	0.30		5.6	4.1
19.	Bromobenzene			1.557	156	1.04	2.7	0.32		5.4	4.1
20.	Iodobenzene						2.8	0.35			
21.	Phenyl ether			1.580	258	3.3	3.4			3.7	3.7
22.	Phenetole			1.505	170	1.14	3.3			4.2	4.9
23.	Ethyl ether (**)	LC	218	1.350	35	0.24	2.8	0.38	1.3	4.3	4.0
24.	Benzene	LC	280	1.498	80	0.60	2.7	0.32	0.058	2.3	3.6
25.	Tricresyl phosphate										
26.	Ethyl iodide			1.510	72	0.57	2.2			7.8	4.2
27.	n-Octanol		205	1.427	195	7.3	3.4	0.5	3.9	10.3	5.8
28.	Fluorobenzene			1.46	85	0.55	3.1			5.4	4.6
29.	Benzylether			1.538	288	4.5	4.1				
30.	Methylene chloride (**)	LC	233	1.421	40	0.41	3.1	0.42	0.17	8.9	5.6
31.	Anisole			1.514	154	0.9	3.8			4.3	4.6
32.	i-Pentanol			1.405	130	3.5	3.7	0.61	9.2	14.7	7.3
33.	1,2-Dichloroethane	LC	228	1.442	83	0.78	3.5	0.44	0.16	10.4	6.3
34.	t-Butanol			1.385	82	3.6	4.1	0.7	miscible	12.5	
35.	n-Butanol	LC	210	1.397	118	2.6	3.9	0.7	20.1	17.5	8.3
36.	n-Propanol	LC	240	1.385	97	1.9	4.0	0.82	miscible	20.3	
37.	Tetrahydrofuran (*)	LC	212	1.405	66	0.46	4.0	0.57	miscible	7.6	
38.	Propylamine (*)			1.385	48	0.35	4.2		miscible	5.3	
39.	Ethylacetate (*)	LC	256	1.370	77	0.43	4.4	0.58	8.8	6.0	5.8
40.	i-Propanol	LC	205	1.384	82	1.9	3.9	0.82	miscible	20.3	
41.	Chloroform (*)	LC	245	1.443	61	0.53	4.1	0.40	0.072	4.8	5.6
42.	Acetophenone			1.532	202	1.64	4.8			17.4	8.7
43.	Methylethyl ketone (*)	LC	329	1.376	80	0.38	4.7	0.51	23.4	18.3	9.1
44.	Cyclohexanone			1.450	156	2.0	4.7			18.3	9.1

(1) Solvent	*7<.5cP,>45° **7<.5cP,<45°	(2) Source	(3) UV Cutoff	(4) R.I. _{25°}	Boiling Point (°C)	Viscosity (cP,25°C)	(5) p'	(6) e _a	(7) Water Solubility %* in 20°C Solvent	(8) Dielectric Constant e ²⁰	(9) P'+ 0.25 e
45. Nitrobenzene				1.550	211	1.8	4.4			34.8	13.2
46. Benzonitrile				1.536	191	1.2	4.8			25.2	10.9
47. Dioxane		LC	215	1.420	101	1.2	4.8		miscible	2.2	
48. Tetramethyl urea		LC	265	1.449	175		6.0	0.56		23.0	10.7
49. Quinoline				1.625	237	3.4	5.0			9.0	7.4
50. Pyridine				1.507	115	0.88	5.3		miscible	12.4	
51. Nitroethane			380	1.390	114	0.64	5.2		0.9		
52. Acetone (*)		LC	330	1.356	56	0.30	5.1	0.71	miscible		
53. Benzyl alcohol				1.538	205	5.5	5.7			13.1	8.8
54. Tetramethyl guanidine							6.1	0.6			
55. Methoxyethanol		LC	210	1.400	125	1.60	5.5		miscible	19.9	
56. Tris (cyanoethoxy) propane		GC					6.6	0.56			
57. Propylene carbonate		LC					6.1				
58. Ethanol		LC	210	1.359	78	1.08	4.3		miscible	24.6	
59. Oxydipropionitrile		GC					6.8				
60. Aniline				1.584	184	3.77	6.3			6.9	8.1
61. Acetic acid				1.370	118	1.1	6.0		miscible	6.2	
62. Acetonitrile (*)		LC	190	1.341	82	0.34	5.8		miscible	37.5	
63. N, N-dimethylaceta- mide		LC	268	1.436	166	0.78	6.5	0.88		37.8	
64. Dimethylformamide		LC	268	1.428	153	0.80	6.4			36.7	
65. Dimethylsulfoxide		LC	268	1.477	189	2.00	7.2	0.62	miscible	4.7	
66. N-methyl-2-pyrrolidone		LC	285	1.468	202	1.67	6.7			32	
67. Hexamethyl phosphoric acid triamide				1.457	233	3	7.4	Q.65		30	
68. Methanol (*)		LC	205	1.326	65	0.54	5.1		miscible	32.7	
69. Nitromethane			380	1.380	101	0.61	6.0		2.1		
70. m-Cresol				1.540	202	14	7.4			11.8	10.0
71. N-methylformamide				1.447	182	1.65	6.0		miscible	182	
72. Ethylene glycol				1.431	182	16.5	6.9		miscible	37.7	
73. Formamide				1.447	210	3.3	9.6		miscible	111	
74. Water		LC		1.333	100	0.89	10.2			80	

- (1) (*) indicates solvents most suitable for LC, having convenient boiling points ($>45^{\circ}\text{C}$) and low viscosity ($\leq 0.5\text{cp}$).
(**) indicates solvents with very low viscosity and boiling point.
- (2) "LC" indicates that the solvents are commercially available specifically for LC from the following companies: Burdick & Jackson, Baker Chemical, Mallinkrodt Chemical, Fischer Scientific, Waters Associates and Manufacturing Chemists, Inc.
(Note: In Japan, they are also commercially available from the following companies: Wako Pharmaceutical Co., Ltd., Nakarai Pharmaceutical Co., Ltd. and Kanto Chemical Co., Ltd.)
"GC" indicates that the solvent is used as a stationary phase for gas chromatography, and can be purchased from companies selling GC columns and stationary phases. (These solvents are used as a stationary phase in the liquid-to-liquid LC.)
- (3) The wavelength shorter than which the solvent becomes opaque.
- (4) Refractive index at 25°C
- (5) Polarity parameter of solvent
- (6) Solvent strength parameter of liquid-to-solid adsorption on alumina
- (7) Water solubility (W%) at 20°C for solvent used in liquid-to-solid adsorption
- (8) Value at 20°C
- (9) Function where P' is proportional to solvent strength and dielectric constant in ion pair chromatography.