# **Barnstead Mega-Pure Glass Stills**

Barnstead Mega-Pure Glass Stills effectively remove inorganic solids, organics with boiling points higher than water, bacteria, and pyrogens. The stills are constructed of non-leaching components to ensure high purity water.

#### **Product features**

- Contact with only glass, quartz and Teflon<sup>™</sup> components ensures ultimate purity and eliminates cross-contamination
- · Choice of three production capacities
- Units are easily wall-mounted or bench-mounted
- · Vertical condenser design provides maximum purity
- High temperature cutoff shuts down the still if the temperature is too high, preventing heating element burnout
- Quick release cover provides easy cleaning access
- Connection to Automatic Collection System (ACS) allows 24-hour automatic system operation
- Push-button system start with feed water solenoid valve

#### Barnstead MegaPure-1 Glass Stills

 Compact, easy-to-use 1 L/hr glass still available with or without a storage bottle

#### Barnstead MegaPure-3A and 6A Glass Stills

• Ideal for larger volume demands with storage and options



#### Quick look comparison

	MegaPure-1	MegaPure-3A	MegaPure-6A
System features	Compact, easy-to-use	Stills for larger videmands	volume
Product water capacity, L/hr	1	3	6
Cooling water capacity, L/hr	11.3	30.3	53
Number of heaters	1	1	2
Resistivity at 25°C, MΩ.cm	1–10	1–10	1–10

#### **Applications include**

#### Routine lab work

- Rinsing lab glassware
- Supplying autoclaves and ultrapure water systems
- · Preparing and diluting buffers, reagents, and media

# Glass stills ideal for biological applications

### Product accessories

#### **Automatic Collection System (ACS)**

- All glass system designed to collect water from the Barnstead Mega-Pure Glass Stills and control their operation
- Generous 45 L capacity
- Easy wall- or bench-mount
- Preset to shut off the power to the heaters and the water supply to the still after the bottle has collected approximately 45 L
- The system will automatically restart the still and refill the collection system when the purified water supply has been depleted to approximately 38 L

#### Pretreatment for hard water

- · Allows for pretreating feed water to Barnstead Mega-Pure Glass Stills
- B-Pure single- or dual-cartridge deionizer reduces scale build up and increases distillate purity
- Greater flexibility in meeting varying feed water qualities
- Purity light measures water quality and indicates when cartridges are exhausted
- · Large cartridge capacity for extended cartridge use
- · Unit must be wall mounted



Plastic bottle (413964)



**Automatic Collection System** 



Barnstead B-Pure Cartridge Holder



13 L bottle (413934)

# Barnstead Mega-Pure Glass Stills

#### **Product specifications**

Inlet water temperature				
°C (°F)	Vent temperature °C (°F)	Auto drain	Inlet pressure (psig)	Auto start/stop controls
4-37 (39.2-98.6)	85-96.7 (185-206)	No	20-100	Yes

### System options ordering information

	Volume	Overall dimensions	Electrical (5	60/60 Hz)		Cat. No.	
	output						6 L bottle
	(L/hr)*	W x H x D in (cm)	Volts	Amps	Phase	Unit	included
MegaPure-1**	1.4	18 x 34 x 9.75 (45.7 x 86.4 x 24.8)	120	9	1	A440266	A440267
MegaPure-3A**	3.4	23 x 45 x 12	240	11	1	A440367	_
Megarure-SA	3.4	(58.4 x 114.3 x 30.5)	208	13	1	A440696	_
MegaPure-6A**	6	23 x 45 x 12	240	21	1	A440518	_
iviegarure-bA 6	· · · · · · · · · · · · · · · · · · ·	(58.4 x 114.3 x 30.5)	208	25	1	A440697	_

<sup>\*</sup> Volume Output for 208 V units can be reduced up to 25 % less than specified. \*\*

### Required accessories ordering information

Description	Cat. No.
Storage solutions See Storage Options	
MP-1 Flow Meter (Required for MP-1 units only)	440092

### Optional accessories ordering information

Description			Cat. No.	
			D4511 Single B-Pure system;	
Single-cartridge deionizer Single-cartridge deionizer removes inorganic contaminants. Use with adapter kit to connect			• 331040 Adapter Kit;	
to Megapure still. Wall mount only.	7 x 24 x 7 in (17.8 x 61	x 17.8 cm)	D0803 High Capacity Cartridge (Recommended)	
Double-cartridge deionizer			<ul> <li>D4524 Double B-Pure system with 50 Kohm light;</li> </ul>	
Provides increased capacity to remove inorganic contaminates. Suggested combination— (2) High Capacity cartridges. Use with adapter kit to connect to MegaPure Glass-Still.			• 331041 Adapter kit;	
Wall mount only. 15 x 27 x 7 in (38.1	x 68.65 x 17.8 cm)		D0803 High-capacity cartridge (Qty 2 recommended)	
Deionization cartridges	High capacity cartridge Used in B-Pure system; removes ionized solids for extended capacity		D0803	
	Organic removal coursed in B-Pure for	artridge rorganic and chlorine removal	D0813	
Dual solenoid Required if using service DI or RO as boiler feed water and untreated (tap) cooling water		MP-1 (240 V), MP-3, MP-6A	440236	
		MP-1 (120 V)	RY798X2A	



## Ordering information

Storage options		Cat. No.
Storage bottles	Plastic bottle, 6 L capacity For use with MegaPure-1 and MegaPure-3A Glass Stills	
or manual operation	Glass bottle, 13 L capacity For use with any Mega-Pure Glass Still	413934
Automatic Collection System (ACS)	Automatic Collection System (ACS) unit	B440704
For use with any Mega-Pure Glass Still Collects up to 45 L of water and controls the still Dimensions: W x H x D 16.25 x 38 x 16.25 in. (41.3 x 96.5 x 41.3 cm)	Flexible tubing adapter kit  Needed when the supplied glass tubing is either too long or too short to reach still; flexible tubing is not included. 0.5 in. Teflon tubing able to handle 80-90°C water is recommended	440138
	ACS wall mounting bracket	400634
	Adapter for ACS washer hook-up for glassware washers	440241

