




Communication board

User Guide

REF

160140, 161535, 161537, 161635, 161990, 161520, 161636, 160184,
160185, 10076964, 10076965

624728/04 | 2020-09-15

© 2020 Hamilton Medical AG. All rights reserved. Printed in Switzerland.

No part of this publication may be reproduced, stored in a database or retrieval system, or transmitted in any form or by any means, electronic, mechanical, or by photocopying, recording, or otherwise, without prior written permission of Hamilton Medical AG.

This document may be revised, replaced, or made obsolete by other documents by Hamilton Medical AG at any time and without notice. Ensure that you have the most current applicable version of this document; if in doubt, contact the technical support department of Hamilton Medical AG, Switzerland. While the information set forth herein is believed to be accurate, it is not a substitute for the exercise of professional judgment.

Nothing in this document shall limit or restrict in any way Hamilton Medical AG's right to revise or otherwise change or modify the equipment (including its software) described herein, without notice. In the absence of an express, written agreement to the contrary, Hamilton Medical AG has no obligation to furnish any such revisions, changes, or modifications to the owner or user of the equipment (including software) described herein.

The equipment must be operated, serviced, or upgraded only by trained professionals. Hamilton Medical AG's sole responsibility with respect to the equipment and its use is as stated in the limited warranty provided in the device *Operator's Manual*.

Hamilton Medical AG shall not be liable for any loss, cost, expense, inconvenience, or damage that may arise out of misuse of the product, or if non-Hamilton Medical AG parts were used when replacing parts, or if serial numbers were amended, deleted, or removed.

If returning parts to Hamilton Medical AG, be sure to use the standard Hamilton Medical returned goods authorization (RGA) procedure. Disposal of parts shall follow all local, state, and federal regulation with respect to environmental protection.

For all proprietary as well as third-party trademarks used by Hamilton Medical AG, see www.hamilton-medical.com/trademarks. Product and/or company names marked with a [§] symbol may be the trademarks and/or registered trademarks of their respective owners, including but not limited to Aerogen[§], Nihon Kohden[§], Masimo[§], Capnostat[§], Bluetooth[§], and Wi-Fi[§].

Manufacturer

Hamilton Medical AG
Via Crusch 8, CH-7402 Bonaduz
Switzerland
Phone: (+41) 58 610 10 20
Fax: (+41) 58 610 00 20
info@hamilton-medical.com
www.hamilton-medical.com

Preface

Conventions used in this guide



WARNING

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.



CAUTION

Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.

NOTICE

Emphasizes information of particular importance.

1 Overview

Hamilton Medical ventilators offer extended communication options using the communication boards described in the following tables.

Table 1. HAMILTON-C1/C2/C3/C6/T1 communication boards

Communication board	REF	Description
HAMILTON-C2, HAMILTON-C3		
CO2, SpO2, Nurse Call	160140	Support for CO2 monitoring and remote alarm capability
HAMILTON-T1		
CO2, COM1, Nurse Call	161535	Support for CO2 monitoring, data transfer, and remote alarm capability
CO2	161537	Support for CO2 monitoring
CO2, SpO2,  /COM1 ^{1, 2}	10076965	Support for CO2 and SpO2 monitoring, interface with HAMILTON-H900 humidifier, and data transfer ¹
CO2, SpO2, COM1 ² (red)	161635	Support for CO2 and SpO2 monitoring, and data transfer
CO2, SpO2, COM1 ² (gray)	161990	Support for CO2 and SpO2 monitoring, and data transfer
HAMILTON-C1		
CO2, COM1, Nurse Call	161520	Support for CO2 monitoring, data transfer, and remote alarm capability
CO2, SpO2,  /COM1 ^{1, 2}	10076964	Support for CO2 and SpO2 monitoring, interface with HAMILTON-H900 humidifier, and data transfer ¹
CO2, SpO2, COM1 ²	161636	Support for CO2 and SpO2 monitoring, and data transfer

¹ RS-232 connection over the /COM1 port is only available when using the communication Y-cable (PN 10077038).
² Available for HAMILTON-T1 ventilators with SN 3000 and higher, HAMILTON-C1 ventilators with SN 6000 and higher.

Communication board	REF	Description
HAMILTON-C6		
CO ₂ , SpO ₂	160185	Support for CO ₂ and SpO ₂ monitoring
CO ₂ , SpO ₂ , Aerogen	160184	Support for CO ₂ , SpO ₂ and Aerogen nebulization data monitoring

The communication boards are intended to be used for the following tasks:

Table 2. Supported communication options

Using the connector...	The ventilator can...
CO ₂	Monitor the patient's exhaled CO ₂ with an attached CO ₂ sensor (<i>not included</i>).
SpO ₂ ^{3,4}	<i>HAMILTON-C1/T1/C3/C6 only.</i> Monitor the patient's SpO ₂ , pulse rate, and optional related pulse-oximetry data with an attached SpO ₂ adapter and sensor (<i>not included</i>).
COM	<i>HAMILTON-C1/T1 only.</i> Send monitored data, ventilator settings, and alarms to a patient monitor, patient data management system (PDMS), or other computer system. On the HAMILTON-C6 ventilator, the COM ports are included as extended communication ports; the communication board is not required.
 COM1	<i>HAMILTON-C1/T1 only.</i> Support remote access to status and controls of a connected HAMILTON-H900 humidifier (<i>not included</i>), and send monitored data, ventilator settings, and alarms to a patient monitor, patient data management system (PDMS), or other computer system over an RS-232 connection. ⁵

³ Currently only available for HAMILTON-T1 ventilators with SN 3000 and higher, and HAMILTON-C1 ventilators with SN 6000 and higher.
⁴ SpO₂ monitoring and the COM connection are NOT supported on HAMILTON-C2 ventilators.
⁵ RS-232 connection over the /COM1 port is only available when using the communication Y-cable (PN 10077038).

Using the connector...	The ventilator can...
Nurse Call	<i>HAMILTON-C1/T1/C2/C3 only.</i> Send alarm and timing signals to external stations. On the HAMILTON-C6 ventilator, the Nurse call port is included as an extended communication port; the communication board is not required.
Aerogen	<i>HAMILTON-C6 only.</i> Support Aerogen nebulizer integration (<i>not included</i>).

Table 3. HAMILTON-C 1/T1 communication board configurations





CO2, SpO2, /COM1	CO2, SpO2, COM1
	
CO2, COM1, Nurse Call	CO2 <i>HAMILTON-T1 only</i>
	

Table 4. HAMILTON-C2/C3 communication board configurations



CO2, SpO2, Nurse Call^d


Table 5. HAMILTON-C6 communication board configurations

CO2, SpO2	CO2, SpO2, Aerogen
	

2 Installing a communication board

WARNING

- Additional equipment connected to medical electrical equipment must comply with the respective IEC or ISO standards. All configurations must comply with the requirements for medical electrical systems, IEC 60601-1, clause 16.
- Anybody connecting additional equipment to medical electrical equipment configures a medical system and is responsible for ensuring that the system complies with the requirements for medical electrical systems. Local laws take priority over the above-specified requirements.

CAUTION

To prevent damage to sensitive electronic circuits, do not touch the circuit board, and be sure to use ESD-protective gear.

NOTICE

- The communication board must be installed by a Hamilton Medical authorized service engineer.
- Use only the original screws provided with the ventilator.
- The communication board includes protective covers for the connector ports. When a port is not in use, make sure the cover is in place, sealing the port.


- If you remove the communication board, be sure to reattach the cover plate to protect the connection components.
- All of the listed options are not supported on every ventilator model. See Table 1 to determine what is supported on your device.

Before you begin

- Ensure the ventilator is turned off and disconnected from all power sources:
 - Disconnect the ventilator from external power.
 - Remove the batteries. For details, see your ventilator *Operator's Manual*.
- Have an appropriate screwdriver available.

To install the communication board in a HAMILTON-C1/T1 ventilator


1. Perform the electrical safety text as described in the ventilator *Service Manual*.
2. Remove the cover plate from the ventilator by removing the two flat-head screws (1).



3. Carefully insert the communication board into the open slot and secure the board using the same screws. Do not overtighten the screws.
4. Be sure to store the cover plate after removal for potential future use.

To install the communication board in a HAMILTON-C2/C3 ventilator

1. Perform the electrical safety text as described in the ventilator *Service Manual*.
2. Remove the cover plate from the ventilator by removing the two Phillips screws (1).



3. Carefully insert the communication board into the open slot and secure the board using the same screws. Do not overtighten the screws.
4. Be sure to store the cover plate after removal for potential future use.

To install the communication board in a HAMILTON-C6 ventilator

Installation of the communication board is only permitted by Hamilton Medical technical service personnel. Contact your Hamilton Medical technical representative for additional information.

Next steps

- For additional required components, see Section 3.
- Options (CO₂, SpO₂) and sensors must be enabled on the ventilator before they are operational. For details, see the ventilator *Operator's Manual*.
- For more information, and when applicable, connector diagrams and pin assignments for the various cables (COM, Nurse Call), see the associated user guides and the ventilator *Operator's Manual*.
 - *Nurse Call Cable User Guide* PN 624344
 - *COM1 Cable User Note* PN 624442

3 Additional required components

The communication boards do not include any required cables, sensors, or other components required to enable the supported functionality.

This section describes the additional components required for each connection option, available from Hamilton Medical.

Table 6. Required components for use with communication boards

Communication option	Additional required components
CO2 monitoring	Order the appropriate Hamilton Medical CO2 sensor: • Capnostat 5 CO2 sensor (PN 281718) • Capnostat 5 CO2 sensor with angled connector (PN 282157) • LoFlow CO2 sidestream sensor (PN 281928)
SpO2	<i>HAMILTON-C1/T1/C3/C6 only.</i> • Masimo: Masimo SpO2 Starter Kit (PN 950210, 950201, or 10104745) and the appropriate SpO2 sensors. • Nihon Kohden SpO2 sensor adapter (PN 282010) and the appropriate SpO2 sensors. For details, see the SpO2-measurement-related catalog.
 /COM1	<i>HAMILTON-C1/T1 only.</i> Two communication cables are supported for HAMILTON-H900 humidifier integration: • Communication cable to HAMILTON-H900 (PN 950473) • Communication Y-cable to HAMILTON-H900 and with RS232 interface (PN 10077038) RS-232 connection over the /COM1 port is <i>only</i> available when using the communication Y-cable (PN 10077038).
COM	<i>HAMILTON-C1/T1 only.</i> Order the COM cable (PN 161545, 161650). Additional third-party components may be required to connect to the target device. For interface and connection details, see the ventilator <i>Operator's Manual</i>.

Communication option	Additional required components
Nurse Call	Order the appropriate Nurse Call cable for your requirements (PN 160166, PN 161761, or PN 161196). For additional details, see <i>Nurse Call Cable User Guide</i> .
Aerogen	<i>HAMILTON-C6 only</i> . Order an Aerogen nebulizer system. See the <i>Hamilton Medical e-catalog</i> for details.

4 Enabling options on the ventilator


After installing the communication board, you must also enable the CO2 and/or SpO2 hardware, define SpO2 sensor settings (if applicable), and enable the sensors on the ventilator.

- In **Configuration**:
 - Enable the **CO2** and/or **SpO2** option.
 - Configure SpO2 sensor data settings, if needed.
- In the **System > Sensors On/off** window, enable CO2 and/or SpO2 monitoring.

These steps are described in detail in your ventilator *Operator’s Manual* and the *Pulse Oximetry Instructions for Use*.

HAMILTON-C1/T1, HAMILTON-H900 integration

By default, the humidifier integration interface is enabled. When integration is enabled, the humidifier droplets are displayed in the lower right corner of the ventilator main display.



If the droplets are absent, integration is disabled. In this case, the **Configuration > SW options** window lists the entry, **Humidifier GUI disabled**.

To enable the humidifier integration interface, see the *Software Installation Guide* and contact your Hamilton Medical technical representative for the appropriate option key.

5 Disposal

Dispose of all parts removed from the device according to your institution’s protocol. Follow all local, state, and federal regulations with respect to environmental protection.



For more information:
www.hamilton-medical.com



HAMILTON
MEDICAL

Intelligent Ventilation since 1983



Hamilton Medical AG

Via Crusch 8, 7402 Bonaduz, Switzerland

+41 (0)58 610 10 20

info@hamilton-medical.com

www.hamilton-medical.com