
Philips M1032A VueLink

**External Device
Service Booklet**

PHILIPS

Part Number M1032-9000Y

4512 610 18581

Published in October 2006

Edition 22

Y.00.00

Notice

Philips makes no warranty of any kind with regard to this material, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose. Philips shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this material.

This document contains proprietary information that is protected by copyright. All rights are reserved. No part of this document may be photocopied, reproduced or translated to another language without prior written consent of Philips Medical Systems. The information contained in this document is subject to change without notice. Philips assumes no responsibility for the use or reliability of its software on equipment that is not furnished by Philips.

Important

United States federal law restricts these devices to sale by or on the order of a physician.

Intended Use

Description

The Philips M1032A VueLink module is a plug-in module for use with a Philips IntelliVue, V24 or V26 Patient Monitor or a Philips CMS Patient Monitoring System.

The Philips M1032A VueLink module is powered via the patient monitor.

Purpose

The Philips M1032A VueLink module transmits information from a connected external device to a Philips IntelliVue, V24 or V26 Patient Monitor or a Philips CMS Patient Monitoring System. **The Philips M1032A VueLink module is not a therapeutic device.**

Environment

The Philips M1032A VueLink module is intended to be used in a clinical environment by trained healthcare professionals. It is not intended for home use.

Indications for Use

Condition	The use of a Philips M1032A VueLink module is generally indicated when the clinician decides there is a need to include the physiological and device parameters and waves and alarms generated by external devices into the Philips IntelliVue, V24 or V26 Patient Monitor or a Philips CMS Patient Monitoring System.
Part of Body or Type of Tissue with Which the Device Interacts	The Philips M1032A VueLink module does not contact the body or tissue of the patient. Signals are obtained from external devices.
Frequency of Use	The Philips M1032A VueLink module is indicated for use when prescribed by a clinician.
Physiological Purpose	The Philips M1032A VueLink module is indicated when the purpose is to gain information for treatment, to assess adequacy of treatment, or to rule out causes of symptoms. The Philips M1032A VueLink module is well suited for patient monitoring.
Patient Population	Please refer to the documentation describing the external device.
Prescription Versus Over-the-Counter	The Philips M1032A VueLink module is a prescription device.

Warnings, Cautions and Notes

In this guide:

Warning

- *A warning alerts you to a potential serious outcome, adverse event or safety hazard.
Failure to observe a warning may result in death or serious injury to the user or patient.*
-
-

Caution

A caution alerts you where special care is necessary for the safe and effective use of the product.
Failure to observe a caution may result in minor or moderate personal injury or damage to the product or other property, and possibly in a remote risk of more serious injury.

Note

A note gives special instructions to highlight an operating procedure or practice. Notes may precede or follow the applicable text.

Responsibility of the Manufacturer

Philips Medical Systems only considers itself responsible for any effects on safety, reliability and performance of the equipment if:

- assembly operations, extensions, re-adjustments, modifications or repairs are carried out by persons authorized by Philips, and
- the electrical installation of the relevant room complies with national standards, and
- the Philips M1032A VueLink module is used in accordance with the instructions for use contained in this document, and in the relevant chapter of the user documentation for the patient monitor into which the module is plugged (Philips IntelliVue, V24 or V26 Patient Monitor or a Philips CMS Patient Monitoring System).

Warning

- *Failure on the part of the responsible individual hospital or institution employing the this equipment to implement a satisfactory maintenance schedule may cause undue equipment failure and possible health hazards.*

The maintenance schedule is defined in the Testing and Maintenance chapter of the IntelliVue Patient Monitor Service Guide, the Philips CMS Patient Monitoring System Service Guide, or the Philips V24/V26 Component Monitoring System Service Guide, corresponding to your monitor.

Note

Not all the device drivers are available in all countries.

The Instructions for Use for all drivers are translated into the following languages: Czech, Danish, Dutch, Estonian, Finnish, French, German, Greek, Hungarian, Italian, Japanese, Norwegian, Polish, Portuguese, Romanian, Russian, Simplified Chinese (China), Slovak, Spanish, Swedish, Traditional Chinese (Taiwan), Turkish. The driver software may not be available in all of these languages. Please refer to the CMS-PM softserver (<http://pww.softserver.anr.ms.philips.com/navigation.asp>,

subfolder Vuelink, subfolder Additional Information) for a list of the translated drivers that are available.

Please also check your local regulatory and business requirements for permission to sell.

Philips authorized technical personnel are responsible for the setup, configuration and repair of the VueLink module and all related Philips equipment. The configuration, setup and repair of external devices furnished by manufacturers other than Philips Medical Systems must be carried out by local site personnel, or by representatives of the external device manufacturer.

Philips Medical Systems is not responsible for any problems arising from inaccurate or erroneous data displayed on the Philips patient monitor that is received from external devices furnished by manufacturers other than Philips Medical Systems.

Philips Medical Systems makes every effort to ensure that signal names are maintained between the Philips patient monitor and any interfaced external devices. However, signal names within the Philips patient monitor also need to remain consistent. This means that in some cases the wave and numeric labels used on the Philips monitor may be different from those used on the interfaced external device.

Philips Medical Systems will not be held responsible for any errors in the configuration of the VueLink module for connection to free-analog devices.

If manufacturers other than Philips Medical Systems make changes to any VueLink supported device, Philips Medical Systems will try to ensure that interfacing to that device remains possible. However, Philips Medical Systems reserves the right to discontinue interfacing to any device.

Introduction

This Service Guide contains technical details for the connection of external devices to the VueLink module. It is not a comprehensive, in-depth explanation of the product architecture or technical implementation.

Who Should Use This Guide

This guide is for biomedical engineers or technicians responsible for connecting external devices to Philips patient monitoring systems using the VueLink module.

How to Use This Booklet

This booklet provides information on the worldwide external devices that can be connected to the VueLink module. It should be printed and inserted into Appendix A of the VueLink Module Handbook, completely replacing the existing worldwide driver information.

The information contained in this document is subject to change without notice.

Philips Medizin Systeme Böblingen GmbH
Hewlett-Packard Str. 2
71034 Böblingen
Germany

© 2006 Philips Medizin Systeme Böblingen GmbH
All rights are reserved.

Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Connecting your Device to VueLink

There are two different ways to connect your Device to VueLink:

- Connection via a Supported Device driver.
This uses one of the following device outputs:
 - Analog
 - RS-232 Digital
 - Analog/Digital combination
- Connection via a Free Analog Device driver.
This only uses device analog outputs.

Supported Device Drivers

Supported Device Drivers are preconfigured in the VueLink module. A dedicated ready-prepared VueLink cable is available for connecting the external device to the module. (Please refer to page A-4 for cable part numbers).

Free Analog Device Drivers

Free Analog Device Drivers must be manually configured into the module at installation. A Free Analog VueLink cable is available for connecting the Free Analog external device to the VueLink module. One end of the cable is unterminated and must be prepared for connection to the external device on-site. (Please refer to page B-4 for cable option part numbers).

Connecting your Device to VueLink

Contents

Connecting your Device to VueLink	9
Supported Device Information	15
Driver Part Numbers for Supported Devices	16
VueLink Cable Overview	17
Auxiliary Devices	23
Abbott Oximetry 3 - SO2 Monitor	23
Critikon Dinamap 1846 / 1846 SX - NBP Monitor	25
Danmeter AAI/AEP Monitor	27
Nellcor® N-200 - SpO2 Monitor	29
Nellcor® N-100C - SpO2 Monitor	31
Ventilators	33
Nellcor Puritan-Bennett 740 / 760 / 840 Ventilators	33
Puritan-Bennett 7200a / 7200ae Ventilator	36
Dräger Babylog 8000 Ventilator	38
Dräger Evita Ventilator	40
Dräger Evita 2 Ventilator	42
Dräger Evita 2 Dura, Dräger Evita 4 & Dräger Evita XL Ventilators	44
Dräger Graphic Screen for Babylog 8000	46
Dräger Graphic Screen for Savina	48
Dräger Savina	50
GE Engström Carestation Ventilator	52
Hamilton Veolar, Veolar FT, and Amadeus Ventilators	54
BEAR 1000 Ventilator	56
Infrasonics® Infant Star / ISV Ventilator	58
Siemens SCM 990 Ventilator	61
Siemens 900 C/D/E Ventilator	64
Maquet Servo 300/300A Ventilator	66
Maquet SERVO-i Ventilator	68
Maquet SERVO-s Ventilator	70
Gas Analyzers	73
HP M1025A/B Gas Analyzer	73
Datex Capnomac (II and Ultima) Gas Analyzer	75
Dräger Vamos Variable Anaesthetic Gas Monitor	78

Ohmeda RGM 5250 Respiratory Gas Monitor	80
Ohmeda Rascal II Anesthetic Gas Monitor	82

Anesthesia Machines85

Dräger Cato Anesthesia Device	85
Dräger Cicero Anesthesia Machine	87
Dräger Cicero EM Anesthesia Machine	89
Dräger Fabius GS/Tiro Anesthesia Machine	91
Dräger Julian Anesthesia Workstation	93
Dräger Pallas Anesthesia Workstation	95
Dräger PM 8050 Anesthesia Machine	96
Dräger Primus/Apollo Anesthesia Workstation	98
Dräger Zeus Anesthesia Device	100
GE Aestiva and Avance	102
GE Aisys	105
North American Dräger Narkomed 2B and Narkomed 2C Anesthesia Systems	107
North American Dräger Narkomed 4 Anesthesia System	112
North American Dräger Narkomed GS Anesthesia System	114
North American Dräger Narkomed 6000 Anesthesia System	116
Ohmeda 7800/7810	118
Ohmeda 7900	120
Ohmeda Modulus CD Anesthesia Machine	122
Taema Alys Anesthesia Machine	125

Auxiliary Plus Devices127

Aspect Medical Systems A-2000™ BIS™ Monitor	127
Baxter Explorer Cardiopulmonary Monitor	129
B Braun SpaceCom	131
Braun FM Device Information	134
Diametrics Medical IRMA SL Series 2000	136
Diametrics Medical Trendcare Monitor	138
Edwards CardiacOutput-Oximetry Monitor Family	140
Fresenius Vial Base A/DPS/MVP	144
GE DINAMAP ProCare 100 Series	146
i-STAT® I Analyzer	148
Mortara ELI 100/STM 12-Lead Monitor	150
Organon Teknika TOF-Watch SX	152

Free Analog Device Information155

Free Analog Device Information	155
Configuring the Free Analog Device.	156
Dräger Evita Ventilator	158
Engström ELSA/EAS Anesthesia System	166
Nellcor N1000 and N2500 Multi-Function Monitor	174
Radiometer Transcutaneous pO ₂ /pCO ₂ Monitor.	182

AI

Supported Device Information

This appendix contains information about the supported devices that can be connected to the VueLink module. The following information is provided for each supported device:

- The part number of the corresponding device driver.
- The configuration settings required on the external device.
- The placement of the output connectors on the external device.

A cable matrix has also been compiled, detailing the cable part numbers for each supported device.

Caution

Only use cables specified by Philips Medical Systems when connecting supported devices to the VueLink module. Use of unauthorized or incorrect cables may result in damage to supported devices and may cause incorrect data to be displayed on the Philips patient monitor.


The configuration settings required on the external devices are detailed on the following pages. Procedures on how to configure the external devices are not provided. Configuration procedures can be found in the documentation supplied with the external devices. Philips personnel are only responsible for configuring Philips manufactured external devices.

Note

The signal labels used on the Philips patient monitor may be different from those given on the external device. The labels used on the Philips patient monitor are listed in the Philips M1032A VueLink External Device Instructions for Use.

Driver Part Numbers for Supported Devices

Each driver is allocated a part number, using the driver key detailed below:



VueLink Cable Overview

3 Separate part numbers are assigned to each VueLink Cable:

M1032A KXX option	To order a cable with a new module.
M1182A KXX option	To order a cable for an existing module or for upgrades (Sales).
M1032-XXXXXX option	To order a cable for an existing module from SMO/SLI. (Support or Immediate Requirement)

The Supported Device cable options are detailed in the matrix opposite:

VueLink Cable Overview

Device	New Module Option (M1032A #)		Patient Monitor Cable Option (M1182A #)		SMO/SLI Option (M1032-)	
	2m	4m	2m	4m	2m	4m
Auxiliary Devices						
Abbott Oximetrix 3	K73	-	K73	-	61619	61620
Critikon 1846/1846 SX	K72	-	K72	-	61613	61614
Danmeter AAI/AEP Monitor	K56	-	K56	-	61693	-
Nellcor N-200	K71	-	K71	-	61617	61618
Nellcor N-100C	K70	-	K70	-	61615	61616
Ventilators						
Dräger Babylog 8000	-	K2H	-	K2H	-	61644
Dräger Evita	K24	K2E	K24	K2E	61629	61630
Dräger Evita 2	-	K2G	-	K2G	-	61643

VueLink Cable Overview

Device	New Module Option (M1032A #)		Patient Monitor Cable Option (M1182A #)		SMO/SLI Option (M1032-)	
	2m	4m	2m	4m	2m	4m
Dräger Evita 2 dura Dräger Evita 4 Dräger Evita XL Ventilator	-	K2M	-	K2M	-	61680
Dräger Graphic Screen for Babylog 8000	-	K29	-	K29	-	61604
Dräger Graphic Screen for Savina	-	K28	-	K28	-	61604
Dräger Savina	-	K27	-	K27	-	61644
GE Engström Carestation	-	K2N	-	K2N	-	61607
Hamilton Amadeus	-	271	-	271	-	61673
Hamilton Veolar	-	271	-	271	-	61673
Maquet Servo 300	-	K2F	-	K2F	-	61642
Maquet SERVO-i	-	K23	-	K23	-	61696
Maquet SERVO-s	-	K26	-	K26	-	61696
Nellcor Puritan Bennett 740	-	K25	-	K25	-	61665
Nellcor Puritan Bennett 760	-	K25	-	K25	-	61665
Nellcor Puritan Bennett 840	-	K25	-	K25	-	61665
Puritan-Bennett 7200a/ae	K20	K2A	K20	K2A	61621	61622
Siemens 900 C/D/E	K21	K2B	K21	K2B	61623	61624
Siemens SCM 990	K22	K2C	K22	K2C	61625	61626
BEAR 1000	-	K2K	-	K2K	-	61657

VueLink Cable Overview

Device	New Module Option (M1032A #)		Patient Monitor Cable Option (M1182A #)		SMO/SLI Option (M1032-)	
	2m	4m	2m	4m	2m	4m
Infrasonics Infant Star / ISV	-	K2J	-	K2J	-	61645
Gas Analyzers						
Datex Capnomac (II & Ultima)	K41	-	K41	-	61633	-
Dräger Vamos	K4A	-	K4A	-	61695	-
Philips M1025A/B	K40	-	K40	-	61631	-
Ohmeda RGM 5250 (old) ¹	-	-	-	-	61635	-
Ohmeda RGM 5250 (new)	K43	-	K43	-	61636	-
Ohmeda Rascal II	K45	-	K45	-	61664	-
Anesthesia Machines						
Dräger Cato	K02	-	K02	-	61602	-
Dräger Cicero	K01	-	K01	-	61601	-
Dräger Cicero EM mono	K10	-	K10	-	61675	-
Dräger Cicero EM color	-	K12		K12	-	61685
Dräger Fabius GS/Tiro	K49	-	K49	-	61700	-
Dräger Julian	K11	-	K11	-	61681	-
Dräger Pallas	K4C	-	K4C	-	61694	-
Dräger PM 8050	270	-	270	-	61676	-
Dräger Primus/Apollo	K48	-	K48	-	61694	-
Dräger Zeus	-	K4B	-	K4B	-	61666

VueLink Cable Overview

Device	New Module Option (M1032A #)		Patient Monitor Cable Option (M1182A #)		SMO/SLI Option (M1032-)	
	2m	4m	2m	4m	2m	4m
GE Aestiva	K15	-	K15	-	61684	-
GE Aisys	K16	-	K16	-	61684	-
GE Avance	K15	-	K15	-	61684	-
NAD ² Narkomed 2B/2C	K05	-	K05	-	61605	-
NAD ² Narkomed 4	K06	-	K06	-	61605	-
NAD ² Narkomed GS	K47	-	K47	-	61605	-
NAD ² Narkomed 6000	K46	-	K46	-	61605	-
Ohmeda Modulus CD	K03	-	K03	-	61603	-
Ohmeda 7800	K07	-	K07	-	61658	-
Ohmeda 7810	K08	-	K08	-	61659	-
Ohmeda 7900	K09	-	K09	-	61678	-
Taema Alys	274	-	274	-	61674	-
Auxiliary Plus Devices						
Aspect Medical Systems A-2000 BIS Monitor	-	K52	-	K52	-	61687
Baxter Explorer	K74	-	K74	-	61651	-
Baxter Vigilance	-	-	-	-	61652	-
B Braun SpaceCom	-	K59	-	K59	-	61649
Braun FM	-	275	-	275	-	61691

VueLink Cable Overview

Device	New Module Option (M1032A #)		Patient Monitor Cable Option (M1182A #)		SMO/SLI Option (M1032-)	
	2m	4m	2m	4m	2m	4m
Braun FM software revision < 3.0	-	-	-	-	-	61648
Diametrics Medical IRMA SL Series 2000	-	K54	-	K54	-	61688
Diametrics Medical Trendcare monitor	-	K55	-	K55	-	61689
Edwards CardiacOutput-Oximetry Monitor Family	K75	-	K75	-	61692	-
Fresenius Vial Base A/DPS/MVP	-	277	-	277	-	61682
GE DINAMAP ProCare 100 Series Patient Monitor	-	K78	-	K78	-	61606
i-STAT [®] 1 Analyzer	-	K58	-	K58	-	61667
Mortara ELI 100/STM	-	K76	-	K76	-	61653
Organon Teknika TOF-Watch SX	K53	-	K53	-	61686	-
Free Analog Devices						
Free Analog	K60	K6A	K60	K6A	61611	61612
Open Interface Cable (25-pin digital/analog)	-	K6B	-	K6B	-	61654
Open Interface Cable (9-pin digital)	-	K6C	-	K6C	-	61699


1. Old cable only available for support in case an existing cable fails.
2. NAD = North American Dräger

A2

Auxiliary Devices

Abbott Oximetrix 3 - SO₂ Monitor

Device Driver Name:	Abbott Oximetrix 3
Device Driver P/N:	M1032-A07rl
	where: r = revision l = language
Supported Devices:	Oximetrix (R)3 SO ₂ Monitor (version 104 and 105)
Connection:	RS-232 Digital (Fixed configuration)
	Baud Rate: 1200
	Word Length: 8 bits
	Stop Bits: 1
	Start Bits: 1
	Parity: None
Switch Settings:	Factory default



RS232 DIGITAL
CONNECTOR


Abbott Oximetrix 3 - SO₂ Monitor Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Critikon Dinamap I846 / I846 SX - NBP Monitor

Device Driver Name:	Critikon 1846 / 1846 SX
Device Driver P/N:	M1032-A06rl
	where: r = revision l = language
Supported Devices:	Critikon 1846 (with software revision 1846RCM and 1846PRBG) Critikon 1846 SX (with software revision SXRCH and SXPRDH)
Connection:	RS-232 Digital (Fixed configuration)
	Baud Rate: 600
	Word Length: 8 bits
	Stop Bits: 1
	Parity: None.
Switch Settings:	Factory default




Critikon Dinamap 1846 / 1846 SX - NBP Monitor Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Danmeter AAI/AEP Monitor

Device Driver Name:	AAI/AEP Monitor
Device Driver P/N:	M1032-A08rl
	where: r = revision. l = language.
Supported Devices:	Danmeter AAI/AEP Monitor (SBC(software) Version 1.5)
Connection:	RS-232 Digital (at the rear of the device)
	Baud Rate: 9600 (Low ¹)
	Word Length: 8 bits
	Stop Bits: 1
	Start Bits: 1
	Parity: None
	Protocol: Device Link (ON ¹)



Danmeter AAI/AEP Monitor Rear Panel

1. Must be set and verified according to the document "AAI/AEP Monitor Directions for Use".

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Nellcor® N-200 - SpO₂ Monitor

Device Driver Name: Nellcor® N-200

Device Driver P/N: M1032-A05rl

where:
 r = revision
 l = language


Supported Devices: Nellcor® N-200

Connection: Analog / RS-232 Digital combination

Baud Rate: 9600
 Word Length: 8 bits
 Stop Bits: 2
 Parity: None
 Analog: 1 alarm, 1 wave.

Switch Settings: See table below:

Name	Number	Position
SCALE and OUTPUT VOLT Switches	1	UP
	2	DOWN
DIP Switches	1	Either
	2	DOWN
	3	DOWN
	4	UP
	5	DOWN
	6	Either
	7	UP
	8	DOWN



Nellcor® N-200 - SpO₂ Monitor Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Nellcor® N-100C - SpO₂ Monitor

Device Driver Name: Nellcor® N-100

Device Driver P/N: M1032-A04rl

where:


r = revision

l = language

Supported Devices: Nellcor® N-100

Connection: Analog: 1 wave, 2 numerics, no alarms

Switch Settings: Set range switch to 0-10 V.



Nellcor® N-100C - SpO₂ Monitor Rear Panel

Note

The VueLink cable can get stuck in the outer cabinet without providing proper contact to the connectors. Apply moderate force to overcome this problem.

Note


Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

A3 Ventilators


Nellcor Puritan-Bennett 740 / 760 / 840 Ventilators

Device Driver Name:	NPB740/760/840
Device Driver P/N:	M1032-B48rl
	where: r = revision l = language
Supported Devices:	NPB 740 (sw rev. J). NPB 760 (sw rev. J). NPB 840 (sw rev. J).
Connection:	RS-232 Digital
	Baud Rate: 9600
	Word Length: 8 bits
	Stop Bits: 1
	Start Bits: 1
	Parity: None

Note RS-232 parameters must be set and verified according to the Nellcor Puritan-Bennett Service Manual



NPB 740 / 760 Ventilator Side Panel




NPB 840 Ventilator Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Puritan-Bennett 7200a / 7200ae Ventilator

Device Driver Name:	Bennett 7200a/ae
Device Driver P/N:	M1032-B02rl where: r = revision l = language
Supported Devices:	Puritan-Bennett 7200a (sw rev. V English only). Puritan-Bennett 7200ae (sw rev. D for Eng, Sp, Fr, Ger and Ital).
Connection:	Analog / RS-232 Digital combination Baud Rate: 9600 Word Length: 7 bits Stop Bits: 1 Parity: Even Analog: 2 waves
Switch Setting:	Internal; must be set by customer (see Puritan-Bennett User Guide).



Puritan-Bennett 7200a / 7200ae Ventilator Rear Panel


Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Dräger Babylog 8000 Ventilator

Device Driver Name:	Dräger Babylog 8000 Ventilator
Device Driver P/N:	M1032-B24rl where: r = revision l = language
Supported Devices:	Dräger Babylog 8000 with: <ul style="list-style-type: none">• SW Revision: 5.01• Digital Interface (optional): Babylink on MEDIBUS protocol 3.00.
Supported Device Options:	None
Connection:	RS-232 Digital Baud Rate: 9600 ¹ Word Length: 8 bits (fixed) Stop Bits: 1 (fixed) Parity: None ¹ Handshake: Software (fixed) Protocol: MEDIBUS ¹ (BabyLink)
Switch Settings:	Factory default

1. Must be set and verified in accordance with Dräger Babylog 8000 User's Manual.




Dräger Babylog 8000 Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Dräger Evita Ventilator

Device Driver Name:	Dräger Evita Ventilator
Device Driver P/N:	M1032-B11rl where: r = revision l = language
Supported Devices:	Dräger Evita Typ 8410614 with <ul style="list-style-type: none">• DW-Bus analog interface option number 8303940.• digital interface option number 8305327• SW Revision 11.01, 13.01 and 14.01, 13.02 or 14.02 for the digital interface and the Dräger Evita Ventilator.
Connection:	Analog / RS-232 Digital combination exchange Baud Rate: 9600 Word Length 7 bits (fixed) Stop Bits: 2 (fixed) Parity: Even (fixed) Handshake: Software (fixed) Analog: 2 waves
Switch Settings:	Factory default



Dräger Evita Ventilator Rear Panel

Note

Please ensure that the first two digits of the software version number shown on the Dräger Evita display after switching on the device are the same as the first two digits of the version number on the Evita's interface cards.


Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Dräger Evita 2 Ventilator

Device Driver Name:	Dräger Evita 2 Ventilator
Device Driver P/N:	M1032-B23rl where: r = revision l = language
Supported Devices:	Dräger Evita 2 with: <ul style="list-style-type: none">• SW Revision 2.00.• Digital interface (optional): EvitaLink 2.00 on MEDIBUS protocol 3.00.
Supported Device Options:	None
Connection:	RS-232 Digital Channel A. (Channel B can also be used as an alternative when set up as follows.) Baud Rate: 19200 ¹ Word Length: 8 bits (fixed) Stop Bits: 1 ¹ Parity: Even ¹ Handshake: Software (fixed) Protocol: MEDIBUS ¹

1. These settings are fixed for digital channel A. If channel B is used, they must be set and verified according to the Dräger Evita 2 User's Manual.




Dräger Evita 2 Ventilator Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Dräger Evita 2 Dura, Dräger Evita 4 & Dräger Evita XL Ventilators

Device Driver Name:	Dräger Evita 2 DuraVentilator Dräger Evita 4 Dräger Evita XL
Device Driver P/N:	M1032-B41rl where: r = revision l = language
Supported Devices:	Dräger Evita 2 Dura with: <ul style="list-style-type: none">• SW Revision 4.10.• MEDIBUS protocol 4.00. Dräger Evita 4 with: <ul style="list-style-type: none">• SW Revision 4.10.• MEDIBUS protocol 4.00. Dräger Evita XL with: <ul style="list-style-type: none">• SW Revision 6.00.• MEDIBUS protocol 4.00.
Supported Device Options:	None
Connection:	RS-232 C Baud Rate: 19200 Word Length: 8 bits (fixed) Stop Bits: 1 Parity: Even Handshake: None



Dräger Evita 2 Dura, Evita 4 and Evita XL Ventilator Rear Panel


Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Dräger Graphic Screen for Babylog 8000

Device Driver Name:	Dräger Graphic Screen for Babylog 8000
Device Driver P/N:	M1032-B64rl where: r = revision l = language
Supported Devices:	Dräger Babylog 8000 Software 5.01 MEDIBUS Version 3.00. Dräger Graphic Screen version 2.n MEDIBUS Version 4.00.
Supported Device Options:	None
Connection:	RS-232 port COM3 on rear side of the Graphic Screen Baud Rate: 9600 Word Length: 8 Sart Bits: 1 Stop Bits: 1 Parity: None Protocol: MEDIBUS

Note All RS 232 port settings are fixed.



Dräger Graphic Screen Rear Panel


Note

Please refer to the Philips M1032A VueLink External Device User's Information for Dräger for Graphic Screen for Babylog 8000 for exact details of numerics, INOPs and alarms available from the external device via Philips patient monitoring network.

Dräger Graphic Screen for Savina

Device Driver Name:	Dräger Graphic Screen for Savina
Device Driver P/N:	M1032-B63rl where: r = revision l = language
Supported Devices:	Dräger Savina Software 3.01 MEDIBUS Version 4.00 Dräger Graphic Screen version 2.n MEDIBUS Version 4.00.
Connection:	RS-232 port on rear side Baud Rate: 19200 Data Bits: 8 Start Bits: 1 Stop Bits: 1 Parity: None Protocol: MEDIBUS

Note All RS 232 port settings are fixed.




Dräger Graphic Screen Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device User's Information for Dräger for Graphic Screen for Savina for exact details of numerics, INOPs and alarms available from the external device via Philips patient monitoring network.

Dräger Savina

Device Driver Name:	Dräger Savina
Device Driver P/N:	M1032-B66rl where: r = revision l = language
Supported Devices:	Dräger Savina Software 3.01 MEDIBUS Version 4.00
Connection:	RS-232 port on rear side Baud Rate: 19200 Data Bits: 8 Start Bits: 1 Stop Bits: 1 Parity: None Protocol: MEDIBUS




Dräger Savina Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device User's Information for Dräger for Savina for exact details of waves, numerics, INOPs and alarms available from the external device via Philips patient monitoring network.

GE Engström Carestation Ventilator

Device Driver Name:	GE Engström
Device Driver P/N:	M1032-B67rl where: r = revision l = language
Supported Devices:	GE Engström Software 3.x Ohmeda Com 1.3 Serial Protocol
Connection:	RS-232 port on rear side (Port 4) Baud Rate: 19200 Data Bits: 7 Start Bits: 1 Stop Bits: 1 Parity: Odd



GE Engström Carestation Ventilator Rear Panel


Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs, and alarms available from the external device, and via Philips patient monitoring network.


Hamilton Veolar, Veolar FT, and Amadeus Ventilators

Device Driver Name:	Veolar ¹ Amadeus
Device Driver P/N:	M1032-B19rl where: r = revision l = language
Supported Devices:	<ul style="list-style-type: none">• Veolar and Veolar FT Ventilators Control Processor and Front Panel sw rev: V30 and V33 Interface sw rev: NIK01• Amadeus Ventilator Control Processor and Front Panel sw rev: A31 Interface sw rev: NIK01
Connection:	RS-232 Digital Baud Rate:9600 ² Word Length: 7 bits (fixed) Stop Bits: 2 bits (fixed) Parity: EVEN (fixed) Handshake: Flag control ² Program Selector: Position 9 ²

1. The product name Veolar is used throughout to describe both the Hamilton Veolar Ventilator, and the Hamilton Veolar FT Ventilator.
2. Must be set by Hamilton Service Engineer before connecting to VueLink.



Hamilton Veolar Ventilator Rear Panel



Hamilton Amadeus Ventilator Rear Panel

Note


Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs, and alarms available from the external device, and via Philips patient monitoring network.

BEAR 1000 Ventilator

Device Driver Name:	BEAR 1000
Device Driver P/N:	M1032-B28rl where: r = revision l = language
Supported Devices:	Bear 1000 (Basic and Comprehensive) SW Release 9.7 Bear GD 1000 Graphics Display (optional) Software P/N 51001-02007
Connection:	Analog (2 waves) / RS-232 Digital combination. No graphics display: Baud Rate: 9600 (Bear 1000 ventilator) Graphics display: Baud Rate: 19200 (Bear 1000 ventilator) 9600 (graphics display) Word Length: 8-bits (fixed) Stop Bits: 1 (fixed) Parity: None (fixed)

The 15-pin connector should be connected to the analog port on the ventilator. The 9-pin connector should either be connected to the RS-232 port on the ventilator, if

no Graphics Display is present, or to the Digital Out port on the Graphics Display. In the latter case, the RS-232 gender adapter must also be used.



BEAR 1000 Ventilator Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Infrasonics® Infant Star / ISV Ventilator

Device Driver Name: Infrasonics® Infant Star / ISV

Device Driver P/N: M1032-B30rl

where:

r = revision

l = language

Supported Devices:

- Original Infant Star Ventilator (yellow casing):
 - Standard Model TÜV Approved (software version 49, 51)
 - High Frequency Model TÜV Approved (software version 46, 50, 52)
 - Enh. High Frequency Model TÜV Approved (software version 82, 83)
- New Infant Star ISV Ventilator (grey casing):
 - Standard Model ISV 500 TÜV Approved (software version 105, 107)
 - High Frequency Model ISV 950 (software version 107)

Connection:


The device driver supports the following configuration for the serial communications port:

Baud Rate:	1200
Word Length:	8 bits
Stop Bits:	2
Parity:	None
Analog:	1 wave


These settings are fixed and cannot be changed.

Switch Settings:

None



Original Infrasonics Infant Star Ventilator Rear Panel




New Infrasonics Infant Star Ventilator Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Siemens SCM 990 Ventilator

Device Driver Name:	Siemens SCM 990 AD (Adult) Siemens SCM 990 NEO (Neonatal)
Device Driver P/N:	M1032-B08rl (Adult) M1032-B14rl (Neonatal) where: r = revision l = language
Supported Devices:	Siemens 900 C/D/E with Servo Computer Module 990 (serial number 713 or greater, sw rev 2.0).
Connection:	Analog / RS-232 Digital combination Baud Rate: 9600 Word length: 8 bits Stop Bits: 1 Parity: Even. Analog: 2 waves, 1 alarm, 1 numeric.
Switch Settings:	See following table.




Transmission Switch Settings

Transmission Switch Settings

Name	Switch Number	Position
ADDRESS Code Switches	1	0-9 (ANY)
	2	0-9 (ANY)
PARITY Switches (for RS-232 configuration)	1	ON
	2	ON
	3	OFF
	4	OFF
BAUD RATE switch	1	5

Note

The **restart** switch is to the right of the BAUD RATE switch. If any of the above switches are changed during power-on; the device power must be switched off and on.



Siemens SCM990 Ventilator Rear Panel

Note


Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Siemens 900 C/D/E Ventilator

Device Driver Name:	Siemens 900 C/D/E AD (Adult) Siemens 900 C/D/E NEO (Neonatal)
Device Driver P/N:	M1032-B03rl (Adult) M1032-B07rl (Neonatal) where: r = revision l = language
Supported Devices:	Siemens 900 C/D/E
Connection:	Analog: 2 waves, 5 numerics, 1 alarm.
Switch Settings:	Factory default

Note

Please ensure that the infant / adult switch on the Siemens Ventilator is in the correct position.



Siemens 900 C/D/E Ventilator Rear Panel


Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Maquet Servo 300/300A Ventilator

Device Driver Name:	Maquet Servo 300/300A
Device Driver P/N:	M1032-B20rl
	where: r = revision l = language
Supported Devices:	Maquet Servo 300/300A Software revision 2.01
Connection:	Port N82 on rear side
	Baud Rate: 9600 ¹
	Word length: 8 bits
	Stop Bits: 1
	Parity: Even.
	Protocol: SV300CI EXTENDED

1. Must be set and verified according to the Maquet Servo 300/300A User Manual



Maquet Servo 300 Ventilator Rear Panel

Note

The analog cable previously connecting the VueLink module to port N81 is no longer supported.


All waves and numerics are now transferred to the VueLink module using the digital cable connected to port N82.

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Maquet SERVO-i Ventilator

Device Driver Name:	Maquet SERVO-i
Device Driver P/N:	M1032-B55rl
	where: r = revision l = language
Supported Devices:	SERVO-i System version 2.0, 3.0
Connection:	Computer Interface Emulator
	Baud Rate: 9600 Word length: 8 bits Start Bits: 1 Stop Bits: 1 Parity: Even




Maquet SERVO-i Ventilator Right Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Maquet SERVO-s Ventilator

Device Driver Name:	Maquet SERVO-s
Device Driver P/N:	M1032-B65rl
	where: r = revision l = language
Supported Devices:	SERVO-s System version 3.0
Connection:	Computer Interface Emulator
	Baud Rate: 9600
	Word length: 8 bits
	Start Bits: 1
	Stop Bits: 1
	Parity: Even



Maquet SERVO-s Ventilator Back Panel


Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

A4 Gas Analyzers

HP M1025A/B Gas Analyzer

Device Driver Name:	Philips M1025A/B
Device Driver P/N:	M1032-B06rl where: r = revision l = language
Supported Devices:	Philips M1025A (Protocol Rev. A, serial Number 31xx and lower) and Philips M1025B.
Connection:	RS-232 Digital Baud Rate: 9600 Parity: Even Stop Bit: 1 Data Bit: 8 Handshake: Hardwired Hardwire Mode: Switched Lines Text Line: Ignore Monitor Message: Ignore
Switch Settings:	Use Philips M1025A/B screen configurations.



Philips M1025A/B Gas Analyser Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Datex Capnomac (II and Ultima) Gas Analyzer

Start Up Procedure

The driver automatically identifies the connected device at start-up. For this reason, the following start-up procedure must be followed whenever the user changes the Datex Device, even if it is the same model:

1. Turn off the Datex Device.
2. Plug in and Set up the VueLink Module and connect the cable.
3. Turn on the Datex Device.

Note

To ensure proper operation of the VueLink module and accurate data, the start up procedure described above must be carried out:-

- If the Datex device is turned off.
- If the VueLink module or cable is unplugged.
- Whenever the Datex device is selected.
- After re-entering monitoring mode.

Warning

- *In order to avoid damaging the Datex device, the VueLink cable should only be connected to, or disconnected from, the Datex device when it is switched off.*
-
-

Device Driver Name: DATEX
Device Driver P/N: M1032-B13rl
where:
r = revision

l = language

Supported Devices:

- Datex Capnomac II
- Datex Capnomac II four channel analog output option
- Datex Ultima

All units shipped after Feb 21st 1994

Connection:

Analog / RS-232 Digital combination

Baud Rate: 1200

Word Length: 8 bits

Stop Bits: 1

Parity: None

Analog: 1 wave on Capnomac II

3 waves on Capnomac II

four channel analog output option


4 waves on Capnomac

Ultima

Switch Settings:

Factory default

**Digital and
Analog Connector**




Datex Capnomac (II and Ultima) Gas Analyzer rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Dräger Vamos Variable Anaesthetic Gas Monitor

Device Driver Name:	Dräger Vamos
Device Driver P/N:	M1032-B54 <i>r/l</i>
	where: <i>r</i> = revision <i>l</i> = language
Supported Devices:	Dräger Vamos Software 2.0 MEDIBUS Version 4.03
Connection:	MEDIBUS RS-232 port on rear side
	Baud Rate: 9600 ¹
	Data Bits: 8
	Start Bits: 1
	Stop Bits: 1
	Parity: even



Dräger Vamos Rear Panel

1. Must be set and verified according to the Dräger Vamos User's Manual.

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Ohmeda RGM 5250 Respiratory Gas Monitor

Device Driver Name: Ohmeda RGM 5250 Respiratory Gas Monitor

Device Driver P/N: M1032-B22Bl

where:

r = revision

l = language

Supported Devices: Ohmeda RGM 5250 with SW Revision 6.0

Supported Device Options:

- Anesthetic Agent

- SpO₂

Connection: Analog / RS-232 Digital

Baud Rate: 1200

Word Length: 7 bits


Stop Bits: 1

Parity: Odd

Handshake: None

Hardwire mode: Uni-directional

Print Period: 10 sec



Ohmeda RGM 5250 Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Ohmeda Rascal II Anesthetic Gas Monitor

Device Driver Name: Ohmeda Rascal II

Device Driver P/N: M1032-B44rl

where:

r = revision

l = language

Supported Devices: Ohmeda Rascal II Anesthetic Gas Monitor with SW Revision 1.12


Connection: RS-232 Digital and Analog

Baud Rate: 9600 (fixed)

Word Length: 8 bits (fixed)

Stop Bits: 1 (fixed)

Parity: None (fixed)



Ohmeda Rascal II Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.


A5

Anesthesia Machines

Dräger Cato Anesthesia Device

Device Driver Name:	Dräger Cato Anesthesia Device
Device Driver P/N:	M1032-B25rl where: r = revision l = language
Supported Devices:	Dräger Cato - Screen Type PM 8050 CD with <ul style="list-style-type: none">• SW Revision 2.01, 2.02• Digital interface: MEDIBUS protocol 3.00
Supported Device Options:	None
Connection:	RS-232 Digital Baud Rate: 9600 ¹ Word Length: 8 bits (fixed) Stop Bits: 1 (fixed) Parity: Even (fixed) Handshake: Software (fixed) Protocol: MEDIBUS ¹

1. Must be set and verified according to the Dräger Cato User's Manual.



Dräger Cato Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Dräger Cicero Anesthesia Machine

Device Driver Name: Dräger Cicero Anesthesia Machine

Device Driver P/N: M1032-B12rl

where:

r = revision

l = language

Supported Devices: PM8020 datamanager
Software revision 4.6 only

Connection: RS-232 Digital


Baud Rate: 9600 (fixed)

Word Length: 8 bits fixed)

Stop Bits: 1 (fixed)

Parity: Even (fixed)

Switch Settings: A jumper must be set inside the Cicero in order to get wave output. This must be done by your local Dräger Representative.



RS232 Connector

Dräger Cicero Anesthesia Machine Rear Panel

Note

VueLink supports Isoflurane, Enflurane and Halothane gases only. If the Cicero is configured with any other unsupported gas, Philips cannot accept any responsibility for information displayed on the Philips patient monitor via VueLink.


Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Dräger Cicero EM Anesthesia Machine

Device Driver Name:	Dräger Cicero EM Anesthesia Machine
Device Driver P/N:	M1032-B36rl
	where:
	r = revision
	l = language
Supported Devices:	PM 8060 datamanager Software revision 2.01 only.
Connection:	RS-232 Digital
	Baud Rate: 9600 ¹
	Word Length: 8 bits (fixed)
	Stop Bits: 1 (fixed)
	Parity: Even (fixed)
	Handshake: Software (fixed)
	Protocol: MEDIBUS ¹

1. Must be set and verified according to the Dräger Cicero EM User's Manual.



Dräger Cicero EM Rear Panel

Note To access the connectors, open the Cicero EM housing.


Note VueLink supports Isoflurane, Enflurane, Halothane, Sevoflurane, and Desflurane gases only. If the Cicero EM is configured with any other unsupported gas, Philips cannot accept any responsibility for information displayed on the Philips patient monitor via VueLink.

Note Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Dräger Fabius GS/Tiro Anesthesia Machine

Device Driver Name:	Dräger Fabius GS/Tiro
Device Driver P/N:	M1032-B53rl
	where: r = revision l = language
Supported Devices:	Dräger Fabius GS/Tiro Software 3.01
Connection:	RS-232 port COM1
	Baud Rate: 9600 ¹
	Word Length: 8 ¹
	Start Bits: 1
	Stop Bits: 1 ¹
	Parity: Odd ¹
	Protocol: Medibus ¹

1. Must be set and verified according to the Dräger Fabius GS/Tiro Technical Service Manual.



Dräger Fabius GS/Tiro Rear Panel showing the location of COM1


Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Dräger Julian Anesthesia Workstation

Device Driver Name:	Dräger Julian
Device Driver P/N:	M1032-B42rl
	where:
	r = revision
	l = language
Supported Devices:	PM 8055 datamanager
	<ul style="list-style-type: none"> • Software Revision 3.02. • Digital Interface: MEDIBUS Protocol 3.00.
Connection:	RS-232 Digital
	Baud Rate: 9600 ¹
	Word Length: 8 bits (fixed)
	Stop Bits: 1 (fixed)
	Parity: Even (fixed)
	Handshake: Software (fixed)
	Protocol: MEDIBUS ¹

1. Must be set and verified according to the Dräger Julian Anesthetic Workstation User's Manual.



Dräger Julian Rear Panel

Note

VueLink supports Isoflurane, Enflurane, Halothane, Sevoflurane and Desflurane gases only. If the Dräger Julian is configured with any other unsupported gas, Philips cannot accept any responsibility for information displayed on the Philips patient monitor via VueLink.

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Dräger Pallas Anesthesia Workstation

Device Driver Name: Dräger Pallas

Device Driver P/N: M1032-B62*r/l*

where:
r = revision
l = language

Supported Devices: Dräger Pallas Software 3.01
 MEDIBUS Version 4.03

Connection: RS-232 port COM1 on rear side

Baud Rate: 9600¹
 Data Bits: 8
 Start Bits: 1
 Stop Bits: 1
 Parity: even
 Protocol: MEDIBUS

COM 1
(Medibus)



Dräger Pallas Rear Panel

Note


Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

1. Must be set and verified according to the Dräger Pallas User's Manual.

Dräger PM 8050 Anesthesia Machine

Device Driver Name:	Dräger PM 8050
Device Driver P/N:	M1032-B32rl where: r = revision l = language
Supported Devices:	Dräger PM 8050 with: <ul style="list-style-type: none">• Software Revision 2.02, 2.03.• Digital Interface: MEDIBUS Protocol 3.00.
Connection:	RS-232 Digital Baud Rate: 9600 ¹ Word Length: 8 bits (fixed) Stop Bits: 1 (fixed) Parity: Even (fixed) Handshake: Software (fixed) Protocol: MEDIBUS ¹

1. Must be set and verified according to the Dräger PM 8050 User's Manual.



Dräger PM 8050 Rear Panel

Note

VueLink supports Isoflurane, Enflurane, Halothane, Sevoflurane and Desflurane gases only. If the Dräger PM 8050 is configured with any other unsupported gas, Philips cannot accept any responsibility for information displayed on the Philips patient monitor via VueLink.


Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Dräger Primus/Apollo Anesthesia Workstation

Device Driver Name:	Dräger Primus/Apollo
Device Driver P/N:	M1032-B52 <i>r/l</i>
	where: <i>r</i> = revision <i>l</i> = language
Supported Devices:	Dräger Primus Software 2.02 MEDIBUS Version 4.03 Dräger Apollo Software 3.0 MEDIBUS Version 4.03
Connection:	RS-232 port COM1 on rear side
	Baud Rate: 9600 ¹
	Data Bits: 8
	Start Bits: 1
	Stop Bits: 1
	Parity: even
	Protocol: MEDIBUS

1. Must be set and verified according to the Dräger Primus/Apollo User's Manual.



Dräger Primus/Apollo Rear Panel


Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Dräger Zeus Anesthesia Device

Device Driver Name:	Dräger Zeus
Device Driver P/N:	M1032-B58rl where: r = revision l = language
Supported Devices:	Dräger Zeus Software 3.n <ul style="list-style-type: none">• MEDIBUS Version 4.03
Connection:	RS-232 port COM1 or COM2 on rear side Baud Rate: 9600 ¹ Word Length: 8 Stop Bits: 1 Parity: Even Protocol: MEDIBUS

1. Must be set and verified according to the Dräger Zeus User's Manual.




Dräger Zeus Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for Dräger Zeus for exact details of waves, numerics, INOPs and alarms available from the external device via Philips patient monitoring network.

GE Aestiva and Avance

Device Driver Name:	GE Avance and Aestiva
Device Driver P/N:	M1032-B39rl where: r = revision l = language
Supported Devices:	GE Avance Software 3.0 Ohmeda Com 1.2 Serial Protocol GE Aestiva Software 4.5 Ohmeda Com 1.0 Serial Protocol
Connection:	RS-232 Digital Baud Rate: 19200 Word Length: 7 bits Stop Bits: 1 Parity: Odd



RS-232

Rear of GE Aestiva



Rear of GE Avance

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use Information for GE Avance and Aestiva for exact details of waves, numerics, INOPs and alarms available from the external device via Philips patient monitoring network.

GE Aisys

Device Driver Name:	GE Aisys
Device Driver P/N:	M1032-B60rl
	where: r = revision l = language
Supported Devices:	GE Aisys Software 1.x Ohmeda Com 1.2 Serial Protocol
Connection:	RS-232 port on rear
	Baud Rate: 19200
	Word Length: 7 bits
	Start Bits: 1
	Stop Bits: 1
	Parity: Odd

RS-232



Rear of GE Aisys

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device via Philips patient monitoring network.


North American Dräger Narkomed 2B and Narkomed 2C Anesthesia Systems

Device Driver Name:	NARKOMED 2B/2C/3200
Device Driver P/N:	M1032-B31rl where: r = revision l = language
Supported Devices:	Narkomed 2B Anesthesia System Software Revision 2.05 Narkomed 2C Anesthesia System Software Revision 1.11 Vitalert 3200 Monitor Software Revision 1.08 NAD Communications HubSoftware Revision 1.07
Connection:	RS-232 Digital Protocol: Vitalink Word Length: 8 bits Stop Bits: 1 Parity: None Baud Rate: Narkomed 2B 1200 Narkomed 2C 19200 Vitalert 3200 19200 Communications Hub19200


Note

- Connect VueLink to Port A only (Narkomed 2B, 2C and Vitalert 3200). Otherwise, incorrect data may be displayed.


- Waves are only available directly from the device (Narkomed 2C or Vitalert 3200); therefore, connect VueLink to the device whose waves are desired. (The Communications Hub does not provide waves.)
 - If the patient monitor RS-232 MECIF interface is being used, then configure MECIF only for the highest (in letter, e.g., C is higher than B) used port on the Communications Hub or Vitalert 3200, or only to Port B on the Narkomed 2C. Otherwise, incorrect data may be displayed.
 - If all of the following conditions are true:
 - The patient monitor RS-232 MECIF interface is connected to the Narkomed 2C or Communications Hub.
 - Another Type B VueLink module (Ventilator, Gas Analyzer, Anesthesia Machine or Auxiliary Plus) is being used.
 - The above VueLink module is measuring one or more numerics which are also output by the Narkomed VueLink interface. (For example, the Philips M1025A/B Anesthetic Gas Analyzer numerics include !imCO2, !etCO2, !inN2O, etN2O, !FIO2, !inAGT, !etAGT, !PULSE, !SpO2 and !AWRR. These correspond to the Narkomed VueLink interface numerics: IMCO2, ETCO2, inN2O, etN2O, FIO2, inAGT, etAGT, PULSE, SpO2 and AWRR.)
 - Then the following must be performed to avoid the display of incorrect data: Each duplicate numeric (from condition (c)) should NOT be selected on the Narkomed VueLink interface. For example, if the Philips M1025A/B Anesthetic Gas Analyzer VueLink interface is also being used, then the IMCO2, ETCO2, inN2O, etN2O, FIO2, inAGT, etAGT, PULSE, SpO2 and AWRR numerics must not be selected on the Narkomed VueLink interface.
-




Narkomed 2B: Rear View



Narkomed 2C: Rear View



Vitalert 3200: Rear View



Communications Hub: Rear View

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.


North American Dräger Narkomed 4 Anesthesia System

Device Driver Name:	NAD NARKOMED 4
Device Driver P/N:	M1032-B38rl
	where: r = revision l = language
Supported Devices:	Narkomed 4 Anesthesia System Software Revision 1.30
Connection:	RS-232 Digital
	Protocol: Vitalink
	Word Length: 8 bits
	Stop Bits: 1
	Parity: None
	Baud Rate: 19200

Note

- Connect VueLink to Port A only on Narkomed 4. Otherwise, incorrect data may be displayed.
- If the patient monitor RS-232 MECIF interface is being used, then configure MECIF only for the highest (in letter, e.g., C is higher than B) used port on the Narkomed 4. Otherwise, incorrect data may be displayed.
- If all of the following conditions are true:
 - The patient monitor RS-232 MECIF interface is connected to the Narkomed 4.
 - Another Type B VueLink module (Ventilator, Gas Analyzer, Anesthesia Machine or Auxiliary Plus) is being used.
 - The above VueLink module is measuring one or more numerics which are also output by the Narkomed VueLink interface. (For example, the Philips M1025A/B Anesthetic Gas Analyzer numerics include !imCO2, !etCO2, !inN2O, etN2O, !FIO2, !inAGT, !etAGT, !PULSE, !SpO2 and !AWRR. These correspond to the Narkomed VueLink interface numerics: IMCO2, ETCO2, inN2O, etN2O, FIO2, inAGT, etAGT, PULSE, SpO2 and AWRR.)
- Then the following must be performed to avoid the display of incorrect data: Each duplicate numeric (from the third condition) should NOT be selected on

the Narkomed VueLink interface. For example, if the Philips M1025A/B Anesthetic Gas Analyzer VueLink interface is also being used, then the IMCO₂, ETCO₂, inN₂O, etN₂O, FIO₂, inAGT, etAGT, PULSE, SpO₂ and AWRR numerics must not be selected on the Narkomed VueLink interface.



Narkomed 4: Rear View

Note


Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

North American Dräger Narkomed GS Anesthesia System

Device Driver Name:	NAD NARKOMED GS
Device Driver P/N:	M1032-B45rl
	where: r = revision l = language
Supported Devices:	Narkomed GS Anesthesia System Software Revision 1.10
Connection:	RS-232 Digital
	Protocol: Vitalink
	Word Length: 8 bits
	Stop Bits: 1
	Parity: None
	Baud Rate: 19200

Note

- Connect VueLink to Port A only on Narkomed GS. Otherwise, incorrect data may be displayed.
 - If the patient monitor RS-232 MECIF interface is being used, then connect VueLink to Port B only on Narkomed GS.
-



Narkomed GS: Rear View

Note


Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

North American Dräger Narkomed 6000 Anesthesia System

Device Driver Name:	NAD NARKOMED 6000
Device Driver P/N:	M1032-B49rl
	where: r = revision l = language
Supported Devices:	North American Dräger Narkomed 6000 Anesthesia System Software Revision 1.06
Connection:	RS-232 Digital (use one of the ports: com5 - com8)
	Protocol: Vitalink (fixed)
	Word Length: 8 bits
	Start Bits: 1
	Stop Bits: 1
	Parity: None
	Baud Rate: 19200

The word length, start and stop bits, parity and baud rate must be set and verified according to the NAD Narkomed 6000 Technical Service Manual

1. Connect VueLink to the one of the RS-232 ports COM5 to COM8 on Narkomed 6000.



NAD Narkomed 6000: Rear View

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Ohmeda 7800/7810

Device Driver Name: OHMEDA 7800/7810

Device Driver P/N: M1032-B37rl

where:


r = revision

l = language

Supported Devices: Ohmeda 7800 in Excel Anesthesia Machine
Software Revision 4.22


Ohmeda 7810 in Modulus II Anesthesia Machine
Software Revision 4.22

Connection: Baud Rate: 1200
Word Length: 7 bits
Stop Bits: 1
Parity: Odd



RS232 and Analog Connector

Rear Panel of Ohmeda 7800



Rear Panel of Ohmeda 7810

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Ohmeda 7900

Device Driver Name: OHMEDA 7900

Device Driver P/N: M1032-B39rl

where:

r = revision

l = language

Supported Devices: Software Revision 3.5

Connection: RS-232 Digital

Baud Rate: 19200

Word Length: 7 bits

Stop Bits: 1

Parity: Odd

RS-232



Underside of Ohmeda 7900

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.


Ohmeda Modulus CD Anesthesia Machine

Device Driver Name:	Ohmeda Modulus CD
Device Driver P/N:	M1032-B27rl where: r = revision l = language
Supported Devices:	Ohmeda Modulus CD (Base Type) Software versions 3.50 and 4.00 are supported.
Connection:	Connect the interface cable to the A1 serial port of the <i>Ohmeda Modulus CD</i> , which is located on the back panel. The device driver supports the following configuration for the serial communication port. Baud Rate: 1200 Word Length: 7 bits Stop Bits: 1 Parity: Odd These settings are fixed and cannot be changed.
Switch Settings:	None

Set up the Ohmeda Modulus CD as follows:

- **For software version 3.50:**
 1. Press MANUAL key twice on display pod.
 2. Select EXTERNAL PORT DEVICE and change it to ARKIVE.
 3. Change the Printer Logging Interval to 20 seconds.

- **For software version 4.00:**
 1. Press MANUAL key on display pod.
 2. Select EXTERNAL PORT DEVICE and change it to COMPUTER LOGIN.



Rear Panel of the Ohmeda Modulus CD Anesthesia Machine

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Taema Alys Anesthesia Machine

Device Driver Name: Taema Alys

Device Driver P/N: M1032-B29rl

where:

r = revision

l = language

Supported Devices: Taema Alys, all versions with software driver version V02.40 and upwards.

Connection: Connect the device driver to the serial port of the Taema Alys. The device driver supports the following RS-232 configuration for the serial communication port.

Baud Rate: 19200

Word Length: 8 bits

Stop Bits: 1


Parity: Even

Handshake: Software.

These settings are fixed and cannot be changed.

Switch Settings: None

The Taema Alys is an anesthetic device that is intended for use on adult and pediatric patients (not neonates). The Taema Alys is equipped with a standard RS 232 interface, through which data such as waves, vital parameters, machine settings, alarms, and inops can be retrieved.



Rear Panel of the Taema Alys Ventilator

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

A6


Auxiliary Plus Devices

Aspect Medical Systems A-2000™ BIS™ Monitor

Device Driver Name:	Aspect A-2000 BIS
Device Driver P/N:	M1032-B47rl
	where:
	r = revision
	l = language
Supported Devices:	Aspect Medical Systems A-2000 BIS Serial Protocol Revision 1.07
Connection:	RS-232 Digital
	Protocol: ASCII
	Word Length: 8 bits
	Start Bits: 1
	Stop Bits: 1
	Parity: None
	Baud Rate: 9600

Note

- Connect VueLink to the RS-232 port on A-2000 BIS Monitor.
 - Make sure that the serial port protocol is set to ASCII on the A-2000 BIS Monitor (see the documentation supplied with the BIS monitor for information on how to do this).
-



A-2000 BIS Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Baxter Explorer Cardiopulmonary Monitor


Device Driver Name: Baxter Explorer
Device Driver P/N: M1032-B15rl
where:
r = revision
l = language

Supported Devices: Baxter Explorer (sw rev. 06.37)

Connection: RS-232 Digital (Fixed configuration)

Baud Rate: 9600
Word Length: 8 bits
Stop Bits: 1
Parity: None.

Switch Settings: Factory default



Baxter Explorer Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

B|Braun SpaceCom

Device Driver Name:	B Braun SpaceCom
Device Driver P/N:	M1032-B61 <i>r</i> / where: <i>r</i> = revision <i>/</i> = language
Supported Devices:	B Braun SpaceCom with Protocol Version V3.25 Infusomat Space Version 686 Infusomat Space P Version 687 Perfusor Space Version 688
Connection:	RS232 port on rear side Baud Rate: 9600 ¹ Data Bits: 8 Start Bits: 1 Stop Bits: 1 Parity: None Protocol: BCC V3.25

Note The B|Braun SpaceCom VueLink driver supports 24 pumps. A maximum of 12 pumps can be placed in each of the first two columns. The pumps must be placed in the first three segments of each pillar.

Note If the VueLink Cable is disconnected during normal operation, the communication between the BBraun SpaceCom (Communication SW for the BBraun Space System)

1. Must be set and verified according to the B|Braun SpaceCom User Manual.

and the VueLink module is lost. You have to restart the BBraun SpaceCom to reestablish the communication.

Caution

Please upgrade the IntelliVue Software to Rev. C.00.90 (patch) or higher or IntelliVue Rev. D.0 to avoid any loss of VueLink data from the resting display.

Please refer to the respective IntelliVue FCO and Service Bulletin.

Explanation: If you change the VueLink labels on the main screen, labels and values get mixed up. The values appear with the wrong label. The values in the VueLink task window are correct.



The location of the communication ports of the B | Braun SpaceCom.

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device via Philips patient monitoring network.

Braun FM Device Information


Device Driver Name:	Braun fm-system
Device Driver P/N:	M1032-B33rl
	where: r = revision l = language
Supported Devices:	Braun fm system with fm computer Software Revision 3.10
Connection:	RS-232
	Baud Rate: 9600
	Word Length: 8 bits
	Stop Bits: 1
	Parity: None

Caution

Please upgrade the IntelliVue Software to **Rev. C.00.90 (patch)** or higher or IntelliVue Rev. D.0 to avoid any loss of VueLink data from the resting display.

Please refer to the respective IntelliVue FCO and Service Bulletin.

Explanation: If you change the VueLink labels on the main screen, labels and values get mixed up. The values appear with the wrong label. The values in the VueLink task window are correct.



Braun FM Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device via Philips patient monitoring network.

Diametrics Medical IRMA SL Series 2000

Device Driver Name: IRMA SL Series 2000

Device Driver P/N: M1032-B50rl

where:
r = revision
l = language


Supported Devices: Diametrics Medical IRMA SL Series 2000

Software version 5.1.x

Connection: RS-232 Digital

Baud Rate: 9600
Stop Bits: 1
Data Bits: 8
Parity: None.
Handshake: None

Switch Settings: None



Diametrics Medical IRMA SL Series 2000: Rear View

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device via Philips patient monitoring network.

Diametrics Medical Trendcare Monitor

Device Driver Name: Trendcare Monitor

Device Driver P/N: M1032-B51rl

where:

r = revision

l = language

Supported Devices: Diametrics Medical Trendcare Monitor

Software version 2.70.01

Connection: RS-232 Digital (Fixed configuration)

Baud Rate: 9600, 19200¹


Data Bits: 8

Start Bits: 1

Stop Bits: 1


Parity: None.

Switch Settings: The data types 1 to 4 (DT1, DT2, DT3, DT4) and all numeric values (PH, PCO2, PO2, ETC) must be selected for output at the external device.



1. The baud rate is set in the VueLink configuration (config mode) on the monitor.

Satellite Monitor



Senior Monitor

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device via Philips patient monitoring network.

Edwards CardiacOutput-Oximetry Monitor Family

Device Driver Name:	Edwards CardiacOutput-Oximetry Monitor Family
Device Driver P/N:	M1032-B16rl where: r = revision l = language
Supported Devices:	Edwards Vigilance Edwards Vigilance II Edwards Vigileo
Connection:	RS-232 Digital Baud Rate: 19200 ¹ Parity: None Stop Bits: 1 Data Bits: 8 Flow Control: 2 seconds.

1. Must be set and verified according to the Edwards CardiacOutput-Oximetry Monitor Family User's Manual



Fig. 1: COM1 position on Vigilance (backside)

Edwards Vigilance Rear Panel



Edwards Vigilance II Rear Panel




Edwards Vigileo Rear Panel


Connecting to devices with male connectors

Some releases of the *Edwards CardiacOutput-Oximetry Monitor Family* devices have a male COM-Port socket. The Philips cable provides a “gender-changer” to connect these devices to the *VueLink* module. The following instructions describe, how to assemble the “gender-changer” with the *VueLink* cable.


1. Unscrew the nuts on the female side of the gender-changer..



2. Screw the *VueLink* cable to the gender-changer



3. Connect the assembled *VueLink*-cable and gender-changer to COM1 of the *Edwards CardiacOutput-Oximetry Monitor Family* device.



Note


Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device via Philips patient monitoring network.

Fresenius Vial Base A/DPS/MVP

Warning

*The Fresenius Vial Base A can be configured to support either six or eight DPS or MVP modules.
If the base is configured to support six modules, and more than six modules are connected, VueLink provides only numerics, alarms, and inops from the first six modules.*

Device Driver Name:	F-V Base A/DPS/MVP
Device Driver P/N:	M1032-B43rl
	where: r = revision l = language
Supported Devices:	Fresenius-Vial Base A with Software Revision 2.2 Fresenius Vial DPS Module with Software Revision from 2.2 to 6.6 Fresenius Vial MVP Module with Software Revision 3.1
Connection:	RS-232 C
	Baud Rate: 19200
	Word Length: 7 bits (fixed)
	Stop Bits: 1 (fixed)
	Parity: Even (fixed)
	Handshake: None (fixed)



Fresenius Vial Base A Module Rear Panel

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for Fresenius Vial Base A/DPS/MVP for exact details of numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

GE DINAMAP ProCare 100 Series

Device Driver Name:	GE DINAMAP ProCare 100 Series
Device Driver P/N:	M1032-B72rl where: r = revision l = language
Supported Devices:	GE DINAMAP ProCare 100 Series
Connection:	RS-232 port on rear side Baud Rate: 9600 Data Bits: 8 Start Bits: 1 Stop Bits: 1 Parity: None



The location of the RS 232 port on rear side of the GE DINAMAP ProCare 100 Series

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of numerics and INOPs available from the external device.

i-STAT[®] 1 Analyzer

Device Driver Name: i-STAT[®] 1 Analyzer

Device Driver P/N: M1032-B59*r/l*
where:
r = revision
l = language


Supported Devices: i-STAT[®] 1 Analyzer with
i-STAT[®] 1 RIBS Downloader/Recharger

Cartridge Types:

6+	CG4+	CG8+
Chem8+	CK-MB	Crea
cTnI	E3+	EC4+
EC8+	EG6+	EG7+
G	G3+	PT/INR
Celite _{ACT}	Kaolin _{ACT}	

Connection: RJ-12 port on rear side of the i-STAT[®] 1 RIBS Downloader/Recharger

Baud Rate:	9600
Data Bits:	8
Start Bits:	1
Stop Bits:	1
Parity:	none




The location of the RJ-12 port on rear side of the i-STAT[®] 1 RIBS Downloader/Recharger.

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of numerics and INOPs available from the external device.

Mortara ELI 100/STM 12-Lead Monitor

- Device Driver Name:** Mortara ELI 100/STM
- Device Driver P/N:** M1032-B26rl
- where:
r = revision
l = language
- Supported Devices:** Mortara ELI 100/STM with:
- Software part no: 11001-001-20
 - Software revision: 1.95, 11.95
- Connection:** RS-232 Digital (fixed configuration)
- Baud Rate: 19200
Word Length: 8 bits
Stop Bits: 1
Parity: None
Handshake: Software
- Switch Settings:** None



Mortara ELI 100/STM Rear Panel

Note


Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Organon Teknika TOF-Watch SX

Note This information only applies to the special Philips version of the TOF-Watch SX.

Device Driver Name:	TOF-Watch SX
Device Driver P/N:	M1032-B46rl
	where:
	r = revision
	l = language
Supported Devices:	Organon Teknika TOF-Watch SX (Philips version) Software Revision 1.2
Connection:	RS-232 Digital
	Word Length: 8 bits
	Start Bits: 1
	Stop Bits: 1
	Parity: None
	Baud Rate: 19200

Connect VueLink to the TOF-Watch SX via the TOF-Watch SX Interface.



TOF-Watch SX: Connecting

Note

Please refer to the Philips M1032A VueLink External Device Instructions for Use for exact details of waves, numerics, INOPs and alarms available from the external device and via Philips patient monitoring network.

Free Analog Device Information

Free Analog Device Information

This appendix contains information about Free Analog Devices that can be connected to the VueLink module. The following information is provided:

- General configuration and cable information about Free Analog Devices.
- Step by step instructions on how to configure the following Free Analog Device into the VueLink module.
 - Dräger Evita Ventilator
 - Engström EAS/ELSA Anesthesia System
 - Nellcor N1000 and N2500 Multi-Function Monitor
 - Radiometer Transcutaneous pO₂/pCO₂ Monitor

Configuring the Free Analog Device

If an external device name does not appear in the Preselection Task Window on the Philips patient monitor, but has an analog interface, it can be connected to the VueLink Module using the *Free Analog Device driver*. The device information must however be manually configured into the module using the Task Windows in the System's Configuration Mode. The following must be specified for every signal produced by the device:

- signal labels
- units
- voltage range
- value range
- grid settings.

These settings are then stored permanently within the VueLink module, in a *Free Analog Device Driver*. A name can also be assigned to the device driver using the handheld keypad.


Warning

- *Voltages greater than $\pm 12V$ must not be applied across any two pins of the VueLink module front connector.*
 - *Only one external device must be connected to a single VueLink module at any one time.*
 - *Patient safety checks must always be performed after repair and installation of the module, as detailed in the patient monitor Service Manual.*
-
-

Preparing the Free Analog Cable

A 2m and 4m moulded cable with an unconnected 25-pin (male) connector is available for connecting the Free Analog Device. If the external device has a 25-pin female plug, the connector should be soldered onto the VueLink cable. If a different connector is required, this must be provided by local technical personnel or by representatives of external device manufacturers. This must then be soldered onto the unterminated VueLink cable as before.

Configuring the Free Analog Device



The Free Analog cable options are detailed in the matrix below:

Devicecol	New Module Option (MI032A-)		Patient Monitor Cable Option (MI182A-)		SMO/SLI Option (MI032-)	
	2m	4m	2m	4m	2m	4m
Free Analog Devices	K60	K6A	K60	K6A	61611	61612

Dräger Evita Ventilator

This note provides information on configuring the Evita Free Analog Device into the VueLink module. This procedure should take approximately 2 hours and consists of 6 steps:-

	Estimated Time
Step 1 Ensure that all the necessary hardware is available	10 minutes
Step 2 Prepare the free analog cable	45 minutes
Step 3 Enter a device name, device-id and channel labels	10 minutes
Step 4 Configure each wave and numeric channel	30 minutes
Step 5 Connect the free analog cable to the Evita analog output	10 minutes
Step 6 Verify configuration with the user	15 minutes
Total Time:	120 minutes

General information on configuring the free analog devices has not been included in this technical note. Please refer to Chapter 3 of the VueLink Handbook, if you are not familiar with this procedure.

Step 1: Ensure that all the necessary hardware is available

(10 Minutes)

You will need:-

- An M1032A #A02 VueLink Module
- 4m Free Analog (open end) cable M1032A #K6A or M1182A #K6A
or
2m Free Analog cable M1032A #K60 or M1182A #K60.
- A connector to fit the Evita Analog Interface connector:
Standard 37 PIN DSUB connector (male)
Standard hood with snap adapter

Note

This connector is provided locally by the hospital or by the manufacturer; it will not be set up at SMO/SLI.

- The Evita Analog Interface - This is an optional feature. If the Evita is not equipped with this interface it needs to be ordered by the customer.

Step 2: Prepare the free analog cable

(45 Minutes)

1. Connect the black and white wires that are twisted together to the Evita Analog ground (Evita pin 1). Do NOT separate these wires.
2. Connect the thick black wire, or the copper shield, to the case resp. ground on the device (Evita pin 32).



Note

If there is a conductive connection from the copper shield, via the connector, to the case resp. ground, then the black wire does not need to be connected.

3. Connect the colored wires to the analog output connector on the external device. The respective channels for each wire color are given in the following table.

Free Analog Cable		Designation	EVITA PIN
Channel	Wire Colours		
---	Black and White (twisted)	Analog ground	1
---	Black (thick)	Shield	32
Wave 1	Brown	Airway Pressure	7
Wave 2	Red	Expiratory Flow	23
Numeric 1	Orange	Minute Volume	2
Numeric 2	Yellow	Airway Resp. Rate	5
Numeric 3	Green	Tidal Volume	20
Numeric 4	Blue	Inspired O ₂ Conc.	21
Numeric 5 ¹	Purple	Compliance	15
		Resistance	16

1. Only one available at a time, depending on how the jumpers are set, as shown in the diagram below.



Pins used for the Evita free analog interface

Step 3: Enter a device name, device-id and channel labels

(10 Minutes)

Only the Expiratory Flow Label (EXPFL) needs to be entered with the handheld keypad. All the other Evita waves and numerics are offered in the configuration Task Window of the Ventilator Free Analog Device Driver. The Evita Analog Interface provides the following waves and numerics:

Name	Recommended Label
Waves	
Airway Pressure Paw	AWP
Expiratory Flow	!EXPFL
Numerics	
Minute Volume V	!MIVOL
Ventilation frequency f^1	!AWRR

Name	Recommended Label
Expired Tidal Volume V_T	!TV
Inspired O ₂ Concentration	!FIO2

1. Referred to as Airway Respiration Rate in patient monitor

Alarms

No alarms are activated for the Evita Free Analog Device Driver.

Step 4: Configure each wave and numeric channels

(30 Minutes)


The configurations for each wave and numeric channel are provided in the following table.

Channel	Label	Unit	Value Range = Voltage Range	Decimal ¹ Represent.
Wave 1	AWP	mbar	-100 to 0 mbar-10 to 0 V 0 to 100 mbar0 to 10 V	x
Wave 2	!expF ²	l/min	0 to 150 l/min0 to 10 V	x
Numeric 1	!MIVOL	l/min	0 to 35 l/min0 to 10 V	x.x
Numeric 2	!AWRR	l/min	0 to 60 l/min0 to 10 V	x
Numeric 3	!TV	l	0 to 2 l0 to 10 V	x.xx
Numeric 4	!FIO ₂	none	0.00 to 1.000 to 10 V	x.xx
Numeric 5 ³	!COMP	ml/mbar	0 to 255 ml/mbar0 to 10V	x
	!Rstat	mbar/L/s	0 to 100 mbar/L/s0 to 10V	x

1. You can configure the number of decimal places by moving the decimal point on the input value range, using the left/right arrow keys. The decimal representation in the table above refers to the Value Range. Please ensure that they are configured as shown.
2. Evita only provides positive voltages for expFl (Expiratory Flow).
3. Only 1 available at a time, depending on configuration.

Step 5: Connect the free analog cable to the Evita analog output

(10 Minutes)



Step 6: Verify configuration with the user

(15 Minutes)

Ensure that the correct waves and numerics are being interfaced from the Evita.

Engström ELSA/EAS Anesthesia System

(This device can be connected in Germany - SUV approval is available). This note provides information on configuring the ELSA/EAS Free Analog device into the VueLink module. This procedure should take approximately 130 minutes and consists of 7 steps:-

	Estimated Time
Step 1 Ensure that all the necessary hardware is available	10 minutes
Step 2 Selection of two waveforms	10 minutes
Step 3 Prepare the free analog cable	45 minutes
Step 4 Enter a device name, device-id and channel labels	10 minutes
Step 5 Configure each wave and numeric channel	30 minutes
Step 6 Connect the free analog cable to the EAS analog output	10 minutes
Step 7 Verify configuration with the user	15 minutes
Total Time:	130 minutes

General information on configuring the free analog devices has not been included in this technical note. Please refer to Chapter 3 of the VueLink Handbook, if you are not familiar with this procedure.

Step 1: Ensure that all the necessary hardware is available

(10 Minutes)

You will need:-

- An M1032A #A04 VueLink Module
 - 4m Free Analog (open end) cable M1032A #K6A or M1182A #K6A
- or
- 2m Free Analog cable M1032A #K60 or M1182A #K60.
- A connector to fit the ELSA/EAS Analog Interface connector:
ITT Cannon 25-pole, female
Standard hood with snap adapter

Note

This connector is provided locally by the hospital or by the manufacturer; it will not be set up at SMO/SLI. (The connector provided with the free analog cable cannot be used).

Step 2: Select two waveforms

(10 Minutes)

The following parameters are provided as a waveform at the Analog Output port of ELSA/EAS:

- Expiration Flow (ELSA/EAS pin 2)
- Airway Pressure (ELSA/EAS pin 15)
- CO₂ Concentration (ELSA/EAS pin 5)
- Inspiration Flow (ELSA/EAS pin 14)
- Anesthetic Agent (Elsa pin 16)

Let the customer decide which two of the provided waveforms should be interfaced to VueLink. The following configuration is described using a typical selection of two waveforms. These are:

- CO₂ Concentration
- Airway Pressure

Step 3: Prepare the free analog cable

(45 Minutes)

1. Connect the black and white wires that are twisted together to the ELSA/EAS Analog ground (pin 18). Do NOT separate these wires.
2. Connect the thick black wire, or the copper shield, to the case resp. ground on the device (ELSA/EAS pin 1).


Note

If there is a conductive connection from the copper shield, via the connector, to the case resp. ground, then the black wire does not need to be connected.

3. Connect the colored wires to the analog output connector on the external device. The respective channels for each wire color are given in the following table.

VueLink Free Analog Cable		Designation	ELSA/ EAS PIN
Channel	Wire Colours		
---	Black and White (twisted)	Analog ground	18
---	Black (thick)	Shield	1
Wave 1	Brown	CO ₂ Concentration ¹	5
Wave 2	Red	Airway Pressure	15
Numeric 1	Orange	O ₂ Concentration	3
Numeric 2	Yellow	Expired Minute Volume	4
Numeric 3	Green	Total Fresh Gas Flow	17
Numeric 4	Blue	-	-

1. This wave is optional for ELSA. Please select another waveform, if CO₂ Concentration is not available on your device.



Pins used for the ELSA/EAS free analog interface

Step 4: Enter a device name, device-id and channel labels

(10 Minutes)

Only the AWP and the MIVOL labels are offered in the configuration window of the Ventilator Free Analog Driver. All the other ELSA/EAS parameter labels need to be entered with the handheld keypad. The following parameter abbreviations are recommended as labels for these, but the customer can decide on his/her own abbreviations.

The ELSA/EAS Analog Interface provides the following waves and numerics (two of the waveforms were selected when the free analog cable was prepared in Step 3):

Name	Label
Waves	
CO ₂ Concentration ¹	!CO2CC
Airway Pressure	AWP
Expiration Flow	!EXPFL
Inspiration Flow	!INSFL
Anesthetic Agent	!AGT 2
Numerics	
O ₂ Concentration	!O2CC
Exp. Minute Volume	!MIVOL
Total Fresh Gas Flow	!FRESH

1. Optional for ELSA.

Alarms

No alarms are activated for the ELSA/EAS Free Analog Device Driver.

Step 5: Configure each wave and numeric channels²

(30 Minutes)


The configurations for each wave and numeric channel are provided in the following table.

Channel	Label	Unit	Value Range = Voltage Range	Decimal¹ Represent.
Wave 1	CO2CC	mmHg	0 to 76 mmHg = 0 to 1 V	x
Wave 2	AWP	cmH ₂ O	0 to 100 cmH ₂ O = 0 to 1 V	x
Numeric 1	O2CC	%	0 to 100% = 0 to 1 V	x
Numeric 2	MIVOL	l/min	0 to 50 l/min = 0 to 1 V	x.x
Numeric 3	FRESH	l/min	0 to 50 l/min = 0 to 1 V	x.x
Numeric 4	-	-	-	-

1. You can configure the number of decimal places by moving the decimal point, using the left/right arrow keys. The decimal representation in the table above refers to the Value Range. Please ensure that they are configured as shown.

Step 6: Connect the free analog cable to the ELSA/EAS analog output

(10 Minutes)



Step 7: Verify configuration with the user

(15 Minutes)

Ensure that the correct waves and numerics are being interfaced from the ELSA/EAS.

Nellcor NI000 and N2500 Multi-Function Monitor

This note provides information on configuring the Nellcor Free Analog Device into the VueLink module. This procedure should take approximately 3 hours and consists of 7 steps:-

	Estimated Time
Step 1 Ensure that all the necessary hardware is available	10 minutes
Step 2 Prepare the free analog cable	75 minutes
Step 3 Enter a device name, device-id and channel labels	10 minutes
Step 4 Configure each wave and numeric channel	30 minutes
Step 5 Configure each Nellcor analog output	30 minutes
Step 6 Connect the free analog cable to the Nellcor analog output ports	10 minutes
Step 7 Verify configuration with the user	15 minutes
Total Time:	180 minutes

Note

- Nellcor N1000/N2500 is a Class I, Type CF instrument and requires a SUV approval in Germany.
- Differences in the displayed numerics can occur between the Nellcor and the Philips patient monitor, if exceptionally unstable monitoring occurs (such as crying babies). This is due to the different numeric update rates between both devices.

General information on configuring the free analog devices has not been included in this technical note. Please refer to Chapter 3 of the VueLink Handbook, if you are not familiar with this procedure.

Step 1: Ensure that all the necessary hardware is available

(10 Minutes)

You will need:-

- An M1032A #A05 VueLink Module
- 4m Free Analog (open end) cable M1032A #K6A or M1182A #K6A
- or
- 2m Free Analog cable M1032A #K60 or M1182A #K60.
- Seven plugs to fit the Nellcor analog output ports:
7 x 3/32 inch subminiature phone plugs (single tip-mono)
- Seven cables to connect the plugs to the open end of the Free Analog cable:
7 x shielded cable (length: 10 - 15cm) (diameter: 2.5mm, RG145 or similar)

Note

These cables and plugs are provided locally by the hospital or by the manufacturer; they will not be set up at SMO/SLI.


- The connector hood, the divider (5 holes) and five sleeves to connect the seven cables to the Free Analog cable are included with the free analog cable.

Step 2: Prepare the free analog cable

(75 Minutes)

1. Connect each of the seven short cables (cable 1 to 7) to a phone plug.
2. Put all of the included sleeves through the holes of the divider (A₁ to A₅).

3. Put Cable 1 through the sleeve in A₁ and crimp the sleeve on the inside of the connector hood, so that the cable cannot be torn out. Repeat this procedure for all the other cables. (As there are only 5 holes, two cables must be fitted through single holes in two cases).
4. Connect the black and white wires that are twisted together to the shield wires of all the cables. Do NOT separate these wires.
5. Connect the brown wire of the Free Analog cable to the signal wire of Cable 1 for Numeric 1. Continue for all cables as shown in the table opposite.
6. Label each cable (1 to 7)



VueLink Free Analog Cable		Designation	Cable
Channel	Wire Colours		
---	Black and White (twisted)	Analog ground	1-7
---	Black (thick)	Shield	Black & White
Wave 1	Brown	Pleth	1
Wave 2	Red	CO ₂	2
Numeric 1	Orange	SpO ₂	3

VueLink Free Analog Cable		Designation	Cable
Channel	Wire Colours		
Numeric 2	Yellow	EtCO ₂	4
Numeric 3	Green	Resp	5
Numeric 4	Blue	EtN ₂ O	6
Numeric 5	Purple	Pulse	7

Step 3: Enter a device name, device-id and channel labels

(10 Minutes)

Only the Pleth, Resp and Pulse labels need to be entered with the handheld keypad. All the other Nellcor waves and numerics are offered in the configuration window of the Nellcor Free Analog Device Driver. The Nellcor Analog Interface provides the following waves and numerics:

Name	Label
Waves	
Plethysmogram	!PLETH
Carbon Dioxide	CO ₂
Numerics	
Oxygen Saturation	!SpO ₂
End-Tidal CO ₂	!EtCO ₂
Respiration Rate	!RESP
End-Tidal N ₂ O	!EtN ₂ O
Pulse Rate	PULSE
Inspired CO ₂	!InCO ₂
Inspired N ₂ O	!InN ₂ O
Pulse Amplitude	-

Alarms

No alarms are activated for the Nellcor Free Analog Device Driver.

Step 4: Configure each wave and numeric channel

(30 Minutes)

The configurations for each numeric are provided in the following table.

Channel	Label	Unit	Value Range = Voltage Range	Decimal Represent.
Wave 1	Pleth	-	0 to 100 = 0 to 10 V	x
Wave 2	CO2	mmHg	0 to 76 mmHg = 0 to 10 V	x
Numeric 1	SpO2	%	0 to 100% = 0 to 10 V	x
Numeric 2	EtCO2	mmHg	0 to 76 mmHg = 0 to 10 V	x
Numeric 3	Resp	rpm	0 to 75 rpm = 0 to 10 V	x
Numeric 4	EtN2O	mmHg	0 to 760 mmHg = 0 to 10 V	x
Numeric 5	Pulse	rpm	0 to 250 bpm = 0 to 10 V	

Note

The low and high grid values have to meet the restrictions outlined in the VueLink Handbook.

Step 5: Configure each Nellcor analog output


(30 Minutes)

The following settings should be stored in slot 10 (defaults) of the Nellcor device:

Analog Output Port	Signal	Full scale output voltage
1	Pulse wave	10.00V
2	CO ₂ wave	10.00V = 76 mmHg
3	SaO ₂ %	10.00V = 100%
4	End-Tidal CO ₂	10.00V = 76 mmHg
5	Resp. Rate	10.00V = 75/min
6	End-Tidal N ₂ O	10.00V = 760 mmHg
7	Pulse Rate	10.00V = 250/min

Step 6: Connect the free analog cable to the Nellcor analog output

(10 Minutes)



Step 7: Verify configuration with the user

(15 Minutes)

Ensure that the correct waves and numerics are being interfaced from the Nellcor monitor.

Radiometer Transcutaneous pO₂/pCO₂ Monitor

This note provides information on configuring the TCM3 Free Analog device into the VueLink module. This procedure should take approximately 140 minutes and consists of 8 steps:-

	Estimated Time
Step 1 Ensure that all the necessary hardware is available	10 minutes
Step 2 Select two numerics	10 minutes
Step 3 Prepare the free analog cable	45 minutes
Step 4 Enter a device name, device-id and channel labels	10 minutes
Step 5 configure the alarm channel	10 minutes
Step 6 Configure and numeric channel	30 minutes
Step 7 Connect the free analog cable to the TCM3 analog output	10 minutes
Step 8 Verify configuration with the user	15 minutes
Total Time:	140 minutes

General information on configuring the free analog devices has not been included in this technical note. Please refer to Chapter 3 of the VueLink Handbook, if you are not familiar with this procedure.

Step 1: Ensure that all the necessary hardware is available

(10 Minutes)

You will need:-

- An M1032A #A01 VueLink Module
 - 4m Free Analog (open end) cable M1032A #K6A or M1182A #K6A
- or
- 2m Free Analog cable M1032A #K60 or M1182A #K60.
- A connector to fit the TCM3 Analog Interface connector:
Connector D-sub 15 pole, female.
Standard hood with M3 securing screws.

Note

This connector is provided locally by the hospital or by the manufacturer; it will not be set up at SMO/SLI. (The connector provided with the free analog cable cannot be used).

Step 2: Select two Numerics

(10 Minutes)

The following parameters are provided at the Analog Output port of TCM3:

- Alarm (TCM3 pin 5) *Transducer Temperature (TCM3 pin 8)
- Transcutaneous partial oxygen pressure (TCM3 pin 9)
- Transcutaneous partial CO₂ pressure (TCM3 pin 10)

In the following, the configuration is described using a typical selection of the two numerics and the alarm:

- Transcutaneous partial oxygen pressure
- Transcutaneous partial CO₂ pressure


Step 3: Prepare the free analog cable

(45 Minutes)

1. Connect the black and white wires that are twisted together and the thick black wire to the TCM3 Analog ground (pin 15). Do NOT separate these wires. Wire pin 13 and pin 15 together.
2. Connect the colored wires to the analog output connector on the external device. The respective channels for each wire color are given in the following table.

VueLink Free Analog Cable		Designation	TCM3 PIN
Channel	Wire Colours		
---	Black and White (twisted)	Analog ground	15
---	Black (thick)	Shield	153
Bridge			13 to 15

VueLink Free Analog Cable		Designation	TCM3 PIN
Channel	Wire Colours		
Numeric 1	Red	TCPO ₂	9
Numeric 2	Orange	CPCO ₂	10
Alarm	Yellow	Alarm	5



Pins used for the TCM3 free analog interface

Step 4: Enter a device name, device-id and channel labels

(10 Minutes)

All the TCM3 parameter labels need to be entered with the handheld keypad. The following parameter abbreviations are recommended as labels for these, but the customer can decide on his/her own abbreviations.

The TCM3 Analog Interface provides the following 2 numerics:

Name	Label
Numerics	
Partial Oxygen Pressure	!TCPO ₂
Partial CO ₂ Pressure	!CPCO ₂ ¹

1. Label restricted to 5 characters.

Step 5: Configure the alarm channels

(10 Minutes)

The configurations for the alarm channel are provided below:

Input Voltage Range	0...1 V
Active	Highlablists

Step 6: Configure each numeric channels

(30 Minutes)


The configurations for each numeric channel are provided in the following table.

Channel	Label	Unit	Value Range = Voltage Range	Decimal ¹ Represent.
Numeric 1	TCPO ₂	mmHg	0 to 200 mmHg = 0 to 1 V	x
Numeric 2	CPCO ₂	mmHg	0 to 200 mmHg = 0 to 2 V	x

1. You can configure the number of decimal places by moving the decimal point, using the left/right arrow keys. The decimal representation in the table above refers to the Value Range. Please ensure that they are configured as shown.

Step 7: Connect the free analog cable to the TCM3 analog outputs

(10 Minutes)



Step 8: Verify configuration with the users

(15 Minutes)

Ensure that the correct numerics are being interfaced from the TCM3. Simulate an alarm to ensure that the alarming capability functions correctly.

Note VueLink only provides a visual alarm, not an audible alarm.
