Klausen Process Machinery Pty. Ltd.

20 Bessemer Street
BLACKTOWN NSW 2148
AUSTRALIA

Tel: (02) 9672-3222
Fax: (02) 9672-8733

Mobile Vacuum Loading System

for

Powders & Granules

Customer	: Roche Products Pty. Ltd.
Machine Type	: KM38-310SS
Manufactured	: February, 2001
Serial No.	: M010227SS/M010228SS

Table of **Contents**

Operating Instructions	Page 3
Section 1.0 Receiver.....	Page 4
Section 2.0 Suction Wand	Page 10
Section 3.0 Angle Valves	Page 12
Section 4.0 Vacuum Pump	Page 16
Section 5.0 Lifting Column	Page 18
Section 6.0 Control System.....	Page 19
Section 7.0 Wiring Diagrams.....	Page 23
Section 8.0 Parts List	Page 29

Operating Instructions

The Klausen KM-Series mobile vacuum transfer system is designed for use with a minimal amount of installation. Proceed as follows:

- After unpacking, if vacuum transfer receiver is not mounted on the support arm on the column, the receiver should be placed on the arm and secured with the thumb screws or knobs.
- plug the connectors into the corresponding colour coded sockets on the junction box which is mounted on the column. The sockets are:
 - yellow - solenoid for filter cleaning air blasts
 - black - level sensor on receiver
 - green - solenoid for discharge butterfly valve (if applicable)
 - blue - vibrator on receiver (if applicable)
- if an external sensor is to be used to start/stop the receiver, the sensor should be mounted into the receiving hopper and the plug connected to the socket located under the main control box.
- connect the power (415 V, 3 Ph) by plugging into suitable power outlet
- connect compressed air to the connection box located on the mobile trolley
- turn main power switch on at control panel
- press start button and allow to run briefly before pressing stop button
- check direction of rotation of vacuum pump
- if direction is incorrect, cables connection must be exchanged by an electrician to reverse the direction
- connect the hose and suction wand to the inlet of the receiver
- raise receiver above filling hopper by turning column lift switch
- push the mobile trolley into position so that the receiver will discharge into the receiving hopper or container and lock the wheels into position
- press start button on control panel and check that vacuum is available at the end of the suction wand
- the system is now ready for operation

The optimal conditions for operating the vacuum transfer system will depend on the product being transferred, the method of pick-up of the product (through suction wand, dump hopper etc) and how it is discharged from the receiver.

We would recommend that the system be first operated with:

- operation with sensor
- vibrator off
- adjustment collar on suction wand closed or air inlet valve on dump hopper closed if dump hopper is used

If after operation, transfer time is short and the sensor triggers often, powder may be building up on the sensor causing it to trigger. This generally happens with fine powders on granules that can sit on top of the sensor. In this case, it would be better to operate on time (See control system instructions).

If product does not discharge freely from the receiver, the use of the vibrator may be required.

If product does not flow freely into the wand or from the dump hoppers, adjusting the collar of the suction wand or the air inlet valve of the dump hopper can increase the amount of air into the system.

1.0 Receiver

The receiver can be used for feeding processing and packaging machines using a vacuum conveying operation. Feeding is done batchwise.

The receiver is suitable for fluid powders as well as granular bulk materials.

General Description

The receiver consists of the following main parts:

- Stainless steel cover plate on which the compressed air tank is mounted. Included is a solenoid valve for cleaning the filters and also the clean air connection.
- Filter socks mounted on a filter plate to form the filter assembly.
- Main body of the receiver consisting of a cylindrical filter housing and conical receiver. The product inlet is welded to the cylindrical part of the receiver. The proximity switch is also mounted in this section.
- Butterfly discharge valve is mounted on the bottom of the receiver.
- The cover is mounted onto the body using three hook clamps.
- Compressed air connection to be installed from air up to the air pressure regulator. Hose $\text{Ø } 8 \times \text{Ø } 6$. Compressed air has to be free from oil and water, at a pressure of 4 - 5 bar over pressure.

Maintenance Instructions

- For products which are injurious to health, please consult the corresponding safety rules.
- Filters, cartridge filters or hose filters should be checked for tightness and damages:
 - once a week, and
 - if dust escapes at the clean air exhaust outlet

In order to check filter insert, remove top cover by loosening quick release clamp and lift filter insert out of the intermediate housing.

- Cleaning before disassembly of the conveying unit, switch off control, interrupt powder supply and compressed air supply. Release pressure of air tank by activating the diaphragm valve. Remove product conveying line and if applicable clean air exhaust line. By means of the quick release clamps, all separate parts can be disassembled for cleaning and inspection.

Important

The filter on the vacuum pump must be checked and cleaned on a daily basis. If the filter is dirty or clogged with material, then it is a sign that there is a leak in the filter socks or plate, or the product may be too fine for the filter. The fault must be investigated before proceeding with any further material transfer.


Functional Description

The control system should work as follows:

- On starting conveying to a receiver, the angle valve to that receiver should open and the other angle valve closed.
- If the receiver butterfly valve is closed, the vacuum pump is started
 - ♦ The vacuum is built up via the receiver and the conveying line to the product inlet wand
 - ♦ The product is sucked via the wand and the conveying line into the receiver. It collects in the receiver
 - ♦ The conveying air is cleaned in the filter and passes through to the vacuum pump where it escapes at the blow-off outlet
- Once conveying starts, a conveying timer starts and conveying continues until product level builds up to the proximity switch or the timer times out depending on the mode selected
- Upon completion of conveying, the vacuum pump stops and the butterfly valve is opened. Two timers
 - butterfly valve open timer, and
 - delay timerare started
 - ♦ Due to the weight of the material, the hopper discharges. (For non- free flowing materials an optional vibrator can be installed. When connected, the vibrator would operate whenever the valve is open)
- When the delay time expires, filter cleaning is started by allowing compressed air in the air tank to flow to the filters
 - ♦ The pulse of compressed air cleans the filters and aerates the receiver
 - ♦ The delay timer allows some material to discharge from the hopper before filter cleaning to reduce the amount of dust in the receiver, caused by the air jet
- When the butterfly valve timer expires, the valve closes and the system starts on the next cycle

Technical Data

Receiver	:	KM38-310 SS
Normal Width	:	38 mm
Capacity	:	13 Litres
Permissible Ambient Temperature	:	50°C
Filter Area	:	2 x 0.23 m ²
Filter Material	:	Polycster
Compressed Air	:	Water and oil free
Compressed Air Pressure	:	4 - 5 bar
Compressed Air Consumption	:	Approx. 5 litres / impulse

WITH COLUMN RAISED

WITH COLUMN LOWERED


TITLE: DURA-VAC DV310-38 TRANSFER SYSTEM
 OUTLINE DRAWING SERIAL No.M010227SS

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
REV	DESCRIPTION	DATE	DRN	APP'D
A	FIRST ISSUE	26.02.01	JAL	

WITH COLUMN RAISED

WITH COLUMN LOWERED

TITLE: DURA-VAC MOBILE VTS - ROCHE
 OUTLINE DRAWING SERIAL No. M010228SS
 REFERENCE DRAWING MON016-004


DATE: 26.02.01	JOB No.: MON016
DATE: 26.02.01	SCALE: 1:20
DATE: 26.02.01	SHEET: 1 OF 1
DATE: 26.02.01	REV: A

DRAWN: J LEECH
 CHECKED: B. LEE
 DWG No. MON016-002

REV	DESCRIPTION	DATE	DRN	APP'D
A	FIRST ISSUE	26/02/01	JAL	

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TITLE: DURA-VAC VACUUM RECEIVER
ROCHE - OUTLINE DIMENSION DRAWING

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
1 2 3 4 5 6 6 7 8

A B D

**Specification Sheet
for
Klausen Vacuum Transfer System Filter Cartridges**

Filter Model No. : 135 - 200 Pharma.
Main Body : Stainless Steel
End Cap : Stainless Steel
Filter Material : Spun bonded, heat pleated Axtar Polyester (G2260-IS)

Material Specifications

Weight (g/m ²)	:	260
Thickness (mm)	:	0.61
Tensile Strength (kg/15 mm)	Warp :	27.3
	Weft :	16.2
Tearing Strength (kg)	Warp :	10.8
	Weft :	10.3
Bursting Strength (kg)	:	27.3
Resistance to Bending (at 50000 times)	Warp :	Not Cut Off
	Weft :	Not Cut Off
Air permeability (cc/cm ² /sec. at 1/2° WG)	:	11.5
Pressure Loss (mmAq)	:	3.8
Efficiency (%)	0.5 u :	60.8
	1.0 u :	72.2
Pore Size (u)	min. :	23
	Peak :	29
	Max. :	54

2.0 Suction Wand

The suction wand is used for product input into the feed lines. This wand is designed for product input out of barrels, containers etc.

General Description


The suction wand consists of an outer tube (Item 1) and inner tube with air slit (Item 2), together with a breather ring (Item 3) for air supply.

The two tubes are joined together with spacer bushes (Item 4) and blind rivets (Item 5). A space stirrup (Item 6) is also mounted on the bottom of the wand to prevent the wand from attaching itself to the container bottom and blocking the product line.




Functional Description

The vacuum in the feed draws through the wand into the feed line. The throughput can be modified by adjusting the breather plug to change the amount of air that is sucked in.

The suction wand should be inserted into the product to such a depth so that the air holes in the breather plug remain clear.

ITEM No.	QTY	DESCRIPTION	MATERIAL	REFERENCE
1	1	CENTER TUBE	SS316	DV-STD-026
2	1	OUTER TUBE	SS316	DV-STD-027
3	1	AIR INTAKE COLLAR	SS316	DV-STD-028
4	1	STIRRUP	SS316	DV-STD-029
5	1	HEX. HD SCREW M6 x 15	GALV.	
6	4	SOLID RIVOT	SS316	
7	4	SPACER	SS316	REFER DETAIL

NOTES:

- WELDING:**
ALL WELDS WIRE BRUSHED (SS WIRE BRUSH).
ALL WELD DISCOLOURATION REMOVED.
ALL IRREGULAR WELDS DRESSED.
- SURFACE FINISH:**
PHARMACEUTICAL GRADE POLISH

KLAUSEN PROCESS MACHINERY PTY LTD
20 BERSBAMER ST BLAOKETOWN NSW 2148
TELE: (02)9672 3522 FAX: (02)9672 8733


CONTRACTS IN WHICH THIS DRAWING IS USED ARE THE PROPERTY OF KLAUSEN PROCESS MACHINERY PTY LTD AND IT IS NOT TO BE REPRODUCED BY ANY MEANS OR GIVEN TO ANY PARTY WITHOUT THE WRITTEN PERMISSION OF KLAUSEN PROCESS MACHINERY PTY LTD.

REV	DESCRIPTION	DATE	BY	APP'D
A	FIRST ISSUE	03/24	JAL	

--	--	--


DRAWN: J. LEECH	DATE: 13.02.01	JOB No.: STANDARD
CHECKED: S. LEL	DATE: 26.02.01	SCALE: 1:2.5
APP'D: A3	DWG No.: DV-STD-025	REV'T: A

TITLE:
DURA-VAC 2" SUCTION WAND
GENERAL ARRANGEMENT

3.0 Angle Valves

Angle valves in clean air lines
Angle valves in conveying air lines

Application

Angle valves are used for shutting off clean air respect. conveying air lines.

General Description


The angle valves may be delivered with 24 V or 220 V solenoid valves. Mounting in the pipe line is done with pipe clamps or with the welded, oval flanges.

Functional Description

- By moving the compressed air cylinder with disc out and in, the pipe line is opened and closed.
- With angle valves with aeration openings in the cover, one part of the pipe line is aerated when shutting off the other part. This is used for example in clean air lines coming from a separator. By closing the part of the pipe line to the underpressure generator, the holes in the cover are opened. The fresh air coming in aerates the separator. The pressure is balanced quicker.

Assembly instructions

The assembly depends on the use:

with or without aeration

- The version with oval flange is fastened to the corresponding counter flange. Do not forget the gasket.
- The version without oval flange is fastened in the pipe line with clamps. Pay attention to the fact, that the gasket lies in the slot provided for this in the halves of the pipe clamp.
- Connect the compressed air hose to the compressed air supply. Mount a shut-off cock into the line. This allows shutting off the for inspection and maintenance purposes.
- Compressed air free from water and oil, pressure approx. 6 bar absolute.
- Connect solenoid valve.
- With version for conveying air lines, the proximity switches for supervising the end positions must be connected (if available).
- Check, if compressed air connection from the solenoid valve to the compressed air cylinder is connected according to the function.

Technical Data

Nominal widths	:	38, 50 and 65 mm
Connection	:	220V/24V according to type specification plate, 50/60 Hz
Compressed air	:	free from water and oil
Comp. Air pressure	:	approx. 6 bar absolute
Comp. Air consumption	:	EV 38... and EV 50... 1 Ltr / switching impulse EV 65... and EVG.... 2 Ltr / switching impulse

Maintenance Instructions

- Check gaskets
 - every three months
 - with reduced conveying capacity (by - padded air is sucked in).
- Renew if necessary
 - Check the function of the solenoid valve and, if applicable, the proximity switches for supervising the end positions.

Instructions for remedying faults

Fault	Cause	Remedying
Conveying capacity reduces	Wear of gasket: by padded air is sucked in	Check gasket, renew if necessary
Compressed air cylinder is not activated	Electrical connection from solenoid valve to compressed air cylinder interrupted	Connect
	Compressed air pressure too low	Increase to approximately 6 bar absolute
	Compressed air supply interrupted	Connect
	Compressed air hose from solenoid valve to compressed air hose interrupted	Renew compressed air hose

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Further pump descriptions state the nominal displacement and the construction level:
Example:

RA 0100 E
RA = Ultimate pressure 0,5 mbar
0100 = Nominal displacement = 100 m³/h
E = Construction level

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The RA version, if used constantly, requires an oil-return suction to be installed into the B-endplate. Operating pressure must be below 300 mbar. When removing condensable vapours, a gas ballast valve must be installed. Both versions are air-cooled.

In case of questions about the application and versions, please contact the manufacturer.

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Transport and Packing

R 5 vacuum pumps pass a rigorous operating test in the factory and are packed carefully to avoid transit damage.

The inlet flange and as the case may be, the exhaust are sealed with a plug, so that no dirt can enter the pump during transport. Before start-up these plugs have to be removed. Please check packing on delivery for transport damage. The pump can be lifted from the packing with a suitable lifting device using the lifting bracket on the pump. (Fig. 7.1)

Do not lift the pump via the fan covers (fig. 2; 8, 9)!

Packing materials should be disposed of according to environmental laws or re-used.

These operating instructions are part of the consignment.

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**Pumps are generally shipped without oil.
Operating the vacuum pump without oil will damage the pump!**

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Start-up

It is essential to observe the following instructions step by step to ensure a safe start-up.

Start-up may only be conducted by trained specialists.

1. Setting-up

The pump must be set up or mounted horizontally on a flat surface. Special mounting is not required, as the pump can be mounted with screws via the threads of the rubber feet of the pump. The following ambient operating environment must be observed:

Ambient temperature: 12 to 30° C
Ambient pressure = Atmosphere

In order to avoid over-heating of the pump, an undisturbed fresh-air-flow to the pump is necessary.

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2. Inlet connection

The inlet flange can be connected with a vacuum-tight flexible hose or pipe.

The mounted pipes should cause no tension on the pump's flanges. If necessary, compensators must be installed.

Restriction of the pipes must be avoided in order not to decrease the displacement of the pump. The nominal diameter of the pipes has to be at least the same as the diameter of pump's inlet flange.

No foreign particles (e.g. solder soot) or liquids may enter the inlet line, as they could destroy the vacuum pump.

Restricting devices should not be installed in the exhaust line. Always connect the exhaust pipe in a manner, so that no condensate can enter the pump (slope, siphon).

3. Oil filling (Fig. 9.1)

The vacuum pump always leaves the factory without oil. Operation without oil destroys the pump! Therefore the pump must be charged with oil before start-up.

This is accomplished by charging the oil through the oil fill plug of the service unit, until the oil level reaches the MAX-mark on the oil sight glass.

Type and quantity of oil see "Maintenance" and page 9.

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Electrical connection

Electrical installation may only be conducted by a specialist. Regulations following EMV Directive 89/ 336 EWG, low Voltage Directive 73/ 23 EWG, and the

appropriate EN Standards have to be applied as well as VDE/ EVU regulations and local or national regulations. The operator of the vacuum pump has to inform the manufacturer, if electric or electromagnetic interference from his mains is to be expected.

1. Voltage and frequency on the nameplate must agree with the supply voltage.

2. The drive motor must be protected against overloads according to VDE 0113. In the case of portable installation of the vacuum pump, the electrical connection must be equipped with cable guides that have the function of traction relief.

3. To check the direction of rotation of the pump, flick the ON/OFF switch. In case of incorrect direction reverse the polarity of any two of the electrical phases.

Looking at the motor fan cover, the direction of rotation is counter-clockwise (fig. 6.1).

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Operation advice

1. These vacuum pumps can be used to evacuate air or dry gases, which are not aggressive, poisonous or explosive. Other agents may not be transported. In case of doubt, please contact the manufacturer.

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Attention! Not to be used with aggressive and explosive gases or gas mixtures and fluids.

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2. To pump condensable vapours, a gas ballast valve should be installed. The vacuum pump should run for 30 minutes prior to operation with the inlet connection closed, in order to reach the operating temperature of 75°C. Only at this operating temperature can condensable vapours be transported. After use, the pump should be left running for an additional 30 minutes to clear the oil of condensate.

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Maintenance

The vacuum pump must be switched off and secured against accidental switch-on for all maintenance.

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1. The oil level should be checked at least once a day.

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2. If the oil level is below the MAX-mark on the oil sight glass (3), more oil should be added.

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Fill with oil (5) until the level reaches the MAX-mark on the oil sight glass.

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3. Oil changing and oil filter changing
 Oil must be changed after the first 100 hours of operation.

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Further oil changes depend on operating conditions. The oil must be changed after 500 - 2000 hours of operation, but at least semi-annually.

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If there is considerable pollution it could be necessary to change the oil more frequently.

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To change the oil, the warm pump must be switched off and ventilated to reach atmospheric pressure. Drain the oil through the oil plug (4).

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Safety advice

In this operating instruction safety measures are advised before each step. It is imperative that these safety precautions are observed.

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Principle of operation

These pumps work according to the rotary vane principle. An eccentrically installed rotor (8) rotates in the cylinder. The centrifugal force of the rotation pushes the vanes (7), which glide in slots in the rotor, towards the wall of the cylinder. The vanes separate the sickle-shaped space between rotor and cylinder into chambers. When the chambers are connected with the inlet channel, gas is sucked in, compressed by the next rotation and pushed into the oil separator (1). The differential pressure constantly causes oil to be passed into the compression chambers. The oil and the medium are then discharged into the oil separator and there separated from the exhaust air by gravity, a demister (12) and the exhaust filter (9). The oil collects on the bottom of the oil separator and is then pushed into the compression chamber again (oil circulation). The oil-free medium is discharged through the exhaust cover plate (10) into atmosphere.

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Versions

The difference between the versions RA and RC is the ultimate pressure that can be reached:

- RA reaches an ultimate pressure of 0,5 mbar
- RC reaches an ultimate pressure of 20 mbar

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When oil stops running, close the plug and start up the pump again for a few seconds. Reopen the oil drain plug and discharge the remaining oil. Refasten the oil drain plug, remove the oil filter and replace it. Fill with fresh oil through the oil fill plug (5).

For oil quantity see nameplate.
 Used oil and used oil filters are to be disposed of according to environmental laws.

4. Types of oil

Oils according to DIN 51506, lubricating oil group VC must be used.
 We recommend original Busch oils of VM series, which comply with this DIN. If you need further information, request our leaflet "Special Oils for Vacuum Pumps".

5. Monitoring and change of exhaust filter
 The exhaust filter is best monitored using a filter pressure gauge (see "Accessories" on reverse side). This filter pressure gauge can be screwed into the thread of the oil fill plug. When the measured pressure reaches > 0,6 bar (Overpressure), the exhaust filter (9) has to be changed.

Increased energy intake by the motor could also be the result of soiled exhaust filters.

If oil mist escapes from the exhaust during the operation of the vacuum pump the exhaust filter must be changed. When installing a new exhaust filter, a new seal for the exhaust cover plate must also be installed.

The RA version may discharge oil in case of continuous operation (see "Versions", page 3 and 4).

6. Cleaning of the gas ballast valve

In case of visible dirt in the gas ballast valve, the sinter must be unscrewed and cleaned. Clean with blast air.

7. Cleaning of inlet flange

To clean the inlet flange screen (2), disconnect the four screws and remove the inlet flange. Take out the screen and clean with blast air.

8. Cleaning of fan covers

Both fan covers should be inspected regularly for dirt. Soiling of the fan covers prevents cool air intake and may lead to overheating of the vacuum pump.

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Information

We would be pleased to supply further information if required.
 Available are:

- Type sheet R 5 0063/ 0100 E
- Leaflet "Special Oils for Vacuum Pumps"
- R 5 Trouble shooting
- R 5 Storage Instructions

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
Spare parts and accessories

To guarantee safe operation of the vacuum pump, only original spare parts and accessories should be used. When ordering spare parts and accessories, always state pump specification and serial number.
 Please find the part number in the spare parts list.

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Order sample for spare parts:
 Pump: R 5 0100 E
 Spare part: Axial fan
 = Part number **0524108653**

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
- RA 0063 E
- RC 0063 E
- RA 0100 E
- RC 0100 E

These Installation and Operating Instructions are valid for the following pumps:

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It is mandatory that these operating instructions be read and understood prior to the vacuum pump installation and start-up.

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
Index	page
Safety	2
- Application	2
- Safety advice	3
Principle of operation	3
Versions	3-4
Transport and packing	4
Start-up	4-5
- Setting-up	4
- Inlet connection	5
- Oil filling	5
Electrical connection	5-6
Operating advice	6
Maintenance	6-7
- Oil level	6
- Oil changing and changing of oil filter	6-7
- Types of oil	7
- Quantity of oil required	9
- Monitoring and change of exhaust filter	7
- Cleaning of gas ballast valve	7
- Cleaning of inlet flange	7
Technical data	9
Spare parts	9-10
Exploded view drawing	8
Accessories	reverse side

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Safety

This vacuum pump has been manufactured according to the latest technical standards and safety regulations. If not installed properly or not used as directed, dangerous situations or damage might occur.

Application

This vacuum pump is designed for use in the fields of coarse vacuum and can be used to suck off air or dry gases, which are not aggressive, poisonous or explosive.

Other agents must not be transported. In case of doubt, please contact your local Busch Agency.

**Liquid and solid particles must not enter the pump.
 In case of doubt consult the manufacturer.**

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5.0 Trolley / Lifting Column

The vacuum transfer receiver is mounted onto a lifting column. The lifting column is electrically powered and has a lifting range of approximately 1 meter. The upper and lower limit of the column movement is protected by limit switches.

The lifting column, vacuum pump, control and ancillary items are mounted onto a trolley. The trolley can be wheeled into position so that the vacuum receiver discharge is positioned above the receiving hopper. The wheels on the trolley can then be locked.

Yearly Column Maintenance

- ♦ Turn on mains power and lower column to automatic stop
- ♦ Isolate transfer system from all power
- ♦ Remove limit switch cover on the lower end of the outer column tube
- ♦ Lift limit switches (rollers) so no roller will be able to touch the plastic lining
- ♦ Unscrew three screws in column below the motor plate
- ♦ Remove motor by lifting up motor mounting plate
- ♦ Unscrew and remove three M12 screws on the outer column tube and let the tube slide down
- ♦ Undo grub screw on motor couple sleeve. Remove sleeve and bearing block
- ♦ Secure ACME screw (crane, hoist etc) from falling into inner tube
- ♦ Use 19 mm socket and 19 mm open-end spanner to loosen nut in inner tube
- ♦ Use pliers to hold nut (to prevent nut from falling into tube) while removing screw
- ♦ Lift ACME screw with screw nut and screw nut mounting block out of inner tube
- ♦ Clean all parts with kerosene
- ♦ Inspect all parts for wear and proper functioning
- ♦ Grease ACME screw and assemble all parts in reversed order

6.0 Control System

Standard Control Panel Operation

The control system is designed for operating one vacuum transfer receiver with a single 1.5 kW vacuum pump. The control panel/control box contains the main switch, start/stop buttons and indication lamps as well as the switch for the lifting column. The box houses all the required contactors, the PLC and the selector switches for selection between timer or proximity switch mode.

In timer mode conveying continues until the set time has expired before discharging. In proximity switch mode, conveying stops when the proximity switch is triggered. Switching is via the selector switches located in the control box however, it is expected that once the preferred mode is determined and selected, you will need to change the mode again.

After turning on the main switch on the main panel, the system is active.

- 1 Angle valve for the selected receiver will open to the vacuum line
- 2 The pump will start
- 3 Conveying will continue until the proximity switch is triggered or the timer times out
- 4 The pump will stop
- 5 The discharge butterfly valve will open for the set time
- 6 Filter cleaning will occur
- 7 After the butterfly valve closes, the system will be ready for the next cycle
- 8 The system will operate until the stop button on the main panel is pushed, at which time, step 4 - 7 will occur.

Whilst transferring, a indicator light will display on the panel indicating transfer is being made to the selected receiver.

- ♦ The safety timer is to ensure that unit does not continue to operate if there is no product and it is unattended. If conveying stops because of the safety timer, conveying will start again on a new cycle however, if there are three consecutive stops because of the safety timer, a fault will appear and the unit must be reset before starting again.
- ♦ In timer modes the system switches off automatically after it reaches the set cycle count. To continue transferring press reset and then start again.

Column Operation

The lifting column is raised/lowered using the two way switch on the side of the control panel. The switch must be held in position to raise or lower the column.

External Operation

The transfer system can be started and stopped with an external switch e.g. lever sensor.

- ♦ External switch has to be connected to the socket on the bottom of the control panel
- ♦ Turn on main power/switch

The system will show start and stop according to the opening and closing of the external switch. The emergency button on the first panel will stop and reset the system with external operation or control panel operation. The start button on the panel will override the external switch.

Operator Interface Menu

The operator interface menu allows the operator easy access to the control settings of the vacuum transfer system to give optimal operating conditions.

On switching on the system:

VACUUM TRANSFER SYSTEM

is displayed

To change parameters press the ↓ key to give the suction/conveying time screen.

SUCTION TIME -- SEC
SET -- SEC

↓ key again to give the empty/discharge time screen.

DISCHARGE -- SEC
SET -- SEC

↓ key to give the filter clean screen.

FILTER CLEAN PLS -- SEC
-- SEC

↓ to give the cycle count screen.

CYCLE COUNT -- SEC
SET -- SEC

mode selection screen.

↓ key to give sensor/conveying time

LEVEL SENSOR ON
TURN OFF PRESS F2

Setting Suction / Conveying Time

The vacuum transfer can operate in sensor or conveying time mode. If the sensor is set to off, the system will operate in conveying time mode. In this mode the suction and therefore the conveying will be applied for the time set as the suction time. The optimal suction time will depend on the product and should be set to a value so that the receiver is filled to a level just below the filters before conveying stops.

To set the suction time, bring up the suction time screen. The currently set suction time (in seconds) is shown. Press 'ENTER' and then enter the new time using the up and down arrow keys and press 'ENTER' to confirm the setting.

Settings between 5 and 60 seconds are valid.

Setting Empty / Discharge Time

The empty or discharge time is the time the valve on the receiver is opened during each vacuum transfer cycle. With granular or free flowing products the time can be reduced to speed up the cycle time. With powdery products, the valve would need to be opened longer to enable the product to completely discharge from the vacuum receiver.

To set the empty time, bring up the empty time screen. The currently set empty time (in seconds) is shown. Press 'ENTER' and then enter the new time using the up and down arrow keys and press 'ENTER' to confirm the setting. Settings between 3 and 35 seconds are valid.

Setting Filter Cleaning Blasts

At the completion of each transfer cycle, there is a filter cleaning blast of compressed air. This cleans the filters and can also help dislodge product in the vacuum receiver. With granular products only one filter clean blast is necessary. With powdery products, more filter blasts may be required. Up to five blasts per transfer cycle can be set.

To set the number of filter cleaning air blasts, bring up the filter cleaning screen. The currently set number of cleaning air blasts is shown. Press 'ENTER' and then enter the required number using the up and down arrow keys and press 'ENTER' to confirm the setting.

Setting Maximum Cycle Count

The vacuum transfer system can be set so that it will stop after a set number of cycles. To set the maximum number of cycles, bring up the cycle count screen. Press 'ENTER' and then enter the new cycle count using the up and down arrow keys and press 'ENTER' again to confirm the setting.


Settings between 1 and 1000 are valid. If you do not typically want the system to stop after a set number of cycles, then set the maximum cycle count to 1000.

Selecting Sensor / Conveying Time Mode

To change between the sensor and conveying time mode, bring up the sensor/conveying time selection mode screen. The currently set position of the timer is displayed. Sensor on indicates sensor mode and sensor off indicates timer mode. Turn to sensor on and sensor off by pressing 'F1' or 'F2' keys as indicated. On selecting, the main starting screen 'Vacuum Transfer System' is shown.

7.0 Wiring Diagrams

1.0	LIFTING COLUMN
1.1	COLUMN LOW LEVEL LIMIT SWITCH
1.2	COLUMN HIGH LEVEL LIMIT SWITCH
1.3	COLUMN LIFTER MOTOR 0.18kW
2.0	VACUUM TRANSFER RECEIVER
2.1	RECEIVER HIGH LEVEL PROBE
2.2	DISCHARGE BUTTERFLY VALVE
2.3	BUTTERFLY VALVE ACTUATOR
2.4	BUTTERFLY VALVE SOLENOID
2.5	DISCHARGE VIBRATOR
2.6	PRODUCT LINE
2.7	AIR REGULATOR / GAUGE
3.0	FILTER CLEANING AIR TANK
3.1	FILTER CLEANING SOLENOID VALVE
3.2	COMPRESSED AIR DISCHARGE VALVE
3.3	AIR REGULATOR / GAUGE
3.4	NON RETURN VALVE
4.0	VACUUM PUMP 15kW
4.1	VACUUM RELIEF VALVE
4.2	VACUUM GAUGE
4.3	VACUUM REGULATING GATE VALVE
4.4	VACUUM LINE


TITLE:
DURA-VAC MOBILE VTS - ROCHE
CONTROL P&I DIAGRAM

REFERENCE DRAWING	A
DRAWN:	J. LEECH
CHECKED:	J. LEECH
DATE:	24.02.01
JOB No.:	MON016
SCALE:	NTS
DATE:	24.02.01
DWG No.:	MON016-005
3rd. ANGLE PROJECTION	A3
REV.	A

CONFIDENTIAL:
THIS DRAWING IN DESIGN AND DETAIL IS THE
PROPERTY OF KLAUSEN PROCESS MACHINERY
PTY LTD AND IS NOT TO BE REPRODUCED BY
ANY MEANS OR GIVEN TO A THIRD PARTY
WITHOUT THE SPECIFIC PERMISSION OF
KLAUSEN PROCESS MACHINERY PTY LTD.

KLAUSEN PROCESS
MACHINERY PTY LTD
20 BESSEMER ST BLACKTOWN NSW 2148
TEL: (02)9672 3222 FAX: (02)9672 8733


1 2 3 4 5 6 7 8

A B C D

8 7 6 5 4 3 2 1

D C B A

D C B A


TITLE: ROCHE - DURA-VAC DV -310-38 VTS
CIRCUIT DIAGRAM

REFERENCE DRAWING: VTSEST.DWG


DRAWN: AUHANG	DATE: 12.03.98	JOB No.: MON016
CHECKED: Lee	DATE: 26.02.98	SCALE: NONE
3rd. ANGLE PROJECTION	DWG No. MON016-006	SHEET: 1 OF 4
		REV.: B

CONFIDENTIAL:
THIS DRAWING IN DESIGN AND DETAIL IS THE PROPERTY OF KLAUSEN PROCESS MACHINERY PTY LTD AND IS NOT TO BE REPRODUCED BY ANY MEANS OR GIVEN TO A THIRD PARTY WITHOUT THE SPECIFIC PERMISSION OF KLAUSEN PROCESS MACHINERY PTY LTD.

KLAUSEN PROCESS MACHINERY PTY LTD
20 BESSEMER ST BLACKTOWN NSW 2148
TEL:(02)9672 3222 FAX:(02)9672 8733

8 7 6 5 4 3 2 1

TITLE: ROCHE - DURA-VAC DV-310-38 VTS
 CIRCUIT DIAGRAM



REFERENCE DRAWING: VTSESTDWG

DATE:	12/03/98	JOB No.:	MON016
DRAWN:	A. UHANG	CHECKED:	B. LAY
SCALE:	NONE	SHEET:	2 OF 4
DWG No.:	MON016-006	REV.:	B


3rd. ANGLE PROJECTION

CONFIDENTIAL: THIS DRAWING IN DESIGN AND DETAIL IS THE PROPERTY OF KLAUSEN PROCESS MACHINERY PTY LTD AND IS NOT TO BE REPRODUCED BY ANY MEANS OR GIVEN TO A THIRD PARTY WITHOUT THE SPECIFIC PERMISSION OF KLAUSEN PROCESS MACHINERY PTY LTD.

KLAUSEN PROCESS MACHINERY PTY LTD
 20 BESSEMER ST BLACKTOWN NSW 2148
 TEL: (02)9672 3222 FAX: (02)9672 8733

8 7 6 5 4 3 2 1


TITLE: ROCHE - DURA-VAC DV-310-38 VTS
CIRCUIT DIAGRAM

REFERENCE DRAWING: VISEST/DWG

DRAWN: A. UHANG	DATE: 12.03.98	JOB No.: MON016
CHECKED: J. JEC	DATE: 26.02.01	SCALE: NONE
DRAWN: A3		REV: 3 OF 4
CHECKED: A3		REV: B
3rd ANGLE PROJECTION		

CONFIDENTIAL: THIS DRAWING IN DESIGN AND DETAIL IS THE PROPERTY OF KLAUSEN PROCESS MACHINERY PTY LTD AND IS NOT TO BE REPRODUCED BY ANY MEANS OR GIVEN TO A THIRD PARTY WITHOUT THE SPECIFIC PERMISSION OF KLAUSEN PROCESS MACHINERY PTY LTD.


KLAUSEN PROCESS MACHINERY PTY LTD
20 BESSEMER ST BLACKTOWN NSW 2148
TEL: (02) 9672 3222 FAX: (02) 9672 8733

1 2 3 4 5 6 7 8

1 2 3 4 5 6 7 8

D C B A


TITLE:
ROCHE - DURA-VAC DV-310-38 VTS
PANEL LAYOUT

REFERENCE DRAWING: VTSEST.DWG

DATE:	12.03.98	JOB No.:	MON016
DRAWN:	A.UHANG	CHECKED:	3. LSC
DATE:	24.02.01	SCALE:	NONE
3rd. ANGLE PROJECTION		DWG No.	MON016-006
		APP'D	

CONFIDENTIAL:
 THIS DRAWING IN DESIGN AND DETAIL IS THE PROPERTY OF KLAUSEN PROCESS MACHINERY PTY LTD AND IS NOT TO BE REPRODUCED BY ANY MEANS OR GIVEN TO A THIRD PARTY WITHOUT THE SPECIFIC PERMISSION OF KLAUSEN PROCESS MACHINERY PTY LTD.

KLAUSEN PROCESS MACHINERY PTY LTD
 20 BESSEMER ST. BLACKTOWN NSW 2146
 TEL:(02)9672 3222 FAX:(02)9672 8733



1 2 3 4 5 6 7 8

D C B A

8.0 Parts List


ROCHE-VISMANUAL

ITEM NO.	QTY	DESCRIPTION	MATERIAL	REFERENCE
1	1	VACUUM RECEIVER ASSEMBLY		DV-310-38-41
2	1	RECEIVER SUPPORT BRACKET		DV-STD-07
3	1	LIFTING COLUMN ASSEMBLY		DV-STD-054
4	1	CONTROL PANEL		
5	1	CONTROL PANEL SUPPORT COLUMN		
6	1	PUSH BAR		DV-STD-06
7	1	VACUUM PUMP		BUSCH 0063
8	1	VACUUM PUMP GUARD		DV-STD-013
9	1	VACUUM LINE ASSEMBLY		DV-STD-041
10	1	BASE PLATE		DV-STD-005
11	2	LOCKABLE SWIVEL CASTORS		
12	2	FIXED CASTORS		
13	2	JUNCTION BOX		
14	1	DISCHARGE VALVE REGULATOR BOX	SS 304	
15	6	HEX HEAD SCREW M5 x 10	SS 304	
16	17	BUTTON HEAD SOCKET SCREW M5 x 10	SS 304	
17	2	HEX HEAD SCREW M5 x 20	SS 304	
18	25	FLAT WASHER M5	SS 304	
19	4	HEX HEAD SCREW M6 x 10	SS 304	
20	4	FLAT WASHER M6	SS 304	
21	7	HEX HEAD SCREW M8 x 20	SS 304	
22	4	HEX HEAD SCREW M8 x 25	SS 304	
23	16	C/SUNK SOCKET HEAD SCREW M8 x 25	SS 304	
24	3	PLASTIC HANDLE, THREAD M8 x 20	SS 304	
25	31	FLAT WASHER M8	SS 304	
26	19	SPRING WASHER M8	SS 304	
27	4	LODGE NUT M8	SS 304	
28	6	SPACER M8	SS 304	
29	16	HEX NUT M8	SS 304	
30	2	HEX HEAD SCREW M10 x 20	SS 304	
31	4	HEX HEAD SCREW M10 x 25	SS 304	
32	6	FLAT WASHER M10	SS 304	
33	2	DOME NUT M5	SS 304	
34	1	15" PIPE REST		DV-STD-056
35	1	15" PIPE SADDLE		DV-STD-057

TITLE:
**DURA-VAC MOBILE VTS
 GENERAL ARRANGEMENT
 DV-310-38 RECEIVER - SERIAL No. M0102755**

REFERENCE DRAWING


DATE:	16.02.01
CHECKED:	B. V.S.
DATE:	26.02.01
SCALE:	1:1
DWG No.	MON016-003
3rd. ANGLE PROJECTION	



REV	DESCRIPTION	DATE	DRN	APP'D
A	FIRST ISSUE	5.02.01	JAL	

CONFIDENTIAL:
 THIS DRAWING IN DESIGN AND DETAIL IS THE
 PROPERTY OF KLAUSEN PROCESS MACHINERY
 PTY LTD AND IS NOT TO BE REPRODUCED BY
 ANY MEANS OR GIVEN TO A THIRD PARTY
 WITHOUT THE SPECIFIC PERMISSION OF
 KLAUSEN PROCESS MACHINERY PTY LTD.

**KLAUSEN PROCESS
 MACHINERY PTY LTD**
 20 BESSEMER ST BLACKTOWN NSW 2148
 TEL:(02)9672 3222 FAX:(02)9672 6733

ITEM No.	QTY	DESCRIPTION	MATERIAL	REFERENCE
1	1	VACUUM RECEIVER ASSEMBLY		DV-310-38-01
2	1	RECEIVER SUPPORT BRACKET		DV-STD-007
3	1	LIFTING COLUMN ASSEMBLY		DV-STD-004
4	1	CONTROL PANEL		
5	1	CONTROL PANEL SUPPORT COLUMN		
6	1	PUSH BAR		DV-STD-018
7	1	VACUUM PUMP		DV-STD-042
8	1	VACUUM PUMP GUARD		BUSCH 0063
9	1	VACUUM LINE ASSEMBLY		DV-STD-013
10	1	BASE PLATE		DV-STD-041
11	2	LOCKABLE SWIVEL CASTORS		DV-STD-005
12	2	FIXED CASTORS		
13	2	JUNCTION BOX		
14	1	DISCHARGE VALVE REGULATOR BOX		DV-STD-005
15	6	HEX. HEAD SCREW M5 x 10	SS 304	
16	18	BUTTON HEAD SOCKET SCREW M5 x 10	SS 304	
17	2	HEX. HEAD SCREW M5 x 20	SS 304	
18	25	FLAT WASHER M5	SS 304	
19	4	HEX. HEAD SCREW M6 x 10	SS 304	
20	4	FLAT WASHER M6	SS 304	
21	7	HEX. HEAD SCREW M8 x 20	SS 304	
22	4	HEX. HEAD SCREW M8 x 25	SS 304	
23	16	L/SUNK SOCKET HEAD SCREW M8 x 25	SS 304	
24	3	PLASTIC HANDLE THREAD M8 x 20		
25	31	FLAT WASHER M8	SS 304	
26	19	SPRING WASHER M8	SS 304	
27	4	DOVE NUT M8	SS 304	
28	6	SPACER M8	SS 304	
29	16	HEX. NUT M8	SS 304	
30	2	HEX. HEAD SCREW M10 x 20	SS 304	
31	4	HEX. HEAD SCREW M10 x 25	SS 304	
32	6	FLAT WASHER M10	SS 304	
33	2	DOVE NUT M5	SS 304	
34	1	15° PIPE REST		DV-STD-066
35	1	15° PIPE SADDLE		DV-STD-067


TITLE:
**DURA-VAC MOBILE VTS
 GENERAL ARRANGEMENT
 DV-310-38 RECEIVER - SERIAL No. M0102285S**


DATE: 16.02.01	JOB No. MON016	SHEET: 1 OF 1
DATE: 26.02.01	SCALE: 1:10	REV: A
DRAWN: JLEECH	CHECKED: B. IRE	DWG No. MON016-004
3rd. ANGLE PROJECTION		

CONFIDENTIAL: THIS DRAWING IN DESIGN AND DETAIL IS THE PROPERTY OF KLAUSEN PROCESS MACHINERY PTY LTD AND IS NOT TO BE REPRODUCED BY ANY MEANS OR GIVEN TO A THIRD PARTY WITHOUT THE SPECIFIC PERMISSION OF KLAUSEN PROCESS MACHINERY PTY LTD.

KLAUSEN PROCESS MACHINERY PTY LTD
 20 BESSEMER ST BLACKTOWN NSW 2146.
 TEL:(02)9672 3522 FAX:(02)9672 8733

REV	DATE	DRN	APP'D	DESCRIPTION
A	02/01	JAL		FIRST ISSUE

ITEM No.	QTY	DESCRIPTION	MATERIAL	REFERENCE
1	1	VACUUM PUMP - BUSCH RC 0063 E		
2	1	VACUUM GUAGE		
3	1	BALL VALVE		
4	1	VACUUM RELIEF VALVE		
5	1	ELBOW		
6	1	PIPE CLAMP		
7	1	1.5" ANGLE VALVE		
8	1	SOCKET HEAD SCREW M8 x 25	SS 304	
9	4	NYLOC NUT M8	SS 304	
10	4			

TITLE: DURA-VAC MOBILE VTS
VACUUM LINE ASSEMBLY

REFERENCE DRAWING:


DRAWN: JLEECH
CHECKED: G. UEL
DATE: 22.02.01
DATE: 8.11.02.01

DWG No. DV-STD-041

JOB No. STANDARD
SCALE: 1:5
SHEET: 1 OF 1
REV: A

CONFIDENTIAL:
THIS DRAWING IN DESIGN AND DETAIL IS THE
PROPERTY OF KLAUSEN PROCESS MACHINERY
PTY LTD AND IS NOT TO BE REPRODUCED BY
ANY MEANS OR GIVEN TO A THIRD PARTY
WITHOUT THE SPECIFIC PERMISSION OF
KLAUSEN PROCESS MACHINERY PTY LTD.


**KLAUSEN PROCESS
MACHINERY PTY LTD**
20 BESSEMER ST BLACKTOWN NSW 2148
TEL: (02)9672 3222 FAX: (02)9672 8733

REV	DESCRIPTION	DATE	DRN	APP'D
A	FIRST ISSUE	22/2/01		

1 2 3 4 5 6 7 8

ITEM NO.	QTY	DESCRIPTION	MATERIAL	REFERENCE
1	1	VACUUM CONNECTION ELBOW	RUBBER	DV-310-38-012
2	1	OVAL GASKET		
3	1	REGULATOR / GAUGE		
4	1	NON-RETURN VALVE		
5	A	AIR BLEED DIAPHRAGM VALVE		
6	1	COMPRESSED AIR TANK		
7	1	FILTER CLEANING SOLINOID VALVE		DV-310-38-013
8	1	VACUUM RECEIVER LID		DV-310-38-014
9	1	FILTER CLAMP PLATE		DV-310-38-017
10	1	FILTER SEAL		DV-310-38-018
11	1	FILTER CARTRIDGE		
12	1	FILTER PLATE		
13	1	FILTER PLATE SEAL		
14	3	RECEIVER LID CLAMP		
15	1	RECEIVER BODY		
16	1	LEVEL SENSOR		
17	1	VIBRATOR MOUNTING PLATE		DV-310-38-016
18	1	VIBRATOR		
19	1	BUTTERFLY VALVE TOP FLANGE		
20	1	BUTTERFLY VALVE SEAL		
21	1	BUTTERFLY VALVE BLADE		
22	1	BUTTERFLY VALVE BOTTOM FLANGE		
23	1	DISCHARGE SPIGOT		DV-310-38-005
24	1	BUTTERFLY VALVE ACTUATOR MOUNTING PLATE		
25	1	BUTTERFLY VALVE ACTUATOR		
26	8	WANG STUD M8 x 40	SS 304	
27	4	CSK SCREW M8 x 20	SS 304	
28	8	BUTTON HD SOCKET SCREW M6 x 6	SS 304	
29	1	HEX. HD. SCREW M6 x 6	SS 304	
30	5	FLAT WASHER M6	SS 304	
31	4	DOME NUT M6	SS 304	
32	2	HEX. HD. SCREW M8 x 25	SS 304	
33	2	HEX. HD. SCREW M8 x 20	SS 304	
34	2	HEX. HD. SCREW M8 x 40	SS 304	
35	2	HEX. HD. SCREW M8 x 45	SS 304	
36	4	HEX. HD. SCREW M8 x 16	SS 304	
37	2	FLAT WASHER M8	SS 304	
38	1	DOME NUT M8	SS 304	
39	1	OVAL GASKET	RUBBER	


TITLE: DURA-VAC DV-310-38
VACUUM RECEIVER
GENERAL ARRANGEMENT

DRAWN: JLEECH	DATE: 26.02.01	JOB No.: STANDARD
CHECKED: S. LEE	DATE: 26.02.01	SCALE: 1/5
DWG No. DV-310-38-011		SHEET: 1 OF 1
REV: A		REV.: A


REV	DESCRIPTION	DATE	APP'D
A	FIRST ISSUE	26.02.01	JAL

CONFIDENTIAL IN DESIGN AND DETAIL IS THE
PROPERTY OF KLAUSEN PROCESS MACHINERY
PTY LTD AND IS NOT TO BE REPRODUCED BY
ANY MEANS OR GIVEN TO A THIRD PARTY
WITHOUT THE SPECIFIC PERMISSION OF
KLAUSEN PROCESS MACHINERY PTY LTD.

**KLAUSEN PROCESS
MACHINERY PTY LTD**
20 BESSEMER ST BLACKTOWN NSW 2148
TEL:(02)9672 3222 FAX:(02)9672 8733

ITEM No.	QTY	DESCRIPTION	MATERIAL	REFERENCE
1	1	MOTOR COVER		DV-STD-002
2	1	MOTOR XXOXY		
3	1	LIFT COLUMN OUTER TUBE		DV-STD-003
4	1	LIFT COLUMN INNER TUBE		DV-STD-004
5	1	MOTOR MOUNTING PLATE		DV-STD-001
6	2	DRIVE KEY - 4.8 (3/16") SQR x 25	KEY STEEL	
7	1	DRIVE COUPLING		DV-STD-049
8	1	LIFT COLUMN BEARING BLOCK		DV-STD-047
9	1	LIFT COLUMN DRIVE SCREW M6 x 1066		DV-STD-050
10	1	SAFETY NUT - HEX. T ACME THREAD	M.S	
11	2	SAFETY NUT RETAINING CLIPS		
12	1	LIFT COLUMN DRIVE NUT	BRONZE	DV-STD-051
13	1	LIFT COLUMN DRIVE NUT RETAINING BLOCK	M.S	DV-STD-048
14	2	LIFT COLUMN DRIVE NUT LOCATING PIN		DV-STD-053
15	1	LIFT COLUMN LIMIT SWITCH COVER		DV-STD-066
16	2	LIMIT SWITCH		
17	2	LIFT COLUMN LIMIT SWITCH BRACKET		DV-STD-046
18	1	BEARING		
19	1	HEX. HEAD SCREW M4 x 35	SS 304	
20	4	PAN PHILLIPS HEAD SCREW M4 x 4	SS 304	
21	4	PAN PHILLIPS HEAD SCREW M4 x 6	SS 304	
22	6	PAN PHILLIPS HEAD SCREW M4 x 6	SS 304	
23	18	FLAT WASHER M4	SS 304	
24	4	HEX. NUT M4	SS 304	
25	8	HEX. HEAD SCREW M5 x 10	SS 304	
26	8	FLAT WASHER	SS 304	
27	3	HEX. HEAD SCREW M6 x 10	SS 304	
28	4	SOCKET HEAD SCREW M6 x 20	SS 304	
29	3	FLAT WASHER M6	SS 304	
30	2	HEX. NUT M2	SS 304	
31	3	HEX. HEAD SCREW M2 x 25	SS 304	
32	3	FLAT WASHER M2	SS 304	

TITLE:
**DURA-VAC MOBILE VTS
 LIFT COLUMN WITH 760mm TRAVEL
 GENERAL ARRANGEMENT**

DATE:	22.02.01	JOB No.:	STANDARD
CHECKED:	B. LEECH	SCALE:	1:1
DATE:	22.02.01	REV.:	A
3rd. ANGLE PROJECTION			

REV	DESCRIPTION	DATE	APP'D
A	FIRST ISSUE	22.02.01	JAL

CONFIDENTIAL:
 THIS DRAWING IN DESIGN AND DETAIL IS THE
 PROPERTY OF KLAUSEN PROCESS MACHINERY
 PTY LTD AND IS NOT TO BE REPRODUCED BY
 ANY MEANS OR GIVEN TO A THIRD PARTY
 WITHOUT THE WRITTEN PERMISSION OF
 KLAUSEN PROCESS MACHINERY PTY LTD.

**KLAUSEN PROCESS
 MACHINERY PTY LTD**
 20 BESSEMER ST BLACKTOWN NSW 2146
 TEL:(02)9672 3222 FAX:(02)9672 8733