

OLYMPUS

Your Vision, Our Future

ENDOSCOPE OVERVIEW 2017/2018



IMAGINE...

**Your Thumb Could
Maximize ADR**

Explore the evidence at:

 www.olympus.co.uk/proven

EVIS LUCERA
ELITE

MEDICAL ENDOSCOPES

Medical Video Endoscopes										
	Insertion Tube OD (mm)	Distal End OD (mm)	Channel ID (mm)	Bending (Up/Down)	Bending (Left/ Right)	Working Length (mm)	Field of View	Depth of Field (mm)	NBI	Comments
Video Gastrosopes										
GIF-HQ290*	9.9	10.2	2.8	210°/90°	100°/100°	1030	140° (near) 140° (normal)	3 – 7 (near) 7 – 100 (normal)	✓	DF, OT, WJ, HDTV
GIF-H290Z*	9.6	9.9	2.8	210°/90°	100°/100°	1030	95° (tele) 140° (wide)	1.5 – 3 (tele) 7 – 100 (wide)	✓	Zoom, OT, WJ, HDTV
GIF-H290*	8.9	8.9	2.8	210°/90°	100°/100°	1030	140°	3 – 100	✓	OT, HDTV
GIF-XP290N*	5.8	5.4	2.2	210°/90°	100°/100°	1100	140°	3 – 100	✓	Transnasal, OT, HF compatible
GIF-H260	9.5	9.8	2.8	210°/90°	100°/100°	1030	140°	3 – 100	✓	HDTV
GIF-PQ260	7.7	7.9	2.0	210°/90°	100°/100°	1030	140°	3 – 100	✓	–
GIF-XP260NS	5.8	5.4	2.0	210°/90°	100°/100°	1100	120°	3 – 100	✓	Transnasal
GIF-FQ260Z**	10.5	11.0	2.8	210°/90°	100°/100°	1030	60° (tele) 140° (wide) 140° (AFI)	2 – 3.5 (tele) 7 – 100 (wide) 7 – 100 (AFI)	–	AFI, Zoom
GIF-2TQ260M	11.7	11.7	2x 3.2	210°/90° Sec.: 70°/70°	100°/100°	1030	140°	3 – 100	✓	WJ, Multi- bending

* with CV-290 only ** with CV-290/-260SL only

OT = One-touch Connector DF = Dual Focus CF = Close Focus
WJ = WaterJet RIT = Responsive Insertion Technology SG = ScopeGuide
VS = Variable Stiffness HFT = High Force Transmission PB = Passive Bending
AFI = Auto-Fluorescence Imaging

Subject to change without notice. Liability for errors and omissions excepted.

MEDICAL ENDOSCOPES

	Insertion Tube OD (mm)	Distal End OD (mm)	Channel ID (mm)	Bending (Up/Down)	Bending (Left/ Right)	Working Length (mm)	Field of View	Depth of Field (mm)	NBI	Comments
Video Colonoscopes										
CF-HQ290L/I*	12.8	13.2	3.7	180°/180°	160°/160°	1680/1330	160° (near) 170° (normal)	4 – 9 (near) 9 – 100 (normal)	✓	DF, OT, RIT, SG, WJ, HDTV
CF-HQ290ZL/I*	12.8	13.2	3.7	180°/180°	160°/160°	1680/1330	90° (tele) 170° (wide)	2 – 3 (tele) 7 – 100 (wide)	✓	OT, RIT, WJ, HDTV, Zoom
CF-H290L/I*	12.0	12.2	3.2	180°/180°	160°/160°	1680/1330	170°	5 – 100	✓	OT, RIT, WJ, HDTV
PCF-H290DL/I*	11.8	11.7	3.2	180°/180°	160°/160°	1680/1330	170°	5 – 100	✓	OT, RIT, SG, WJ, HDTV
PCF-H290ZL/I*	11.8	11.7	3.2	180°/180°	160°/160°	1680/1330	85° (tele) 170° (wide)	1 – 2 (tele) 7 – 100 (wide)	✓	HDTV, OT, RIT, WJ, Zoom
PCF-H290L/I*	11.5	11.7	3.2	180°/180°	160°/160°	1680/1330	170°	5 – 100	✓	OT, RIT, WJ, HDTV
CF-H260DL	13.2	13.6	3.7	180°/180°	160°/160°	1680	140°	5 – 100	✓	HDTV, SG, VS
CF-H260AZL/I	12.9	13.6	3.2	180°/180°	160°/160°	1680/1330	80° (tele) 140° (wide)	2 – 3 (tele) 7 – 100 (wide)	✓	HDTV, Zoom, VS
CF-FH260AZL/I**	13.2	14.8	3.2	180°/180°	160°/160°	1680/1330	80° (tele) 140° (wide) 140° (AFI)	2 – 3 (tele) 7 – 100 (wide) 5 – 100 (AFI)	✓	HDTV, AFI, WJ, Zoom, VS
CF-Q260DL	12.4	12.2	3.2	180°/180°	160°/160°	1680	140°	5 – 100	✓	SG, VS
PCF-Q260AZI	11.8	11.7	3.2	180°/180°	160°/160°	1330	60° (tele) 140° (wide)	2 – 3.5 (tele) 7 – 100 (wide)	✓	WJ, Zoom, VS
PCF-Q260JL	10.5	10.5	3.2	180°/180°	160°/160°	1680	140°	5 – 100	✓	WJ
PCF-PQ260L/I	9.2	9.2	2.8	180°/180°	160°/160°	1680/1330	140°	5 – 100	✓	RIT (HFT, PB)

* with CV-290 only ** with CV-290/-260SL only

OT = One-touch Connector DF = Dual Focus CF = Close Focus
WJ = WaterJet RIT = Responsive Insertion Technology SG = ScopeGuide
VS = Variable Stiffness HFT = High Force Transmission PB = Passive Bending
AFI = Auto-Fluorescence Imaging

Subject to change without notice. Liability for errors and omissions excepted.

MEDICAL ENDOSCOPES

	Insertion Tube OD (mm)	Distal End OD (mm)	Channel ID (mm)	Bending (Up/Down)	Bending (Left/ Right)	Working Length (mm)	Field of View	Depth of Field (mm)	NBI	Comments
Video Duodenoscopes										
TJF-260V	11.3	13.5	4.2	120°/90°	120°/90°	1240	100°	5 – 60	✓	V-System
JF-260V	11.3	12.6	3.7	120°/90°	120°/90°	1240	100°	5 – 60	✓	V-System
Video Choledochoscope										
CHF-V	5.2	4.9	2.0	160°/130°	–	380	120°	3 – 50	✓	–
Video Enteroscope										
SIF-Q260	9.2	9.2	2.8	180°/180°	160°/160°	2000	140°	5–100	✓	–
Capsule Endoscope										
EC-S10***	11.0	11.0	–	–	–	26	160°	0–20	–	–

*** with EC-10 ENDOCAPSULE system only

OT = One-touch Connector DF = Dual Focus CF = Close Focus
WJ = WaterJet RIT = Responsive Insertion Technology SG = ScopeGuide
VS = Variable Stiffness HFT = High Force Transmission PB = Passive Bending
AFI = Auto-Fluorescence Imaging

Subject to change without notice. Liability for errors and omissions excepted.

MEDICAL ENDOSCOPES

	Insertion Tube OD (mm)	Distal End OD (mm)	Channel ID (mm)	Bending (Up/Down)	Bending (Left/ Right)	Working Length (mm)	Field of View	Depth of Field (mm)	NBI	Comments
Video Bronchoscopes										
BF-H290*	5.7	6.0	2.0	210°/130°	–	600	120°	3 – 100	✓	HDTV, OT
BF-Q290*	4.9	4.8	2.0	210°/130°	–	600	120°	2 – 100	✓	OT
BF-P290*	4.1	4.2	2.0	210°/130°	–	600	110°	2 – 50	✓	OT
BF-XP290*	2.8	3.1*	1.2	210°/130°	–	600	110°	2 – 50	✓	OT
BF-1TQ290*	6.0	5.9	3.0	180°/130°	–	600	120°	2 – 100	✓	OT
BF-1T260	6.0	5.9	2.8	180°/130°	–	600	120°	2 – 100	–	–
BF-F260**	5.4	5.5	2.0	180°/130°	–	600	120°	3 – 100	–	AFI
BF-P260F	4.4	4.0	2.0	180°/130°	–	600	120°	3 – 50	–	Fibervideoscope
BF-XP260F	2.8	2.8	1.2	180°/130°	–	600	120°	2 – 50	–	Fibervideoscope
Video Pleuroscope (Medical Thoracoscopy)										
LTF-160	7.0	7.0	2.8	160°/130°	–	270	120°	3 – 100	–	Autoclavable

*with CV-290 only **with CV-290/-260SL only
♦tapered part of distal end 2.9mm

OT = One-touch Connector DF = Dual Focus CF = Close Focus
WJ = WaterJet RIT = Responsive Insertion Technology SG = ScopeGuide
VS = Variable Stiffness HFT = High Force Transmission PB = Passive Bending
AFI = Auto-Fluorescence Imaging

Subject to change without notice. Liability for errors and omissions excepted.

MEDICAL ENDOSCOPES

Medical Fiber Endoscopes									
	Insertion Tube OD (mm)	Distal End OD (mm)	Channel ID (mm)	Bending (Up/Down)	Bending (Left/Right)	Working Length (mm)	Field of View	Depth of Field (mm)	Comments
Gastrointestinal Fiberscopes									
CHF-BP30	3.4	3.1	1.2	160°	130°	1870	90°	1 – 150	Baby scope
Fiber Bronchoscopes									
BF-1T60	6.0	5.9	3.0	180°/130°	–	600	120°	3 – 50	–
BF-P60	5.0	4.9	2.2	180°/130°	–	600	120°	3 – 50	–
BF-XP60	2.8	2.8	1.2	180°/130°	–	600	90°	2 – 50	–
BF-3C40	3.6	3.3	1.2	180°/130°	–	550	120°	3 – 50	–

ENDOSCOPY INTEGRATION

Olympus' systems offer a complete platform of solutions for endoscopy, consisting of modular components that are combined to meet the individual needs of the hospital. Utmost support throughout the workflow from patient reception to the examination, via the reprocessing and recovery areas, ending in the doctor's office and training room facilities.



Examination Room

A high level of ergonomics, efficiency, and communication in a modern and friendly ambience.

- Optimized workspace and patient comfort
- Integrated HD documentation system
- Saving space, examination, and cleaning time



Reprocessing Area

Hygiene supports and secures the endoscope's workflow with full traceability and documentation of all relevant reprocessing steps.

- ETD Double and ETD4 washer disinfectors
- EDC Plus Drying and storage cabinets



Recovery Area

Patient comfort and safety in perfection – a central overview of the status of the rooms, support smooth processes and a stressless atmosphere.

- Access to status of endoscopy room
- Patient-care documentation for recovery nurse



Doctor's Office

In the office but with an optimal audio/video connection to the examination rooms. Access from the PC to all data and images and reporting system.

- Access to centralized patient/procedure information
- Audio/video connection to examination rooms
- Tools for activity and cost analysis



Training / Meeting Room

The media center for training, courses, and seminars with advanced video management, including live video and bidirectional communication from examination rooms.

- Access/display of departmental activity information
- IP audio/video routing and multi-display functionality

ENDOSCOPIC ULTRASOUND

Endoscopic Ultrasound									
	Insertion Tube OD (mm)	Channel ID (mm)	Bending (Up/Down Left/Right)	Working Length (mm)	Ultrasound Field of View	Endoscopic Field of View (Forward/Oblique)	Depth of Field (mm)	Frequency (MHz)	Scanning Mode
Ultrasonic Videoscopes for EUS									
GF-UE260-AL5	11.8	2.2	130°/90° 90°/90°	1250	360°	100°/55°	3 – 100	5/6/7.5/10*	Electronic 360° radial array
TGF-UC260J	12.6	3.7	180°/90° 90°/90°	1245	90°	Forward viewing	3 – 100	5/6/7.5/10/12*	Electronic curved linear array
GF-UCT260	12.6	3.7	180°/90° 90°/90°	1250	180°	100°/55°	3 – 100	5/6/7.5/10*	Electronic curved linear array
Ultrasonic Endoscope for EBUS									
BF-UC260FW	6.3	2.2	120°/90° —	600	60°	80°/35°	2 – 50	5/6/7.5/10/12*	Electronic curved linear array

*with EU-ME2. Please contact your Olympus representative for compatibility with other systems ** Equipped with a 0.035" guidewire port at the distal end

ENDOSCOPIC ULTRASOUND

	Frequency (MHz)	Working Length (mm)	Insertion Tube OD (mm)	Comp. Biopsy Channel (mm)	Model	Comp. Biopsy Channel (mm)	Scanning Mode
Ultrasonic Miniature Probes							
UM-2R/3R	12/20	2050	2.5	2.8	MH-246R Balloon Sheath	3.7	Mechanical radial
UM-S20-17S	20	2150	max. 1.8	2.0	K-201 / K-202 GuideSheathKit	2.0	Mechanical radial
UM-S20-20R	20	2050	max. 2.0	2.2•	K-203 / K-204 GuideSheathKit	2.6	Mechanical radial
UM-BS20-26R	20	2050	Use only with balloon sheath		MAJ-643R Balloon Sheath	2.8	Mechanical radial
UM-G20-29R**	20	2050	max. 2.9	3.2	Not used with balloon sheath		Mechanical radial
UM-DG20-31R**	20	2050	max. 3.1	3.7	Not used with balloon sheath		Mechanical helical for Dual Plane Reconstruction (DPR)
UM-DP12-25R	12	2050	2.5	2.8	Not used with balloon sheath		Mechanical helical for Dual Plane Reconstruction (DPR)
UM-DP20-25R	20	2050	2.5	2.8	Not used with balloon sheath		Mechanical helical for Dual Plane Reconstruction (DPR)

* with EU-ME2. Please contact your Olympus representative for compatibility with other systems ** Equipped with a 0.035" guidewire port at the distal end
• 2.0 mm in bronchoscopes


ADVANCING THE ART OF BRONCHOSCOPY

Finding and Sampling of Peripheral Pulmonary Lesions

Peripheral Pulmonary Lesion

Real-time confirmation of the lesion through endobronchial ultrasound using an EBUS miniature probe.

Combining EBUS with a GuideSheath enables easy and repeated access to the lesion via sampling devices.



7023



NEW

EZDILATE

Multi-Stage Endoscopic Balloon Dilator

EZDilate balloons were specifically designed to meet your needs while setting a higher standard for performance and precision.

- Smooth insertion and navigation to the stricture while protecting sensitive anatomy
- Precise placement and confidence in positioning
- Optimal clarity
- Controlled accurate dilation

The EZDilate 3-stage balloon dilatation catheters are available in a full range of fixed wire and wire-guided size options for use in esophageal, pyloric, colonic, and biliary applications.

 www.olympus.co.uk/etcatalogue



Subject to change without notice. Liability for errors and omissions excepted.

OLYMPUS

OLYMPUS KEYMED

KeyMed House, Stock Road, Southend-on-Sea,
Essex, SS2 5QH, UK

Telephone: +44 (0)1702 616333

info@olympus.co.uk www.olympus.co.uk