

FP Double Head (34-2 & 46-2)

IWK Verpackungstechnik GmbH Lorenzstraße 6 D – 76297 Stutensee Phone: +49 7244 968-0

www.iwk.de




#### Demand




"An efficient and economical tube filling system at an attractive price"






# Concept FP 34-2









## Concept FP 46-2






# Concept FP 64-2









# Station overview






# FP 46-2 (with robot tube inserting)





# Footprint FP 46-2





#### Characteristics



- Performance
- Tube diameter (D) plastic/laminate tube
- Tube diameter (D) metal tube
- Tube length with cap (L)
- Tube coat length (L1) plastic/laminate tube
- Tube coat length (L1) metal tube (double fold)
- Tube coat length (L1) metal tube (triple fold)
- Tube coat length (L1) metal tube (saddle fold)
- Filling volume

up to 220 tubes/min

Ø 10-52 mm

Ø 10-52 mm

60-250 mm


50-230 mm

50-230 mm

56-230 mm

62-230 mm

1,5-350 ml



#### Transport system





- no stretching of the belts
- enables a high degree of positioning precision



Tube holders depend on tube length and diameter

For plastic tubes up to Ø52 mm, which clear dents by the tubes' elastic force. With tilting segments



For metal and laminate tubes up to Ø52 mm. With 4-part metal-clamping segments.





- cast chain links
- mounted in precision needle bearings
- servo-motor drive with high-precision planetary gear drive

# Tube inserting








Size parts



Insertion rail\*





Pusher cone

\* also depends on tube length

## Print registration










Centering cone print registration depends on tube diameter




# Tube dedusting









#### Dosing & filling




servo driven lifter and dosing piston

Stable design of filling frame

high filling accuracy

almost tool-free changing of wet parts changing of wet parts on block (option) optimal accessibility




stream-optimized low shearing

easily accessible

easy to clean

- FDA-approved
- Limited dead legs





#### Product-contact format parts



Dosing piston depends on filling volume



Filling pipe depends on tube diameter



mass = 3.13 kg filling pipe suspension



mass = 1.64 kg rotary valve




mass = 6.8 kg rotary valve housing

Change of rotary valve housing and of filling pipe suspension is only required when changing the volume </> 15 ml

# Lifting device application with trolley of product-contact parts






FP Double Head I IWK Verpackungstechnik GmbH

# Multicolor filling





**FP Double & Triple Head** 

#### Fold closure



#### **Double-/ triple-fold**




Saddle-fold



We do not use lubricated guides above the tube (cGMP!)


#### **Folding tools**






# Folds folding width 5mm






double-fold



triple-fold




saddle fold


# Hot Air-closing



robust design

high performance heating element
energy-efficient pre-Heating (without water)
servo drive as an option for contour welds
change over without tools







## Format Parts Hot-Air tool



Hot-air tool depends on tube diameter



maximum sealing height for contour sealing is 40 mm



DISTURIZING LOTIO






Contour seal - ARC Contour seal - Wav








Hot-air tool with slot



Hot-air tool with holes










pre-forming tool as size part depend on tube diameter and length





change over without tools robust mechanical design cam driven compact design pressure adjustment: continuously variable

Chilled pressing jaws depend on tube fold embossing image





adjustable pre-forming tool

## Trimming station for straight cut



double-sided driven shear blades


robust mechanics

compact design

cam driven

shear blades made of powder-

metallurgical structure





#### Punching station





double-sided driven punching tool

tool-free exchange of punching tools



robust mechanics

cam driven

compact design

#### Tube discharge



**FP Double Head** 

Direct discharge into

SC5

SC4



- high reproducibility
- high positioning accuracy
- high repetition accuracy

**Size parts** 



**Tube flipper** 



**Brace** 



**Suction cups holder** 

see simulation

#### Drive concept



#### Main functions servo driven

- Print registration
- Lifter filling and dosing piston
- Rotary valve (option)
- Lowering Hot Air (option)
- Transport, main shaft

#### Height adjustment of transport system is motor-driven

- high processing reliability
- high accuracy
- highest flexibility, because of
- fast change over
- fast running in




# Drive concept – FP Double Head



Function	Basic machine	Standard	Option
Tube Transport	X	Servo Drive	
Main Shaft	X	Servo Drive	
Tube Infeed	X	Servo Drive	
Print Registration – Tube Turning	X	Servo Drive	
Print Registration – Tube Lifting	X	Closed Cam	
Tube Dedusting		Servo Drive	
Dosing Piston	X	Servo Drive	
Dosing – Tube Lifting	X	Servo Drive	
Rotary Valve	X	Closed Cam	Servo Drive
Shut-off Valve	X	Pneumatic	
Folding Tools		Closed Cam	
Hot-Air – Diving		Closed Cam	Servo Drive
Press	X	Closed Cam	
Trimm / Punch		Closed Cam	
Tube Discharge	X	Closed Cam	
Tube Lifting – Cartoner Transfer		Servo Drive	

#### Machine control







#### Intuitive and ergonomic user interface



Display of all relevant process parameters

Fast & comfortable change over

Counting of good/faulty tubes

**Tele-service via Internet** 

**User profiles** 

# HMI concept (Optional) 15,4" Touch Panel



# Industrial PC with diagnose/report functionality



All standard functions

History / statistic of messages

**Charge protocol** 

21 CFR Part 11 conform (audit trail)

All functionality for user administration

# Storage boxes







- Storage boxes for size parts
- Storage boxes for product-contact format parts



#### Further options



- Code scanning at print registration
- Monitoring by camera:
  - Tube edge
  - Particles
  - Tube cap presence
- Different hoppers, Accumulation cylinder (Level monitoring)
- Product contact format parts change on block
- Rotary lobe pump for CIP
- Gas purging, Droplet dosing
- Heated pressing jaws for metal tubes
- Labeling (Embossing, Laser-, Ink-jet-printing)
- Automatic central lubrication
- and more

#### Conclusion



- Outstanding transparency & accessibility
- Innovative solutions
- Customized solutions
- High flexibility
- Tool-free change over
- Ergonomic operating concept

#### Pictograms Overview





tube infeed



saddle fold



tube discharge



print registration/code scanning



heated pressing jaw



**IPC-discharge** 



dedusting



hot-air



bulk ware discharge



drip dosing



press



particle monitoring



gas purging



trimming cutters punching



roundness monitoring



filling



\*

laser labeling Inkjet labeling



cap monitoring



tube fill level monitoring



code monitoring



code reading





folding 1



cap tightener



air expeller





folding 2



faulty tube discharge



Robot tube infeed system

# Thank you for your attention!



