




Ethicon **Product Catalog**

Products by Category

At Ethicon, we develop meaningful solutions that enhance surgical outcomes and improve the lives of patients. For more than 80 years, we have listened to the voice of our customers, working in partnership with clinicians, hospitals, payors, and other healthcare stakeholders to gain deep experience and perspective into ways to redefine surgery. Our global scope, along with our broad surgical portfolio, let us provide market-appropriate, value-driven solutions that help patients around the world live longer, healthier lives.

Access

Ethicon offers you flexibility across surgical specialties, so you can choose the right trocar, with the right level of performance, when and where it is needed. The ENDOPATH® Trocar Portfolio brings value to every procedure.



3

Surgical Stapling

Our surgical stapling portfolio addresses the challenges of tissue thickness, tissue movement during firing, staple line security, and tissue tension to provide you with greater stapling precision in surgery.



48

Adjunctive Hemostats

The Ethicon portfolio of Adjunctive Hemostats is an area dedicated to aid in surgery by focusing on an array of technologies to minimize bleeding challenges. Solutions comprise innovative, biologically based products for surgical conditions that are often difficult or expensive to manage.



9

Uterine & Pelvic Surgery

Ethicon has worked with gynecologists to create innovative new solutions for a greater choice of procedure. Enhanced technique combined with innovative technologies helps to improve surgical performance.



57

Energy

With the broadest advanced energy portfolio, Ethicon provides you with unmatched choice for your laparoscopic and open surgical needs.



13

Wound Closure

Ethicon has pioneered advances in skin closure that enable you to minimize tissue trauma, optimize healing, and improve patient outcomes. We offer a portfolio of innovative wound closure solutions designed for the specific needs of skin, fascia, and organ tissues.



61

Hernia Repair & Fixation

Ethicon has been pioneering advances in hernia surgery to help restore patients back to their pre-hernia life. Our innovative value-led solutions are meant to provide strong and safe hernia repair and are designed for good long-term patient outcomes.



37

Vascular Access & Infection Control

BIOPATCH® is the ONLY protective disk with a cleared indication and PROVEN to reduce the incidence of catheter-related bloodstream infections (CRBSIs), local infections and skin colonization in patients with central venous and arterial catheters.



63

Ligation

The Ethicon portfolio of endoscopic and open ligation devices are designed to provide optimal patient surgical experiences. Our deep understanding of the interaction between surgical tools and living tissue has enabled us to design solutions for confident occlusion.



44

Other Devices

Whatever your needs during the course of surgery, Ethicon is committed to providing you a solution. From laparoscopic hand instruments to gastric banding, drains, and reservoirs, you can count us to provide you with outstanding value, excellent service, and a broad range of product offerings.



70

Products by Surgical Specialties

Bariatric Surgery

Colorectal Surgery

Gynecologic Surgery

Hernia Repair

Hepato-Biliary Surgery

Thoracic Surgery

Access

Ethicon offers you flexibility across surgical specialties, so you can choose the right trocar, with the right level of performance, when and where it is needed. The ENDOPATH® Trocar Portfolio brings value to every procedure.

Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)



Access

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENDOPATH XCEL® Bladeless Trocars with OPTIVIEW® Technology
5mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
2B5ST	5mm	75mm	Stability	6
2B5LT	5mm	100mm	Stability	6
2B5XT	5mm	150mm	Stability	6
2B12LT	12mm	100mm	Stability	6
2B12XT	12mm	150mm	Stability	6

ENDOPATH XCEL® Blunt Tip Trocar with OPTIVIEW® Technology
12mm




Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
2H12LP	12mm	100mm	Stability	6

ENDOPATH XCEL® Dilating Tip Trocar with OPTIVIEW® Technology
5mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
2D5ST	5mm	75mm	Stability	6
2D5LT	5mm	100mm	Stability	6
2D12LT	12mm	100mm	Stability	6
2D12XT	12mm	150mm	Stability	6

ENDOPATH XCEL® Universal Sleeves with OPTIVIEW® Technology
5mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
2CB5ST	5mm	75mm	Stability	6
2CB5LT	5mm	100mm	Stability	6
2CB12LT	12mm	100mm	Stability	6

Access

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENDOPATH XCEL® Bladeless Trocars

5mm, 8mm, 11mm, 12mm, 15mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
B5ST	5mm	75mm	Stability	6
B5LT	5mm	100mm	Stability	6
B5XT	5mm	150mm	Stability	6
B8LT	8mm	100mm	Stability	6
B11LT	11mm	100mm	Stability	6
B11LP	11mm	100mm	Smooth	6
B11LTH	11mm	100mm, handled	Stability	6
B11LPH	11mm	100mm, handled	Smooth	6
B12SRT	12mm	75mm	Stability	6
B12LT	12mm	100mm	Stability	6
B12LP	12mm	100mm	Smooth	6
B12LTH	12mm	100mm, handled	Stability	6
B12LPH	12mm	100mm, handled	Smooth	6
B12XT	12mm	150mm	Stability	6
B15LT	15mm	100mm	Stability	6

ENDOPATH XCEL® Dilating Tip Trocars

5mm, 11mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
D5ST	5mm	75mm	Stability	6
D5LT	5mm	100mm	Stability	6
D11LT	11mm	100mm	Stability	6
D12LT	12mm	100mm	Stability	6
D12XT	12mm	150mm	Stability	6

Access

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENDOPATH XCEL® Blunt Tip Trocar
12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
H12LP	5mm	100mm	Smooth	6

ENDOPATH XCEL® Universal Sleeves
5mm, 11mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
CB5ST	5mm	75mm	Stability	6
CB5LT	5mm	100mm	Stability	6
CB11LT	11mm	100mm	Stability	6
CB12LT	12mm	100mm	Stability	6

ENDOPATH BASX® Bladeless Trocars
5mm, 11mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
TB5ST	5mm	75mm	Stability	6
TB5LT	5mm	100mm	Stability	6
TB11LT	11mm	100mm	Stability	6
TB12ST	12mm	75mm	Stability	6
TB12LT	12mm	150mm	Stability	6

ENDOPATH BASX® Universal Sleeves
5mm, 11mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
CTB5ST	5mm	75mm	Stability	6
CTB5LT	5mm	100mm	Stability	6
CTB11LT	11mm	100mm	Stability	6
CTB12LT	12mm	100mm	Stability	6

ENDOPATH® Mini/Micro
Trocars Nonoptical
2mm/3mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
23NBS	2mm/3mm	65mm	Smooth	6
23NBL	2mm/3mm	100mm	Smooth	6

Access

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENDOPATH® Non-Shielded Trocars
5mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
355NS	5mm	65mm	Non-shielded	6
35LNS	5mm	100mm	Non-shielded	6

ENDOPATH® Rigid Thoracic Trocar
10mm/12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
TT012	10mm/12mm	75mm	Rigid	6

FLEXIPATH® Flexible Thoracic Trocars



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
FP007	7mm	80mm	Flexible	6
FP015	15mm	80mm	Flexible	6
FP020	20mm	80mm	Flexible	6

FLEXIPATH® Flexible Thoracic Trocar Pack



Product Code	Description	EA/BX
FPK02	Contains: (1) 15mm obturator and (3) 15mm flexible sleeves	6

Access – HALS/Insufflation Needles

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENDOPATH DEXTRUS® Minimally Invasive Access System



Product Code	Description	Abdominal Wall Thickness	EA/BX
LFRO1	Fixed-length access retractor-small	<4cm	6
FLRO2	Fixed-length access retractor-medium	4-7cm	6
FLRO3	Fixed-length access retractor-larger	>7cm	6
HAPO2	Seal cap assembly with accessories	N/A	3

ENDOPATH® Insufflation Needles Pneumoperitoneum



Product Code	Shaft Length	EA/BX
PN120	120mm	12
PN150	150mm	12

ENDOPATH® Insufflation Needle Ultra Veress



Product Code	Diameter	EA/BX
UV120	120mm	12

Adjunctive Hemostat

The Ethicon portfolio of Adjunctive Hemostats is an area dedicated to aid in surgery by focusing on an array of technologies to minimize bleeding challenges. Solutions comprise innovative, biologically based products for surgical conditions that are often difficult or expensive to manage. Our Hemostasis Optimization Program is a systematic approach that enables improvement of adjunctive hemostat utilization. This program looks at solutions for surgical bleeding situations, and potentially reduces procedural costs and variability in product utilization.

Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)



Adjunctive Hemostat

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

SURGICEL SNoW® Hemostat



Product Code	Description	EA/BX
2081	1in x 2in	10
2082	2in x 4in	10
2083	4in x 4in	10

SURGICEL® FIBRILLAR™ Hemostat



Product Code	Description	EA/BX
1961	1in x 2in	10
1962	2in x 4in	10
1963	4in x 4in	10

SURGICEL® NU-KNIT™ Hemostat



Product Code	Description	EA/BX
1943	3in x 4in	24
1946	6in x 9in	10

SURGICEL® Original Hemostat



Product Code	Description	EA/BX
1951	2in x 14in	24
1952	4in x 8in	24
1953	2in x 3in	24
1955	0.5in x 2in	24

SURGIFOAM® Absorbable Gelatin Powder



Product Code	Description	EA/BX
1978	1.0g	6

Adjunctive Hemostat

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

SURGIFOAM® Absorbable Gelatin Powder Kit



Product Code	Description	EA/BX
1979	1.0g	6

SURGIFOAM® Absorbable Gelatin Sponge



Product Code	Description	EA/BX
4972	2cm x 6cm x 7mm (12-7)	12
1973	8cm x 6.25cm x 10mm (50)	4
1974	8cm x 12.5cm x 10mm (100)	6
1975	8cm x 12.5cm x 2mm (100c)	6

SURGIFOAM® Absorbable Gelatin Sponge Hemorrhoidectomy Sponge



Product Code	Description	EA/BX
1977	8in x 3in	5

INSTAT® MCH Hemostat



Product Code	Description	EA/BX
1985	1.0g	12

Adjunctive Hemostat

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

EVARREST® Fibrin Sealant Patch



Product Code	Description	EA/BX
EVT5024	EVARREST® Fibrin Sealant Patch 2"x4" (5.1cm x 10.2cm)	2

SURGIFLO® Hemostatic Matrix



Product Code	Description	EA/BX
2993	8mL SURGIFLO® Hemostatic Matrix Kit with thrombin	6
2994	SURGIFLO® Hemostatic Matrix Kit with thrombin (fully sterile)	6
2991	8mL SURGIFLO® Hemostatic Matrix	6
MS1995	34cm SURGIFLO® Endoscopic Applicator	6

EVICEL® Fibrin Sealant (Human)



Product Code	Description	EA/BX
3905	EVICEL 5ML KIT	1
3902	EVICEL 2ML KIT	1
3921S	EVICEL AIRLESS PRAY ACC US EA	1
3920	EVICEL 4CM CONTROL TIP (SET OF 10)	1
3909	EVICEL 45CM FLEXIBLE TIP (SET OF 3)	1
3908	EVICEL 35 CM APPLICATOR TIP	1
3905ASA	10 EACH ES OF 3905 AND 1 BOX OF 39215	1
3902ASA	10 EACH ES OF 3902 AND 1 BOX 39215	1
3901ASA	10 EACH ES OF 3901 AND 1 BOX OF 39215	1
3901	EVICEL 1ML KIT	1

Energy

Ethicon provides you with unmatched choice for your Energy needs in laparoscopic and open surgery. HARMONIC® Ultrasonic devices offer precision, multi-functionality and strong vessel sealing. ENSEAL® Advanced Bipolar devices bring confidence and consistency to vessel sealing and deliver versatility in every device. MEGADYNE™ delivers an expansive suite of innovative Core Energy products.

Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)



Advanced Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

HARMONIC® HD 1000i Shears



Product Code	Shaft Diameter	Shaft Length	Tip	EA/BX
HARHD20	5mm	20cm	Curved	6
HARHD36	5mm	36cm	Curved	6

HARMONIC ACE®+7 Shears



Product Code	Shaft Diameter	Shaft Length	Tip	EA/BX
HARH23	5mm	23cm	Curved	6
HARH36	5mm	36cm	Curved	6
HARH45	5mm	45cm	Curved	6

HARMONIC ACE®+ Shears



Product Code	Shaft Diameter	Shaft Length	Tip	EA/BX
HAR23	5mm	23cm	Curved	6
HAR36	5mm	36cm	Curved	6

HARMONIC FOCUS®+ Shears



Product Code	Length	Tip	EA/BX
HAR9F	9cm	Curved	6

HARMONIC FOCUS®+ Long Shears



Product Code	Length	Tip	EA/BX
HAR17F	17cm	Curved	6

Advanced Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

HARMONIC SYNERGY® Blades



Product Code	Shaft Length	Tip	EA/BX
SNGCB	Adjustable 4-9cm	Curved	6
SNGHK	Adjustable 4-9cm	Hook	6
SNGHK2	Adjustable 4-9cm	Combination Hook	6

HARMONIC® Blades



Product Code	Shaft Diameter	Shaft Length	Tip	EA/BX
HC325	5mm	23cm	Curved	6
HDH05	5mm	32cm	Hook	6
HBC05	5mm	31cm	Ball	6

HARMONIC® Hand Piece



Product Code	Compatible with	EA/BX
HP054	HAR23, HAR36, HARH23, HARH36, HARH45, HBC05, HC325, HDH05	1

HARMONIC® BLUE Hand Piece



Product Code	Compatible with	EA/BX
HPBLUE	HAR9F, HAR17F, SNGCB, SNGHK, SNGHK2	1

Advanced Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

HARMONIC® Accessories

Product Code	Description	EA/BX
HAOSSL	HARMONIC ACE® 5mm open sleeve	6
HP54CTA	HARMONIC® Hand Piece cleaning tool	1
HSA08	Disposable hand switch adapter (only use with HPO54)	3
HSTOI	Sterilization tray	1
HSTO2	Test tip for HPO54	1
TLBOI	Torque lock blade wrench	2
TTBLUE	Test tip for HPBLUE	1

ENSEAL® X1 Large Jaw



Product Code	Shaft Diameter	Shaft Length	Jaw	EA/BX
NSLX120L	13mm	20cm	Curved	6

ENSEAL® G2 Articulating Tissue Sealers



Product Code	Shaft Diameter	Shaft Length	Jaw	EA/BX
NSLG2C35A	5mm	35cm	Curved	6
NSLG2C45A	5mm	45cm	Curved	6
NSLG2S35A	5mm	35cm	Straight	6
NSLG2S45A	5mm	45cm	Straight	6

ENSEAL® G2 Curved Tissue Sealers



Product Code	Shaft Diameter	Shaft Length	Jaw	EA/BX
NSLG2C14	5mm	14cm	Curved	6
NSLG2C25	5mm	25cm	Curved	6
NSLG2C35	5mm	35cm	Curved	6
NSLG2C45	5mm	45cm	Curved	6

Advanced Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENSEAL® G2 Straight Tissue Sealers



Product Code	Shaft Diameter	Shaft Length	Jaw	EA/BX
NSLG2S14	5mm	14cm	Straight	6
NSLG2S25	5mm	25cm	Straight	6
NSLG2S35	5mm	35cm	Straight	6
NSLG2S45	5mm	45cm	Straight	6

ENSEAL® Trio Tissue Sealers



Product Code	Shaft Diameter	Shaft Length	Jaw	EA/BX
ETRIO314H	5mm	14cm	Curved	6
ETRIO325H	5mm	25cm	Curved	6
ETRIO335H	5mm	35cm	Curved	6
ETRIO345H	5mm	45cm	Curved	6

ENSEAL® Round Tip Tissue Sealers



Product Code	Shaft Diameter	Shaft Length	Jaw	EA/BX
NSEAL514RH	5mm	14cm	Straight	6
NSEAL525RH	5mm	25cm	Straight	6
NSEAL535RH	5mm	35cm	Straight	6
NSEAL545RH	5mm	45cm	Straight	6

ETHICON GEN11 Generator



Product Code	Compatible with	EA/BX
GEN11	All HARMONIC® and ENSEAL® devices	1

Advanced Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ETHICON GEN11 Generator Accessories



Product Code	Description	Compatible with	EA/BX
HGA11	HARMONIC® connector	HPO54 and HPBLUE hand pieces	1
EGA11	ENSEAL® connector	ENSEAL® Trio Tissue Sealer and ENSEAL® Round Tip Tissue Sealer	1
GEN11VK	Output verification key	GEN11	1
FSW11	Foot switch and cable	GEN11	1
CRT11	Generator cart	GEN11	1
GEN11WAR1	1-Year extended warranty		1
GEN11WAR2	2-Year extended warranty		1

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ MEGA POWER™ Generator



Product Code	Compatible with	Ships with	EA/BX
1000	Monopolar and Bipolar devices	One electrosurgical generator and one power cable	1

MEGADYNE™ MEGA POWER™ Generator accessories



Product Code	Product Description	EA/BX
1300SJ	Mega Cart with Accessory Shelf	1
1450J	Mega Power Bipolar Foot Switch	1
1400JJ	Mega Power Monopolar Foot Switch	1
96007	Monopolar Adapter	1

MEGADYNE™ ACE BLADE™ 700 Soft Tissue Dissector



Product Code	Description	Tip Length	EA/BX
ME7251C	Standard Blade C Connector	2.5"	6
ME7251E	Standard Blade E Connector	2.5"	6
ME725M1C	Modified Blade C Connector	2.5"	6
ME725M1E	Modified Blade E Connector	2.5"	6

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ E-Z CLEAN™ Coated Blade Electrodes



Product Code	Description	Tip Length	EA/BX
0012	Standard Blade Electrode	2.5"	12
0012M	Standard Blade Electrode Modified	2.5"	12
0012A	X-Long Standard Blade Electrode	2.75"	12
0012AM	X-Long Standard Blade Electrode Modified	2.75"	12
0014A	Extended Blade Electrode	4"	12
0014AM	Extended Blade Electrode Modified	4"	12
0014	Extended Blade Electrode	6.5"	12
0014M	Extended Blade Electrode Modified	6.5"	12

MEGADYNE™ E-Z CLEAN™ Coated Needle Electrodes



Product Code	Description	Tip Length	EA/BX
0013	Needle Electrode	2.75"	12
0013M	Needle Electrode Modified	2.75"	12
0016A	Extended Needle Electrode	4"	12
0016AM	Needle Electrode Modified	4"	12
0016	Extended Needle Electrode	6"	12
0016M	Extended Needle Electrode Modified	6"	12

MEGADYNE™ E-Z CLEAN™ Coated Blunt Needle Electrodes



Product Code	Description	Tip Length	EA/BX
0113	Blunt Needle Electrode	2.75"	12
0113M	Blunt Needle Electrode Modified	2.75"	12
0113A	Blunt Needle Electrode	4"	12

MEGADYNE™ E-Z CLEAN™ Coated Ball Electrodes



Product Code	Description	Tip Length	EA/BX
0015	Ball Electrode	2"	12
0009	Ball Electrode	5"	12

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ E-Z CLEAN™ Coated Specialty Electrodes



Product Code	Description	Tip Length	EA/BX
0028	Bayonet Needle	5.8"	12
0028M	Bayonet Needle Modified	5.8"	12
0029M	Bayonet Blade Modified	6.3"	12

MEGADYNE™ MEGA FINE™ Needle Electrodes



Product Code	Description	Tip Length	EA/BX
0118	Needle Electrode	2"	12
0118A	Extended Needle Electrode	2.5"	12
0119	Angled 45 Degree	2"	12
0119A	Angled 45 Degree 3mm	2"	12
0120	Angled 90 Degree	2"	12
0121	Needle Electrode	6.5"	12

MEGADYNE™ Precision Electrodes



Product Code	Description	Tip Length	EA/BX
0012AP	Blade Electrode	2.75"	12
0012AMP	Blade Electrode Modified	2.75"	12
0014AP	Blade Electrode	4"	12
0014AMP	Blade Electrode Modified	4"	12

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ E-Z CLEAN™ Protective Nose Cone Electrodes



Product Code	Description	Tip Length	EA/BX
0012MD	Nose Cone Blade	2.5"	24
0012AMD	Nose Cone Blade	2.75"	24
0014AMD	Nose Cone Blade	4"	24
0014MD	Nose Cone Blade	6.5"	24
0013MD	Nose Cone Needle	2.75"	24

MEGADYNE™ E-Z CLEAN™ Laparoscopic Electrodes



Product Code	Tip Configuration	Compatible With	Length	EA/BX
0017	Standard	Standard Pencils	33cm	6
0018	Spatula	Standard Pencils	33cm	6
0018C	Curved Spatula	Standard Pencils	33cm	6
0018CS	Curved Spatula	Reusable Foot Control Cords	33cm	6
0019	Curved	Standard Pencils	33cm	6
0019S	Curved	Reusable Foot Control Cords	33cm	6
0019L	Curved	Standard Pencils	45cm	6
0019LS	Curved	Reusable Foot Control Cords	45cm	6
0020	L-Hook	Standard Pencils	33cm	6
0020S	L-Hook	Reusable Foot Control Cords	33cm	6
0020L	L-Hook	Standard Pencils	45cm	6
0020LS	L-Hook	Reusable Foot Control Cords	45cm	6
0021	J-Hook	Standard Pencils	33cm	6
0021S	J-Hook	Reusable Foot Control Cords	33cm	6
0021L	J-Hook	Standard Pencils	45cm	6
0021LS	J-Hook	Reusable Foot Control Cords	45cm	6
0022	Needle	Standard Pencils	33cm	6
0022S	Needle	Reusable Foot Control Cords	33cm	6
0024	Ball	Standard Pencils	33cm	6
0026S	O-Wire	Reusable Foot Control Cords	33cm	6

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ E-Z CLEAN™ Laparoscopic Electrodes



Product Code	Tip Configuration	Compatible With	Length	EA/BX
0100	L-Wire	Standard Pencils	33cm	6
0100S	L-Wire	Reusable Foot Control Cords	33cm	6
0100L	L-Wire	Standard Pencils	45cm	6
0100LS	L-Wire	Reusable Foot Control Cords	45cm	6
0105S	J-Wire	Reusable Foot Control Cords	35cm	6

MEGADYNE™ Megatip Electrodes



Product Code	Tip Configuration	EA/BX
0600	L-Wire	6
0600M	L-Wire Modified	6
0605	J-Wire	6
0618	Spatula	6
0619	Curved Blade	6
0620	L-Hook	6
0620M	L-Hook Modified	6
0621	J-Hook	6
0621M	J-Hook Modified	6

MEGADYNE™ Reusable Indicator Shafts



Product Code	Description	Tip Length	EA/BX
0690	Shaft	32cm	2
0690S	Shaft with Split Stem	32cm	2
0695	Shaft	38cm	2
0695S	Shaft with Split Stem	38cm	2

MEGADYNE™ Reusable Foot Control Cord



Product Code	Description	Length	EA/BX
0075	Split stem	3m	1

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ Rocker Switch Pencil



Product Code	Description	Cord Length	EA/BX
0030	Pencil with E-Z CLEAN™ Electrode	10ft	50
0030L	Pencil with E-Z CLEAN™ Electrode	15ft	25
0030H	Pencil with E-Z CLEAN™ Electrode and holster	10ft	50
0030HL	Pencil with E-Z CLEAN™ Electrode and holster	15ft	25
0036	Pencil with E-Z CLEAN™ Electrode Modified	15ft	50
0036H	Pencil with E-Z CLEAN™ Electrode Modified and holster	15ft	50



MEGADYNE™ Button Switch Pencil



Product Code	Description	Cord Length	EA/BX
0035	Pencil with E-Z CLEAN™ Electrode	10ft	50
0035L	Pencil with E-Z CLEAN™ Electrode	15ft	25
0035H	Pencil with E-Z CLEAN™ Electrode and holster	10ft	50
0035HL	Pencil with E-Z CLEAN™ Electrode and holster	15ft	25
0037	Pencil with E-Z CLEAN™ Electrode Modified	15ft	50
0037H	Pencil with E-Z CLEAN™ Electrode Modified and holster	15ft	50



MEGADYNE™ E-Z Pen Reusable Electrosurgical Pencil



Product Code	Description	Cord Length	EA/BX
3030H	Rocker switch pencil with holster	10ft	10
3030HL	Rocker switch pencil with holster	15ft	10
3035H	Button switch pencil with holster	10ft	10
3035HL	Button switch pencil with holster	15ft	10



Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ MEGA SOFT™ Universal
Reusable Patient Return Electrode
91cm x 51cm x 0.3cm



Product Code	Description	Electrodes per box	Cables per box
0845M2K01	MEGA SOFT™ Universal with one M2K01 cable: Standard connector, 2.4m length	1	1
0845M2K02	MEGA SOFT™ Universal with one M2K02 cable: Standard connector, 4.4m length	1	1
0845M2K03	MEGA SOFT™ Universal with one M2K03 cable: Phone plug connector, 2.4m length	1	1
0845M2K04	MEGA SOFT™ Universal with one M2K04 cable: Phone plug connector, 4.4m length	1	1
0845M2K05	MEGA SOFT™ Universal with one M2K05 cable: Extended phone plug connector, 4.4m length	1	1
0845M2K06	MEGA SOFT™ Universal with one M2K06 cable: Argon beam connector, 2.4m length	1	1
0845M2K08	MEGA SOFT™ Universal with one M2K08 cable: Compatibility cable, 2.4m length	1	1
0845M2K09	MEGA SOFT™ Universal with one M2K09 cable: Compatibility cable, 4.4m length	1	1

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ MEGA SOFT™ Universal Dual Reusable Patient Return Electrode
91cm x 50cm x 0.3175cm



Product Code	Description	Electrodes per box	Cables per box
0846M2K01	MEGA SOFT™ Universal Dual with two M2K01 cables: Standard connector, 2.4m length	1	2
0846M2K02	MEGA SOFT™ Universal Dual with two M2K02 cables: Standard connector, 4.4m length	1	2
0846M2K03	MEGA SOFT™ Universal Dual with two M2K03 cables: Phone plug connector, 2.4m length	1	2
0846M2K04	MEGA SOFT™ Universal Dual with two M2K04 cables: Phone plug connector, 4.4m length	1	2
0846M2K05	MEGA SOFT™ Universal Dual with two M2K05 cables: Extended phone plug connector, 4.4m length	1	2
0846M2K06	MEGA SOFT™ Universal Dual with two M2K06 cables: Argon beam connector, 2.4m length	1	2
0846M2K08	MEGA SOFT™ Universal Dual with two M2K08 cables: Compatibility cable, 2.4m length	1	2
0846M2K09	MEGA SOFT™ Universal Dual with two M2K09 cables: Compatibility cable, 4.4m length	1	2

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ MEGA SOFT™ Reusable Patient Return Electrode
117cm x 51cm x 1.25cm



Product Code	Description	Electrodes per box	Cables per box
0830M2K01	MEGA SOFT™ with one M2K01 cable: Standard connector, 2.4m length	1	1
0830M2K02	MEGA SOFT™ with one M2K02 cable: Standard connector, 4.4m length	1	1
0830M2K03	MEGA SOFT™ with one M2K03 cable: Phone plug connector, 2.4m length	1	1
0830M2K04	MEGA SOFT™ with one M2K04 cable: Phone plug connector, 4.4m length	1	1
0830M2K05	MEGA SOFT™ with one M2K05 cable: Extended phone plug connector, 4.4m length	1	1
0830M2K06	MEGA SOFT™ with one M2K06 cable: Argon beam connector, 2.4m length	1	1
0830M2K08	MEGA SOFT™ with one M2K08 cable: Compatibility cable, 2.4m length	1	1
0830M2K09	MEGA SOFT™ with one M2K09 cable: Compatibility cable, 4.4m length	1	1

MEGADYNE™ MEGA SOFT™ Dual Reusable Patient Return Electrode
117cm x 51cm x 1.25cm



Product Code	Description	Electrodes per box	Cables per box
0835M2K01	MEGA SOFT™ Dual with two M2K01 cables: Standard connector, 2.4m length	1	2
0835M2K02	MEGA SOFT™ Dual with two M2K02 cables: Standard connector, 4.4m length	1	2
0835M2K03	MEGA SOFT™ Dual with two M2K03 cables: Phone plug connector, 2.4m length	1	2
0835M2K04	MEGA SOFT™ Dual with two M2K04 cables: Phone plug connector, 4.4m length	1	2
0835M2K05	MEGA SOFT™ Dual with two M2K05 cables: Extended phone plug connector, 4.4m length	1	2
0835M2K06	MEGA SOFT™ Dual with two M2K06 cables: Argon beam connector, 2.4m length	1	2
0835M2K08	MEGA SOFT™ Dual with two M2K08 cables: Compatibility cable, 2.4m length	1	2
0835M2K09	MEGA SOFT™ Dual with two M2K09 cables: Compatibility cable, 4.4m length	1	2

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ MEGA SOFT™ Pediatric
Reusable Patient Return Electrode
66cm x 30.5cm x 1.3cm



Product Code	Description	Electrodes per box	Cables per box
0840M2K01	MEGA SOFT™ Pediatric with one M2K01 cable: Standard connector, 2.4m length	1	1
0840M2K02	MEGA SOFT™ Pediatric with one M2K02 cable: Standard connector, 4.4m length	1	1
0840M2K03	MEGA SOFT™ Pediatric with one M2K03 cable: Phone plug connector, 2.4m length	1	1
0840M2K04	MEGA SOFT™ Pediatric with one M2K04 cable: Phone plug connector, 4.4m length	1	1
0840M2K05	MEGA SOFT™ Pediatric with one M2K05 cable: Extended phone plug connector, 4.4m length	1	1
0840M2K06	MEGA SOFT™ Pediatric with one M2K06 cable: Argon beam connector, 2.4m length	1	1
0840M2K08	MEGA SOFT™ Pediatric with one M2K08 cable: Compatibility cable, 2.4m length	1	1
0840M2K09	MEGA SOFT™ Pediatric with one M2K09 cable: Compatibility cable, 4.4m length	1	1

MEGADYNE™ MEGA SOFT™ 2000
Reusable Patient Return Electrode
91cm x 51cm

Product Code	Description	Electrodes per box	Cables per box
0800M2K01	MEGA SOFT™ 2000 with one M2K01 cable: Standard connector, 2.4m length	1	1
0800M2K02	MEGA SOFT™ 2000 with one M2K02 cable: Standard connector, 4.4m length	1	1
0800M2K03	MEGA SOFT™ 2000 with one M2K03 cable: Phone plug connector, 2.4m length	1	1
0800M2K04	MEGA SOFT™ 2000 with one M2K04 cable: Phone plug connector, 4.4m length	1	1
0800M2K05	MEGA SOFT™ 2000 with one M2K05 cable: Extended phone plug connector, 4.4m length	1	1
0800M2K06	MEGA SOFT™ 2000 with one M2K06 cable: Argon beam connector, 2.4m length	1	1
0800M2K08	MEGA SOFT™ 2000 with one M2K08 cable: Compatibility cable, 2.4m length	1	1
0800M2K09	MEGA SOFT™ 2000 with one M2K09 cable: Compatibility cable, 4.4m length	1	1


Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ MEGA SOFT™ Reusable Patient Return Electrode

Product Code	Description	EA/BX
0825	MEGA 2000™ Sheath	30

MEGADYNE™ MEGA SOFT™ Cables and Accessories



Product Code	Description	Cord Length	EA/BX
M2KO1	Standard Connector	8ft	1
M2KO2	Standard Connector	14.4ft	1
M2KO3	Phone Plug Connector, same as B-205 adapter	8ft	1
M2KO4	Phone Plug Connector, same as B-205 adapter	14.4ft	1
M2KO5	Extended Phone Plug Connector, same as B-210 adapter	14.4ft	1
M2KO6	Argon Beam Connector, same as A-238 adapter	8ft	1
M2KO7	Replacement Pigtail Cable	2ft	1
M2KO8	Compatibility Cable	8ft	1
M2KO9	Compatibility Cable	14.4ft	1
PKIT001	MEGA SOFT™ Repair Kit		1

MEGADYNE™ Mega Vac Smoke Evacuators and Accessories



Product Code	Description	Ships with	EA/BX
2200J	MEGADYNE™ Mega Vac PLUS Smoke Evacuator with Laparoscopic Mode	One smoke evacuator, one power cable, one charcoal filter (2220), and one RF sensor (2250)	1
2100J	MEGADYNE™ Mega Vac Smoke Evacuator	One smoke evacuator, one power cable, one charcoal filter (2220), and one RF sensor (2250)	1
2211J	ULPA Filter with Fluid Trap		10
2210J	ULPA Filter with Fluid Trap		9
2220J	Charcoal Filter		1

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ Mega Vac Smoke Evacuators and Accessories



Product Code	Description	Ships with	EA/BX
2250J	MEGADYNE™ RF Sensor for Mega Vac		1
2160J	MEGADYNE™ Mega Vac Smoke Evacuator Footswitch		1

MEGADYNE™ MINI VAC™ Smoke Evacuators and Accessories




Product Code	Description	Ships with	EA/BX
ECVV120	MINI VAC™ Smoke Evacuator 120v	One smoke evacuator, one power cable, and one MicroSafe Filter (MGVS353)	1
ECVV220	MINI VAC™ Smoke Evacuator 220v	One smoke evacuator, one power cable, and one MicroSafe Filter (MGVS353)	1
MGEZLINK01	MEGADYNE™ RF Sensor for MINI VAC™		
MGVS353	MEGADYNE™ Micro Safe Filter 1 Pack		1
MGVS35302	MEGADYNE™ Micro Safe Filter 2 Pack		2
MGVSFT10	Fluid Trap for MINI VAC™		10

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

Smoke Evacuation Pencils and Accessories



Product Code	Description	EA/BX
252510	MEGADYNE™ ZIP-PEN™ Smoke Evacuation Pencil	20
252510EC	MEGADYNE™ ZIP-PEN™ Smoke Evacuation Pencil with 22mm connector	20
2540J	MEGADYNE™ ZIP-PEN™ Extension Nozzle 2.7"	10
2560J	MEGADYNE™ ZIP-PEN™ Extension Nozzle 5.2"	10
211010J	MEGADYNE™ ULTRA VAC™ Smoke Evacuation Pencil	20
211010EC	MEGADYNE™ ULTRA VAC™ Smoke Evacuation Pencil with 22mm connector	20
212009J	MEGADYNE™ Attacha Vac Smoke Evacuation Shroud	25
212009EC	MEGADYNE™ Attacha Vac Smoke Evacuation Shroud with 22mm connector	25
2190J	Laparoscopic Tubing	25
2195J	Speculum Tubing	25
2140J	MEGADYNE™ Mega Vac Smoke Evacuator Connector	1
2145J	22mm male to 10mm male Connector	1
2150J	Universal Connector	1
2151J	22mm ID Tube Connector	1
2155	22mm female to 10mm male Connector	1

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ Stainless Steel Blade Electrodes



Product Code	Description	Tip Length	EA/BX
0312	Blade Electrode	2.5"	12
0312M	Blade Electrode Modified	2.5"	12
0312A	Blade Electrode	2.75"	12
0312AM	Blade Electrode Modified	2.75"	12
0314A	Blade Electrode	4"	12
0314AM	Blade Electrode Modified	4"	12
0314	Blade Electrode	6.5"	12
0314M	Blade Electrode Modified	6.5"	12

MEGADYNE™ Stainless Steel Needle Electrodes



Product Code	Description	Tip Length	EA/BX
0313	Needle Electrode	2.75"	12
0313M	Needle Electrode Modified	2.75"	12
0316AM	Needle Electrode Modified	4"	12
0316	Needle Electrode	6"	12
0316M	Needle Electrode Modified	6"	12

MEGADYNE™ Stainless Steel Ball Electrodes



Product Code	Description	Tip Length	EA/BX
0300	Ball Electrode	5"	12
0315	Ball Electrode	2"	12

MEGADYNE™ Pencils With Stainless Steel Electrodes



Product Code	Description	EA/BX
0038	Rocker switch pencil	50
0038H	Rocker switch pencil with holster	50
0039	Button switch pencil	50
0039H	Button switch pencil with holster	50
FC39H	Foot control pencil with holster	50
TC100	Tip Cleaner	50

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ Disposable Patient Return Electrodes and Accessories



Product Code	Description	Cord Length	EA/BX
0850C	Single plate return electrode	9ft	50
0855C	Dual plate return electrode	10ft	50
0855CL	Dual plate return electrode	15ft	50
0855CN	Dual plate return electrode	9ft	50
0865C	Pediatric dual plate return electrode	10ft	50
0855	Dual plate return electrode	No Cord	50
0870	Return Electrode Cable, Single Style	10ft	1
0875	Return Electrode Cable, Dual Style	10ft	1

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENDOPATH® Electrosurgery PROBE PLUS® II System



Product Code	Grip	Control	EA/BX
EPHO2	Pistol	Hand	6

EPHO4	Pencil	Hand	6
-----------------------	--------	------	---

ENDOPATH® Electrosurgery PROBE PLUS® II Electrode Shafts



Product Code	Shaft Diameter	Shaft Length	Electrode	EA/BX
EPSO1	5mm	34cm	Hook	6
EPSO2	5mm	34cm	Spatula	6
EPSO3	5mm	34cm	Right angle	6
EPSO4	5mm	34cm	Curved dissector	6
EPSO5	5mm	29cm	Hook	6
EPSO6	5mm	29cm	Spatula	6
EPSO7	5mm	29cm	Right angle	6

ENDOPATH® Electrosurgery PROBE PLUS® II Tubing



Product Code	Description	EA/BX
EPTO1	PROBE PLUS® II, dual-bottle irrigation tube set	6

EPTO3	PROBE PLUS® II, pressurized bag irrigation tube set with spike	6
-----------------------	--	---

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ All-In-One Hand Control Unit



Product Code	Description	EA/BX
0055	All-In-One Hand Control Unit	5

MEGADYNE™ All-In-One Electrodes



Product Code	Tip Configuration	Length	EA/BX
0060	L-Hook	32cm	5
0085	J-Wire Hook	38cm	5
0062	Curved Blade	32cm	5
0063	L-Wire Hook	32cm	5
0080J	L-Wire Hook	38cm	5
0064	Standard Blade	32cm	5

MEGADYNE™ All-In-One Electrosurgical Tip



Product Code	Tip Configuration	Tip Length	EA/BX
0066	Standard Blade	2.5"	5

MEGADYNE™ Suction Coagulators



Product Code	Description	French Size	EA/BX
004010J	Foot control, 6" cannula, 10 ft cord	10 Fr	10
004025	Foot control, 6" cannula, 10 ft cord	10 Fr	25
004125	Hand control, 6" cannula, 10 ft cord	10 Fr	25
004225	Hand control, 6" cannula, 10 ft cord	8 Fr	25
004325	Hand control, 6" cannula, 10 ft cord	12 Fr	25

Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ LLETZ Loop Electrodes Square



Product Code	Description	Color	EA/BX
0410	LLETZ Loop 10 x 7 mm	Orange	12

MEGADYNE™ LLETZ Loop Electrodes Round



Product Code	Description	Color	EA/BX
0420	LLETZ Loop 10 x 5 mm	Purple	12
0430	LLETZ Loop 15 x 6 mm	Gray	12
0440	LLETZ Loop 20 x 10 mm	Red	12
0450	LLETZ Loop 10 x 10 mm	Yellow	12
0460	LLETZ Loop 15 x 12 mm	Green	12
0470	LLETZ Loop 20 x 12 mm	White	12
0480	LLETZ Loop 20 x 15 mm	Blue	12
0490	LLETZ Loop 20 x 8 mm	Tan	12

MEGADYNE™ Bipolar Forceps



Product Code	Description	Tip	Length	EA/BX
4010J	Reusable Scoville-Greenwood Forcep	1.5mm	7"	1
4020J	Reusable Jeweler Forcep	0.4mm	4"	1
4030J	Reusable Cushing Forcep	1.5mm	6"	1
4040J	Reusable Adson Forcep	1mm	4"	1
4050J	Reusable Simkin Forcep	0.5mm	5"	1
4060J	Reusable Gerald Forcep	1mm	7"	1
4000J	Disposable Flying Lead Bipolar Cable		12'	10
4005J	Disposable Fixed Lead Bipolar Cable		12'	10

Hernia Repair & Fixation

Ethicon has been pioneering advances in hernia surgery to help restore patients back to their pre-hernia life. The Ethicon portfolio offers choice to help ensure you have the right solutions for the right procedure. Our innovative value-led solutions are meant to provide strong and safe hernia repair and are designed for good long-term patient outcomes.

Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)



Hernia Repair & Fixation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

PROLENE® Polypropylene Mesh



Product Code	Description	Size	EA/BX
PMH	Square	6" x 6" (15cm x 15cm)	6
PMII	Rectangle	3" x 6" (8cm x 15cm)	6
PML	Square	12" x 12" (30cm x 30cm)	3
PMLK	Large pre-shaped keyhole	2.4" x 5.4" (6cm x 13.7cm)	6
PMSK	Small pre-shaped keyhole	1.8" x 4.0" (4.6cm x 10cm)	6
PMXL	Rectangle	2" x 12" (5cm x 30cm)	6
PMXS	Rectangle	1" x 4" (2.54cm x 10cm)	6

PROLENE® Soft Polypropylene Mesh



Product Code	Description	Size	EA/BX
SPMXS	Rectangle	1" x 4" (2.54cm x 10cm)	6
SPMS	Rectangle	2" x 4" (5cm x 10cm)	6
SPMII	Rectangle	3" x 6" (7.6cm x 15cm)	6
SPMH	Square	6" x 6" (15cm x 15cm)	6
SPMLI	Square	10" x 10" (25cm x 25cm)	3
SPMXXL	Rectangle	14" x 12" (35.6cm x 30cm)	3
SPM3XL	Square	19.6" x 19.6" (50 cm x 50 cm)	3

ULTRAPRO® Partially Absorbable Lightweight Mesh




Product Code	Description	Size	EA/BX
UMS3	Rectangle	2.4" x 4.3" (6cm x 11cm)	3
UMR3	Rectangle	3" x 6" (7.6cm x 15cm)	3
UMM3	Square	6" x 6" (15cm x 15cm)	3
UML1	Square	12" x 12" (30cm x 30cm)	1

Hernia Repair & Fixation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ULTRAPRO ADVANCED™



Product Code	Description	Size	EA/BX
UPA3612	Rectangle	6cm x 12cm	3
UPA37615	Rectangle	7.6cm x 15cm	3
UPA31015	Rectangle	10cm x 15cm	3
UPA31515	Square	15cm x 15cm	3
UPA1530	Rectangle	15cm x 30cm	1
UPA3030	Square	30cm x 30cm	1

ULTRAPRO COMFORT PLUG™ Partially Absorbable Hernia Repair Device



Product Code	Plug Diameter	Onlay Mesh Size	EA/BX
UPLUG401	40mm	7cm x 14cm	1
UPLUG403	40mm	7cm x 14cm	3
UPLUG406	40mm	7cm x 14cm	6
UPLUG551	55mm	7cm x 14cm	1
UPLUG553	55mm	7cm x 14cm	3
UPLUG556	55mm	7cm x 14cm	6

PROLENE® Polypropylene Hernia System



Product Code	Description	Underlay Patch Size	Onlay Patch Size	Connector Size	EA/BX
PHSM	Medium	7cm	4.5cm x 10cm	1.9cm x 1.3cm	1
PSHM6	Medium	7cm	4.5cm x 10cm	1.9cm x 1.3cm	3
PHSL	Large	10cm	4.5cm x 10cm	1.9cm x 1.3cm	6
PHSL6	Large	10cm	4.5cm x 10cm	1.9cm x 1.3cm	1
PHSE	Extended	10cm	5.5cm x 12.8cm	1.9cm x 1.3cm	3
PHSE6	Extended	10cm	5.5cm x 12.8cm	1.9cm x 1.3cm	6

Hernia Repair & Fixation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ULTRAPRO® Hernia System



Product Code	Description	Underlay Patch Size	Onlay Patch Size	EA/BX
UHSL	Large	2.4" x 4.7" (6cm x 12cm)	4.0" (-10cm)	3
UHSL6	Large	2.4" x 4.7" (6cm x 12cm)	4.0" (-10cm)	6
UHSM	Medium	2.4" x 4.7" (6cm x 12cm)	3.0" (7.5cm)	3
UHSM6	Medium	2.4" x 4.7" (6cm x 12cm)	3.0" (7.5cm)	6
UHSOV	Oval	2.4" x 4.7" (6cm x 12cm)	4.0" x 4.7" (-10cm x 12cm)	3

PROCEED® Surgical Mesh



Product Code	Description	Size	EA/BX
PCDB1	Rectangle	2" x 4" (5cm x 10cm)	1
PCDD1	Rectangle	4" x 8" (10cm x 20cm)	1
PCDG1	Oval	6" x 8" (15cm x 20cm)	1
PCDH1	Oval	8" x 10" (20cm x 25cm)	1
PCDJ1	Rectangle	8" x 12" (20cm x 30cm)	1
PCDL1	Square	12" x 12" (30.5cm x 30.5cm)	1
PCDM1	Square	6" x 6" (15cm x 15cm)	1
PCDN1	Oval	4" x 6" (10cm x 15cm)	1
PCDR1	Rectangle	3" x 6" C 7 .. 5cm x 15cm)	1
PCDT1	Oval	10" x 13" (26cm x 34cm)	1
PCDW1	Rectangle	10" x 14" (25cm x 35cm)	1

PROCEED® Ventral Patch




Product Code	Description	Size	EA/BX
PVPS	Circle	1.7" x 1.7" (4.3cm x 4.3cm)	2
PVPM	Circle	2.5" x 2.5" (6.4cm x 6.4cm)	2

Hernia Repair & Fixation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ETH ICON PHYSIOMESH™ Open Flexible Composite Mesh Device



Product Code	Description	Size	EA/BX
OPHY1010	≤3cm	10cM "ROUND"	1
OPHY1215	3cm-5cm	12cmx15cm	1
OPHY1520	3cm-≤6cm	15cmx20cm	1
OPHY1525	<3cm-≤6cm	15cmx25cm	1
OPHY2030	6cm-10cm	20cmx30cm	1
OPHY2536	≥10cm-≤19cm	25cmx36cm	1

FlexHD® Structural Acellular Hydrated Dermis (allograft)



	Tissue Code	Thickness	Dimension
Thick (0.8mm-1.7mm)	471608	0.4mm - 0.8mm	6cm x 8cm
	471812	0.8mm -1.7mm	8cm x 12cm
	471016	0.8mm -1.7mm	10cm x 16cm
	471216	0.8mm -1.7mm	12cm x 16cm
	471620	0.8mm -1.7mm	16cm x 20cm
	471202	0.8mm -1.7mm	20cm x 20cm
	471224	0.8mm -1.7mm	12cm x 24cm
	471225	0.8mm -1.7mm	20cm x 25cm
Ultra Thick (1.8mm-4mm)	471230	0.8mm -1.7mm	20cm x 30cm
	472812	1.8mm-4mm	8cm x 12cm
	472016	1.8mm-4mm	10cm x 16cm
	472216	1.8mm-4mm	12cm x 16cm
	472620	1.8mm-4mm	16cm x 20cm
	472224	1.8mm-4mm	12cm x 24cm
	472225	1.8mm-4mm	20cm x 25cm
	472230	1.8mm-4mm	20cm x 30cm

Hernia Repair & Fixation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

FlexHD® Diamond XL Acellular Hydrated Dermis (allograft)



Product Code	Thickness	Dimensions	EA/BX
4D1331	0.8mm-1.7mm	22cm x 30cm	1
4D1335	0.8mm-1.7mm	24cm x 35cm	1
4D2331	1.8mm-4mm	22cm x 30cm	1
4D2335	1.8mm-4mm	24cm x 35cm	1

XCM BIOLOGIC® Tissue Matrix sterile (hydrated), Thick Average 1.5mm thick= 0.3mm



Product Code	Dimensions	Product Code	EA/BX
XM106.0612S	6cm x 12cm	XM106.1620S	16cm x 20cm
XM106.0616S	6cm x 16cm	XM106.1625S	16cm x 25cm
XM106.0812S	8cm x 12cm	XM106.2020S	20cm x 20cm
XM106.0816S	8cm x 16cm	XM106.2025S	20cm x 25cm
XM106.0820S	8cm x 20cm	XM106.2030S	20cm x 30cm
XM106.1016S	10cm x 16cm	XM106.2040S	20cm x 40cm
XM106.1020S	10cm x 20cm	XM106.2540S	25cm x 40cm
XM106.1220S	12cm x 20cm	XM106.2545S	25cm x45cm

XCM BIOLOGIC® Tissue Matrix, Tailored Coverage™, Thick



Product Code	Dimensions
XM106.2535TCS	25cm x 35cm

Hernia Repair & Fixation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ETHICON SECURESTRAP®
Absorbable Strap Fixation
Device



Product Code	Description	EA/BX
STRAP12	5mm single use device with 12 absorbable straps	6
STRAP25	5mm single use device with 25 absorbable straps	6

ETH ICON SECURESTRAP® Open
Absorbable Strap Fixation Device



Product Code	Description	EA/BX
OPENSTRAP20	Single use device with 20 absorbable straps	6

Ligation

The Ethicon portfolio of endoscopic and open ligation devices are designed to provide optimal patient surgical experiences. Our deep understanding of the interaction between surgical tools and living tissue has enabled us to design solutions for confident occlusion.

Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)



Ligation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

LIGAMAX™ 5 Endoscopic Clip Applier



Product Code	Shaft Diameter	Clip Size	Number of clips	EA/BX
EL5ML	5mm	Medium/Large	15	3

10mm LIGACLIP® Endoscopic Rotating Multiple-Clip Applier



Product Code	Shaft Diameter	Clip Size	Number of clips	EA/BX
ER320	10mm	Medium/Large	20	3

12mm LIGACLIP® Endoscopic Rotating Multiple-Clip Applier



Product Code	Shaft Diameter	Clip Size	Number of clips	EA/BX
ER420	12mm	Large	20	3

LIGACLIP® Endoscopic Clip Applier



Product Code	Clip Size	EA/BX
EL214	Medium	1
EL314	Medium/Large	1
EL414	Large	1

Ligation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

LIGACLIP® Multi-Clip Appliers (MCA)



Product Code	Shaft Length	Clip Size	Number of clips	EA/BX
MCS20	9 3/8"	Small	20	6
MSM2	9 3/8"	Medium	20	6
MCM20	11"	Medium	20	6
MCM30	11"	Medium	30	6
MCL20	13"	Large	20	6

LIGACLIP® Multi-Patient Single-Clip Applier - Blue



Product Code	Length	Handle Color	Clip Size	Specification	EA/BX
LX105	5 3/4"	Blue	Small	-	1
LX107	7 1/2"	Blue	Small	-	1
LX110	10 1/2"	Blue	Small	-	1
LX130	10 1/2"	Blue	Small	Right angle	1

LIGACLIP® Multi-Patient Single-Clip Applier - Silver



Product Code	Length	Handle Color	Clip Size	Specification	EA/BX
LX205	5 3/4"	Silver	Medium	-	1
LX207	7 1/2"	Silver	Medium	-	1
LX210	10 1/2"	Silver	Medium	-	1
LX220	10 1/2"	Silver	Medium	Angled jaw	1
LX230	10 1/2"	Silver	Medium	Right jaw	1

LIGACLIP® Multi-Patient Single-Clip Applier - Green



Product Code	Length	Handle Color	Clip Size	Specification	EA/BX
LC307	7 1/2"	Green	Medium/Large	-	1
LC310	10 1/2"	Green	Medium/Large	-	1
LC320	10 1/2"	Green	Medium/Large	-	1

Ligation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

LIGACLIP® Multi-Patient
Single-Clip Applier - Yellow



Product Code	Length	Handle Color	Clip Size	Specification	EA/BX
LC407	7 1/2"	Yellow	Large	-	1
LC410	10 1/2"	Yellow	Large	-	1
LC420	10 1/2"	Yellow	Large	Angled jaw	1
LC430	10 1/2"	Yellow	Large	Right angle	1

LIGACLIP® Ligating Clip
Cartridge Base



Product Code	Description	EA/BX
LC800	LIGACLIP® Ligating Clip Cartridge Base	1

LIGACLIP® EXTRA Ligating Clips
- For use with open and endoscopic
single clip appliers



Product Code	Length	Handle Color	Clip Size	Specification	EA/BX
LT100	Open single	Blue	Small	6	36
LT102	Open single	Blue	Small	20	15
LT200	Open/endo	Silver	Medium	6	36
LT202	Open/endo	Silver	Medium	20	15
LT300	Open/endo	Green	Medium/ Large	6	18
LT400	Open/endo	Yellow	Large	6	18

Surgical Stapling

Reliable tissue repair is a critical factor for all surgery. That's why Ethicon has worked with surgeons worldwide to develop and refine our surgical staplers and cutters. Our surgical stapling portfolio addresses the issues of tissue thickness, staple line security, and tissue tension to provide you solutions that are designed to enhance control and patient safety.

Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)



Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ECEHLON FLEX™ Powered Plus
60mm Stapler



Product Code	Shaft Length	Maximum Firings	EA/BX
PCEE60A	280mm	12	3
PSEE60A	340mm	12	3
PLEE60A	440mm	12	3

ECHELON FLEX™ Powered
ENDOPATH® 45mm Stapler



Product Code	Shaft Length	Maximum Firings	EA/BX
PCE45A	280mm	12	3
PSE45A	340mm	12	3
PLE45A	440mm	12	3

ECHELON FLEX™ Powered Vascular
Stapler with Advanced Placement Tip



Product Code	Shaft Length	Maximum Firings	EA/BX
PVE35A	320mm	12	3

ECHELON FLEX™ ENDOPATH®
60mm Stapler



Product Code	Shaft Length	Articulation	Maximum Firings	EA/BX
SC60A	280mm	Yes	12	3
EC60A	340mm	Yes	12	3
LONG60A	440mm	Yes	12	3

Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ECHELON FLEX™ ENDOPATH®
45mm Stapler



Product Code	Shaft Length	Articulation	Maximum Firings	EA/ BX
SC45A	280mm	Yes	12	3
EC45A	340mm	Yes	12	3
EC45AL	440mm	Yes	12	3

ECHELON 60mm Reloads
with Gripping Surface Technology



Product Code	Color Cartridge	Number of Staples	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/ BX
GST60W	White	6	Mesentery/ thin	2.6mm	1.0mm	12
GST60B	Blue	6	Regular	3.6mm	1.5mm	12
GST60D	Gold	6	Regular/ thick	3.8mm	1.8mm	12
GST60G	Green	6	Thick	4.1mm	2.0mm	12
GST60T	Black	6	Very thick	4.2mm	2.3mm	12

ECHELON 60mm Reloads
- For use with all ECHELON
60mm Endoscopic Staplers



Product Code	Color Cartridge	Number of Staples	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/ BX
ECR60M	Gray	6	Mesentery/ thin	2.0mm	.75mm	12
ECR60B	Blue	6	Regular	3.5mm	1.5mm	12
ECR60D	Gold	6	Regular/ thick	3.8mm	1.8mm	12
ECR60G	Green	6	Thick	4.1mm	2.0mm	12

Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ECHELON 45mm Reloads
- For use with all ECHELON 4
5mm Endoscopic Staplers



Product Code	Color Cartridge	Number of Staples	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/ BX
ECR45M	Gray	6	Mesentery/ thin	2.0mm	.75mm	12
ECR45W	White	6	Vascular/ thin	2.5mm	1.0mm	12
ECR45B	Blue	6	Regular	3.5mm	1.5mm	12
ECR45D	Gold	6	Regular/ thick	3.8mm	1.8mm	12
ECR45G	Green	6	Thick	4.1mm	2.0mm	12

ENDOPATH ECHELON™ Vascular
White Reload for Advanced
Placement Tip - 35mm



Product Code	Color Cartridge	Number of Staples	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/ BX
VASECR35	White	4	Vascular/thin	2.5mm	1.0mm	12

ENDOPATH® ETS 45mm Articulating
Linear Cutter



Product Code	Shaft Length	Articulation	Reload Included	Maximum Firings	EA/ BX
ATS45	340mm	Yes	No	8	3
ATS45NK	340mm	Yes	No	8	3

Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENDOPATH® ETS 45mm Linear Cutters - For use with all ETS 45mm Endoscopic Cutters



Product Code	Color Cartridge	Number of Staples	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/ BX
6R45M	Gray	6	Mesentery/ thin	2.0mm	.85mm	12
TR45W	White	6	Vascular/ thin	2.5mm	1.0mm	12
6R45B	Blue	6	Regular	3.5mm	1.5mm	12
TR45G	Green	4	Thick	4.1mm	2.0mm	12

PROXIMATE® Curved Intraluminal Stapler



Product Code	Shaft Length	Adjustable Closed Staple Height	Head Diameter	Lumen	Open leg Length	EA/ BX
CDH21A	18cm	1.0-2.5mm	21mm	12.4mm	5.5mm	3
CDH25A	18cm	1.0-2.5mm	25mm	16.4mm	5.5mm	3
CDH29A	18cm	1.0-2.5mm	29mm	20.4mm	5.5mm	3
CDH33A	18cm	1.0-2.5mm	33mm	24.4mm	5.5mm	3

PROXIMATE® Curved Intraluminal Stapler



Product Code	Shaft Length	Adjustable Closed Staple Height	Head Diameter	Lumen	Open leg Length	EA/ BX
ECS21A	28cm	1.0-2.5mm	21mm	12.4mm	5.5mm	3
ECS25A	28cm	1.0-2.5mm	25mm	16.4mm	5.5mm	3
ECS29A	28cm	1.0-2.5mm	29mm	20.4mm	5.5mm	3
ECS33A	28cm	1.0-2.5mm	33mm	24.4mm	5.5mm	3

Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

PROXIMATE® ILS Circular Sizer Set



Product Code	Description	EA/BX
CSS	Bowel sizers for the ENDOPATH® ILS Curved Staplers, the set includes one each of a 25mm, 29mm, and 33mm sizer	1

CONTOUR® Curved Cutter Stapler



Product Code	Cartridge Color	Number of Staple Rows	Tissue Thickness	Open Staple Height	Closed Staple Height	Reload Included	Maximum Firings	EA/BX
CS40B	Blue	Regular	3.5mm	3.5mm	1.44mm	Yes	6	3
CS40G	Green	Thick	4.7mm	4.7mm	2.0mm	Yes	6	3

CONTOUR® Curved Cutter Stapler Cartridge



Reload Code	Cartridge Color	Number of Staple Rows	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/BX
CR40B	Blue	4	Regular	3.5mm	1.44mm	3
CR40G	Green	4	Thick	4.7mm	2.0mm	3

ETHICON ENDO-SURGERY™ Linear Cutter



Product Code	Shaft Length	Articulation	Reload Included	Maximum Firings	EA/BX
NTLC55	55mm	6	1.5mm, 1.8mm, 2.0mm	12	3
NTLC75	75mm	6	1.5mm, 1.8mm, 2.0mm	12	3

ETHICON ENDO-SURGERY™ Linear Cutter Cartridge



Product Code	Number of Stapled Rows	Selectable Closed Staple Height	EA/BX
SR55	6	1.5mm, 1.8mm, 2.0mm	12
SR75	6	1.5mm, 1.8mm, 2.0mm	12

Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

PROXIMATE® Linear Cutter



Product Code	Cartridge Color	Number of Staple Rows	Tissue Thickness	Open Staple Height	Closed Staple Height	Reload Included	Maximum Firings	EA/BX
TVC55	White	4	Vascular/thin	250mm	1.0mm	Yes	8	3
TLC55	Blue	4	Regular	385mm	1.5mm	Yes	8	3
TCT55	Green	4	Thick	450mm	2.0mm	Yes	8	3
TLC75	Blue	4	Regular	385mm	1.5mm	Yes	8	3
TCD75	Gold	4	Regular/thick	420mm	1.8mm	Yes	8	3
TCT75	Green	4	Thick	450mm	2.0mm	Yes	8	3
TLC10	Blue	4	Regular	385mm	1.5mm	Yes	8	3
TCT10	Green	4	Thick	450mm	2.0mm	Yes	8	3

PROXIMATE® Linear Cutter Cartridge



Reload Code	Cartridge Color	Number of Staple Rows	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/BX
TVR55	White	4	Vascular/thin	2.50mm	1.0mm	3
TCR55	Blue	4	Regular	3.85mm	1.5mm	3
TRT55	Green	4	Thick	4.50mm	2.0mm	3
TCR75	Blue	4	Regular	3.85mm	1.5mm	3
TRD75	Gold	4	Regular/thick	4.20mm	1.8mm	3
TRT75	Green	4	Thick	4.50mm	2.0mm	3
TCR10	Blue	4	Regular/thick	3.85mm	1.5mm	3
TRT10	Green	4	Thick	4.50mm	2.0mm	3

Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

PROXIMATE® Reloadable Staplers (TL/TLH/TLV)



Product Code	Tissue Thickness	Maximum Closed Staple Height	Reload Included	Maximum Firings	EA/BX
TL30	1.0mm-2.5mm	1.0mm	Yes	4	3
TL60	1.0mm-2.5mm	1.0mm	Yes	4	3
TL90	1.0mm-2.5mm	1.0mm	Yes	4	3
TLH30	1.0mm-2.5mm	1.0mm	Yes	4	3
TLH60	1.0mm-2.5mm	1.0mm	Yes	4	3
TLH90	1.0mm-2.5mm	1.0mm	Yes	4	3
TLV30	1.0mm	1.0mm	Yes	4	3

TL Linear Stapler with Dial



Product Code	EA/BX
TR30	3
TL6R	3
TR90	3

TLH Linear Stapler with Dial



Product Code	EA/BX
TVR30	12
TRH30	12
TRH60	12
TRH90	12

Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

PROXIMATE® Reloadable Staplers (TX)



Product Code	Cartridge Color	Number of Staple Rows	Tissue Thickness	Open Staple Height	Closed Staple Height	Reload Included	Maximum Firings	EA/ BX
TX30V	White	3	Vascular	2.5mm	1.0mm	Yes	8	3
TX30B	Blue	2	Regular	3.5mm	1.5mm	Yes	8	3
TX30G	Green	2	Thick	4.8mm	2.0mm	Yes	8	3
TX60B	Blue	2	Regular	3.85mm	1.5mm	Yes	8	3
TX60G	Gold	2	Thick	4.8mm	1.8mm	Yes	8	3

PROXIMATE® Reloadable Staplers Cartridge (TX)



Product Code	Cartridge Color	Number of Staple Rows	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/ BX
XR30V	White	3	Vascular	2.5mm	1.0mm	12
XR30B	Blue	2	Regular	3.5mm	1.5mm	12
XR30G	Green	2	Thick	4.8mm	2.0mm	12
XR60B	Blue	2	Regular	3.5mm	1.5mm	12
XR60G	Green	2	Thick	4.8mm	2.0mm	12

PROXIMATE® Reloadable Staplers Cartridge (TX)



Product Code	Number of Rows	Closed Staple Height	Knife Diameter	Number of Staples	Staple Leg Length	EA/ BX
PPHO3	2	0.75mm-1.5mm	24.4mm	28	4mm	3

Uterine & Pelvic Surgery

Ethicon has worked with gynecologists to create innovative new solutions for a greater choice of procedure. Plus, we've aligned our portfolio to the surgical treatment pathway so you can ensure you have the right solutions for the right procedure. Enhanced technique combined with innovative technologies helps to improve surgical performance.

Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)



Uterine & Pelvic Surgery

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

GYNECARE INTERCEED® Absorbable Adhesion Barrier



Product Code	Description	EA/BX
4350	GYNECARE INTERCEED® Absorbable Adhesion Barrier, 3" x 4" (7.6cm x 10.2cm)	10
4350XL	GYNECARE INTERCEED® Absorbable Adhesion Barrier XL, 5" x 6" (12.7cm x 15.2cm)	10

GYNECARE TVT ABBREVO® Continence System



Product Code	Description	EA/BX
TVTOML	GYNECARE TVT ABBREVO®	1

GYNECARE TVT EXACT® Continence System



Product Code	Description	EA/BX
TVTRL	GYNECARE TVT EXACT®	1
810061	Rigid catheter guide (reusable)	1

GYNECARE TVT™ Obturator System tension-free support for incontinence



Product Code	Description	EA/BX
810081	GYNECARE TVT™ Obturator System	1
810081L	GYNECARE TVT™ Laser Cut Obturator System	1

GYNECARE TVT™ Retropubic System tension-free support for incontinence



Product Code	Description	EA/BX
810041A	GYNECARE TVT™ with abdominal guides	1
810041B	GYNECARE TVT™ (mechanical cut mesh)	1
810041BL	GYNECARE TVT™ (laser cut mesh)	1
830041	GYNECARE TVT™ Undyed	3
830041B	GYNECARE TVT™ (mechanical cut mesh)	3
830041BL	GYNECARE TVT™ (laser cut mesh)	3
810051	Introducer handle (reusable)	1
810061	Rigid catheter guide (reusable)	1

Uterine & Pelvic Surgery

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

GYNECARE VERSAPOINT™ Bipolar Electrosurgery System



Product Code	Description	EA/BX
00481	Foot switch	1
00480	Connector cable (Handpiece)	1
04995	Universal Power Cord (U.S. and Canada)	1
01936	Hysteroscope 4.0mm 30 degree wide angle lens	1
01938	Hysteroscope 4.0mm 12 degree high magnification lens	1
01735	Resectoscope hysteroscopic bridge/visual obturator	1
01939	Resectoscope rotating outer sheath 27Fr	1
01934	Resectoscope obturator - standard	1
01931	Resectoscope working element - passive	1
00466	Ball tip electrode	5
00467	Twizzle tip electrode	5
00468	Spring tip electrode	5
01950	Bipolar electrode, resectoscopic - vaporizing tip CO degree)	5
01985	Bipolar electrode, resectoscopic - angled loop tip	5
00463	VERSAPOINT II Hysteroscopic Connector Cable	1
00464	VERSAPOINT II Hysteroscopic Footswitch	1
00465	VERSAPOINT II Hysteroscopic Electrosurgical Generator	1
01942	VERSAPOINT II Resectoscopic System- Resectoscope Electrode 4.0mm Angled Loop	5

Uterine & Pelvic Surgery

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ARTISYN® Y-Shaped Mesh



Product Code	Size	EA/BX
ARTY	27cm x 5cm	1
ARTY5L	27cm x 5cm	5

GYNECAREGYNEMESH® PS
Nonabsorbable PROLENE®
Soft Mesh



Product Code	Size	EA/BX
GPSL	10cm x 15cm	1
GPSXL3	25cm x 25cm	3

Wound Closure

Secure skin closure is an integral part of nearly every surgical procedure. Ethicon has pioneered advances in skin closure that enable you to minimize tissue trauma, optimize healing, and improve patient outcomes. We offer a portfolio of innovative wound closure solutions designed for the specific needs of skin, fascia, and organ tissues.

Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)



Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

STRATAFIX™ Spiral PDO Knotless Tissue Control Device - **Bidirectional dyed**



Cutting Edge Reverse 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 FS-2 19mm 3/8 Circle	5-0	Violet	7cm X 7cm	1	SXPDP2B421
	4-0	Violet	7cm X 7cm	1	SXPDP2B422
	4-0	Violet	14cm X 14cm	1	SXPDP2B423
 FS 26mm 3/8 Circle	3-0	Violet	14cm X 14cm	1	SXPDP2B417
	2-0	Violet	14cm X 14cm	1	SXPDP2B418
	3-0	Violet	14cm X 14cm	1	SXPDP2B419
	2-0	Violet	14cm X 14cm	1	SXPDP2B420
 CP-2 26mm 1/2 Circle	0	Violet	24cm X 24cm	1	SXPDP2B415
	1	Violet	24cm X 24cm	1	SXPDP2B416
 OS-6 36mm 1/2 Circle	1	Violet	24cm X 24cm	1	SXPDP2B200
	1	Violet	30cm X 30cm	1	SXPDP2B201
 OS-8 40mm 1/2 Circle	1	Violet	36cm X 36cm	1	SXPDP2B202

Tapercut 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 V-4 17mm 3/8 Circle	4-0	Violet	24cm X 24cm	1	SXPDP2B426
 V-26 26mm 3/8 Circle	4-0	Violet	14cm X 14cm	1	SXPDP2B425
	4-0	Violet	24cm X 24cm	1	SXPDP2B424

Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

STRATAFIX™ Spiral PDO Knotless
Tissue Control Device -
Bidirectional dyed



Taperpoint 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 CT-3 22mm 1/2 Circle	0	Violet	14cm X 14cm	1	SXPDP2B404
 SH 22mm 1/2 Circle	2-0	Violet	10cm X 10cm	1	SXPDP2B413
	2-0	Violet	14cm X 14cm	1	SXPDP2B414
 MH 36mm 1/2 Circle	2-0	Violet	7cm X 7cm	1	SXPDP2B406
	2-0	Violet	14cm X 14cm	1	SXPDP2B407
	0	Violet	14cm X 14cm	1	SXPDP2B409
	2-0	Violet	24cm X 24cm	1	SXPDP2B408
	0	Violet	24cm X 24cm	1	SXPDP2B410
	0	Violet	30cm X 30cm	1	SXPDP2B411
	2-0	Violet	36cm X 36cm	1	SXPDP2B412
 CT-1 36mm 1/2 Circle	0	Violet	14cm X 14cm	1	SXPDP2B401
	0	Violet	24cm X 24cm	1	SXPDP2B402
	0	Violet	30cm X 30cm	1	SXPDP2B403
 MO-4 36mm 1/2 Circle	0	Violet	36cm X 36cm	1	SXPDP2B400
 CTX 48mm 1/2 Circle	0	Violet	36cm X 36cm	1	SXPDP2B405

Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

STRATAFIX™ Spiral PGA-PCL
Knotless Tissue Control Device
Bidirectional undyed



Cutting Edge Reverse ▾					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 FS-2 19mm 3/8 Circle	4-0	Undyed	7cm X 7cm	1	SXMD2B405
	3-0	Undyed	7cm X 7cm	1	SXMD2B404
	4-0	Undyed	14cm X 14cm	1	SXMD2B407
	3-0	Undyed	14cm X 14cm	1	SXMD2B406
	4-0	Undyed	30cm X 30cm	1	SXMD2B409
	3-0	Undyed	30cm X 30cm	1	SXMD2B408
 FS-1 24mm 3/8 Circle	4-0	Undyed	30cm X 30cm	1	SXMD2B150
	3-0	Undyed	30cm X 30cm	1	SXMD2B410
 FS 26mm 3/8 Circle	3-0	Undyed	14cm X 14cm	1	SXMD2B411
	4-0	Undyed	30cm X 30cm	1	SXMD2B413
	3-0	Undyed	30cm X 30cm	1	SXMD2B412
 CP-2 26mm 1/2 Circle	2-0	Undyed	14cm X 14cm	1	SXMD2B414

Taperpoint ⊙					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 RB-1 17mm 1/2 Circle	4-0	Undyed	16cm X 16cm	1	SXMD2B403
	3-0	Undyed	16cm X 16cm	1	SXMD2B402
 MH 36mm 1/2 Circle	2-0	Undyed	14cm X 14cm	1	SXMD2B400
	2-0	Undyed	36cm X 36cm	1	SXMD2B401

Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

STRATAFIX™ Spiral Polypropylene
Knotless Tissue Control Device
Bidirectional undyed



Taperpoint 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 CT-3 22mm 1/2 Circle	0	Undyed	14cm X 14cm	1	SXPL2B401
 CT-1 36mm 1/2 Circle	1	Undyed	24cm X 24cm	1	SXPL2B400

STRATAFIX™ Spiral PGA-PCL
Knotless Tissue Control Device
Unidirectional undyed



Cutting Edge Reverse 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 PS-2 19mm 3/8 Circle	4-0	Undyed	30cm	1	SXMD1B105
	3-0	Undyed	30cm	1	SXMD1B103
	4-0	Undyed	60cm	1	SXMD1B106
	3-0	Undyed	60cm	1	SXMD1B104
 PS-1 24mm 3/8 Circle	3-0	Undyed	30cm	1	SXMD1B100
	3-0	Undyed	60cm	1	SXMD1B101
 PSL 30mm 3/8 Circle	2-0	Undyed	60cm	1	SXMD1B102

Taperpoint 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 RB-1 17mm 1/2 Circle	4-0	Undyed	20cm	1	SXMD1B402
	3-0	Undyed	20cm	1	SXMD1B403
	2-0	Undyed	20cm	1	SXMD1B404

Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

STRATAFIX™ Spiral PGA-PCL
Knotless Tissue Control Device
Unidirectional undyed



Taperpoint 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 SH 26mm 1/2 Circle	3-0	Undyed	20cm	1	SXMD1B405
	2-0	Undyed	20cm	1	SXMD1B406
 CT-1 36mm 1/2 Circle	3-0	Undyed	45cm	1	SXMD1B407
	2-0	Undyed	45cm	1	SXMD1B408
	0	Undyed	45cm	1	SXMD1B409

STRATAFIX™ Spiral PDO Knotless
Tissue Control Device
Unidirectional dyed



Cutting Edge Reverse 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 PS-2 19mm 3/8 Circle	4-0	Violet	45cm	1	SXPD1B101
	3-0	Violet	45cm	1	SXPD1B100

Taperpoint 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 CT-1 36mm 1/2 Circle	2-0	Violet	20cm	1	SXPD1B400
	2-0	Violet	30cm	1	SXPD1B401

STRATAFIX™ Spiral
Polypropylene Knotless
Tissue Control Device
Unidirectional undyed



Taperpoint 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 SH-1 22mm 1/2 Circle	2-0	Undyed	30cm	1	SXPL1B400

Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.



Product Code	Description	Adhesive Volume	EA/BX
DNX12	Unique formulation with an easy-to-use applicator for precise control for wide and fine line applications	0.7mL	12

DNX6	Unique formulation with an easy-to-use applicator for precise control for wide and fine line applications	0.7mL	6
----------------------	---	-------	---



Product Code	Description	Adhesive Volume	EA/BX
CLR222US	2-octylcyanoacrylate topical skin adhesive applicator. Flexible, self-adhesive polyester mesh	22cm	2

CLR602US	2-octylcyanoacrylate topical skin adhesive applicator. Flexible, self-adhesive polyester mesh tape dispenser	60cm	2
--------------------------	--	------	---




Product Code	Description	Adhesive Volume	EA/BX
DHVM12	2-octylcyanoacrylate high viscosity formulation for small incisions and lacerations	0.36mL	12



Product Code	Description	EA/BX
EC11	EN-3 Needle - O, 42" Suture	12
EC12	EN-3 Needle - 2-O, 42" Suture	12
ESS15	ENDO-HOLDER™ Cannula (for use with 5mm diameter trocar)	1

Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.



Product Code	Description	EA/BX
EXIOG	ETHIBOND® Extra Suture, ST-3 needle, O, 42"	12
JK10G	Coated VICRYL® Suture, ST-3 needle, O, 42"	12
ZK10G	PDS® II Suture, ST-3 needle, O, 42"	12

ENDOLOOP® Ligature- made from PDS (polydioxanone) suture



Product Code	Description	EA/BX
EZIOG	ENDOLOOP® Ligature made from PDS® II O 18" 12 (polydioxanone) Suture	12



Product Code	Description	EA/BX
EJIOG	ENDOLOOP® Ligature made from VICRYL® O 18" 12 (polyglactin 910) Suture	12



Product Code	Description	EA/BX
E705R	ENDOPATH® Needle Holder, reusable, 5mm, 30cm length	1
E705RL	ENDOPATH® Needle Holder, reusable, 5mm, 45cm length	1
SRNHI	Self Righting Needle Holder, reusable, 5mm, 30cm length	1
SRNHIL	Self Righting Needle Holder, reusable, 5mm, 45cm length	1



Product Code	Description	EA/BX
KA200	Sterile, one-piece clips molded from the polymer poly (p-dioxanone).	1
KA200L	Sterile, one-piece clips molded from the polymer poly (p-dioxanone).	1
XC200	Sterile, one-piece clips molded from the polymer poly (p-dioxanone).	6

Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

Suture Assistant



Product Code	Needle Type	EA/BX
SW100	N/A	4
SW110	6" ETHIBOND EXCEL® (size 0) and EEN (canoe) needle	10
SW112	6" ETHIBOND EXCEL® (size 2-0) and EEN (canoe) needle	10
SW120	6" ETHIBOND EXCEL® (size 0) and SH needle	10
SW122	6" ETHIBOND EXCEL® (size 2-0) and SH Needle	10
SW222	6" VICRYL® (size 2-0) and SH Needle	10

PROXIMATE® Rotating Head Skin Staplers



Product Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Size	EA/BX
PRR35	35	3.9	0.53	5.7	Regular	6
PRW35	35	3.9	0.58	6.9	White	6

PROXIMATE® Plus MD Skin Stapler



Product Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Size	EA/BX
PMR35	35	3.9	0.53	5.7	Regular	6
PMW35	35	3.9	0.58	6.9	Wide	6

PROXIMATE® PX Fixed Head Skin Staplers



Product Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Size	EA/BX
PXR35	35	3.9	0.53	5.7	Regular	6
PXW35	35	3.9	0.58	6.9	Wide	6

Staple Extractor




Product Code	Number of Staples	EA/BX
PSX	Squeeze handle staple extractor	12

Vascular Access & Infection Control

Infection Control has become a high priority for hospitals and throughout the continuum of care. Clinicians focus on fewer breaches of the skin and proper site maintenance which improves outcomes, lowers costs and increases patient satisfaction. Ethicon provides solutions for your vascular access and infection control needs.

Learn more about the portfolio on Ethicon.com



Indications for Use

BIOPATCH® containing Chlorhexidine Gluconate is intended for use as a hydrophilic wound dressing that is used to absorb exudate and to cover a wound caused by the use of vascular and non-vascular percutaneous medical devices such as: IV catheters, central venous lines, arterial catheters, dialysis catheters, peripherally inserted coronary catheters, mid-line catheters, drains, chest tubes, externally placed orthopedic pins, and epidural catheters. It is also intended to reduce local infections, catheter-related blood stream infections (CRBSI), and skin colonization of microorganisms commonly related to CRBSI, in patients with central venous or arterial catheters.

Vascular Access & Infection Control

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

BIOPATCH Protective Disk with CHG



Product Code	Size	French Size Range	Common Uses
4150	1" disc (2.5cm) w/4mm center hole	6-12Fr	Central Lines, PICCs

BIOPATCH Protective Disk with CHG



Product Code	Size	French Size Range	Common Uses
4151	3/4" disc (1.9cm) w/1.5mm center hole	<6Fr	Peripheral IVs, Huber Needles (ports), Arterial Lines, Extended Dwell PIVs, Midlines/ PICCS, Pins

BIOPATCH Protective Disk with CHG



Product Code	Size	French Size Range	Common Uses
4152	1" disc (1.9cm) w/2.5mm center hole	13-20Fr	Dialysis, Drians, Sheaths, Cordis, VAD drive Lines

Other Devices

Whatever your needs during the course of surgery, Ethicon is committed to providing you a solution. From laparoscopic hand instruments to gastric banding, drains, and reservoirs, you can count us to provide you with outstanding value, excellent service, and a broad range of product offerings.

Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)



Other Devices

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENDOPATH® Bipolar Forceps



Product Code	Shaft Diameter	Jaw	EA/BX
EBFO1	5mm	Macro	5
EBFO2	5mm	Micro	5

ENDOPATH® Active Cords - for use with bipolar forceps

Product Code	Description	EA/BX
EBCO2	Reusable coaxial pin connector bipolar cord, gray, 1 prong	1
EBCO5	Active return cord for bipolar instruments fixed pin connector	1

ENDOPATH® Scissors



Product Code	Description	EA/BX
5DCS	Scissors, curved, with monopolar cautery, 5mm	6

ENDOPATH® Graspers



Product Code	Description	EA/BX
10AG	Anvil grasper	6
5DSG	Grasper, 5mm	6

Other Devices

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENDOPATH® Bipolar Forceps



Product Code	Description	EA/BX
SDCD	Dissector, curved, with monopolar cautery, 5mm	6
BCD10	Dissector, blunt tip, 10mm, 3 per pouch, 12 pouches per sales unit	36
BTDO5	Dissector, blunt tip, 5mm, 3 per pouch, 12 pouches per sales unit	36

ENDOPATH® Babcocks



Product Code	Description	EA/BX
SBB	Babcock, 5mm	6
10BB	Babcock, 10mm	6

Purse String Clamp - Multi-patient use



Product Code	Description	EA/BX
EH40	Reusable, purse string clamp	1

ENDOPOUCH® Specimen Retrieval Bag



Product Code	EA/BX
POUCH	6

ENDOPATH® Electrosurgery PROBE PLUS® II System



Product Code	Grip	Control	Common Uses
EPHO2	Pistol	Hand	6
EPHO4	Pencil	Hand	6

Other Devices

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENDOPATH® Electrosurgery
PROBE PLUS® II Electrode Shafts



Product Code	Shaft Diameter	Shaft Length	Electrode	EA/BX
EPS01	5mm	34cm	Hook	6
EPS02	5mm	34cm	Spatula	6
EPS03	5mm	34cm	Right angle	6
EPS04	5mm	34cm	Curved dissector	6
EPS05	5mm	29cm	Hook	6
EPS06	5mm	29cm	Spatula	6
EPS07	5mm	29cm	Right angle	6

ENDOPATH® Electrosurgery
PROBE PLUS® II Tubing



Product Code	Description	EA/BX
EPT01	PROBE PLUS® II, dual-bottle irrigation tube set	6
EPT03	PROBE PLUS® II, pressurized bag irrigation tube set with spike	6

ETHICON ENDO-SURGERY®
Gastric Band Dissector



Product Code	Description	EA/BX
RGBRSI	Gastric Band dissector device	1

ETHICON ENDO-SURGERY®
Gastric Calibration Tube



Product Code	Description	EA/BX
RGCT180	Gastric calibration tube, 180 asymmetrical	1
RGCT360	Gastric calibration tube, 360 asymmetrical	1

Other Devices

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

BLAKE® Cardio Connectors

Product Code	Description	EA/BX
BCC1	BLAKE® Cardio Connector 1:1 - compatible with 19 and 24Fr round hubless drains	1
BCC2	BLAKE® Cardio Connector 2:1 - compatible with 19 and 24Fr round hubless drains	1
BCC3	BLAKE® Cardio Connector 3:1 - compatible with 19 and 24Fr round hubless drains	1

BLAKE® Silicone Drains - Flat



Product Code	Material	EA/BX
2210	7mm flat, 3/4 fluted	1
2211	7mm flat, full fluted	1
2212	7mm flat, full fluted, 3/16" trocar	1
2213	10mm flat, 3/4 fluted	1
2214	10mm flat, full fluted	1
2215	10mm flat, full fluted, 3/16" trocar	1
2216	7mm flat, 3/4 fluted, 3/16" trocar	1
2217	10mm flat, 3/4 fluted, 3/16" trocar	1

BLAKE® Silicone Drains - Hubless



Product Code	Material	EA/BX
2226	Hubless, 10Fr round	1
2227	Hubless, 10Fr round, 1/8" trocar	1
2228	Hubless, 15Fr round	1
2229	Hubless, 15Fr round, 3/16" trocar	1
2230	Hubless, 19Fr round	1
2231	Hubless, 19Fr round, 1/4" trocar	1
2232	Hubless, 19Fr round, 1/4" bendable trocar	1
2233	Hubless, 15Fr round, 3/16" bendable trocar	1
2234	Hubless, 24Fr round	1

Other Devices

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

VICRYL® (polyglactin 910) Knitted Mesh
VICRYL® (polyglactin 910) Woven Mesh



Product Code	Description	Size	EA/BX
VKML	Knitted	12" x 12" (30cm x 30cm)	3
VKMM	Knitted	6" x 6" (15cm x 15cm)	3
VWML	Woven	12" x 12" (30cm x 30cm)	3
VWMM	Woven	6" x 6" (15cm x 15cm)	3

VICRYL® (polyglactin 910) Mesh Bag

Product Code	Size	EA/BX
VM106K	10" x 12" (26cm x 30cm)	1
VM210K	17" x 12" (44cm x 30cm)	1
VM220K	20" x 8" (50cm x 20cm)	1

J-VAC® Reservoirs - Sterile

Product Code	Description	EA/BX
2161	150ml closed wound drainage reservoir	1
2162	450ml closed wound drainage reservoir	1
2163	300ml closed wound drainage reservoir	1

J-VAC® Bulb Suction Reservoirs
- Sterile Bulbs

Product Code	Material	EA/BX
2160	100cc bulb suction reservoir	1

J-VAC® Silicone Drains
- Drain Adapters

Product Code	Material	EA/BX
2209	J-VAC Drain Adapters 1/8"	1
2199	J-VAC Drain Adapters 1/4"	1
2298	J-VAC Drain Adapters 3/16"	1

Products by Surgical Specialty

Bariatric Surgery

Ethicon is committed to shaping the future of surgery in the fight against obesity and metabolic diseases and helping to elevate the standard of care through research and evidence, best-in-class education and training, innovating products and expanding patient access to care. Together, we're advancing surgical care, enabling patients to live longer, more fulfilling lives. For more information, [click here](#).

Ethicon offers a comprehensive portfolio of market-leading, innovative devices to meet the changing demands of today's most effective surgical treatments for severe obesity, including:



ECHELON FLEX™ GST System



HARMONIC ACE®+ 7 SHEARS with Advanced Hemostasis



STRATAFIX™ Knotless Tissue Control Devices

Colorectal Surgery

Our commitment to colorectal surgery focuses on reducing risk factors that drive variability in outcomes and costs. By combining innovation with education and research, we enable more predictability in outcomes. For more information, [click here](#).

Ethicon is addressing the most pressing challenges of today's colorectal procedures through innovations such as:



ETHICON ENDO-SURGERY® Endoscopic Curved Intraluminal Circular Stapler



ECHELON FLEX™ GST System



HARMONIC ACE®+ 7 SHEARS with Advanced Hemostasis

Products by Surgical Specialty

Gynecologic Surgery

Our goal is to collaborate with the gynecologic community to help promote positive clinical outcomes. We offer solutions to support a number of gynecological surgical procedures – both open and laparoscopic – as well as a variety of techniques. For more information, [click here](#).

Ethicon provides a broad portfolio of devices and products to help meet your clinical and economic goals, including:



HARMONIC ACE®+ 7 SHEARS with Advanced Hemostasis



SURGICEL® Family of Absorbable Hemostats



STRATAFIX™ Knotless Tissue Control Devices

Hernia Repair

Ethicon’s commitment to advancing hernia repair is reflected in our ongoing investments in evidence, innovative products, and best-in-class professional education. For more information, [click here](#).

Ethicon provides essential solutions inspired by surgeon and patient needs, including:



ETHICON PHYSIOMESH™ Open Flexible Composite Mesh Device



ETHICON SECURESTRAP® Absorbable Strap Fixation Device



ULTRAPRO COMFORT PLUG™ Partially Absorbable Hernia Repair Device

Products by Surgical Specialty

Hepato-Biliary Surgery

We are committed to developing flexible solutions for every aspect of cancer procedures – from supporting education to funding research to delivering innovations that reduce the incidence and cost of managing bleeding and infection-related complications. For more information, [click here](#).

Ethicon offers a broad portfolio that delivers surgical precision and optimal hemostasis at every critical step of complex HPB procedures. Our portfolio includes such key innovations as:



HARMONIC™ HD 1000i Shears



SURGICEL® Family of Absorbable Hemostats



STRATAFIX™ Knotless Tissue Control Devices

Thoracic Surgery

We are dedicated to procedural precision and have a renewed focus on bringing greater precision to every aspect of thoracic surgery. Through training, research, and innovative products, we provide tools to enable better outcomes in thoracic surgery. For more information, [click here](#).

Our portfolio of innovative products is designed to address the most pressing challenges of today’s thoracic procedures through innovation, including:



ECHELON FLEX™ Powered Vascular Stapler



ECHELON FLEX™ GST System



HARMONIC™ HD 1000i Shears

Essential Product Information

SURGICEL Essential Product Information

INDICATIONS

SURGICEL® Absorbable Hemostat (oxidized regenerated cellulose) is used adjunctively in surgical procedures to assist in the control of capillary, venous, and small arterial hemorrhage when ligation or other conventional methods of control are impractical or ineffective. SURGICEL® ORIGINAL, SURGICEL® FIBRILLAR™ and SURGICEL® NU-KNIT® Hemostats can be cut to size for use in endoscopic procedures.

PRECAUTIONS

Use only as much SURGICEL® Absorbable Hemostat as is necessary for hemostasis, holding it firmly in place until bleeding stops. Remove any excess before surgical closure in order to facilitate absorption and minimize the possibility of foreign body reaction.

In urological procedures, minimal amounts of SURGICEL® Absorbable Hemostat should be used and care must be exercised to prevent plugging of the urethra, ureter, or a catheter by dislodged portions of the product.

Since absorption of SURGICEL® Absorbable Hemostat could be prevented in chemically cauterized areas, its use should not be preceded by application of silver nitrate or any other escharotic chemicals.

If SURGICEL® Absorbable Hemostat is used temporarily to line the cavity of large open wounds, it should be placed so as not to overlap the skin edges. It should also be removed from open wounds by forceps or by irrigation with sterile water or saline solution after bleeding has stopped.

Precautions should be taken in otorhinolaryngologic surgery to assure that none of the material is aspirated by the patient. (Examples: controlling hemorrhage after tonsillectomy and controlling epistaxis.)

Care should be taken not to apply SURGICEL® Absorbable Hemostat too tightly when it is used as a wrap during vascular surgery (see Adverse Reactions).

ADVERSE REACTIONS

"Encapsulation" of fluid and foreign body reactions have been reported.

There have been reports of stenotic effect when SURGICEL® Absorbable Hemostat has been applied as a wrap during vascular surgery.

Paralysis and nerve damage have been reported when SURGICEL® Absorbable Hemostat was used around, in, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm.

Blindness has been reported in connection with surgical repair of a lacerated left frontal lobe when SURGICEL® Absorbable Hemostat was placed in the anterior cranial fossa.

Possible prolongation of drainage in cholecystectomies and difficulty passing urine per urethra after prostatectomy have been reported.

For more information, please consult your doctor or for product quality and technical questions, call 1-800-795-0012.

Essential Product Information

INSTAT Essential Product Information

INDICATION

INSTAT MCH is recommended for use in surgical procedures (other than Urological and Ophthalmological surgery) as an adjunct to hemostasis when control of bleeding by ligature or conventional procedures is ineffective or impractical.

PRECAUTIONS

As with other hemostatic agents, it is not recommended that INSTAT MCH be left in an infected or contaminated space.

Safety of this product has not been established in pregnant women; therefore, it should only be used when benefit to risk clearly warrants its use.

Another microfibrillar collagen agent has been reported to cause interference with the healing of skin edges when used in the closure of skin incisions, and to reduce the strength of methyl methacrylate adhesives used to attach prosthetic devices to bone surfaces.

INSTAT MCH is not intended to be used to treat systemic coagulation disorders.

The long term effects of leaving INSTAT MCH in situ are unknown.

ADVERSE EVENTS

Adverse reactions reported with another microfibrillar collagen hemostatic agent that were possibly related to its use were adhesion formation, allergic reaction, foreign body reaction, and subgaleal seroma (report of a single case). The use of microfibrillar collagen in dental extraction sockets has been reported to increase the incidence of alveolalgia.

For more information, please consult your doctor or for product quality and technical questions, call 1-800-795-0012.

Essential Product Information

SURGIFLO® Hemostatic Matrix Essential Product Information (Made from Absorbable Gelatin Sponge, USP)

DESCRIPTION

SURGIFLO® Hemostatic Matrix Kit is intended for hemostatic use by applying to a bleeding surface.

ACTIONS

When used in appropriate amounts SURGIFLO® Hemostatic Matrix is absorbed completely within 4 to 6 weeks.

INTENDED USE/INDICATIONS

SURGIFLO® Hemostatic Matrix, mixed with sterile saline or thrombin solution, is indicated in surgical procedures (other than ophthalmic) as an adjunct to hemostasis when control of bleeding by ligature or other conventional methods is ineffective or impractical.

CONTRAINDICATIONS

Do not use SURGIFLO® Hemostatic Matrix in intravascular compartments because of the risk of embolization.

Do not use SURGIFLO® Hemostatic Matrix in patients with known allergies to porcine gelatin.

Do not use SURGIFLO® Hemostatic Matrix in closure of skin incisions because it may interfere with the healing of skin edges. This interference is due to mechanical interposition of gelatin and is not secondary to intrinsic interference with wound healing.

WARNINGS

SURGIFLO® Hemostatic Matrix should not be used in the presence of infection and should be used with caution in contaminated areas of the body.

SURGIFLO® Hemostatic Matrix should not be used in instances of pumping arterial hemorrhage. SURGIFLO® Hemostatic Matrix will not act as a tampon or plug in a bleeding site.

SURGIFLO® Hemostatic Matrix should be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm because it may swell resulting in nerve damage.

Excess SURGIFLO® Hemostatic Matrix should be removed once hemostasis has been achieved.

The safety and effectiveness of SURGIFLO® Hemostatic Matrix for use in ophthalmic procedures has not been established.

SURGIFLO® Hemostatic Matrix should not be used for controlling post-partum intrauterine bleeding or menorrhagia.

The safety and effectiveness of SURGIFLO® Hemostatic Matrix has not been established in children and pregnant women.

The blue flexible applicator tip should not be trimmed to avoid exposing internal guidewire.

The white straight applicator tip should be trimmed away from the surgical area. Cut a square angle to avoid creating a sharp tip.

PRECAUTIONS

Safe and effective use of SURGIFOAM® Sponge has been reported in a published neurologic retrospective study involving 1700 cases in Europe. Safe and effective use in neurosurgery has not been proven through randomized, controlled clinical studies in the United States.

SURGIFLO® Hemostatic Matrix is supplied as a sterile product and cannot be resterilized.

SURGIFLO® Hemostatic Matrix should not be used for packing unless excess product that is not needed to maintain hemostasis is removed.

SURGIFLO® Hemostatic Matrix may swell up to 20% upon contact with additional fluid.

SURGIFLO® Hemostatic Matrix should not be used in conjunction with autologous blood salvage circuits.

SURGIFLO® Hemostatic Matrix should not be used in conjunction with methylmethacrylate adhesives.

In urological procedures, SURGIFLO® Hemostatic Matrix should not be left in the renal pelvis or ureters to eliminate the potential foci for calculus formation.

Essential Product Information

SURGIFLO® Hemostatic Matrix Essential Product Information (Made from Absorbable Gelatin Sponge, USP)

ADVERSE EVENTS

A total of 142 patients received SURGIFOAM® Sponge during a clinical trial comparing SURGIFOAM® Sponge to another absorbable gelatin sponge. In general, the following adverse events have been reported with the use of absorbable porcine gelatin-based hemostatic agents:

Gelatin-based hemostatic agents may serve as a nidus for infection and abscess formation and have been reported to potentiate bacterial growth.

Giant cell granulomas have been observed at implant sites when used in the brain.

Compression of the brain and spinal cord resulting from the accumulation of sterile fluid have been observed.

Multiple neurologic events were reported when absorbable gelatin-based hemostatic agents were used in laminectomy operations, including cauda equina syndrome, spinal stenosis, meningitis, arachnoiditis, headaches, paresthesias, pain, bladder and bowel dysfunction, and impotence.

The use of absorbable gelatin-based hemostatic agents during the repair of dural defects associated with laminectomy and craniotomy operations, has been associated with fever, infection, leg paresthesias, neck and back pain, bladder and bowel incontinence, cauda equina syndrome, neurogenic bladder, impotence, and paresis.

The use of absorbable gelatin-based hemostatic agents has been associated with paralysis, due to device migration into foramina in the bone around the spinal cord, and blindness, due to device migration in the orbit of the eye, during lobectomy, laminectomy, and repair of a frontal skull fracture and lacerated lobe.

Foreign body reactions, "encapsulation" of fluid, and hematoma have been observed at implant sites.

Excessive fibrosis and prolonged fixation of a tendon have been reported when absorbable gelatin-based sponges were used in severed tendon repair.

Toxic shock syndrome was reported in association with the use of absorbable gelatin-based hemostats in nasal surgery.

Fever, failure of absorption, and hearing loss have been observed when absorbable hemostatic agents were used during tympanoplasty.

Essential Product Information

SURGIFLO® Hemostatic Matrix Kit Essential Product Information (Made from Absorbable Gelatin Sponge, USP) with Thrombin

DESCRIPTION

SURGIFLO® with Thrombin (SURGIFLO® Hemostatic Matrix Kit) is intended for hemostatic use by applying to a bleeding surface.

ACTIONS

When used in appropriate amounts SURGIFLO® is absorbed completely within 4 to 6 weeks.

INTENDED USE/INDICATIONS

SURGIFLO®, mixed with thrombin solution, is indicated in surgical procedures (other than ophthalmic) as an adjunct to hemostasis when control of bleeding by ligature or other conventional methods is ineffective or impractical.

CONTRAINDICATIONS

Do not use SURGIFLO® in intravascular compartments because of the risk of embolization.

Do not use SURGIFLO® in patients with known allergies to porcine gelatin.

Do not use SURGIFLO® in closure of skin incisions because it may interfere with the healing of skin edges. This interference is due to mechanical interposition of gelatin and is not secondary to intrinsic interference with wound healing.

WARNINGS

SURGIFLO® should not be used in the presence of infection and should be used with caution in contaminated areas of the body.

SURGIFLO® should not be used in instances of pumping arterial hemorrhage. SURGIFLO® will not act as a tampon or plug in a bleeding site.

SURGIFLO® should be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm because it may swell resulting in nerve damage.

Excess SURGIFLO® should be removed once hemostasis has been achieved.

The safety and effectiveness of SURGIFLO® for use in ophthalmic procedures has not been established.

SURGIFLO® should not be used for controlling post-partum intrauterine bleeding or menorrhagia.

The safety and effectiveness of SURGIFLO® has not been established in children and pregnant women.

The blue flexible applicator tip should not be trimmed to avoid exposing internal guidewire.

The white straight applicator tip should be trimmed away from the surgical area. Cut a square angle to avoid creating a sharp tip.

PRECAUTIONS

Safe and effective use of SURGIFOAM® Sponge has been reported in a published neurologic retrospective study involving 1700 cases in Europe. Safe and effective use in neurosurgery has not been proven through randomized, controlled clinical studies in the United States.

SURGIFLO® is supplied as a sterile product and cannot be resterilized.

SURGIFLO® should not be used for packing unless excess product that is not needed to maintain hemostasis is removed.

SURGIFLO® may swell up to 20% upon contact with additional fluid.

SURGIFLO® should not be used in conjunction with autologous blood salvage circuits.

SURGIFLO® should not be used in conjunction with methylmethacrylate adhesives.

In urological procedures, SURGIFLO® should not be left in the renal pelvis or ureters to eliminate the potential foci for calculus formation.

Essential Product Information

SURGIFLO® Hemostatic Matrix Kit Essential Product Information (Made from Absorbable Gelatin Sponge, USP) with Thrombin

ADVERSE EVENTS

A total of 142 patients received SURGIFOAM® Sponge during a clinical trial comparing SURGIFOAM® Sponge to another absorbable gelatin sponge. In general, the following adverse events have been reported with the use of absorbable porcine gelatin-based hemostatic agents:

Gelatin-based hemostatic agents may serve as a nidus for infection and abscess formation and have been reported to potentiate bacterial growth.

Giant cell granulomas have been observed at implant sites when used in the brain.

Compression of the brain and spinal cord resulting from the accumulation of sterile fluid have been observed.

Multiple neurologic events were reported when absorbable gelatin-based hemostatic agents were used in laminectomy operations, including cauda equina syndrome, spinal stenosis, meningitis, arachnoiditis, headaches, paresthesias, pain, bladder and bowel dysfunction, and impotence.

The use of absorbable gelatin-based hemostatic agents during the repair of dural defects associated with laminectomy and craniotomy operations, has been associated with fever, infection, leg paresthesias, neck and back pain, bladder and bowel incontinence, cauda equina syndrome, neurogenic bladder, impotence, and paresis.

The use of absorbable gelatin-based hemostatic agents has been associated with paralysis, due to device migration into foramina in the bone around the spinal cord, and blindness, due to device migration in the orbit of the eye, during lobectomy, laminectomy, and repair of a frontal skull fracture and lacerated lobe.

Foreign body reactions, "encapsulation" of fluid, and hematoma have been observed at implant sites.

Excessive fibrosis and prolonged fixation of a tendon have been reported when absorbable gelatin-based sponges were used in severed tendon repair.

Toxic shock syndrome was reported in association with the use of absorbable gelatin-based hemostats in nasal surgery.

Fever, failure of absorption, and hearing loss have been observed when absorbable hemostatic agents were used during tympanoplasty.

Essential Product Information

EVITHROM® Thrombin, Topical (Human) for Topical Use Only Lyophilized Powder for Solution

EVITHROM® is a topical thrombin indicated as an aid to hemostasis whenever oozing blood and minor bleeding from capillaries and small venules is accessible and control of bleeding by standard surgical techniques (such as suture, ligature or cautery) is ineffective or impractical. EVITHROM® may be used in conjunction with an Absorbable Gelatin Sponge, USP.

IMPORTANT SAFETY INFORMATION

- For topical use only.
- Do not inject.
- Apply EVITHROM® on the surface of bleeding tissue only.
- The amount of EVITHROM® required depends upon the area of tissue to be treated and the method of application. In clinical studies, volumes up to 10 ml were used in conjunction with Absorbable Gelatin Sponge.
- Do not use for the treatment of severe or brisk arterial bleeding.
- Do not use in individuals known to have anaphylactic or severe systemic reaction to human blood products. Hypersensitivity reactions, including anaphylaxis, may occur.
- There is a potential risk of thrombosis if absorbed systemically.
- May carry a risk of transmitting infectious agents such as viruses and theoretically, the Creutzfeldt-Jakob disease (CJD) agent, despite manufacturing steps designed to reduce the risk of viral transmission.
- The most common adverse reactions during clinical trial (reported in at least 2% of subjects treated with EVITHROM®) were prolonged activated partial thromboplastin time, increased INR, decreased lymphocyte count, prolonged prothrombin time and increased neutrophil count.
- None of the patients treated with EVITHROM developed antibodies to human thrombin or to human Factor V/Va. The clinical significance of these findings is unknown.

For complete indications, contraindications, warnings, precautions, and adverse reactions, please reference full package insert.

Essential Product Information

EVARREST® Fibrin Sealant Patch **IMPORTANT SAFETY INFORMATION**

INDICATIONS AND USAGE

EVARREST® is a fibrin sealant patch indicated for use with manual compression as an adjunct to hemostasis in adult patients undergoing surgery, when control of bleeding by conventional surgical techniques (such as suture, ligation, and cautery) is ineffective or impractical.

LIMITATIONS FOR USE

- Cannot be used in place of sutures or other forms of mechanical ligation in the treatment of major arterial or venous bleeding.
- Not for use in children under one month of age.
- Laparoscopic and other minimally invasive surgeries where manual compression would be difficult to achieve.

IMPORTANT SAFETY INFORMATION

For topical use only. Apply immediate manual compression over the entire surface of the patch and maintain contact pressure for 3 minutes to control the bleeding.

Do not apply intravascularly. This can result in life threatening thromboembolic events.

Do not use to treat bleeding from large defects in arteries or veins where the injured vascular wall requires conventional surgical repair and maintenance of vessel patency or where there would be persistent exposure of EVARREST® to blood flow and/or pressure during absorption of the product. Thrombosis can occur if absorbed systemically.

Do not use in individuals known to have anaphylactic or severe systemic reaction to human blood products. EVARREST® can cause hypersensitivity reactions including anaphylaxis.

Avoid application to contaminated areas of the body or in the presence of active infection. Infection can occur.

EVARREST® contains oxidized regenerated cellulose which adheres to bleeding surfaces. Inadvertent adhesions can occur.

Avoid use in, around, or in proximity to, foramina in bone or areas of bony confines where swelling may cause compression.

Use the least number of patches required to cover the entire bleeding area. Portions of excess patch material can become dislodged and migrate to other areas of the body.

Do not use more than eight 2x4 inch (5.1 x 10.2 cm) or more than four 4x4 inch (10.2 x 10.2 cm) patches.

Use in patients who have been previously exposed to EVARREST® has not been studied.

May carry a risk of transmitting infectious agents, e.g., viruses, the variant Creutzfeldt-Jakob disease (vCJD) agent and, theoretically, the Creutzfeldt-Jakob disease (CJD) agent.

The adverse reactions reported during clinical trials occurred in less than 1% of all cases and included deep venous thrombosis, pulmonary embolism, blood fibrinogen increase, anastomotic hemorrhage, post procedural and intra-abdominal hemorrhage, abdominal distension, anemia, gastrointestinal hemorrhage, thoracic cavity drainage, pleural effusion, abdominal abscess, ascites, localized intra-abdominal fluid collection, cardiac failure, operative hemorrhage, and ischemic bowel.

Pediatrics: Use in children under the age of one month may be unsafe or ineffective due to small size and limited ability to apply the patch as recommended.

Please see package insert for EVARREST® Full Prescribing Information.

To report SUSPECTED ADVERSE REACTIONS, contact ETHICON Customer Support Center at 1-877-384-4266 or FDA at 1-800-FDA-1088 or www.fda.gov/medwatch.

Essential Product Information

EVICEL® Fibrin Sealant (Human) **IMPORTANT SAFETY INFORMATION**

INDICATIONS

EVICEL® Fibrin Sealant (Human) is indicated as an adjunct to hemostasis for use in patients undergoing surgery, when control of bleeding by standard surgical techniques (such as suture, ligation, or cautery) is ineffective or impractical.

CONTRAINDICATIONS

Do not inject directly into the circulatory system. Intravascular application of EVICEL® may result in life-threatening thromboembolic events.

Do not use in individuals known to have anaphylactic or severe systemic reaction to human blood products.

Do not use for the treatment of severe or brisk arterial bleeding.

Do not use EVICEL® for spraying in endoscopic or laparoscopic procedures where the minimum recommended distance from the applicator tip to the target site cannot be ensured.

WARNINGS AND PRECAUTIONS

Life-threatening gas embolism has occurred with the use of spray devices employing pressure regulator to administer fibrin sealants. These events appear to be related to the use of the spray device at pressures higher than recommended and/or at distances closer than recommended to the surface of the tissue. Follow labeled application instructions regarding pressure range and distance when using a spray device and monitor patients for the possibility of gas embolism.

Monitor changes in blood pressure, pulse oxygen saturation, and end-tidal CO₂ when spraying EVICEL® because of the possibility of gas embolism.

To reduce the risk of potentially life-threatening gas embolism, spray EVICEL® using only pressurized CO₂ gas at the recommended pressures and distances

Use EVICEL® spray application only if it is possible to accurately judge the spray distance, especially during endoscopic or laparoscopic procedures. Apply as a thin layer.

Prior to applying EVICEL®, dry surface areas of the wound by standard techniques (e.g. intermittent application of compresses, swabs, use of suction devices). Prepare and administer EVICEL® according to the instructions and only with devices recommended for this product.

May carry a risk of transmitting infectious agents, e.g. viruses, the variant Creutzfeldt-Jakob disease (vCJD) agent and theoretically, the Creutzfeldt-Jakob disease (CJD) agent

Most common adverse events reported in clinical trials (≥5%) are bradycardia, nausea, hypokalemia, insomnia, hypotension, pyrexia, graft infection, vascular graft occlusion, peripheral edema, and constipation.

ADVERSE REACTIONS

The type and frequency of adverse events reported are consistent with events typically seen following surgery. Postsurgical adhesions may occur in the presence of GYNECARE INTERCEED®.

For complete indications, contraindications, warnings, precautions, and adverse reactions, please reference full package insert.

Essential Product Information

GYNECARE INTERCEED® ABSORBABLE ADHESION BARRIER

INDICATIONS

GYNECARE INTERCEED® Absorbable Adhesion Barrier is indicated as an adjuvant in open Claparotomy) gynecologic pelvic surgery for reducing the incidence of postoperative pelvic adhesions after meticulous hemostasis is achieved consistent with microsurgical principles.

CONTRAINDICATIONS

The use of GYNECARE INTERCEED® is contraindicated in the presence of frank infection. GYNECARE INTERCEED® is not indicated as a hemostatic agent. Appropriate means of achieving hemostasis must be employed.

WARNINGS

The safety and effectiveness of GYNECARE INTERCEED® Adhesion Barrier in laparoscopic surgery or any procedures other than open (laparotomy) gynecologic microsurgical procedures have not been established.

Postoperative adhesions may be induced by GYNECARE INTERCEED® application if adjacent tissues (eg, ovary and tube) and structures are coapted or conjoined by the device, or if GYNECARE INTERCEED® is folded, wadded or layered. Postoperative adhesions may occur in the presence of GYNECARE INTERCEED® if meticulous hemostasis is not achieved prior to application. As with all foreign substances, GYNECARE INTERCEED® should not be placed in a contaminated surgical site.

PRECAUTIONS

Use only a single layer of GYNECARE INTERCEED®, since multiple layers of packing or folding will not enhance the adhesion barrier characteristics and may interfere with the absorption rate of GYNECARE INTERCEED®. Care should be exercised in applying GYNECARE INTERCEED® to a pelvic organ not to constrict or restrict it. If the product comes in contact with blood prior to completing the procedure, it should be discarded, as fibrin deposition cannot be removed by irrigation and may promote adhesions formation. Ectopic pregnancies have been associated with fertility surgery of the female reproductive tract. No data exist to establish the effect, if any, of GYNECARE INTERCEED® on the occurrence of ectopic pregnancies. No adequate studies have been conducted in women who have become pregnant within the first month after exposure to GYNECARE INTERCEED®. No teratogenic studies have been performed. Therefore, avoidance of conception should be considered during the first complete menstrual cycle after use of GYNECARE INTERCEED®. The safety and effectiveness of using GYNECARE INTERCEED® in combination with other adhesion prevention treatments have not been clinically established. GYNECARE INTERCEED® is supplied sterile. As the material is not compatible with autoclaving or ethylene oxide sterilization, GYNECARE INTERCEED® must not be resterilized. Foreign body reactions may occur in some patients. Interactions may occur between GYNECARE INTERCEED® and some drugs used at the surgical site. Pathologists examining sites of GYNECARE INTERCEED® placement should be made aware of its usage and of the normal cellular response to GYNECARE INTERCEED® 'to facilitate proper evaluation of specimens'.

ADVERSE REACTIONS

The type and frequency of adverse events reported are consistent with events typically seen following surgery. Postsurgical adhesions may occur in the presence of GYNECARE INTERCEED®.