



## MIDMARK 75L PROCEDURE TABLE

# FEATURES

- \* Four Positioning Functions: Back up and down, Tilt up and down, Foot up and down, and Table up and down.
- \* Ergonomic hand control or foot control.
- \* Manually-adjustable head rest.
- \* Stirrups: Left and right that may be positioned into any number of ways as well as laterally in three positions.
- \* Two 1" (2.5 cm) tall x 3/8" wide accessory rails.
- \* Two duplex, hospital-grade 115 VAC receptacles, one on each side of table.
- \* Paper roll holder hidden in head end of top.



**CeviMed**  
MEDICAL EQUIPMENT & SUPPLIES

## MIDMARK 75L PROCEDURE TABLE

### SPECIFICATIONS

Patient Weight (Maximum) ..... 325 lbs (147.4 kgs) maximum  
 Weight of Table..... 515 lbs (233.6 kgs)  
 Back Section Range..... 0° (horizontal) to +85° ± 5°  
 Table Top Height Range..... 22.5 ± 0.5 in. to 40.0 ± 1.0 in.  
 (55.9 ± 1.3 cm to 101.6 ± 2.5 cm)  
 Foot Section Range ..... 0° (horizontal) to -90° ± 5°  
 Table Top Tilt Range ..... 0° (horizontal) to +30° ± 5° (foot up)  
 Dimensions:  
 Upholstered top..... 27.0 in. wide x 69.25 in. long  
 (68.6 cm x 175.9 cm)  
 With Headrest Extended ..... 82.75 in. (210.2 cm)  
 Electrical Requirements ..... 115 VAC nominal (126 VAC maximum - 110 VAC minimum), 60 HZ, 12 amps maximum  
 Paper Rolls ..... Can accept two standard paper rolls of  
 18 in. x 3.0 in. (45.7 cm x 7.6 cm)  
 Duty Cycle..... 15 seconds on / 5 minutes off (motor run time)  
 Power Cord..... extends 70 in. (Minimum) from table. 12 AWG / 3 conductor, SJT grey jacketed junior hard service with hospital grade grounding type plug.



1 year warranty.



Brand new upholstery.



Save 20% - 50% on Procedure tables by buying refurbished.

