

Infant Flow[®] SiPAP Plus

The next generation of the Infant Flow[®] SiPAP now comes with a standard auxiliary O₂ outlet that was previously an optional feature.

The active alarms can now be silenced for 60 seconds.

On the rear of the device, we have added an RS-232 port. The port provides an output-only interface to provide patient monitored values to a client such as an interface to an Electronic Medical Record system.



Infant Flow SiPAP with roll stand and humidifier bracket (*front view*)



Infant Flow SiPAP with roll stand and humidifier bracket (*rear view*)



CareFusion

Infant Flow SiPAP Plus


Part no.	Description	Qty.	Comments
27416-001	Infant Flow SiPAP Plus	1	
25648-001	Kit including air/oxygen hoses and power cord	1	
28966-001	SiPAP-PC interface kit	1	

Standard and accessories

Part no.	Description	Qty.	Comments
27158-001	Humidifier base bracket <i>(Fisher & Paykel)</i>	1	These ship as standard accessories with purchase of configured stand <i>(part no. 27353-001)</i>
27194-001	Storage basket	1	
27357-001	Top pole <i>(for water bag)</i>	1	
27353-001	SiPAP stand, configured	1	

Optional/Replacement accessories

27159-001	Humidifier adaptor <i>(Teleflex® Neptune®)</i>	1	
27366-001	Dual-basket and humidifier attachment mount bar	1	<p>Dimensions allow:</p> <ul style="list-style-type: none"> Fisher & Paykel or Teleflex Neptune humidifier to be mounted side by side with the basket Clearance for mounting optional dual-tank holder to rear of stand <i>(part no. 27159-001 is required to mount a Teleflex/Neptune humidifier)</i>
27195-001	Cord wrap	1	
27356-001	Oxygen tank mount	1	Dual-tank holder for D or E cylinders
28290-001	75 mm caster, non-locking	1	Individual twin-wheel replacement caster
28291-001	75 mm caster, locking	1	Individual twin-wheel replacement caster



⚠ WARNING—U.S. Federal Law restricts this device to sale by or on the order of a physician.

CareFusion 800.231.2466 toll-free
 22745 Savi Ranch Parkway 714.283.2228 tel
 Yorba Linda, CA 92887 714.283.8493 fax

carefusion.com

