



Unit comes patient ready and includes:

FEATURES:

- Shock advisory system guides user with voice and visual prompts through each step.
- easy manual mode, including waveform displays.
- Noninvasive pacing, Masimo SET pulse oximetry, ECG wirigin (3- or 5-wire) and synchronized cardioversion.
- Designed for indoor use, compact and lightweight makes it easy to rush to the scene
- Optional docking station enables device to be firmly attached to crash cart for safe transport
- Superb color graphics, easily viewed from many angles, enable clinicians to clearly see color matched waveforms and values for efficient and informed patient care decisions.
- Battery indicator on device display

The LIFEPAK 20e defibrillator/monitor is highly intuitive to use and adapts to various patient environments. It skillfully combines an AED function for the infrequent, BLS-trained responder, with manual capability so that ALS-trained clinicians can quickly and easily deliver advanced diagnostic and therapeutic care. Lifepak 20e defibrillator delivers higher energy that is needed to defibrillate patients who require more escalating biphasic energy to restore normal heart rate and rhythm and induce fibrillation.

- A/C Power cord
- SpO2 Sensor w/ adapter
- Patient cable w/ leads
- Battery
- Handsfree cable