



Midmark 304 EXAM TABLE

FEATURES

Manually adjustable back section controlled by spring assisted release handle, located on left and right side of table.

Extendable footrest shelf include upholstered pad and removable treatment pan.

Right and left stirrups.

Five storage drawers accommodate supplies.

Four leveling feet (one on each corner) to ensure solid, level installation.

Paper roll holder.

Pelvic Lift Bar (optional): Allows seat to be positioned in Lithotomy position.

50 Watt Heater (optional): Warms contents of top, foot end drawer.



CeviMed
MEDICAL EQUIPMENT & SUPPLIES

Midmark 304 EXAM TABLE

SPECIFICATIONS

Maximum Patient Weight: 300 lbs (136 kg).

Weight of Table: 250 lbs (113.4 kg).

Back Section Adjustment: 0 - 70° ($\pm 5^\circ$).

Stirrup Extension (Maximum): 16" (40.6 cm).

Stirrup Lateral Positions: 4.

Footrest Extension (Maximum): 10" (25.4 cm).

Table Top Dimensions:

Upholstered Top: 27" (W) x 56" (L) (68.6 cm x 142.2 cm).

With Footrest Extended: 75.5" (191.8 cm).

Seat Height: 31.75" (80.6 cm) to top of upholsters seat at foot end.

Step: 10" x 17.75" x 10.25" (H) (25.4 cm x 45.1 cm x 26 cm).

Base: 17.5" x 42.375" x 4.125" (H) (44.4 cm x 107.6 cm x 10.5 cm).

Paper Roll: Can accept up to 21" x 3.5" Diameter (53.3 cm x 8.9 cm).

Electrical Requirements: 220 - 240 $\pm 6\%$ VAC, 50/60 HZ, 5 amps.

Replacement Fuse Spec: 250 VAC, 6.3 amp, IEC 127 rated, 5 x 20 mm, Time Lag Fuse.



1 year warranty.



Brand new upholstery.



Save 20% - 50% on Exam tables by buying refurbished.

