



MIDMARK 641 PROCEDURE TABLE

FEATURES

Narrow patient support surface, slim back section, and tuck away side rails provide maximum patient accessibility and supports proper posture.

Powered seat, back, foot, and tilt offer multi-positioning capabilities to improve comfort and reduce likelihood of injury.

Premium Comfort System enhances patient comfort for a chair-like feel.

Optional wireless controls allow for placement of hand and foot controls anywhere in the room, removing chords and eliminating obstacles.

Multiple options and accessories available that improve the functionality of any office.



MIDMARK 641 PROCEDURE TABLE

SPECIFICATIONS

Patient Weight (maximum)..... 450 lbs (204.1 kg)
Weight of Chair641 (w/standard base)
w/packaging & skid641 (w/rotational base)
w/packaging & skid440 lbs (200 kg)
495 lbs (224 kg)
525 lbs (238 kg)
585 lbs (265 kg)
Power Cord 8 ft. (243.8 cm) long

Electrical Requirements: [See Model I.D. / Compliance Charts]
Foot / Hand Control Voltage 10 VAC, SELV (Safety Extra Low Voltage)
Duplex Receptacle (optional)
maximum load 115 VAC, 3A, 50/60 Hz

Fuses:
Located at power cord inlet:
F1 & F2 6.3 A, 250V, Type T, 5 x 20 mm
Duty Cycle

(motor run time): Intermittent Operation
[30 seconds ON / 5 minutes OFF]
Classifications: Class 1, Type B Applied Part
Protection against ingress of fluids:
Foot Control only:
Ordinary Equipment
IPX1
Regulatory Compliance: [See Model I.D. / Compliance Charts]



1 year warranty.



Brand new upholstery.



Save 20% - 50% on
Procedure tables by
buying refurbished.

