



MIDMARK 625 Bariatric Treatment HIGH-LOW TABLE

FEATURES

Powered movement of the seat and back gives you the flexibility of positioning patients where you need them.

Standard hand and foot controls designed to make patient positioning quick, effective, and always within reach.

Premium Comfort System stretches elastomeric material over a seat cutout designed to enhance patient comfort for a chair-like feel on an exam table.

Easily remove upholstery without the use of tools, allowing for a color change in seconds or additional access for a thorough cleaning.



CeviMed
MEDICAL EQUIPMENT & SUPPLIES

MIDMARK 625 Bariatric Treatment HIGH-LOW TABLE

SPECIFICATIONS

Patient Weight (maximum): 650 lbs
(295 kg)
Weight of Table:
w/upholstery
w/packaging & skid (no uph.)
Uph. w/packaging (shipped separately)
470 lbs (213 kg)
487 lbs (221 kg)
45 lbs (20 kg)
Power Cord Length: 8 ft. (244 cm)
Electrical Requirements: (See Regulatory Compliance Chart)
Foot / Hand Control Voltage: 14 VAC,
SELV (Safety Extra Low Voltage)
Duplex Receptacle(s) (maximum
load) 115 VAC, 5A, 50/60 Hz
Fuses:
IEC inlet fuses (2)
F1 & F2 (on main PC board):
F3 (on main PC board):
F4 (on main PC board)
Models w/ drawer heaters



1 year warranty.



Brand new upholstery.



Save 20% - 50% on
High Low Tables by
buying refurbished.

