



## Midmark 646 Barrier PODIATRY CHAIR

# FEATURES

The Midmark 646 Barrier-Free Power Podiatry Procedures Chair Base includes power cord and chair arms.

The 646 Barrier-Free Power Podiatry Procedure Chair is a step forward in foot and ankle care with an exclusive, ergonomic foot section design that puts control right where it's needed - the foot section.

**Control:** Integrated touch controls on both sides of the foot section are not only convenient but also help reduce costly breakage of hand and foot controls.

**Positioning:** Powered back and tilt to simplify positioning. Easy-to-reach release handles on both sides of the foot section allow for smooth, one-handed operation.



**CeviMed**  
MEDICAL EQUIPMENT & SUPPLIES

## Midmark 646 Barrier PODIATRY CHAIR

### SPECIFICATIONS

Patient Weight (max): 450 lbs  
(204 kg)

Weight of Chair: 290 lbs (131 kg)

Power Cord Length: 8 ft (244 cm)

Foot Control Voltage: 10 VAC,  
SELV (Safety Extra Low Voltage)

Chair Receptacle

Maximum Load:

115 VAC, 3 amps, 50 / 60 Hz

Duty Cycle

(Motor Run Time):

Intermittent Operation

(30 seconds ON - 5 minutes OFF)

Protection against  
ingress of fluids:

Ordinary Equipment

Foot control only: IPX1

Classifications: Class 1, Type B,  
Applied Part



1 year warranty.



Brand new upholstery.



Save 20% - 50% on  
Podiatry chairs by  
buying refurbished.

