



CLASSIC MF MOP BLUE KENTUCKY 250 G

Mop system, for use with removable wet mop holder. Suitable for vertical or jaw presses. Ideal for cleaning residual contamination where frequent wringing of the mop is required.



Order no. 9.212-066.0

Technical data

Textile usage		Reusable textiles
Textile attachment		Mop clamp
Material		80% PET / 20% PA
Textile material		microfibre
Washing temperature	°C	max. 90
Washing recommendation	°C	60
Quantity	Piece(s)	1
Weight per product / Weight per m ²	kg	0.3
Package weight	kg	0.3
Dimensions (L × W)	mm	300 × 100
Dimensions, packaged	mm	300 × 100 × 100

Application

Bucket method

For optimal moisture and a single-step wiping process, dip the mop fully into the cleaning water, allow it to absorb as much water as possible and then wring it out in the press.

Properties

- Very easy to use, meaning no training is required

Washing instructions

- Do not wash with other textiles
- Do not use bleach, softener or strongly alkaline detergents
- Washing or drying at high temperatures can damage the fibres
- Thoroughly rinse the textile before washing it
- Please do not dry in the dryer